

ESHRE PGD Consortium Best Practice Guidelines for Amplification-Based Preimplantation Genetic Diagnosis (PGD)

Harton, GL^{1,8}, De Rycke, M², Fiorentino, F³, Moutou, C⁴, SenGupta, S⁵, Traeger-Synodinos, J⁶, and Harper, JC^{5,7}

- 1 Reprogenetics LLC, Livingston, NJ 07039, USA. Formerly Genetics & IVF Institute, Preimplantation Genetic Diagnosis Laboratory, Fairfax, VA 22031, USA
- 2 Centre for Medical Genetics, Universitair Ziekenhuis Brussel, Brussels, Belgium
- 3 Genoma Laboratories, Rome, Italy
- 4 Université de Strasbourg, Hôpitaux Universitaires de Strasbourg, Laboratoire de Biologie de la Reproduction. F-67000 Strasbourg
- 5 UCL Centre for PG & D, Institute for Women's Health, University College London, UK
- 6 Medical Genetics, University of Athens, "Aghia Sophia" Children's Hospital, Athens, Greece
- 7 Centre for Reproductive and Genetic Health, UCLH, London, UK
- 8 To whom correspondence should be addressed

Abstract:

In 2005, the European Society for Human Reproduction and Embryology (ESHRE) PGD Consortium published a set of Guidelines for Best Practice PGD to give information, support, and guidance to potential, existing and fledgling PGD programs (Thornhill et al., 2005). The subsequent years have seen the introduction of a number of new technologies as well as the evolution of current techniques. Additionally, in the light of ESHRE's recent advice (http://www.eshre.com/binarydata.aspx?type=doc/Manual_for_guideline_development_JK_WN_SB_FINAL.pdf) on how practice guidelines should be written and formulated, the Consortium believed it was timely to revise and update the PGD guidelines. Rather than one document that covers all of PGD as in the original publication, these guidelines are separated into four new documents that apply to different aspects of a PGD program i.e. Organization of a PGD centre, FISH-based testing, Amplification-based testing and Polar Body and Embryo Biopsy for Preimplantation Genetic Diagnosis/Screening (PGD/PGS). Here we have updated the sections that pertain to amplification-based PGD.

Topics covered in this guideline include inclusion/exclusion criteria for amplification-based PGD testing, preclinical validation of tests, amplification-based testing methods, tubing of cells for analysis, set-up of local IVF center and transport PGD centers, quality control/assurance and diagnostic confirmation of untransferred embryos.

Introduction

The rapidly changing nature of PGD, specifically the technologies associated with its use and increasing patient access, has necessitated review and revision of the original ESHRE PGD Consortium guidelines (Thornhill et al., 2005). As a result, the ESHRE PGD Consortium has prepared four guidelines (Harton et al., 2010a, b, c, submitted): one relating to the organisation of the PGD Centre and three relating to the methods used: DNA amplification, FISH and biopsy/embryology. The method guidelines should be read in conjunction with the organisation of the PGD Centre guidelines which contains information on personnel, inclusion/exclusion criteria, genetic counselling and informed consent, setting up an IVF centre, transport PGD, QA/QC and accreditation (which is also further discussed in the paper by Harper et al, 2010). In this document, the laboratory performing the diagnosis will be referred to as the PGD/PGS Centre and the centre performing the IVF as the IVF Centre. Topics covered in this guideline include general uses of amplification-based PGD testing, laboratory issues relating to amplification-based testing, pre-examination validation, examination process, and post examination process.

Preimplantation genetic diagnosis (PGD) has been developed for patients at high risk of transmitting a genetic abnormality to their children which includes all monogenic defects (autosomal recessive, autosomal dominant and X-linked disorders). More recently, DNA amplification-based PGD applications have broadened and include sibling-donor selection through HLA-matching (Van de Velde et al., 2009), and the analysis of familial chromosomal rearrangements (Fiorentino et al., in press). PGD is still relatively unregulated and lacks standardization compared to other forms of diagnostic testing, however, more federal, state and local governments are beginning to take an interest in PGD and some have begun accrediting laboratories that offer PGD (Harper et al., 2010). This is a logical step considering the comparative difficulty in achieving the highest levels of accuracy and reliability when analysing single cells as part of PGD versus more routine genetic testing. Many regulations, laws and voluntary networks exist in the mainstream diagnostic community to maintain the highest quality in diagnostic testing. For example, the European Quality Molecular Network has attempted to improve and standardize molecular diagnostic testing across Europe (Dequecker et al., 2001). One step towards higher quality overall and standardisation for PGD is to build consensus

opinion on best practices within the PGD community; a component of the mission of the ESHRE PGD Consortium (hereafter referred to as the Consortium - ESHRE PGD Consortium Steering Committee 1999; 2000; 2002).

The Consortium recognizes that owing to variations in local or national regulations and specific laboratory practices, there will remain differences in the ways in which PGD is practiced (from initial referral through IVF treatment, single cell analysis to follow-up of pregnancies, births and children). However, this does not preclude a series of consensus opinions on best practice based upon experience and available evidence. Indeed, the American Society for Reproductive Medicine published a practice committee report for PGD in 2008 (American Society of Reproductive Medicine and Society for Assisted Reproductive Technology. Practice Committee Report, 2008) essentially reviewing PGD practice in the United States. The PGD International Society (PGDIS) has also drafted guidelines which were recently updated and are more in-depth than the ASRM report. These guidelines are concise and remain so in their recent revised edition (Preimplantation Genetic Diagnosis International Society, 2004, 2007). The consensus opinions provided in this document and the accompanying guidelines, not only reflect current use of PGD but also offer consensus-based specific guidance regarding how best to practice clinical PGD based upon clinical experience, and data, both published and unpublished.

The Consortium hopes that a minimum standard might be achieved across all centres actively providing clinical PGD. Achieving this goal could ultimately have the net effect that patients receive the best care possible regardless of the centre at which they are treated. Rather than a drift towards the lowest common denominator, established and fledgling centres alike can learn from global experiences and be guided by consensus opinion.

These opinions are not intended as rules or fixed protocols that must be followed, nor are they legally binding. The unique needs of individual patients may justify deviation from these opinions, and they must be applied according to individual patient needs using professional judgement. However, guidelines and opinions may be incorporated into laws and regulations and practitioners should check the status of clinical practice guidelines in their own countries to determine the status of this document.

1. General uses of DNA amplification-based tests

- 1.1.** Amplification-based tests can be used for the diagnosis of monogenic defects at the DNA level (Sermon, 2002; Thornhill and Snow, 2002). This includes specific diagnosis for X-linked disease, as well as diagnosis of autosomal recessive and dominant diseases.
- 1.2.** Due to the risk of contamination and allele drop-out (ADO), it is **recommended** that DNA amplification protocols include the use of linked or unlinked markers in addition to the disease locus (Sermon, 2002; Thornhill and Snow, 2002).
- 1.3.** For X-linked diseases, analysis of the mutation and linked markers allows for the transfer of unaffected males as well as the exclusion of carrier females if the patient is so inclined.
- 1.4.** When sexing only is performed for X-linked diseases by DNA amplification, it is recommended that several loci are included to monitor contamination and preclude misdiagnosis due to ADO (Renwick et al., 2006; Renwick and Ogilve, 2007).

2. Laboratory issues relating to single-cell DNA amplification

2.1. Laboratory Infrastructure, Equipment and Materials

The following recommendations are made:

- 2.1.1.** There should be physical separation of pre-amplification (preferentially a positive pressure room), amplification and post-amplification (preferentially a negative pressure room) laboratories and the biopsy laboratory.
- 2.1.2.** When positive and negative pressure rooms are present, they are enclosed by a lock chamber.
- 2.1.3.** It is recommended that primary amplification reactions are set-up in a laminar flow hood in the pre-amplification area.
- 2.1.4.** When performing two rounds of PCR, the area for setting up the first round should still be separated from the area for the second round of PCR.
- 2.1.5.** Each area has its dedicated materials and equipment.
 - 2.1.5.1.** The pre-amplification reagents and materials should be kept from any DNA source.
 - 2.1.5.2.** An appropriate unidirectional workflow should be in place, avoiding any backfire of amplified products to the pre-amplification area.
- 2.1.6.** All clinical equipment should meet the criteria set for the intended application, be appropriately maintained and serviced, with all aspects supported by written standard operating procedures (SOP).

- 2.1.7.** Protective clothing for single cell DNA amplification work should be worn, including full surgical gown (clean, not sterile, and changed after each case), hair cover/hat and face mask (covering nose and mouth) and shoe covers. Gloves should be worn at all times and changed frequently. These should be well-fitting (e.g. latex or nitrile, but not vinyl examination gloves).
- 2.1.8.** All plastic-ware used, including filter tips, should be certified DNA, DNAase free
- 2.1.9.** Work surfaces, equipment, etc. should be cleaned with DNA decontamination solutions or 20 % bleach prior to each case. All tube racks should be either one-use or, if re-used for PGD, be autoclaved/cleaned in 20 % bleach after each use.
- 2.1.10.** All batches of reagents should be recorded so that they may be traceable to specific assays
- 2.1.11.** Whenever possible, all solutions should be purchased 'ready to use' and should be of 'molecular biology' grade or equivalent
- 2.1.12.** All reagent-solutions should be split into single-use aliquots and no aliquot should be re-used for a clinical PGD case
- 2.1.13.** Solutions made in-house may be autoclaved (using a dedicated autoclave) or if this is not possible, passed through 0.2 µM filters under stringent conditions.
- 2.1.14.** A DNA amplification master mix may be pre-treated either by UV irradiation or restriction enzyme digestion for decontamination purposes.

2.2. Work practice controls

The following **recommendations** are made:

- 2.2.1.** It is essential that an adequate labelling system is used to match the cell diagnostic result with the embryo from which that cell was biopsied. Comprehensive and robust labelling is required; printed ticker labelling may be superior to pens as labelling should be legible and indelible.
- 2.2.2.** Labelling and sample identification should be confirmed for critical and high risk steps. It is **recommended** that the unique patient identifier and embryo/cell number should be witnessed and signed by two members of the PGD team at the following steps:
- a) at cell tubing to confirm that the cell identification matches the labelling on the relevant tube
 - b) when diagnostic results are recorded to ensure accuracy and correlation with the correct cell and/or embryo identification

2.2.3. All personnel undertaking single-cell tubing and testing should be adequately trained and shown competent for single-cell procedures, prior to working on clinical specimens. This training should be documented and competencies should be retested regularly.

2.3. Tubing cells

2.3.1. Tubing should be carried out under stringent precautions to minimize contamination. This implies that materials and reagents should be prepared in a dedicated pre-amplification laminar flow (UV light to be turned on before (and after) each cycle and regularly cleaning of working areas with DNA degrading detergents is required).

2.3.2. Specific measures to avoid contamination are:

2.3.2.1. Granulosa cells are meticulously removed from the oocytes to avoid maternal contamination and fertilization is achieved through ICSI to avoid contamination with sperm. The embryo biopsy procedure should preclude operator contamination.

2.3.2.2. Notably, ill-fitting gloves can be a source of carry-over contamination.

2.3.2.3. For each cell assessed (either in the validation phase, or during clinical PGD cases), a blank containing medium collected from the last wash droplet and all amplification components except DNA should be run. However, a contaminated blank is only an indication of a more general problem and does not mean that the corresponding sample is contaminated.

2.3.3. Single cells should be washed at least twice using a sterile transfer pipette before transfer into amplification tubes. The amount of medium co-transferred with the cell in the tube should be minimal.

2.3.4. It is acceptable to transfer embryonic cells to amplification tubes with or without microscopic visualisation.

2.3.5. If the embryonic cells are lysed during washing or transfer, the pipette is possibly contaminated and has to be discarded. For cleavage stage biopsy, another blastomere should be sampled.

2.3.6. If accurate and reliable amplification methods are in place, it is recommended to biopsy only one cell from cleavage-stage embryos as it has been shown that the biopsy of two blastomeres from 8-cell embryos results in a 40% lower live birth rate as compared to 8-cell embryos having lost only one blastomere (Goossens et al., 2008; De Vos et al., 2009)

2.3.7. Efficient cell lysis and adequate DNA denaturation are required. Both alkaline lysis and proteinase K/SDS are accepted as the best cell lysis methods (Thornhill et al., 2005).

3. Single-cell amplification characteristics: specificity, efficiency, ADO and contamination

Many variations in DNA amplification methods have been published. Both PCR-based and non-PCR-based amplification methods can be used, aiming at the amplification of specific region(s) or the whole human genome. Each method has its merits and its limitations. The methods used should have been previously implemented, tested and validated in the PGD Centre.

3.1. PCR-based amplification

Fluorescent PCR amplification followed by post-PCR methods for allele discrimination is regarded as the gold standard. Analysis of fluorescent PCR fragments on an automated sequencer is much more sensitive and yields a more accurate fragment size determination than analysis of conventional PCR fragments on ethidium bromide stained gels (Hattori et al., 1992; Lissens and Sermon, 1997)

- 3.1.1.** Fluorescent PCR is also an efficient way to significantly reduce ADO (Sermon et al., 1998).
- 3.1.2.** Other PCR design factors leading to better specificity, higher PCR efficiency and low ADO rates are smaller amplicon size (< 350 bp) and appropriate primer design using primer software tools together with BLAST (Basic Local Alignment Search Tool, <http://blast.ncbi.nlm.nih.gov/Blast.cgi>) searches and single nucleotide polymorphism (SNP) detection to ensure specificity (Piyamongkol et al., 2003).
- 3.1.3.** DNA degradation and the choice of DNA polymerases also have an impact on efficiency and ADO of single-cell PCR (Piyamongkol et al., 2003).
- 3.1.4.** The co-amplification of polymorphic marker(s), linked or unlinked, with the mutation of interest in a multiplex PCR is recommended as it allows a more accurate diagnosis and simultaneously reveals the presence of ADO and contamination (Pickering et al., 1994).
- 3.1.5.** When no linked markers are available, or the couple is not informative for available markers, or the set-up of a reliable multiplex PCR proves to be too difficult, the biopsy and testing of two cells is an acceptable alternative. Their subsequent independent analysis will help in identifying ADO, which will be seen as a discrepancy in genotype between the cells.
- 3.1.6.** Multiplex PCR in one round likely reduces the chances for contamination and tube transfer errors compared to (hemi)nested PCR protocols. (Hemi)nested PCR protocols are acceptable as long as they are reliable and accurate (Stern et al., 2002).

3.2. Other amplification-based techniques

Different methods have been developed to amplify the whole genome (WGA), including primer extension preamplification (PEP), degenerate oligonucleotide primed PCR (DOP-PCR), and multiple displacement amplification (MDA). The MDA method is non-PCR-based and seems a better WGA method with very good amplification rates and a lower amplification bias than PEP and DOP-PCR at the present time (Lovmar et al., 2003; Panelli et al., 2006). So far, the major problem has been the relatively high preferential amplification and ADO rate (on average 25% which is about 3 – 5 higher than with PCR-based protocols). Analysis of multiple loci near the gene of interest can overcome this problem and allow reliable reconstitution of the parental haplotypes in the embryos, ensuring a diagnostic efficiency similar to PCR-based protocols. Statistically, the probability that two alleles present in a heterozygous cell are amplified at a single marker is $(1 - a)^2$ for ADO rate a . For N such markers, the probability that at least one marker amplifies both alleles is $1 - (1 - (1 - a)^2)^N$ (Renwick et al., 2007). An average ADO of 25% implies that four fully informative markers are required for the detection of both alleles with 95% confidence. Conversely, a lower ADO rate of 5% in case of PCR-based amplification yields a power of 99% for two fully informative markers. Further improvement of ADO rates, amplification yield and genome coverage of the MDA method can be obtained from new primer designs combined with elevated amplification temperature (Alsmadi et al., 2009).

The recently described approach of Preimplantation Genetic Haplotyping (PGH) offers a more generic linkage-based approach to preimplantation diagnosis, and is especially useful for diseases with a wide spectrum of mutations, such as cystic fibrosis (CF) and Duchenne muscular dystrophy (DMD) (Renwick et al., 2006). This generic approach is also useful for HLA haplotyping, which is applied to select embryos that are suitable donors for an affected sibling requiring stem cell transplantation, either alone or in cases of acquired diseases (e.g. leukaemia) or for HLA typing combined with genetic diseases. Finally, a recently described use of DNA amplification-based technology involves the detection of chromosomal imbalances in embryos derived from chromosome rearrangement carriers (Fiorentino et al., in press), which may be done by nested PCR or using whole genome amplification methods as described above.

3.2. Mutation detection strategies

3.2.1. Many mutation detection strategies exist. The following recommendations are made:

3.2.1.1. All methods require validation on suitable genomic DNA controls.

- 3.2.1.2.** Strategies should allow the ability to distinguish between mutant alleles, wild type alleles and amplification failure.
- 3.2.1.2.1.** For instance it is not recommended to only amplify the normal allele in case of a large deletion/insertion as the mutant allele can not be distinguished from amplification failure.
- 3.2.1.3.** It is recommended in case of triplet repeat expansions that strategies based on amplification of the normal alleles are only applied for fully informative couples.
- 3.2.1.4.** For small deletions or insertions, direct fragment size determination can be performed.
- 3.2.1.5.** Detection of single nucleotide substitution via restriction enzyme digestion should include a strategy to control the digestion reaction.
- 3.2.1.6.** Indirect mutation detection is applied for exclusion testing, for HLA typing, in case of unknown mutation or large deletions/insertions with unknown breakpoints, as general strategy for genes with high numbers of (private) mutations or in case direct mutation testing is not successful (presence of pseudogene(s), GC rich sequences refractory to single-cell amplification).
- 3.2.1.7.** Indirect protocols are not applicable in case of *de novo* mutations for couples without previous pregnancies.
- 3.2.2.** It is recommended to include at least two flanking markers in indirect mutation PCR-multiplex protocols.
- 3.2.2.1.** Analysis of at least two loci closely linked to the gene underlying the disease will reduce the risk of unacceptable misdiagnosis (i.e. transfer of an affected embryo) due to ADO (presumed to be around 5%) to a minimum (< 1%).
- 3.2.2.2.** More than two markers will make the test more robust: an assay with just one marker at each side of the mutation will yield “no diagnosis” when one marker fails to amplify. Therefore, two upstream and two downstream markers are preferably applied.
- 3.2.3.** Ideally the polymorphic marker should be intragenic with a high degree of heterozygosity and produce a clearly interpretable peak pattern.
- 3.2.3.1.** The risk of recombination should be considered for every marker and is especially important in case of large genes and genes with recombination hot spots.
- 3.2.3.2.** It is recommended to choose extragenic markers within 1 Mb of the mutation of interest to reduce the risk of recombination events. Even though loci 1 cM apart are

expected to show only approximately 1% recombination, crossover events between markers and a specific mutation cannot be completely ruled out even for very closely linked markers (Altarescu et al., 2008).

3.2.3.3. The recombination risk is important for the risk assessment of accuracy and reliability of the results.

3.2.4. Specifically for indirect HLA typing, a sufficient number of microsatellites should be chosen to be able to detect recombination within the HLA region:

3.2.4.1. A minimum one marker located upstream of HLA-A, minimum one marker between HLA-A and HLA-B, minimum one marker between HLA-B and HLA-DRA, one marker between HLA-DRA and HLA-DQB1 and minimum one marker downstream of HLA-DQB1.

3.2.4.2. The use of two markers will make the test more robust. Around 2% of the embryos with a conclusive diagnosis show recombination within the HLA locus (Rechitsky et al., 2004; Verlinsky et al., 2004; Kuliev et al., 2005; Van de Velde et al., 2009). This is higher than the 0.3–1% described within the locus (Martin et al., 1995; Malfroy et al., 1997).

3.2.5. When linked markers are used, the same minimum haplotyping standards exist for PGD testing as for regular prenatal diagnosis.

3.2.5.1. A molecular biologist experienced in pedigree and linkage analysis should determine which samples are needed as well as the number of samples required for reliable and accurate diagnosis.

3.2.5.2. More than one meiosis is preferentially used in building accurate haplotypes.

3.2.5.3. When no additional family members are available for analysis in cases of paternally linked disorders, single sperm should be used for haplotype construction.

4. Pre-examination Validation

4.1. Confirmatory testing of the mutation on DNA or segregation results (using standard tests) is **recommended**. It is **acceptable** to outsource confirmatory testing to outside laboratories if the PGD Centre is not able to offer appropriate testing of the disease causing mutation.

4.2. Confirmatory testing of the clinical PGD assay on parental and proband blood DNA is **recommended** since the PGD assay could be non-informative due to PCR failure as a result of polymorphisms or DNA sequence variations present in the population (Thornhill et al., 2005).

4.3. Before moving on to single cell validation, it is **recommended** that testing of various genotypes using DNA samples be performed, including affected (autosomal dominant), carrier

(autosomal recessive, X-linked diseases), and unaffected samples for the mutation to be tested or DNA samples with heterozygous markers for tests using a linkage-only approach.

4.4. Single cell analysis should include cells from affected, carrier (autosomal recessive, X-linked diseases) and unaffected individuals (Sermon, 2002, Thornhill and Snow, 2002).

4.4.1. The use of diluted DNA samples to mimic single cells (~6 pg of DNA is equivalent to one cell) is **not recommended** as stochastic variation in pipetting renders the results unreliable.

4.4.2. Acceptable cell types for single cell analysis include lymphocytes, lymphoblastoid cell lines, fibroblasts or buccal cells. Other somatic cell types should be acceptable assuming the laboratory has validated their use internally.

4.4.3. Testing of embryonic cells would be beneficial since they are the target cell type for clinical testing, however their use will depend on availability.

4.4.4. The type of cell used for validation may influence the amplification efficiency and allele drop-out rate (Glentis et al., 2009).

4.4.5. It is **recommended** that at least 50 single cells are used for validation of any amplification-based test and these cells should be run in multiple experiments of 10-20 cells/experiment to assess amplification efficiency, allele drop-out rate and run to run variability.

4.4.6. The same number of cell wash blanks should be run along with the single cells to assess contamination.

4.4.7. If embryonic cells are analyzed, a minimum of 10 samples is **recommended**, and, where possible, several embryonic cells from the same embryo should be tested to determine the expected genotype for the mutation and linked markers.

4.5. Assessing DNA amplification efficiency and ADO rates

4.5.1. Assessing amplification efficiency is an important part of pre-clinical validation. An amplification efficiency of at least 90 % for each marker is **recommended**.

4.5.2. Assessing ADO rates is an important part of pre-clinical validation. ADO rates should be determined using cells that either carry the mutation of interest or, for linkage only tests, wild-type cells that are heterozygous for the markers being used.

4.5.3. ADO rates should be as low as possible, preferably < 10 %. It should be noted that a higher ADO rate can be tolerated when dealing with autosomal recessive diseases as compared to autosomal dominant or compound heterozygous diseases. Higher ADO rates

can also be tolerated when dealing with whole genome amplification-based protocols, however, in this case an increased number of linked markers and/or a 2 cell biopsy strategy can be employed in the protocol.

4.6. Reporting pre-examination validation

A report has to be written detailing the protocol and validation steps of the PGD protocol workup (Harper et al., 2010)

5. Examination Process

5.1. Existence of and adherence to clinical testing protocol

The following **recommendations** are made (Thornhill et al., 2005; Harper et al., 2010)

- 5.1.1.** Clinical testing protocols should include explicit instructions including a summary of results from the validation steps of assay development, scoring criteria, reporting procedures as well as a framework for counseling patients in the presence of diagnostic results.
- 5.1.2.** SOPs are required for all protocols, all equipment, and all processes that take place in the PGD Centre and should include selection and validation of examination procedures, clinical relevance, purpose of exam, specimen requirements and means of identification, equipment and special supplies, reagents, standards or calibrations and internal control materials, instructions for performance of the examination, limitations of the examinations, recording and calculation of results, internal quality control procedures and criteria against which examination processes are judged, reporting reference limits, responsibilities of personnel in authorizing, reporting and monitoring reports, hazards and safety precautions assuring the quality of examinations.
- 5.1.3.** Risk assessments are required for every stage of the PGD process. These assessments should be integrated to the standard operating procedures. Laboratory staff should understand the SOPs clearly as these are the fundamental back bone to the service.
- 5.1.4.** Deviations from protocol should be recorded. If frequent deviations occur, there should be a mechanism in place to change procedures accordingly.
- 5.1.5.** Well structured laboratory forms to report workup, PGD cycle and analysis of untransferred embryo results should also be used.

5.1.6. In all of the critical stages of the PGD process, witnessing and signing is **recommended**.

5.1.7. All steps of all protocols should exist in the laboratory/department/center, including explicit instructions for their execution, evaluation (scoring criteria), reporting procedures and other administrative issues.

5.2. Use of Intra-assay controls

It is **recommended** that positive control genomic DNA is included within each assay. It is **acceptable** to use diluted DNA samples as a positive control but these should never replace the single cell positive controls discussed below.

5.2.1. For dominant diseases, this would include DNA from an affected individual as well as an unaffected individual

5.2.2. For recessive diseases, this would include one heterozygous carrier and one homozygous affected individual

5.2.3. For all diseases, this would include the prospective parents and other relatives to assess the linked markers

5.2.4. Single cells from an appropriate control individual or individuals to be certain that the assay is working at the single cell level on the day of the test.

5.2.5. For each embryonic cell analyzed, a cell-free wash control negative should be analyzed.

5.2.6. At least one reagent negative control (no DNA, no wash) should be run to assess contamination from within the reagents.

5.3. Scoring clinical DNA amplification results

The following **recommendations** are made:

5.3.1. Scoring criteria should be established in written protocol and adhered to for the interpretation of results

5.3.2. Results should be analyzed by two independent observers and discrepancies adjudicated by a third observer (where possible). If no resolution is reached the embryo should not be recommended for transfer, i.e. should be given the diagnosis of uninterpretable or inconclusive.

5.3.3. All data from clinical testing should be kept for Quality Control (QC) purposes and records.

5.3.4. Results should be reviewed and signed or electronically validated by a suitably qualified person

5.3.5. A written or online electronic report should be prepared and double-checked. The checked report is given to the IVF Centre to ensure transfer of the correct embryos.

5.3.6. Reporting of clinical results to the IVF Centre must follow local guidelines or law, or, if nothing local exists, the guidelines in ISO15189 (Harper et al., 2010).

6. Post-examination Process

PGD cycle follow-up

6.1. The following **recommendations** are made (ESHRE PGD Consortium Steering Committee 1999; 2000; 2002, PGDIS, 2008).

6.1.1. Confirmation of the diagnosis should be performed on as many untransferred embryos as possible following diagnosis to provide Quality Assurance (QA) as well as accurate and up to date misdiagnosis rates to prospective PGD patients. As follow-up procedures are costly and time-consuming, alternative methods of QA exist which can be combined with limited follow-up of embryos, including: collection of tissue at prenatal testing, collection of cord blood following birth, and surveys of PGD patients with successful pregnancy.

6.1.2. PGD and IVF Centres should make special efforts to follow-up prenatal testing or birth, especially if confirmatory testing is not possible.

6.1.3. Follow-up of pregnancies (including multiple pregnancy rate and outcome), deliveries, and the health of children at birth and beyond should be attempted and maintained along with the cycle data. These data should be used both for internal QC/QA purposes and sent to the ESHRE PGD Consortium during yearly data collections.

6.1.4. Recently a pediatric follow-up working group has been formed by the ESHRE PGD Consortium for the follow-up of children born after PGD. PGD centers are encouraged to take part in this project.

6.2. Baseline IVF pregnancy rates for PGD

6.2.1. Setting appropriate baseline pregnancy rates should be left up to the individual Centres. However, it is **recommended** that each IVF Centre should compare PGD pregnancy rates and matched non-PGD (routine IVF or ICSI) pregnancy rates within that IVF centre (Thornhill et al. 2005).

6.2.2. Comparison of pregnancy rates with those reported by the annual data collections of the ESHRE PGD Consortium can also be done to set benchmarks for continual improvement of the PGD centre (Goossens et al. 2009)..

Appropriate indications for specific tests

6.2.3. It is **recommended** that specific indications for PGD/PGS should remain within the purview of individual clinics (Thornhill et al., 2005; Harton et al., 2010a)

6.3. Misdiagnosis rate

6.3.1. The ESHRE PGD Consortium makes the following **recommendations**: (Thornhill et al., 2005; PGDIS, 2008; Thornhill and Repping, 2008; Wilton et al., 2009).

6.3.1.1. It is **recommended** that misdiagnosis rates should be calculated for each type of assay and for all assays from a particular Centre (Lewis et al., 2001). Such rates include those clinical cases in which affected pregnancies arose and post-transfer confirmation of diagnosis assays that were discordant with the biopsy result.

6.3.1.2. It is **recommended** that confirmatory testing should be performed at least periodically as a quality assurance.

6.3.1.3. It is **recommended** that the published and in-house estimates of misdiagnosis rates should be available to prospective patients along with pregnancy rates on request to allow informed consent for PGD.

6.3.1.4. Following a misdiagnosis, the laboratory should investigate the possible causes of the misdiagnosis and make changes to protocols to eliminate the risk in the future (Wilton et al., 2009).

6.3.1.5. Misdiagnosis should be reported to the PGD Consortium each year during routine data collection (Thornhill et al., 2005, Goossens et al., 2009).

6.3.1.6. Most of the causes of misdiagnosis are avoidable by taking preventative action and following the principles of quality management present in modern accredited diagnostic testing laboratories. The ESHRE PGD Consortium has recently published an article reviewing the possible causes and adverse outcomes of misdiagnosis (Wilton et al., 2009). It is **recommended** that the suggestions made in this paper for the prevention of specific misdiagnosis causes should be taken into consideration by the PGD center to eliminate the possible causes of misdiagnosis

7. Transport PGD

7.1. For general recommendations on Transport PGD see the Organisation of a PGD centre guideline (Harton et al., 2010a).

References

Alsmadi O, Alkayal F, Monies D, Meyer BF. (2009) Specific and complete human genome amplification with improved yield achieved by phi29 DNA polymerase and a novel primer at elevated temperature. *BMC Res. Notes.*: **24**; 2:48

Altarescu G, Eldar Geva T, Brooks B, Margalioth E, Levy-Lahad E, Renbaum P. (2008) PGD on a recombinant allele: crossover between the TSC2 gene and 'linked' markers impairs accurate diagnosis. *Prenat. Diagn.*: **28**; 929-33.

De Vos A, Staessen C, De Rycke M, Verpoest W, Haentjens P, Devroey P, Liebaers I, Van de Velde H. (2009) Impact of cleavage-stage embryo biopsy in view of PGD on human blastocyst implantation: a prospective cohort of single embryo transfers. *Hum. Reprod.*, in press

ESHRE PGD Consortium Steering Committee (1999) ESHRE Preimplantation Genetic Diagnosis (PGD) Consortium: preliminary assessment of data from January 1997 to September 1998. *Hum. Reprod.* **14**(12); 3138-3148.

ESHRE PGD Consortium Steering Committee (2000) ESHRE Preimplantation Genetic Diagnosis (PGD) Consortium: Data collection II (May 2000). *Hum. Reprod.* **15**(12); 2673-2683.

ESHRE PGD Consortium Steering Committee (2002) ESHRE Preimplantation Genetic Diagnosis (PGD) Consortium: Data collection III (May 2001). *Hum. Reprod.* **17**(1): 233-246.

Fiorentino F, Kokkali G, Biricik A, Stavrou D, Ismailoglu B, De Palma R, Arizzi L, Harton G, Sessa M, and Pantos K. PCR-based detection of chromosomal imbalances on embryos: the evolution of PGD for chromosomal translocations (2010). In press.

Glentis S, SenGupta S, Thornhill A, Wang R, Craft I, Harper JC. (2009) Molecular comparison of single cell MDA products derived from different cell types. *Reprod. Biomed. Online.*: **19**; 89-98.

Goossens V, De Rycke M, De Vos A, Staessen C, Michiels A, Verpoest W, Van Steirteghem A, Bertrand C, Liebaers I, Devroey P, Sermon K. (2008) Diagnostic efficiency, embryonic development

and clinical outcome after the biopsy of one or two blastomeres for preimplantation genetic diagnosis. *Hum. Reprod.*: **23**; 481-92.

V. Goossens, G. Harton, C. Moutou, J. Traeger-Synodinos, M. Van Rij, and J.C. Harper (2009) ESHRE PGD Consortium data collection IX: cycles from January to December 2006 with pregnancy follow-up to October 2007 *Human Reproduction*, doi:10.1093/humrep/dep059

Harton, G, Braude, P, Lashwood, A, Schmutzler, A, Wilton, L, and Harper, JC (2010a) ESHRE PGD Consortium-Best Practice Guidelines for Organization of a PGD Center for Preimplantation Genetic Diagnosis/Screening (PGD/PGS) for Human Reproduction

Harton, G, Coonen, E, Pehlivan, T, Vesela, K, Wilton, L, and Harper, JC (2010b) ESHRE PGD Consortium Best Practice Guidelines for FISH-Based Preimplantation Genetic Diagnosis (PGD) for Human Reproduction

Harton, G, Anick DeVos, Rachel Levy, Kersti Lundin, Cristina Magli, Markus Montag, Monica Parriego, and Harper, JC (2010c) ESHRE PGD Consortium/Embryology Special Interest Group-Best Practice Guidelines for Polar Body and Embryo Biopsy for Preimplantation Genetic Diagnosis/Screening (PGD/PGS). For Human Reproduction

Harper JC, SenGupta S, Vesela K, Thornhill A, Dequeker E, Coonen E, Morris MA (2010) Accreditation of the PGD laboratory. *Hum. Reprod.* (Epub)

Hattori M, Yoshioka K, Sakaki Y. (1992) High-sensitive fluorescent DNA sequencing and its application for detection and mass-screening of point mutations. *Electrophoresis*: **13**; 560–565.

Kuliev A, Rechitsky S, Verlinsky O, Tur-Kaspa I, Kalakoutis G, Angastiniotis M, Verlinsky Y. (2005) Preimplantation diagnosis and HLA typing for haemoglobin disorders. *Reprod. Biomed. Online.*: **11**; 362-70.

Lewis CM, Pinêl T, Whittaker JC, Handyside AH. (2001) Controlling misdiagnosis errors in preimplantation genetic diagnosis: a comprehensive model encompassing extrinsic and intrinsic sources of error. *Hum. Reprod.*: **16**; 43-50.

Lissens, W and Sermon, K. (1997) Preimplantation genetic diagnosis: current status and new developments. *Hum. Reprod.*: **12**; 1756–1761.

Lovmar L, Fredriksson M, Liljedahl U, Sigurdsson S, Syvanen AC. (2003) Quantitative evaluation by minisequencing and microarrays reveals accurate multiplexed SNP genotyping of whole genome amplified DNA. *Nucleic Acids Res.* **31**, e129.

Malfroy L, Roth MP, Carrington M, Borot N, Volz A, Ziegler A, Coppin H. (1997) Heterogeneity in rates of recombination in the 6-Mb region telomeric to the human major histocompatibility complex. *Genomics*: **43**; 226–231.

Martin M, Mann D, Carrington M. (1995) Recombination rates across the HLA complex: use of microsatellites as a rapid screen for recombinant chromosomes. *Hum. Mol. Genet.*: **4**; 423–428.

Panelli S, Damiani G, Espen L, Micheli G, Sgaramella V. (2006) Towards the analysis of the genomes of single cells: further characterisation of the multiple displacement amplification. *Gene.*: **372**; 1–7.

Preimplantation Genetic diagnosis International Society (2004) The preimplantation Genetic diagnosis International Society (PGDIS): Guidelines for good practice in PGD. *Reprod. Biomed. Online* **9**, 430-434

Preimplantation Genetic Diagnosis International Society (PGDIS). (2008) Guidelines for good practice in PGD: programme requirements and laboratory quality assurance. *Reprod. Biomed. Online* **16**; 134-47.

Pickering SJ, McConnell JM, Johnson MH, Braude PR. (1994) Use of a polymorphic dinucleotide repeat sequence to detect non-blastomeric contamination of the polymerase chain reaction in biopsy samples for preimplantation diagnosis. *Hum. Reprod.* : **9**; 1539-45.

Piyamongkol W, Bermudez MG, Harper JC, Wells D. (2003) Detailed investigation of factors influencing amplification efficiency and allele drop-out in single cell PCR: implications for preimplantation genetic diagnosis. *Mol. Hum. Reprod.*: **9**; 411–420.

Renwick PJ, Trussler J, Ostad-Saffari E, Fassihi H, Black C, Braude P, Ogilvie CM, Abbs S. (2006)

Proof of principle and first cases using preimplantation genetic haplotyping- a paradigm shift for embryo diagnosis. *Reprod. Biomed. Online.*; **13**: 110-119.

Renwick PJ, Lewis CM, Abbs S, Ogilvie CM. (2007) Determination of the genetic status of cleavage-stage human embryos by microsatellite marker analysis following multiple displacement amplification. *Prenat. Diagn.*: **27**; 206-15

Renwick P, Ogilvie CM. (2007) Preimplantation genetic diagnosis for monogenic diseases: overview and emerging issues. *Expert Rev. Mol. Diagn.*; **7**: 33-43.

Rechitsky S, Strom, C, Verlinsky, O, Amet T, Ivakhnenko V, Kukharenko V, Kuliev A and Verlinsky Y (1998) Allele dropout in polar bodies and blastomeres. *J. Assist. Reprod. Genet.*; **15**: 253-7

Rechitsky S, Kuliev A, Tur-Kaspa I, Morris R, Verlinsky Y. (2004) Preimplantation genetic diagnosis with HLA matching. *Reprod. Biomed. Online*: **9**; 210–221.

Sermon K, De Vos A, Van de Velde H, Seneca S, Lissens W, Joris H, Vandervorst M, Van Steirteghem A, and Liebaers I. (1998) Fluorescent PCR and automated fragment analysis for the clinical application of preimplantation genetic diagnosis of myotonic dystrophy (Steinert's disease). *Mol. Hum. Reprod.*: **4**; 791–795.

Sermon K, De Rijcke M, Lissens W, De Vos A, Platteau P, Bonduelle M, Devroey P, Van Steirteghem A, Liebaers I. (2002) Preimplantation genetic diagnosis for Huntington's disease with exclusion testing. *Eur. J. Hum. Genet.*; **10**; 591-8

Stern HJ, Harton GL, Sisson ME, Jones SL, Fallon LA, Thorsell LP, Getlinger ME, Black SH, Schulman JD. (2002) Non-disclosing preimplantation genetic diagnosis for Huntington disease. *Prenat Diagn.* Jun;22(6):503-7.

The Practice Committee of the Society for Assisted Reproductive Technology and the Practice Committee of the American Society for Reproductive Medicine (2008) Preimplantation genetic testing: A Practice Committee opinion. *Fertil. Steril.* **90**: 136-43

Thornhill AR and Snow K. (2002) Molecular diagnostics in preimplantation genetic diagnosis. *J. Mol. Diagn.*: **4**; 11-29

Thornhill AR, Guenther AJ, Barbarotto GM, Session DR, Damario MA, Dumesic DA, Snow K. (2002) False-positive Y-microdeletion result for a fertile male caused by an alteration under a PCR primer. *Int J Androl.*; **25**: 352-7.

Thornhill AR, De Die-Smulders CE, Geraedts JP, Harper JC, Harton G L, Lavery SA, Moutou C, Robinson MD, Schmutzler AG, Scriven PN, Sermon KD, Wilton L (2005) ESHRE PGD Consortium "Best practice guidelines for clinical preimplantation genetic diagnosis (PGD) and preimplantation genetic screening (PGS)". *Hum. Reprod.*: **20**; 35–48.

Thornhill AR, Repping S. (2008) Quality control and quality assurance in preimplantation genetic diagnosis. In: Harper J (ed). *Preimplantation Genetic Diagnosis*. UK: Wiley and Sons,

Van de Velde H, De Rycke M, De Man C, De Hauwere K, Fiorentino F, Kahraman S, Pennings G, Verpoest W, Devroey P, Liebaers I. (2009) The experience of two European preimplantation genetic diagnosis centres on human leukocyte antigen typing. *Hum. Reprod.*: **24**; 732-40.

Verlinsky Y, Rechitsky S, Sharapova T, Morris R, M Taranissi M, Kuliev A. (2004) Preimplantation HLA testing. *JAMA*: **291**; 2079–2085.

Wilton L, Thornhill A, Traeger-Synodinos J, Sermon KD, Harper JC. (2009) The causes of misdiagnosis and adverse outcomes in PGD. *Hum. Reprod.*: **24**; 1221-8