

Now sold under the  
Thermo Scientific brand

**Thermo**  
SCIENTIFIC

# Determination of Melamine in Milk Powder by Reversed-Phase HPLC with UV Detection

## INTRODUCTION

Recent investigations of death and health problems of babies in China have revealed that some baby foods (milk powder) have been contaminated by melamine (structure shown in Figure 1). Some manufacturers illegally used melamine as an adulterant to increase the apparent protein content. Standard methods enacted by the Chinese government for determining melamine in raw milk and dairy products included HPLC-UV, LC-MS, and GC-MS methods.<sup>1</sup> However, the high cost of operation and maintenance of GC/LC-MS systems as well as the labor intensive derivatization that GC-MS requires limits their use in the milk product factories. The HPLC-UV method therefore is presently the popular choice for most factories. In this method, melamine is separated on a C<sub>8</sub> or C<sub>18</sub> column using an ion pair buffer (mixture of citric acid and sodium 1-octane sulfonate) and acetonitrile mobile phase.




Figure 1. Structure of melamine.

In this Application Note (AN), we determined melamine in milk powder samples following the regulated method.<sup>1</sup> Melamine was separated from other components in the powdered milk samples using an Acclaim® 120 C18 column and an UltiMate® 3000 HPLC system with UV detection. The results of method detection limits (MDL), recovery, and permitted detection deviation match the requirements in the regulated method.

## EQUIPMENT

Dionex UltiMate 3000 HPLC system consisting of:

HPG 3400A pump

WPS 3000TSL autosampler

TCC-3000 thermostatted column compartment

VWD-3400 UV-vis detector

Chromeleon® 6.80 SP5 Chromatography Data System

Kudos® SK3200LH ultrasonic generator, Kudos

Ultrasonic Instrumental Co., Shanghai, China

Mettler Toledo AL-204 Electrolab balance, Mettler

(Shanghai) Co., Shanghai, China

Anke® TGL-16B centrifuge, Anting Scientific

Instrumental Factory, Shanghai, China

IKA® MS1 Minishaker, IKA Works, Guangzhou, China

Strata™-x-c SCX SPE (Phenomenex) column

SE-506 Nitrogen Purge Instrument, Shine Tech., Beijing, China

## REAGENTS AND STANDARDS

Water, from Milli-Q® Gradient A 10

Methanol (CH<sub>3</sub>OH), HPLC grade, Fisher

Acetonitrile (CH<sub>3</sub>CN), HPLC grade, Fisher

Trichloroacetic acid, analytical grade, SCRC, China

Citric acid (C<sub>6</sub>H<sub>8</sub>O<sub>7</sub>·H<sub>2</sub>O), analytical grade, SCRC, China

Sodium 1-octane sulfonate (98%), Baker Analyzed @  
HPLC Reagent, USA

Melamine (99.0%), HPLC grade, Fluka

Nitrogen (N<sub>2</sub>, 99.999%), Lumin Gas Works, Shanghai,  
China

Ammonia solution (25%–28%), analytical grade, SCRC,  
China

## CHROMATOGRAPHIC CONDITIONS

Guard Column: Acclaim 120 C18, 5 μm,  
4.3 × 10 mm, P/N 059446,  
with guard column holder,  
P/N 59526

Analytical Column: Acclaim 120 C18, 5 μm,  
4.6 × 250 mm, P/N 059149

Mobile Phase: Buffer (dissolve 2.10 g citric acid  
and 2.16 g sodium 1-octane sulfonate  
in 980 mL H<sub>2</sub>O, adjust pH value to  
3.0 with 1 M NaOH solution, add  
water to the mark of 1000-mL  
volumetric flask)—  
CH<sub>3</sub>CN (92 : 8, v/v)

Column Temp.: 40 °C

Flow Rate: 1.0 mL/min

Inj. Volume: 20 μL

UV Detection: Absorbance at 240 nm

## PREPARATION OF STANDARDS

### Stock Standard Solution

Accurately weigh ~100 mg of melamine, dissolve  
in a 100-mL volumetric flask with aqueous methanol  
(50%, v/v). The concentration of melamine is  
1000 μg/mL.

## Working Standard Solutions

Prepare seven working standard solutions for  
calibration by adding defined volumes of the stock  
standard solution and diluting with mobile phase. The  
concentrations of melamine are 0.2, 0.5, 2.0, 20, 25, 50  
and 100 μg/mL, respectively.

## SAMPLE PREPARATION

### Sample Extraction

Put an accurately weighed ~2 g of dried sample  
to a 50-mL centrifuge tube, and then add 15 mL  
aqueous trichloroacetic acid (1%, v/v) and 5 mL  
acetonitrile. After 1 min of vortex shaking, put in an  
ultrasonic bath for 30 min, and then shake for 10 min.  
After 10 min of centrifugation (setting = rpm ≥ 10,000),  
move the supernatant to a 25-mL volumetric flask while  
passing through filter paper, and add the 1% aqueous  
trichloroacetic acid to the mark.

### Cleaning Sample Extract on an SCX SPE Column

Prior to use, the SPE column should be activated by  
passing 3 mL CH<sub>3</sub>OH and 5 mL H<sub>2</sub>O in turn.

Mix 5 mL of the sample extract and 5 mL water,  
move them to the activated SCX SPE column. Wash  
the SPE column with 3 mL methanol and 3 mL water,  
respectively, then elute with 6 mL of aminated methanol  
solution (mixture of 5 mL ammonia solution and  
95 mL methanol). Dry the collected eluent with N<sub>2</sub> at  
50 °C, dissolve the residue in 1-mL mobile phase, and  
then vortex the solution for 1 min. Prior to injection, filter  
the solution through a 0.2-μm filter (Millex®-HV).

### Spiked Milk Powder Samples

Add 20 μL of stock standard solution of melamine to  
the 50-mL centrifuge tubes together with the accurately  
weighed ~2 g of dried sample. The remainder of the  
sample preparation procedure is the same as that for the  
milk powder sample.

## RESULTS

### Method Reproducibility

Figure 2 shows an overlay of chromatograms of  
seven melamine standards with different concentrations.  
The RSD for retention time is 0.143%.




Figure 2. Overlay of chromatograms of seven melamine standards with different concentrations.

### Linearity and Detection Limit

Calibration linearity for melamine was investigated by making three replicate injections of each standard prepared at seven different concentrations. The external standard method was used to establish the calibration curve and to quantify melamine in samples. As shown in Figure 3, excellent linearity was achieved throughout the range from 0.2 to 100  $\mu$ g/mL. The linearity equation of melamine is as follows with the curve forced the through origin.

$$A = 1.2788 c$$

Where,  $A$  stands peak area, and  $c$  stands melamine concentration ( $\mu$ g/mL). The correlation coefficient ( $r$ ) is 0.999961.




Figure 3. Calibration curve for melamine.

Method detection limit (MDL) of melamine was calculated by using  $S/N = 3$ , where  $S$  = signal,  $N$  = noise. The calculated value of MDL is 0.17  $\mu$ g/mL.

### Sample Analysis

Five milk powder samples (1–5) obtained from a manufacturer were analyzed. Melamine was found in samples #1 through #4. An amount of melamine less than the MDL was detected (0.12  $\mu$ g/mL) in sample 5 and therefore can not be reliably identified as melamine. The results are summarized in Table 1. Figure 4 is an overlay chromatograms of these samples.

**Table 1. Sample Analysis Data**

Sample #	1				2	3	4	5
Melamine	Detected (µg/mL)	Added (µg/mL)	Found (µg/mL)	Recovery (%)	Detected (µg/mL)	Detected (µg/mL)	Detected (µg/mL)	Detected (µg/mL)
	0.83	4.0	3.2	80	0.58	0.61	3.6	0.12

Notes: 1. 3 injections were made for each sample

2. Spiked sample was prepared according to the description in the sample section.

3. The detected amount in the list is for the prepared sample (diluted sample).

4. The amount in original sample should be calculated by multiplying "Diluted fold", which is 2.5 in this Application Note.

5. Found = Measured Value of spiked sample - Measured Value of sample.




Figure 4. Chromatograms of five milk powder samples.

## CONCLUSION

The determination of melamine in milk powder was performed on the Acclaim 120 C18 column and UltiMate 3000 HPLC system with UV detection following the regulated method.<sup>1</sup> Table 2 lists the comparison of chromatographic performance data obtained in this experiment to the requirements in the standard method, which demonstrates that the results match the requirements in the regulated method.

**Table 2. Method Performance**

	Method in this Application Note	Requirement in the Regulated Method
Correlation coefficient (r)	0.999961	\
Linearity range (µg/mL)	0.2 ~ 100	0.8 ~ 80
Recovery (%)	80	80 ~ 110

Note: Regulated method, see Ref. 1.

\ Not specified.

## REFERENCE

1. Determination of Melamine in Raw Milk and Dairy Products, *GB/T 22388 2008*.

Acclaim, Chromeleon, and UltiMate are registered trademarks of Dionex Corporation.

Kudos is a registered trademark of Kudos Ultrasonic Instrumental Company.

Anke is a registered trademark of Anting Scientific Instrumental Factory.

IKA is a registered trademark of IKA Works.

Milli-Q and Millex are registered trademarks of Millipore Corporation.

Strata is a trademark of Shine Technology.

Passion. Power. Productivity.



### Dionex Corporation

1228 Titan Way  
 P.O. Box 3603  
 Sunnyvale, CA  
 94088-3603  
 (408) 737-0700

### North America

U.S./Canada (847) 295-7500  
**South America**  
 Brazil (55) 11 3731 5140

### Europe

Austria (43) 1 616 51 25 Benelux (31) 20 683 9768 (32) 3 353 4294  
 Denmark (45) 36 36 90 90 France (33) 1 39 30 01 10 Germany (49) 6126 991 0  
 Ireland (353) 1 644 0064 Italy (39) 02 51 62 1267 Sweden (46) 8 473 3380  
 Switzerland (41) 62 205 9966 United Kingdom (44) 1276 691722

### Asia Pacific

Australia (61) 2 9420 5233 China (852) 2428 3282 India (91) 22 2764 2735  
 Japan (81) 6 6885 1213 Korea (82) 2 2653 2580 Singapore (65) 6289 1190  
 Taiwan (886) 2 8751 6655

[www.dionex.com](http://www.dionex.com)



LPN 2184 PDF 3/09  
 ©2009 Dionex Corporation