

## Thayer Martin Agar (Modified) Procedure

### Principle

Thayer Martin Agar (Modified) is a solid medium used commonly for the primary isolation of *Neisseria gonorrhoeae* from mixed specimens. The agar can also be utilized for primary isolation of *Neisseria meningitidis* from mixed specimens. The agar is classified as a selective enrichment agar. Enrichments added to this medium include both X and V factors. The modified formulation of the Thayer Martin agar includes more agar to help prevent swarming *Proteus*. The agar contains antibiotics to inhibit the growth of normal flora, non-pathogenic *Neisseria* species and most other organisms. *Neisseria gonorrhoeae*, *Neisseria meningitidis*, and *Neisseria lactamica* will grow on the agar. *Neisseria lactamica* is usually non-pathogenic.

The antibiotics in the agar include: vancomycin to inhibit gram-positive organisms, nystatin to inhibit the growth of fungi, colistin to inhibit most gram-negative rods, and trimethoprim helps to prevent *Proteus* from swarming.

### Specimen Collection and Preparation

The original swab specimen should be inoculated at bedside to the plate by rolling the swab in a large "Z" pattern to sufficiently transfer specimen. The plate should be received in the lab within 2 hours. If direct inoculation at bedside is not available the lab should receive the specimen within 2 hours for inoculation to media.

If transport will be delayed a Carbon Dioxide transport system should be used. An agar version called Jembec is available that includes a Carbon Dioxide generations system within the agar. If using a Jembec plate: Using forceps, remove the CO<sub>2</sub> tablet from the foil pouch and plate in the well. Place plate in environmental maintenance pouch, seal and secure. Transport to lab.

### Reagents

Thayer Martin Agar (Modified)	Inoculating loop
Incinerator	Aerobic swab collection system

### Storage

1. Store Thayer Martin agar (Modified) at 2-8°C and bring to room temperature before use.

### Quality Control

Quality control should be performed per lot/shipment date.

<i>Neisseria gonorrhoeae</i> ATCC 43069	Expected results:	Growth
<i>Proteus mirabilis</i> ATCC 43071	Expected results:	Inhibition (partial)
<i>Staphylococcus epidermidis</i> ATCC 12228	Expected results:	Inhibition (partial)

### Procedure

1. Inoculate the original swab specimen to the plate by rolling the swab in a large "Z" pattern. This will sufficiently transfer the specimen.
2. Cross-streak the plate using a sterile wire loop.
3. Incubate in 3-7% CO<sub>2</sub> at 35-37°C and examine at 24 hours.
4. Isolates should be gram stained. If gram-negative diplococci are present colonies should be tested for oxidase production. Suspicious colonies should be identified.
5. If no growth is observed, re-incubate plate for up to 72 hours.
6. **Or, inoculate isolated colonies from culture to plate and observed for growth at 24 hours after incubation in 3-7% CO<sub>2</sub> at 35-37°C.**

### References:

1. Remel package insert, IFU 1880, Lenexa, KS. September 2003
2. Mahon, C.R. and Manuselis, G., Textbook of Diagnostic Microbiology, 3<sup>rd</sup> ED., W.B. Saunders, 2007.