

AsB Detector

Showing Material Contrast and Topographical Information with the AsB Detector



AsB Detector

Showing Material Contrast and Topographical Information with the AsB Detector

Contact: Dr. Ulrich Kohl-Roscher

Date: December 2013

Visualize crystal orientation and compositional contrast.

Introduction

With an integrated GEMINI lens detector, you are able to select and separate the BSE signal by both Z contrast and angle contrast. The high angle BSE electrons are detected in the unique, in-column, energy selective, EsB detector. Large and very large angles, coming from different scattering processes, are collected in the AsB detector or the angular selective BSE detection system.

Availability


The AsB detector is available for the following Field Emission Scanning Electron Microscope series:

- Σ -series
- 15xx (from serial no. 2000)
- SUPRA series
- ULTRA series
- MERLIN series (Smart AsB)

Benefits

- Contrast and topography
The AsB detector allows you to show channelling contrast (crystal orientation), as well as compositional contrast.
- Rutherford (RBS) and Channelling (Mott scattered) BSE electrons can be separated due to the different multiple inelastic (RBS = Z-contrast) and single elastic (MBS = crystallographic contrast) scattering mechanisms. The different scattering processes can be detected as there is almost no magnetic field (<0.3mG) outside of the GEMINI column. This means Z-contrast and crystallographic (orientation or deformation contrast) information can be used to characterize crystalline samples. As multiple inelastic scattering is minimised, it also means the images appear very sharp.
- Ease of use
The AsB detector is completely integrated into the pole piece of the GEMINI lens. This allows you to image BSEs with ultra short working distances without aligning the AsB detector to the optical axis.

Figure 1




Schematics of the AsB detector
 When the primary electron beam hits the specimen, backscattered electrons (BSEs) are released from the specimen. These BSEs are detected by the diodes and generate an image. 1: integrated Detector, 2: e-beam and LL-BSE (only with EsB detectable), 3: Rutherford BSE, 4: Channelling BSE

Figure 2


Operation

The AsB detector is fully software-controlled. It is equipped with four diodes that can be controlled independently via a menu. Compositional mode (COMPO) produces images showing the atomic contrast of the specimen. Topography mode (TOPO) shows surface details.



BSD panel, which controls the AsB detector

Figure 3



A Rutherford BSE image of an aluminium specimen shows barely visible „scratches“. See Figure 4 for an image of the specimen in MBS mode, detecting channelling BSEs.

Upgrade path

Requires SmartSEM V05.01.04 or higher.

On all instruments except MERLIN: requires a free port on the scan generator. Alternatively, a multiplexer board can be retrofitted to allow an additional detector to be operated.


On MERLIN: If there are already four detectors installed on the MERLIN system, an additional ADC4 detector board must be retrofitted (since the standard configuration is meant to handle only four detectors).

ADC4 detector board (MERLIN)	346001-9220-000
------------------------------	-----------------

The retrofit must be done by an authorised Carl Zeiss Microscopy service engineer. For further information, contact: microscopy@zeiss.com


Part	Ordering no.
AsB Upgrade Kit, complete kit including all needed hardware.	347806-9014-990
Smart AsB Upgrade Kit (for MERLIN only)	346006-9010-990

Figure 4



With an AsB image in MBS mode it becomes clear that the „scratches“ are dislocations. Vertical broad shadows are orthogonal to the rolling direction.

Figure 5



Strain deformation in steel with razor sharp dislocations



Carl Zeiss Microscopy GmbH

07745 Jena, Germany

microscopy@zeiss.com

www.zeiss.com/microscopy



We make it visible.