


Imaging of Biological Systems


Size (m) 10^{-8} 10^{-7} 10^{-6} 10^{-5} 10^{-4} 10^{-3} 10^{-2} 10^{-1} 1




Molecules




Complexes




Viruses




Bacteria



Cellular compartments



Cells



Tissues



Whole organisms

Cryo-electron microscopy

Light microscopy

X-ray crystallography

- atomic resolution
- needs protein crystals


MRI

X-ray microscopy

- 20 nm resolution
- thick objects (compared to EM)

X-ray Crystallography

- **Fixed energy beam**
 - rotate crystal through small angle while recording data
 - repeat N times through large angle given by crystal symmetry




Crystal


Diffraction data

Electron density

Model




X-ray Crystallography



Crystal in loop





Loop on pin



X-ray Crystallography

- determine electron density (ie. recover phase information)
- fit known sequence information into electron density



X-ray Crystallography




PS. Just skipped over 40 years of developments!

it used to be more difficult!




STRUCTURE OF THE AcrB PUMP

Bacterial Drug Transporter Reveals Secrets of Drug Resistance




Even ciprofloxacin, an antibiotic used to treat a variety of bacterial infections including inhaled anthrax, is no match for AcrB. In this image, the green-colored drug is firmly ensnared in the protein's cavity.

Transcription Initiation and Elongation by RNA Polymerase II*RNA polymerase II pre-initiation complex model.*

Transcription Initiation and Elongation by RNA Polymerase II

Close-up of strand separation. RNA (red) and DNA (blue) strands are separated by a network of interactions with three protein loops. The loops limit the extent of RNA separation (orange), stabilize the separated single strand, and form part of the RNA exit pore (green and pink).




Imaging of Biological Systems


Size (m) 10^{-8} 10^{-7} 10^{-6} 10^{-5} 10^{-4} 10^{-3} 10^{-2} 10^{-1} 1




Molecules




Complexes




Viruses




Bacteria



Cellular compartments



Cells



Tissues



Whole organisms

Cryo-electron microscopy

Light microscopy

X-ray crystallography


- atomic resolution
- needs protein crystals


MRI

X-ray microscopy

- 20 nm resolution
- thick objects (compared to EM)

X-ray Microscopes






Beamline 6.1.2
Full-field Transmission X-ray
Microscope (XM-1)


Wavelength, 2.4 nm
Photon energy, 517 eV

Center for X-ray Optics
Erik Anderson
David Attwood





A Fresnel Zone Plate Lens





Professor David Attwood
ART PROCESSING
Univ. California, Berkeley

$D_r = 25 \text{ nm}$
 $D = 63 \text{ }\mu\text{m}$
 $N = 618 \text{ zones}$
 $f = 650 \text{ }\mu\text{m}$
 $NA = 0.05$
 $@ 2.4 \text{ nm}\lambda$



Zone Plate Lenses 20 nm Resolution



E. Anderson, D. Attwood, A. Liddle, D. Olynick, B. Harteneck, W. Chao; Center for X-ray Optics






Cryo X-ray Tomography

C.A. Larabell & M. A. Le Gros (2004).
Molecular Biology of the Cell, 15(3), 956-962






C.A. Larabell & M. A. Le Gros (2004).
Molecular Biology of the Cell, 15(3), 956-962



C.A. Larabell & M. A. Le Gros (2004).
Molecular Biology of the Cell, 15(3), 956-962

Saccharomyces cerevisiae





Reconstructed data using different volume analysis algorithms.


- A) Combination of translucent outer surface & opaque surfaces highlight internal organelles; arrow points to nucleus that has been color-coded blue
- B) Surface views combined with volume rendering
- lipid droplets white
 - surface of large vacuole pink
- C) Computer section that was volume-rendered according to x-ray absorption
- lipid droplets white
 - internal structures of vacuole and nucleus red
 - other cytoplasmic structures appear green and orange.

Angle Resolved Photo-Electron Spectroscopy (ARPES)


The Electronic Structure Factory



"2-Dimensional" (E vs k_x)
Bandstructure Measurement
Built up by tiling detector images





Angle Resolved Photo-Electron Spectroscopy (ARPES)



Eli Rotenberg et al, BL 7.0

Angle Resolved Photo-Electron Spectroscopy (ARPES)


- 3d and 4d datasets built up from angle scanning the sample manipulator




Angle Resolved Photo-Electron Spectroscopy (ARPES)

- 2d k-space mapping and building of a 3d map (Chromium)

model



expt




Eli Rotenberg et al, BL 7.0

expt, E. Rotenberg, ALS, and S. D. Kevan, Univ. Oregon


Angle Resolved Photo-Electron Spectroscopy (ARPES)

- 3-D Fermi Surface (extracted from 4-D data)
 - For example, rare earth ferromagnet Terbium

Paramagnetic
 $T=240^{\circ}\text{K}$




Ferromagnetic
 $T=30^{\circ}\text{K}$



exchange-split Fermi Surface

Fermi surface data courtesy of K. Starke, FU Berlin

Angle Resolved Photo-Electron Spectroscopy (ARPES)



- **Count Rate is the most significant detector limitation**
 - Could be improved with faster phosphors, lower readout noise of CCD, faster data collection and processing
- **Angle Resolution limited approx. equally by**
 - sizes of electron spot on phosphor, # pixels on CCD camera, photon beam on sample