


**METADATA AND NUMERICAL DATA CAPTURE:
AZEOTROPIC pressure, p
(2 components)**

***Guided Data*
Capture (GDC)**



This tutorial describes
**METADATA AND NUMERICAL DATA CAPTURE:
AZEOTROPIC pressure p
(2 Components)**
with the Guided Data Capture (GDC) software.

NOTE:

The tutorials proceed sequentially to ease the descriptions. **It is not necessary to enter *all* compounds before entering *all* samples, etc.**

Compounds, samples, properties, etc., can be added or modified at any time.

However, the hierarchy must be maintained (i.e., a property cannot be entered, if there is no associated sample or compound.)

The experimental data used in this example are from:

J. Chem. Eng. Data 2001, 46, 535–540

535

Isothermal Vapor–Liquid Equilibrium of 1-Chlorobutane with Ethanol or 1-Hexanol at Ten Temperatures between 278.15 K and 323.15 K

Santiago Martínez, Rosa Carriga, Pascual Pérez, and Mariano Gracia*

Departamento de Química Orgánica y Química Física (Área de Química Física), Facultad de Ciencias, Universidad de Zaragoza, 50009-Zaragoza, Spain

Vapor pressures of (1-chlorobutane + ethanol or 1-hexanol) at 10 temperatures between 278.15 and 323.15 K were measured by a static method. The reduction of the vapor pressures to obtain activity coefficients and excess molar Gibbs energies was carried out by fitting the vapor pressure data to the Wilson equation according to Barker's method. In the 1-chlorobutane + ethanol system, azeotropic mixtures with a minimum boiling point temperature were observed over the whole temperature range.

AZEOTROPIC pressure for 1-chlorobutane + ethanol

Table 4. Azeotropic Pressures and Mole Fractions for the System $\{(1 - z)$ 1-Chlorobutane + z Ethanol $\}$

T/K	$z(\text{exptl})$	$P_z(\text{exptl})/\text{kPa}$	$z(\text{caled from eq 11})$	$P_z(\text{caled from eq 13})/\text{kPa}$
278.15	0.245	6.172	0.244	6.213
283.15	0.259	8.166	0.259	8.167
288.15	0.274	10.640	0.274	10.633
293.15	0.288	13.756	0.289	13.720
298.15	0.304	17.612	0.304	17.553
303.15	0.318	22.383	0.319	22.274
308.15	0.334	28.128	0.334	28.047
313.15	0.349	35.086	0.349	35.058
318.15	0.364	43.410	0.364	43.516
323.15	0.379	53.328	0.379	53.654

This data set is considered here.

Method Information:

For (1-chlorobutane + ethanol), azeotropic mixtures with a minimum boiling temperature were observed over the whole range of temperature. Azeotropic mole fractions z were graphically calculated, assuming ideal behavior of the vapor, from the well-known equation, $\gamma_1/\gamma_2 = P_2^s/P_1^s$.

NOTE: Although the focus of the data-capture is experimentally determined values, derived azeotropic properties were considered of adequate importance to be included.

Guided Data Capture - Thermophysical and Thermochemical Data

File Edit Tools Help

Reference

Compound

Sample

Mixture

Reaction

Property

Data Tables

- 2001 mar gar 0
 - 1-chlorobutane
 - Sample 1 (cm,99.5m%,nc;mv;99.9m%,glc)
 - ^1: P (L), Set 1, B Method:CCELL dPV=.01 dT=.01
 - ethanol
 - Sample 1 (cm,99.8m%,nc;mv;99.9m%,glc)
 - ^1: P (L), Set 1, B Method:CCELL dPV=.01 dT=.01
 - 1-hexanol
 - Sample 1 (cm,99.0m%,nc;mv;99.4m%,glc)
 - ^1: P (L), Set 1, B Method:CCELL dPV=.01 dT=.01
 - ethanol + 1-chlorobutane**
 - ^2: vle, P (Set 1), B Method:CCELL dP=.01 dT=.01 dX1=0.0003
 - 1-chlorobutane + 1-hexanol
 - ^2: vle, P (Set 1), B Method:CCELL dP=.01 dT=.01 dX2=0.0003

2. CLICK
Property

1. SELECT the *mixture* for which
the data are to be captured.

NOTE: The **bibliographic information, compound identities, sample descriptions, and mixture** were entered previously. (There are separate tutorials, which describe capture of this information, if needed.)

Property and experimental method for ethanol + 1-chlorobutane

Help

Property group: Vapor pressure; Boiling temperature; and Azeotropic T & P

Property: Azeotropic pressure

Units:

Method

Experimental

1. SELECT the **Property group**: *Vapor pressure; Boiling temperature; and Azeotropic T& p* from the menu.

2. SELECT the **Property**: *Azeotropic pressure*.

Comment (optional)

OK

Cancel

Property and experimental method for ethanol + 1-chlorobutane

Help

Property group: Vapor pressure; Boiling temperature; and Azeotropic T & P

Property: Azeotropic pressure

Units: kPa

Method of measurement:

Experimental

Comment (optional)

OK Cancel

SELECT the **UNITS; *kPa* here.**

If necessary, SELECT *ALL OTHER UNITS* to specify an unusual unit.

1. SELECT Method of Measurement from the list provided. **NOTE:** *Other* can be a valid selection and should include a brief description in the **Comment** field.

Units: kPa

Method of measurement: By P(X) extreme

Experimental purpose: Principal objective of the work

2. SELECT the Experimental Purpose from the list provided.

Comment
(optional)

3. CLICK OK

OK

Cancel

Specification of # of Phases in Equilibrium and # of Constraints

Azeotropic pressure (kPa) as function of 1 variable(s)

Mixture: ethanol + 1-chlorobutane

Phases in equilibrium: 2 Constraints: 1 Independent variables: 1 Property set #: 1

Phase of the Property Value(s) Liquid

Phase 2 Gas

Constraint 1 (Fixed value of) Azeotrope

Independent variable 1

NOTE: The number of **Phases in equilibrium** (2, *Liquid* and *gas*) and number of **Constraints** (1, *azeotrope*) are entered automatically by the software.

Azeotropic pressure (kPa) as function of 1 variable(s)

Mixture: ethanol + 1-chlorobutane

Phases in equilibrium: 2 Constraints: 1 Independent variables: 1 Property set #: 1 Sample #: 1 Sample #: 1

Phase of the Property Value(s): Liquid

Phase 2: Gas

Constraint 1 (Fixed value of): Azeotrope

Independent variable: []

Definition of Mixture: []

Data presentation: Experimental values

Comments (Optional): []

Property and method Numerical Data Cancel

Uncertainty: [] %

Position of the Property Value(s): kPa %

NOTE: Multiple samples for a given component can be accommodated, but this is rarely needed.

Specification of variables, units

1. SELECT the **Independent Variable** (*Temperature*) from the lists.

Phase of the Property Value(s) Liquid

Phase 2 Gas

Constraint 1 (Fixed value of) Azeotrope

Independent variable 1 Temperature

Units: K Uncertainty: %

Definition of Measurement Results (Absolute vs Relative) Direct value

Data presentation Experimental values

Comments (Optional):

Property and method Numerical Data Cancel

2. SELECT **Units** for the Variable. Include approximate **Uncertainties**, if known.

Measurement definition and Data presentation

Azeotropic pressure (kPa) as function of 1 variable(s)

Mixture: ethanol + 1-chlorobutane

Phases in equilibrium: 2 Constraints: 1 Independent variable

Phase of the Property Value(s) Liquid

Phase 2 Gas

Constraint 1 (Fixed value of) Azeotrope

Independent variable 1 Temperature of Liquid Units: K Uncertainty: %

Definition of Measurement Results (Absolute vs Relative)
Direct value

Data presentation
Experimental values

Comments (Optional):

Property and method Numerical Data Cancel

1. SELECT *Direct Value* (as compared with a Relative Value) from the list defining the **Measurement Results**

2. SELECT the **Data presentation** method. *Experimental values* here.

3. CLICK *Numerical Data*

Azeotropic pressure (kPa) as function of 1 variable(s)

File Edit Action Help

	Var 1	Property
1		

TYPE, or much preferably, PASTE the variable and property values into the table.

See next page...

Table 4. Azeotropic Pressures and Mole Fractions for the System $\{(1 - z)$ 1-Chlorobutane + z Ethanol $\}$

T/K	$P_z(\text{exptl})/\text{kPa}$	$z(\text{calcd from eq 11})$	$P_z(\text{calcd from eq 13})/\text{kPa}$
278.15	6.172	0.244	6.213
283.15	8.166	0.259	8.167
288.15	10.640	0.274	10.633
293.15	13.756	0.289	13.720
298.15	17.612	0.304	17.553
303.15	22.383	0.319	22.274
308.15	28.128	0.334	28.047
313.15	35.086	0.349	35.058
318.15	43.410	0.364	43.516
323.15	53.328	0.379	53.654

Clear Cancel

Azeotropic pressure (kPa) as function of 1 variable(s)

File Edit Action Help


	Var 1	Property
1	278.15	6.172
2	283.15	8.166
3	288.15	10.640
4	293.15	13.756
5	298.15	17.612
6	303.15	22.383
7	308.15	28.128
8	313.15	35.086
9	318.15	43.410
10	323.15	53.328

Table 4. Azeotropic Pressures and Mole Fractions for the System $\{(1 - z) \text{ 1-Chlorobutane} + z \text{ Ethanol}\}$

T/K		$P_d(\text{exptl})/\text{kPa}$	$z(\text{calcd from eq 11})$	$P_d(\text{calcd from eq 13})/\text{kPa}$
278.15		6.172	0.244	6.213
283.15		8.166	0.259	8.167
288.15		10.640	0.274	10.633
293.15		13.756	0.289	13.720
298.15		17.612	0.304	17.553
303.15		22.383	0.319	22.274
308.15		28.128	0.334	28.047
313.15		35.086	0.349	35.058
318.15		43.410	0.364	43.516
323.15		53.328	0.379	53.654

Clear the Table

NOTE: Simple CUT/PASTE procedures can be used within the table to convert the original table into the required number of columns. (This can also be done externally in spreadsheet software, e.g., EXCEL.)



1. CLICK *View plot* to see an automatic graphical representation of the data.

2. Check for typographical errors, and CLICK *OK*, when done.

Azeotropic pressure (kPa) as function of 1 variable(s)

File Edit Action Help

	Var 1	Property
1	278.15	6.172
2	283.15	8.166
3	288.15	10.640
4	293.15	13.756
5	298.15	17.612
6	303.15	22.383
7	308.15	28.128
8	313.15	35.086
9	318.15	43.410
10	323.15	53.328

CLICK *Accept*, when checking is complete.

Clear the Table

View plot

Accept

Cancel

Guided Data Capture - Thermophysical and Thermochemical Data

File Edit Tools Help

Reference

Compound

Sample

Mixture

Reaction

2001 mar gar 0

1-chlorobutane

Sample 1 (cm,99.5m%,nc;mv;99.9m%,glc)

^1: P (L), Set 1, B Method:CCELL dPV=.01 dT=.01

ethanol

Sample 1 (cm,99.8m%,nc;mv;99.9m%,glc)

^1: P (L), Set 1, B Method:CCELL dPV=.01 dT=.01

1-hexanol

Sample 1 (cm,99.0m%,nc;mv;99.4m%,glc)

^1: P (L), Set 1, B Method:CCELL dPV=.01 dT=.01

ethanol + 1-chlorobutane

^2: vle, P (Set 1), B Method:CCELL dP=.01 dT=.01 dX1=0.0003

^1: PA (Set 1), B Method:PX

1-chlorobutane + 1-hexanol

^2: vle, P (Set 1), B Method:CCELL dP=.01 dT=.01 dX2=0.0003

NOTE: The new data set appears in the tree under the appropriate *mixture*.

NOTE: DOUBLE CLICKING on the *data set* allows editing of all entered information.

END

**Continue with other compounds,
samples, properties, reactions, etc...**

or save your file and exit the program.