


METADATA AND NUMERICAL DATA CAPTURE:

Speed of Sound

(2 – Components)

Guided Data **Capture (GDC)**



This tutorial describes
METADATA AND NUMERICAL DATA CAPTURE:
for **2-components**
Speed of Sound
with the Guided Data Capture (GDC) software.

NOTE:

The tutorials proceed sequentially to ease the descriptions. **It is not necessary to enter *all* compounds before entering *all* samples, etc.**

Compounds, samples, properties, etc., can be added or modified at any time.

However, the hierarchy must be maintained (i.e., a property cannot be entered, if there is no associated sample or compound.)

The experimental data used in this example is from:

312

J. Chem. Eng. Data 2001, 46, 312–316

Temperature Dependence of Densities and Speeds of Sound of Nitromethane + Butanol Isomers in the Range (288.15–308.15) K

Jacobo Troncoso, Clara A. Tovar, Claudio A. Cerdeiriña, Enrique Carballo, and Luís Romani*

Departamento de Física Aplicada, Universidad de Vigo, Facultad de Ciencias del Campus de Ourense, 32004 Ourense, Spain

Densities and speeds of sound of the systems nitromethane + 1-butanol, nitromethane + 2-methyl-propan-1-ol, and nitromethane + 2-butanol were measured in the temperature range (288.15–308.15) K. Excess molar volumes are discussed and compared with those of other systems polar fluid + alcohol founded in the literature.

**Speeds of Sound for the binary system
2,2,2-trifluoroethanol + quinoline
at $p = 101.3$ kPa and various temperatures**

Table 3. Selected Data of Densities ρ and Speeds of Sound u for the Studied Mixtures

x	$\rho/\text{g}\cdot\text{cm}^{-3}$					$u/\text{m}\cdot\text{s}^{-1}$				
	288.15 K	293.15 K	298.15 K	303.15 K	308.15 K	288.15 K	293.15 K	298.15 K	303.15 K	308.15 K
	x Nitromethane + (1 - x) 1-Butanol									
0.056 70	0.82000	0.816 09	0.812 08	0.808 05		1255.86	1238.64	1221.53	1205.03	
0.103 78	0.82925	0.825 17	0.821 04	0.816 88		1255.12	1238.01	1220.76	1203.96	
0.197 35	0.84872	0.844 44	0.840 02	0.835 66		1254.38	1236.86	1218.87	1202.62	
0.290 57	0.87012	0.865 62	0.860 99	0.856 34		1253.87	1236.36	1219.03	1201.74	
0.346 86	0.88404	0.879 43	0.874 78	0.870 13		1254.30	1237.41	1218.59	1202.14	
0.425 91	0.90500	0.900 05	0.895 04	0.890 00		1253.82	1236.96	1219.93	1203.02	
0.482 52	0.92115	0.916 01	0.910 81	0.905 62		1254.03	1237.24	1221.17	1204.63	
0.524 56	0.93390	0.928 61	0.923 30	0.917 98		1255.68	1238.16	1222.92	1206.35	
0.559 49	0.94502	0.939 63	0.934 24	0.928 78		1256.76	1240.56	1224.36	1208.18	
0.587 64	0.95421	0.948 73	0.943 20	0.937 71		1257.09	1241.01	1225.21	1209.97	
0.689 09	0.98996	0.984 27	0.978 52	0.972 77		1267.10	1251.14	1235.51	1219.22	
0.750 78	1.01394	1.008 07	1.002 18	0.996 29		1277.35	1260.39	1244.39	1227.41	
0.852 46	1.05815	1.051 89	1.045 59	1.039 26		1300.15	1282.17	1264.12	1246.13	
0.946 35	1.10580	1.099 35	1.092 86	1.086 02		1326.53	1307.10	1287.76	1268.97	

**This data set is
considered here.**

Experimental Method Info:

Densities ρ and speeds of sound u of the pure components and of the binary mixtures were measured using a DSA-48 densimeter and sound analyzer. This apparatus, con-

The temperature was controlled through a solid-state thermostat that uses the Peltier effect, and the precision in the temperature was ± 0.005 K. The precision is estimated to be about $\pm 2 \times 10^{-5}$ for mole fraction x , $\pm 1 \times 10^{-5} \text{ g}\cdot\text{cm}^{-3}$ for ρ , and $\pm 0.02 \text{ m}\cdot\text{s}^{-1}$ for u .

Guided Data Capture - Thermophysical and Thermochemical Data

File Edit Tools Help

Reference Compound Sample Mixture Reaction **Property** Data Tables

2001 tro tov 0

- 1-butanol
 - Sample 1 (cm,99m%,nc;db,mv;99m%,glc)
- nitromethane
 - Sample 1 (cm,99.8m%,nc;db,mv;99.8m%,glc)
- 1-butanol + nitromethane**

2. CLICK *Property*

1. SELECT the *mixture* for which the data are to be captured.

NOTE: The **bibliographic information, compound identities, sample descriptions, and mixture** were entered previously. (There are separate tutorials, which describe capture of this information, if needed.)

Property and experimental method for 1-butanol + nitromethane

Help

Property group: Refraction; Surface tension; and Speed of sound

Property: Speed of sound

Units: m/s

Method of measurement:

Experimental purpose:

1. SELECT the **Property Group**:
*Refraction; Surface tension; and
Speed of sound* from the menu.

2. SELECT the **Property**:
Speed of sound.

3. SELECT the **Units**: m/s, here.
SELECT *ALL OTHER UNITS* if
another multiplier is needed.

OK

Cancel

1. SELECT Method of Measurement from the list provided. **NOTE:** *Other* can be a valid selection and should include a brief description in the **Comment** field.

Units: m/s

Method of measurement: Sing-around technique in a fixed-path interferometer

Experimental purpose: Principal objective of the work

2. SELECT the Experimental Purpose from the list provided.

Comment (optional) DSA-48.

3. CLICK OK

OK

Cancel

SELECTION of # of Phases in Equilibrium and # of Constraints

Speed of sound (m/s) as function of 2 variable(s)

Mixture: 1-butanol + nitromethane

Phases in equilibrium: 1

Constraints: 1

Independent variables: 2

Phase of the Property Value(s)

SELECT the # of **Phases in equilibrium**. There is **1** phases; *liquid*.

NOTE: For “saturation conditions, this value would be 2; liquid & gas)

SELECT the # of **Constraints**. There is **1** constraint in the example; *pressure = 101.3 kPa*.

Speed of sound (m/s) as function of 2 variable(s)

Mixture: 1-butanol + nitromethane

Phases in equilibrium: 1 Constraints: 1 Independent variables: 2 Property set #: 1

Sample #: 1 Sample #: 1

Phase of the Property Value(s):

Position of the Property Value(s): m/s

Definition of Measurement Results (Absolute vs Relative):

Data presentation: Experimental values

Comments (Optional): DSA-48

Property and method Numerical Data Cancel

Multiple *samples* for a given component can be accommodated, but this is rarely needed.

Speed of sound (m/s) as function of 2 variable(s)

Mixture: 1-butanol + nitromethane

Phases in equilibrium: 1 Constraints: 1 Independent variables: 2 Property set # 1 Sample # 1 Sample # 1

Phase of the Property Value(s) Liquid

Constraint 1 (Fixed value of) of Liquid

Independent variable 1 of Liquid

Independent variable 2 of Liquid

Definition of Measurement Results (Absolute Relative)

NOTE: Constraint and Independent Variable field(s) appear automatically based on the Gibbs Phase Rule.

1) SELECT *Liquid* from the list provided for the **Phase of the Property Value**

Specification of constraints, constraint values, and constraint units

1. SELECT the **Constraint(s)** (p here) and the **Independent Variable(s)** (T and x_1 , here) from the lists provided.

Mixture: 1-butanol + nitromethane

Phases in equilibrium: 1 Constraints: 1 Independent variables: 2 Property set # 1 Sample # 1 Sample # 1

Phase of the Property Value(s) Liquid Precision of the Property Value(s) 0.02 m/s %

Constraint 1 (Fixed value of)	Phase	Value	Units	Uncertainty
Pressure	Liquid	101.3	kPa	
Independent variable 1	Liquid		K	0.005
Independent variable 2	Liquid		Dimensionless	0.00001

Definition of Measurement Results (Absolute vs Relative)

Data pres: Experiment

Comm

Cancel

2. TYPE the Constraint **Value** (*if required*) and SELECT **Units** for the Variable(s) and Constraint(s). Include **Uncertainties**, if known.

Measurement definition and Data presentation

Speed of sound (m/s) as function of 2 variable(s)

Mixture: 1-butanol + nitromethane

Phases in equilibrium: 1 Constraints: 1 Independent variable

Phase of the Property Value(s) Liquid

Constraint 1 (Fixed value of) Pressure of Liquid

Independent variable 1 Temperature of Liquid

Independent variable 2 Mole fraction of nitromethane of Liquid

Units: Dimensionless Uncertainty: 0.00001 %

Definition of Measurement Results (Absolute vs Relative)
Direct value

Data presentation
Experimental values

Comments (Optional): DSA-48

Property and method Numerical Data Cancel

1. SELECT *Direct Value* (as compared with *Relative Value*) from the list defining the **Measurement Results**

2. SELECT the appropriate **Data presentation** method. *Experimental values* here.

3. CLICK *Numerical Data*

Speed of sound (m/s) as function of 2 variable(s)

File Edit Action Help

	Var 1	Var 2	Property
1			

TYPE, or much preferably, PASTE the variable and property values into the table. See next page...

Table 3. Selected Data of Densities ρ and Speeds of Sound u for the Studied Mixtures

x	$\rho/\text{g}\cdot\text{cm}^{-3}$					$u/\text{m}\cdot\text{s}^{-1}$				
	288.15 K	293.15 K	298.15 K	303.15 K	308.15 K	288.15 K	293.15 K	298.15 K	303.15 K	308.15 K
	x Nitromethane + (1 - x) 1-Butanol									
0.056 70		0.82000	0.816 00	0.812 08	0.808 05	1255.86	1238.64	1221.53	1205.03	
0.103 78		0.82925	0.825 17	0.821 04	0.816 88	1255.12	1238.01	1220.76	1203.96	
0.197 35		0.84872	0.844 44	0.840 02	0.835 66	1254.38	1236.86	1218.87	1202.62	
0.290 57		0.87012	0.865 62	0.860 99	0.856 34	1253.87	1236.36	1219.03	1201.74	
0.346 86		0.88404	0.879 43	0.874 78	0.870 10	1254.30	1237.41	1218.59	1202.14	
0.425 91		0.90500	0.900 05	0.895 04	0.890 00	1253.82	1236.96	1219.93	1203.02	
0.482 52		0.92115	0.916 01	0.910 81	0.905 62	1254.03	1237.24	1221.17	1204.63	
0.524 56		0.93390	0.928 61	0.923 30	0.917 98	1255.68	1238.16	1222.92	1206.35	
0.559 49		0.94502	0.939 63	0.934 24	0.928 78	1256.76	1240.56	1224.36	1208.18	
0.587 64		0.95421	0.948 73	0.943 20	0.937 71	1257.09	1241.01	1225.21	1209.97	
0.689 09		0.98996	0.984 27	0.978 52	0.972 77	1267.10	1251.14	1235.51	1219.22	
0.750 78		1.01394	1.008 07	1.002 18	0.996 29	1277.35	1260.39	1244.39	1227.41	
0.852 46		1.05815	1.051 89	1.045 59	1.039 26	1300.15	1282.17	1264.12	1246.13	
0.946 35		1.10580	1.099 35	1.092 86	1.086 02	1326.53	1307.10	1287.76	1268.97	

Speed of sound (m/s) as function of 2 variable(s)

File Edit Action Help

	Var 1	Var 2	Property
1	293.15	0.05670	1255.86
2	293.15	0.10378	1255.12
3	293.15	0.19735	1254.38
4	293.15	0.29057	1253.87
5	293.15	0.34686	1254.30
6	293.15	0.42591	1253.82
7	293.15	0.48252	1254.03
8	293.15	0.52456	1255.68
9	293.15	0.55949	1256.76
10	293.15	0.58764	1257.09
11	293.15	0.68909	1267.10
12	293.15	0.75078	1277.35
13	293.15	0.85246	1300.15
14	293.15	0.94635	1307.10
15	298.15	0.05670	1255.86
16	298.15	0.10378	1255.12
17	298.15	0.19735	1254.38
18	298.15	0.29057	1253.87
19	298.15	0.34686	1254.30
20	298.15	0.42591	1253.82
21	298.15	0.48252	1254.03
22	298.15	0.52456	1255.68
23	298.15	0.55949	1256.76
24	298.15	0.58764	1257.09
25	298.15	0.68909	1267.10

Table 3. Selected Data of Densities ρ and Speeds of Sound u for Studied Mixtures

x	$\rho/\text{g}\cdot\text{cm}^{-3}$					$u/\text{m}\cdot\text{s}^{-1}$			
	288.15 K	293.15 K	298.15 K	303.15 K	308.15 K	293.15 K	298.15 K	303.15 K	308.15 K
	x Nitromethane + (1 - x) 1-Butanol								
0.056 70		0.82000	0.816 09	0.812 08	0.808 05	1255.86	1238.64	1221.53	1205.03
0.103 78		0.82925	0.825 17	0.821 04	0.816 88	1255.12	1238.01	1220.76	1203.96
0.197 35		0.84872	0.844 44	0.840 02	0.835 66	1254.38	1236.86	1218.87	1202.62
0.290 57		0.87012	0.865 62	0.860 99	0.856 34	1253.87	1236.36	1219.03	1201.74
0.346 86		0.88404	0.879 43	0.874 78	0.870 13	1254.30	1237.41	1218.59	1202.14
0.425 91		0.90500	0.900 05	0.895 04	0.890 00	1253.82	1236.96	1219.93	1203.02
0.482 52		0.92115	0.916 01	0.910 81	0.905 62	1254.03	1237.24	1221.17	1204.63
0.524 56		0.93390	0.928 61	0.923 30	0.917 98	1255.68	1238.16	1222.92	1206.35
0.559 49		0.94502	0.939 63	0.934 24	0.928 78	1256.76	1240.56	1224.36	1208.18
0.587 64		0.95421	0.948 73	0.943 20	0.937 71	1257.09	1241.01	1225.21	1209.97
0.689 09		0.98996	0.984 27	0.978 52	0.972 77	1267.10	1251.14	1235.51	1219.22
0.750 78		1.01394	1.008 07	1.002 18	0.996 29	1277.35	1260.39	1244.39	1227.41
0.852 46		1.05815	1.051 89	1.045 59	1.039 26	1300.15	1282.17	1264.12	1246.13
0.946 35		1.10580	1.099 35	1.092 86	1.086 02	1326.53	1307.10	1287.76	1268.97

Clear the Table View plot Accept Cancel

NOTE: Simple CUT/PASTE procedures can be used within the table to convert the original table into the required number of columns. (This can also be done externally in spreadsheet software, e.g., EXCEL.)

Speed of sound (m/s) as function of 2 variable(s)

File Edit Action Help


	Var 1	Var 2	Property
1	293.15	0.05670	1255.86
2	293.15	0.10378	1255.12
3	293.15	0.19735	1254.38
4	293.15	0.29057	1253.87
5	293.15	0.34686	1254.30
6	293.15	0.42591	1253.82
7	293.15	0.48252	1254.03
8	293.15	0.52456	1255.68
9	293.15	0.55949	1256.76
10	293.15	0.58764	1257.09
11	293.15	0.68909	1267.10
12	293.15	0.75078	1277.35
13	293.15	0.85246	1300.15
14	293.15	0.94635	1326.53
15	298.15	0.05670	1238.64
16	298.15	0.10378	1238.01
17	298.15	0.19735	1236.86
18	298.15	0.29057	1236.36
19	298.15	0.34686	1237.41
20	298.15	0.42591	1236.96
21	298.15	0.48252	1237.24
22	298.15	0.52456	1238.16
23	298.15	0.55949	1240.56
24	298.15	0.58764	1241.01
25	298.15	0.68909	1251.14

Clear the Table View plot

1. CLICK *View plot* to see a graphical representation of the data.

2. Check for typographical errors, and CLICK *OK*, when done.

NOTE: The variable associated with the x-axis can be selected for best display.



Speed of sound (m/s) as function of 2 variable(s)

File Edit Action Help

	Var 1	Var 2	Property
1	293.15	0.05670	1255.86
2	293.15	0.10378	1255.12
3	293.15	0.19735	1254.38
4	293.15	0.29057	1253.87
5	293.15	0.34686	1254.30
6	293.15	0.42591	1253.82
7	293.15	0.48252	1254.03
8	293.15	0.52456	1255.68
9	293.15	0.55949	1256.76
10	293.15	0.58764	1257.09
11	293.15	0.68909	1267.10
12	293.15	0.75078	1277.35
13	293.15	0.85246	1300.15
14	293.15	0.94635	1326.53
15	298.15	0.05670	1238.64
16	298.15	0.10378	1238.01
17	298.15	0.19735	1236.86
18	298.15	0.29057	1236.36
19	298.15	0.34686	1237.41
20	298.15	0.42591	1236.96
21	298.15	0.48252	1237.24
22	298.15	0.52456	1238.16
23	298.15	0.55949	1240.56
24	298.15	0.58764	1241.01
25	298.15	0.68909	1251.14

CLICK *Accept*

Clear the Table View plot Accept Cancel

Guided Data Capture - Thermophysical and Thermochemical Data

File Edit Tools Help

Reference

Compound

Sample

- [-] 2001 tro tov 0
 - [-] 1-butanol
 - Sample 1 (cm,99m%,nc;db,mv;99m%,glc)
 - [-] nitromethane
 - Sample 1 (cm,99.8m%,nc;db,mv;9.8m%,dlc)
 - [-] 1-butanol + nitromethane

^2: RSS (Set 1), B Method:FPINT dRSS=0.02 dT=0.005 dX2=0.00001

NOTE: The new data set now appears in the tree under the appropriate *mixture*.

NOTE: DOUBLE CLICKING on the *data set* allows editing of all entered information.

END

**Continue with other compounds,
samples, properties, reactions, etc...**

or save your file and exit the program.