

Adapting the IMPROVE_A Protocol for Multiwavelength Organic and Elemental Carbon Measurements

Judith C. Chow¹, Xiaoliang Wang¹, Benjamin J.
Sumlin¹, Steven B. Gronstal¹, L.-W. Antony Chen^{1,2},
and John G. Watson¹

¹Desert Research Institute, Reno, NV

²University of Nevada, Las Vegas, NV

Presented at:

11th International Conference on Carbonaceous Particles in the Atmosphere
Berkeley, CA

August 11, 2015


Objectives

- Demonstrate equivalence between single- and multiwavelength (λ) systems for organic carbon (OC) and elemental carbon (EC).
- Introduce a calibration procedure for quantifying OC, EC, and brown carbon (BrC).
- Relate relative reflectance (R) and transmittance (T) values to different sources.

Motivation



- Single λ reflectance (R) and transmittance (T) have only been used for pyrolysis adjustment (i.e., reference to initial R or T). It can also be normalized to final R or T to approximate filter attenuation (ATN).
- The multiwavelength light source/detector combination yields different intensities within and among instruments.
- Absolute R and T (in %) can be used to calculate ATN on filters and/or applied to radiative transfer models (e.g., Beer's Law, Kubelka-Munk Theory, or Monte Carlo Ray Tracing, etc.).

Multiwavelength Thermal/Optical Analyzer reports both reflectance and transmittance at 405, 445, 532, 635, 780, 808, and 980 nm




DRI Model 2015
configuration

The multiwavelength optical configurations allow for absolute calibration and wavelength-dependent OC/EC/BrC splits




Similar thermograms are obtained for Models 2001 and 2015

Model 2001 (single wavelength)




Model 2015 (multiwavelength)



Model 2001




Model 2015




Lower R before adding O₂, implying better sensitivity in detecting OC charring

Equivalent OC and EC are obtained for single- and multi-wavelength systems (633 nm vs 635 nm)

Model 2001 (A) vs Model 2001 (B)



Model 2015 vs Model 2001



OCR and ECR are OC and EC by reflectance.

Filter transfer standards with variable deposits can be standardized against Spectralon* standards (diffusive reflective)



(Lambda 35 UV/VIS Spectrometer, Perkin Elmer, Waltham, MA; an Integrating-Sphere Spectrometer; measures R and T at 0 and 100%, 200-1100 nm)




* NIST Certified Labsphere Spectralon® Diffusive Reflectance Standards


R and T are lower for shorter wavelengths

(Vertical lines designate the seven wavelengths in Model 2015)



R and T in diesel samples are lower with increasing loadings.



R and T in smoldering samples show minor changes with loading and clustered at high wavelengths (less useful than the blacker standard).




Calibration curves using transfer standards show linear responses independent of sample type or loading



Transmittance spectral attenuation varies by sample type


$$ATN_T(\lambda) = -\ln\left(\frac{FT_{\lambda,i}}{FT_{\lambda,f}}\right)$$



- Spectral absorption averaged by sample type.
- Smoldering samples acquired in a wood stove connected to DRI dilution chamber.


Reflectance spectral attenuation also varies by sample type

$$ATN_R(\lambda) = -\ln\left(\frac{FR_{\lambda,i}}{FR_{\lambda,f}}\right)$$




- Reflectance usually has lower signal to noise ratios than transmittance.
- R and T can be combined for better quantification of light absorption as indicated by Petzold and Schonlinner (2004).

AAE* can be used to decouple BC and BrC



Ambient Fresno, CA



Smoldering Biomass Burning

Simplified two-component model:





- $\tau_{a,\lambda} = q_{BC} \times \lambda^{-\alpha_{BC}} + q_{BrC} \times \lambda^{-\alpha_{BrC}}$

Assuming $\alpha_{BC}=1$:

- $\tau_{a,\lambda} \times \lambda = q_{BC} + q_{BrC} \times \lambda^{-(\alpha_{BrC}-1)}$ (q_{BC} and q_{BrC} are fitting coefficients; α_{BC} and α_{BrC} are AAEs)

*AAE: Absorption Angstrom Exponent; ** $\tau_{a,\lambda}$ is absorption optical depth, $\tau_{a,\lambda} = ATN_{\lambda}$ if there is no filter effect

BC and BrC contributions to light attenuation (ATN_{405 nm}) vary by sample type



- Assuming only BC absorbs at 980 nm and an AAE_{BC} of 1 to extrapolate BC absorption to 405 nm.
- Samples sorted by BrC fraction (0 to 100%) in ATN_{405 nm}.

Potential future uses of calibrated multiwavelength R and T on thousands of samples

- Identifying light absorbing compounds.
- Separating artifact OC from aerosol OC.
- Ground-truthing remotely-sensed BrC.
- Improving radiation transfer estimates.
- Conducting source apportionment for BC and BrC.

More Information

(Chen et al., 2015; Chow et al., 2015)

Aerosol and Air Quality Research, x: 1–15, xxxx
Copyright © Taiwan Association for Aerosol Research
ISSN: 1680-8584 print / 2071-1409 online
doi: 10.4209/aaqr.2015.02.0106



Optical Calibration and Equivalence of a Multiwavelength Thermal/Optical Carbon Analyzer

Judith C. Chow^{1,2,3*}, Xiaoliang Wang^{1,3}, Benjamin J. Sumlin¹, Steven B. Gronstal¹, L.-W. Antony Chen^{1,4}, Dana L. Trimble¹, Steven D. Kohl¹, Sierra R. Mayorga¹, Gustavo Riggio¹, Patrick R. Hurbain¹, Megan Johnson¹, Ralf Zimmermann⁵, John G. Watson^{1,2,3}

¹Division of Atmospheric Sciences, Desert Research Institute, Reno, Nevada 89512, USA

²The State Key Laboratory of Loess and Quaternary Geology, Institute of Earth Environment, Chinese Academy of Sciences, Xi'an, Shaanxi, 710075, China

³Graduate Faculty, University of Nevada, Reno, Nevada 89503, USA

⁴Department of Environmental and Occupational Health, University of Nevada, Las Vegas, Nevada 89154, USA

⁵Joint Mass Spectrometry Centre, Chair of Analytical Chemistry, Institute of Chemistry, University of Rostock, Rostock, Germany

Atmos. Meas. Tech., 8, 451–461, 2015
www.atmos-meas-tech.net/8/451/2015/
doi:10.5194/amt-8-451-2015
© Author(s) 2015. CC Attribution 3.0 License



Atmospheric
Measurement
Techniques

Multi-wavelength optical measurement to enhance thermal/optical analysis for carbonaceous aerosol

L.-W. A. Chen^{1,2,3}, J. C. Chow^{1,3}, X. L. Wang³, J. A. Robles², B. J. Sumlin², D. H. Lowenthal², R. Zimmermann⁴, and J. G. Watson^{2,3}

¹Department of Environmental and Occupational Health, University of Nevada, Las Vegas, Nevada 89154, USA

²Division of Atmospheric Sciences, Desert Research Institute, Reno, Nevada 89512, USA

³Key Laboratory of Aerosol Science & Technology, SKL-LQG, Institute of Earth Environment, Chinese Academy of Sciences, Xi'an, China

⁴Joint Mass Spectrometry Centre, Chair of Analytical Chemistry, Institute of Chemistry, University of Rostock, Rostock, Germany

Chen, L.-W.A.; Chow, J.C.; Wang, X.L.; Robles, J.A.; Sumlin, B.J.; Lowenthal, D.H.; Watson, J.G. (2015). Multi-wavelength optical measurement to enhance thermal/optical analysis for carbonaceous aerosol. *Atmos. Meas. Tech.*, **8**:451-461. <http://www.atmos-meas-tech.net/8/451/2015/amt-8-451-2015.html>.

Chow, J.C.; Wang, X.L.; Sumlin, B.J.; Gronstal, S.B.; Chen, L.-W.A.; Trimble, D.L.; Kohl, S.D.; Mayorga, S.R.; Riggio, G.; Hurbain, P.R.; Johnson, M.; Zimmermann, R.; Watson, J.G. (2015). Optical calibration and equivalence of a multiwavelength thermal/optical carbon analyzer. *AAQR*, online. doi:10.4209/aaqr.2015.02.0106. http://aaqr.org/JustAcceptedManuscripts/AAQR-15-02-OA-0106_accepted.pdf.

Conclusions

- Reflectance (R) and Transmittance (T) can be traceable to primary standards and made consistent among wavelengths and instruments.
- The detailed absorption spectrum can be approximated by the seven wavelengths.
- Brown carbon (BrC) can be separated from black carbon (BC) by a two-component model.

Acknowledgements

- U.S. National Science Foundation
(CHE 1214163)
- National Park Service IMPROVE Carbon
Analysis Contract (C2350000894)

References

- Chen, L.-W.A.; Chow, J.C.; Wang, X.L.; Robles, J.A.; Sumlin, B.J.; Lowenthal, D.H.; Watson, J.G. (2015). Multi-wavelength optical measurement to enhance thermal/optical analysis for carbonaceous aerosol. *Atmos. Meas. Tech.* , **8**:451-461. <http://www.atmos-meas-tech.net/8/451/2015/amt-8-451-2015.html>.
- Chow, J.C.; Wang, X.L.; Sumlin, B.J.; Gronstal, S.B.; Chen, L.-W.A.; Trimble, D.L.; Kohl, S.D.; Mayorga, S.R.; Riggio, G.; Hurbain, P.R.; Johnson, M.; Zimmermann, R.; Watson, J.G. (2015). Optical calibration and equivalence of a multiwavelength thermal/optical carbon analyzer. *AAQR*, online. doi:10.4209/aaqr.2015.02.0106.
- Petzold, A.; Schonlinner, M. 2004. Multi-angle absorption photometry – A new method for the measurement of aerosol light absorption and atmospheric black carbon. *J. Aerosol Sci.*, **35**:421-441.