# Q4B Evaluation and Recommendation of Pharmacopoeial Texts for Use in the ICH Regions

## **Annex 7: Dissolution Test General Chapter**

This draft guidance, when finalized, will represent the Food and Drug Administration's (FDA's) current thinking on this topic. It does not create or confer any rights for or on any person and does not operate to bind FDA or the public. You can use an alternative approach if the approach satisfies the requirements of the applicable statutes and regulations. If you want to discuss an alternative approach, contact the FDA staff responsible for implementing this guidance. If you cannot identify the appropriate FDA staff, call the appropriate number listed on the title page of this guidance.

For questions regarding this draft document contact (CDER) Robert King 301-796-1242, or (CBER) Christopher Joneckis 301-827-0373.

2	Final for STEP 2 signoff - Annex 7 Dissolution Test [for Brussels November 2008] Corrected 12-18-08  Nov. 13, 2008
3 4 5	INTERNATIONAL CONFERENCE ON HARMONISATION OF TECHNICAL REQUIREMENTS FOR REGISTRATION OF PHARMACEUTICALS FOR HUMAN USE
6	
7	
8	DRAFT CONSENUS GUIDELINE
9	
10	
11	
12	
13	
	Transcript Dragon to the control of
14	EVALUATION AND RECOMMENDATION OF
15	PHARMACOPOEIAL TEXTS FOR USE IN THE ICH REGIONS
16	ON
17	DISSOLUTION TEST GENERAL CHAPTER
18	
19	Q4B ANNEX 7
20	Q 12 11 (1211 )
21	
22	
23	
24	Current Step 2 Version
25	dated 13 November 2008
26	
27	
28	
29	
30	
31	
32	
33	
34 35	At Step 2 of the ICH Process, a consensus draft text or guideline, agreed by the appropriate ICH
36	Expert Working Group, is transmitted by the ICH Steering Committee to the regulatory
37	authorities of the three ICH regions (the European Union, Japan and the USA) for internal and
38	external consultation, according to national or regional procedures.

	Final for STEP 2 signoff - Annex 7 Dissolution Test  [for Brussels November 2008] Corrected 12-18-08  Nov. 13, 2008
39	
40	
41	
42	
43	
44	
45	
46	
47	
48	
49	
50	Q4B Annex 7
51	Document History
52	
53	Current Step 2 version

Code	History	Date
Q4B Annex 7	Approval by the Steering Committee under <i>Step 2</i> and release for public consultation.	13 November 2008

Final for STEP 2 signoff - Ann	nex 7 Dissolution Test
[for Brussels November 2008]	Corrected 12-18-08

Nov. 13, 2008

ON DISSOLUTION TEST GENERAL CHAPTER	IONS
DISSOLUTION TEST GENERAL CHAPTER	
of Dissolution I Lat Oliverial Child I Lin	
62	
63 <b>Q4B Annex 7</b>	
64	
65 Draft ICH Consensus Guideline	
Released for Consultation on 13 November 2008, at <i>Step 2</i> of the ICH Page 100.	rocess
67	
68	
69 TABLE OF CONTENTS	
70	
71	
72 <b>1. INTRODUCTION</b>	1
73	
74 2. Q4B OUTCOME	1
75 76 2.1 Analytical Procedures	1
76 2.1 Analytical Procedures	1
78 2.2 Acceptance Criteria	1
79	
3. TIMING OF ANNEX IMPLEMENTATION	1
81	
4. CONSIDERATIONS FOR IMPLEMENTATION	2
83 84 4.1 General Consideration	9
85	
86 4.2 FDA Consideration	2
87	
88 4.3 EU Consideration.	2
89	
90 4.4 MHLW Consideration	3
91 92 <b>5. REFERENCES USED FOR THE Q4B EVALUATION</b>	વ

Final for STEP 2 signoff - Annex 7 Dissolution Test
[for Brussels November 2008] Corrected 12-18-08

Nov. 13, 2008

95

141

142 143

144

145

#### EVALUATION AND RECOMMENDATION OF 96 PHARMACOPOEIAL TEXTS FOR USE IN THE ICH REGIONS 97 98 ON DISSOLUTION TEST GENERAL CHAPTER 99 100 **Q4B** Annex 7 184 1. INTRODUCTION 106 This annex is the result of the Q4B process for Dissolution Test. 107 108 The proposed texts were submitted by the Pharmacopoeial Discussion Group (PDG). 109 110 2. **Q4B OUTCOME** 111 2.1 **Analytical Procedures** 114 The ICH Steering Committee, based on the evaluation by the Q4B Expert Working 115 Group (EWG), recommends that the official pharmacopoeial texts, Ph.Eur. 2.9.3. 116 Dissolution, JP 6.10 Dissolution Test, and USP <711> Dissolution, can be used as 117 interchangeable in the ICH regions subject to the following conditions: 118 119 **2.1.1** The Dissolution Test is not considered to be interchangeable in the ICH regions 120 when enzymes are used in the media. 121 122 **2.1.2** The dissolution apparatus should be appropriately calibrated to ensure 123 compliance with regional good manufacturing practice (GMP) requirements. 124 125 **2.1.3** Except for Apparatus 1 and 2, apparatus numbers are not consistent in the three 126 pharmacopoeias. Accordingly, other apparatus should be referred to in the 127 dossier by an unambiguous descriptive title or compendial reference. 128 129 The Dissolution Test is not considered to be interchangeable in the ICH regions 130 for dosage forms referred to in the regional compendia as delayed-release, 131 gastro-resistant, or enteric-coated. 132 133 **2.1.5** Validation studies should be conducted to demonstrate that the test results are 134 not adversely affected if the thermometer is to remain in the dissolution vessel 135 per regional good manufacturing practice (GMP). 136 137 **2.1.6** The Dissolution Test is not considered to be interchangeable in the ICH regions 138 for JP Interpretation 2. 139 140 **2.1.7** The Dissolution Test is not considered to be interchangeable in the ICH regions

use and type of sinkers should be specified and justified in the application dossier.

**2.1.8** Product-specific parameters such as media, stirring rate, sampling time, and the

for use of *large* vessels (greater than 1 liter).

Dissolution Test General Chapter

 **2.1.9** When using the small cell tablet holder with the flow-through cell apparatus, only the dimensions described in the PDG harmonised text Figure 5 are considered interchangeable.

### 2.2 Acceptance Criteria

Acceptance criteria should be specified in the application dossier.

#### 3. TIMING OF ANNEX IMPLEMENTATION

When this annex is implemented (incorporated into the regulatory process at ICH Step 5) in a region, it can be used in that region. Timing might differ for each region.

#### 4. CONSIDERATIONS FOR IMPLEMENTATION

#### 4.1 General Consideration

When sponsors or manufacturers change their existing methods to the implemented Q4B-evaluated pharmacopoeial texts that are referenced in Section 2.1 of this annex, any change notification, variation, and/or prior approval procedures should be handled in accordance with established regional regulatory mechanisms pertaining to compendial changes.

#### 4.2 FDA Consideration

Based on the recommendation above, and with reference to the conditions set forth in this annex, the pharmacopoeial texts referenced in Section 2.1 of this annex can be considered interchangeable. However, FDA might request that a company demonstrate that the chosen method is acceptable and suitable for a specific material or product, irrespective of the origin of the method.

An appropriately rigorous mechanical calibration method (such as ASTM International's ASTM E2503-07, Standard Practice for Qualification of Basket and Paddle Dissolution Apparatus, or the procedures for Mechanical Qualification of Dissolution Apparatus 1 and 2, DPA-LOP.002, on the FDA Web site), when properly executed, will satisfy the current good manufacturing practice (CGMP) requirement for dissolution apparatus calibration under § 211.160(b)(4) of Title 21 of the Code of Federal Regulations.

#### 4.3 EU Consideration

For the European Union, the monographs of the Ph. Eur. have mandatory applicability. Regulatory authorities can accept the reference in a marketing authorisation application, renewal or variation application citing the use of the corresponding text from another pharmacopoeia as referenced in Section 2.1, in accordance with the conditions set out in this annex, as fulfilling the requirements for compliance with the Ph. Eur. Chapter 2.9.3. on the basis of the declaration of interchangeability made above.

EU considers that it could accept the approach to the dissolution test for delayedrelease products as published in the USP as meeting the criteria of the Ph. Eur. The

237

238

203 28<del>4</del> 285 Dissolution Test General Chapter 206 207 validation studies referred to in Section 2.1.5 of this annex would normally be 208 submitted in the marketing authorisation dossier. 209 210 4.4 **MHLW Consideration** 213 The pharmacopoeial texts referenced in Section 2.1 of this annex can be used as 214 interchangeable in accordance with the conditions set out in this annex. Details of 215 implementation requirements will be provided in the notification by MHLW when this 216 annex is implemented. 217 218 MHLW considers that it could accept the approach to the dissolution test for 219 reciprocating cylinder apparatus as published in Ph. Eur. and USP, if the validation 220 studies have been submitted in the marketing authorization dossier. 221 222 223 REFERENCES USED FOR THE Q4B EVALUATION 5. 224 5.1 The PDG Stage 5B sign-off document (Rev. 1): Japanese Pharmacopoeial Forum, 225 226 Volume 14, number 4 (December 2005). 5.2 The pharmacopoeial references for Dissolution Test for this annex are: 229 **5.2.1** *European Pharmacopoeia* (Ph. Eur.): 230 6th Edition (official on January 2008) Dissolution Test (reference 01/2008: 231 232 20903). 233 **5.2.2** *Japanese Pharmacopoeia* (JP): 234 6.10 Dissolution Test as it appears in Supplement I to the JP Fifteenth Edition 235 236

(September 28, 2007, The Ministerial Notification No. 316).

<711> Dissolution Test USP 28, 2<sup>nd</sup> Supplement, official August 1, 2005.

**5.2.3** *United States Pharmacopeia* (USP):