SHOR COPY

The Human Airpuck

Ken Altshuler


978-897-4996

We and High School, Wayland, MA 01778

MAYNARO)

Have you ever wanted to have your students experience what life would be like in a frictionless frame of reference? Would Newton's laws of motion and gravity be more apparent if friction weren't around to hide their effects?

During the summer of 1987, I was fortunate to take part in a course for physics teachers at the Virginia Military In: tute. During this course, Barbara Saur (from Kentwood High School in Kent, Washington) described how to make an aircart upon which a person could sit. It proved to be a very simple device to build and an extremely useful demonstration tool for several different topics.



Fi. 1. Design of the airpuck.



Fig. 2. Underside of the airpuck.

_ SCALE=1"=1"

Materials required to build this cart are:

a wooden board (0.5- to 1.0-in thick) that can be shaped into a disk 3.0 ft in diameter

a plastic coffee-can lid

a (1.0- to 1.5-in long) bolt, washer, and nut

a 12.5-sq-ft sheet of plastic at least 4-mil thick

a staple gun with large staples

a bicycle tire

hammer and nails to be used to attach the tire to the wooden wheel, and

an air supply

Cut the board into a 3-ft circle with a hole in its center that is just large enough for the bolt to slip through. The rim of the wheel should be sanded into rounded edges. In addition, a hole should be drilled into the wheel about 10 in from the center. This hole should be slightly larger than the diameter of the nozzle end of the air-supply hose.




Fig. 3. A student moves across the room on the airpuck.




Fig. 4. Resolving vectors. (Photos courtesy of Joseph Damplo.)

Now stretch the sheet of plastic tightly across the bottom of the wooden disk and staple it to the top of the disk along the wheel's rim. Put a hole in the center of the coffee-can lid big enough for the bolt to slide through. This lid is going to function as an extended washer for the bolt as it goes through the plastic. Now put the bolt through the washer, the lid, plastic sheet, and wooden disk. Screw the nut tightly to the bolt on the top of the aircart, as shown in Figs. 1 and 2.

Cut six holes in the plastic at radial distance of about 15 in from the center of the wheel and evenly spaced around the circumference of the wheel. Each of these holes should have a diameter of about 5- to 6-in. Be sure that none of these holes is in direct alignment with the air supply hole.

Cut the bicycle tire and wrap it around the rim of the wooden wheel, nailing it to the wheel on the top side. This will act as a bumper for the wheel so that the plastic along the edge of the wheel won't rip when the aircart hits other objects.

Finally, put the air pump on the cart with the hose nozzle stuck into the hole in the wheel. If the nozzle has a loose fit, use duct tape to hold it in place. You can also place a chair on the disk to make it easier for the operator to maintain his or her balance.

This aircart has proven to be a most useful tool for a large series of demonstrations. For instance, we have used

it to demonstrate that a body in motion remains in motion when not acted upon by an outside force. To show this, a student is told to maintain her balance while she rides the cart across the floor, making sure the tire does not come into contact with the floor. It quickly becomes obvious that a rider has little control over the cart's motion (Fig. 3).

When we talk about vector addition using forces applied to an object along nonparallel lines, we consider the problem generated by having two students pull in known directions with known forces on the aircart. Then we mathematically solve for what equilibrant is required to balance these forces. Once we have mathematically solved this problem, the students test their theory. To do this, force scales are tied to the ends of a rope. These scales are pulled on by two class members. The student on the cart holds the center of this rope against a giant protractor. He can thus assure that the students pull on him at the angles described in the original problem. A third student now pulls on a third scale that has been tied to the back of the chair on the cart. When the cart is returned to a state of equilibrium, we see if the third student is applying a force equal in magnitude and direction to that predicted during our calculation (Fig. 4).

Finally, when we study Newton's third law, students on two carts play catch with a medicine ball (Fig. 5). The class is able to see the student's response both when throwing and catching the ball. In this way, they can see how the ball



Fig. 5. Students demonstrating Newton's third law.

"pushes" one student backward when the other student throws it forward. One problem with which you may have to deal is that during the throwing and catching students may lean so that the tire comes in contact with the floor. Be careful to remind them to maintain their balance throughout this demonstration. These are only a few ideas for the use of the aircart. It is a popular device that permits us to demonstrate topics usually relegated to airtrack experiments. This is a much less expensive piece of equipment than the airtrack. In addition, owing to the personalized nature of these demonstrations, the students feel more a part of the process. This is especially true for the student riding on the cart. For this reason, we try to repeat each experiment with different students riding the cart so that many students get a firsthand experience with the law of physics being studied.

Significantly, and in contrast to airtracks, this airpuck device can easily be built in the typical high school industrial arts room by students. In this way, the students also gain experience in designing and/or building equipment to be used in class.

Oh, one last point. We have found that our carts can support weights of up to 220 lbs easily and, during the demonstration at VMI, people much heavier

than that were able to coast around the room on a cart. So, if you want to take a firsthand role in the demonstration, go ahead and have fun. •

AAPT Apparatus Competition 1998 August 3

Crispin O. Butler and Thomas D. Bishop Physics Department

Middlebury College Middlebury VT 05753

Personal hovercraft with accessories

Hovercraft are vehicles that are supported by a cushion of compressed air or other gasses. Typically a fan or turbine is used to force air under the vehicle and some sort of

INTRODUCTION

skirt around the bottom perimeter of the vehicle serves to hold the air there. Since hovercraft are lifted by a cushion of air, they experience very little friction with the surface under them and can move easily over water as well as land. Commercial hovercraft range

in size from the huge craft that are used to carry both passengers and cars across the

English Channel (25 British pounds for a person, one way, 35 minutes Dover to Calais) to single person inflatable hovercraft built for pleasure use. Through the internet many kits are available for building your own hovercraft and there are hovercraft clubs in many parts of the world devoted to designing hovercraft, building hovercraft, or racing hovercraft. A major deterrent to the general use of hovercraft in our transportation system is the difficulty in controlling their motion. However, it is precisely this difficulty that makes a hovercraft such a useful tool in a physics classroom. Hovercraft have been used by teachers of physics for many years, though there are

few references to their construction or use in the literature. An article in the November

1989 issue of The Physics Teacher by Ken Altshuler, "The Human Airpuck", describes how to make a basic person sized hovercraft out of easily obtainable materials and details several classroom experiments as done by his students2. In this article Altshuler credits Barbara Sauer from Kentwood High School in Kent, Washington with providing him a description of how to make such a device during a course for physics teachers at the Virginia Military Institute in 1987. This hovercraft was particularly suited for use in the classroom as it was an indoor hovercraft. Three feet in diameter it could carry a two hundred fifty pound person and was designed for relatively quiet hovering over smooth

floors. An airtrack blower provided the pressurized air for lift, and the skirt around the bottom perimeter took the form of a sheet-plastic plenum with large holes near the center of the bottom. Many such hovercraft must have been built and used by physics teachers since

that time, but we find very little reference to modifications or uses of Altshuler's human airpuck.

THE PERSONAL HOVERCRAFT The personal hovercraft described here (see Fig. 1.), is a direct descendent of Altshuler's human airpuck. The approximate dimensions and the circular form as well as

many details of the plastic plenum and bumper were taken directly from his article. Our prototype (see Fig. 2.) consisted of just these parts. And, like his airpuck, our prototype was powered by an airtrack blower. Our modifications of his basic design include the addition of two braces to reduce the flexing of the plywood base; weighting the plastic plenum to minimize sounds produced by the hovercraft; the addition of a seat for the

passenger that serves to keep the center of gravity of the rider-hovercraft low; the selection

undergraduate physics students. The youngest students took turns riding and enjoyed pushing the hovercraft around whether a person was riding or not. They tended to wonder how far and how fast it would travel and were able to measure distance traveled and travel time in order to calculate the speed. After a brief introduction the undergraduate students working on a basketball court in groups of five to seven did two preliminary experiments in a semi-quantitative fashion. First, they applied two perpendicular forces to someone on the hovercraft and monitored the resulting motion. Then they made the hovercraft and rider move in a circle with approximately constant speed; measuring the period, the centripetal

force (with scales) and the radius of the circle. With these data they estimated the total mass of the hovercraft and rider and then measured the mass of the rider and the hovercraft on a balance. Finally, the students designed a quantitative lab on the spot and took data for

analysis later. They chose to do experiments such as measuring the frictional forces between the hovercraft and the floor as a function of weight or as a function of speed or measuring the stretch of a spring and the resulting speed of the hovercraft after being launched by that spring. All students at all levels were enthusiastic and involved

themselves in the experiments. The faculty member working with the younger students said, "The hovercraft stold the show." The professor working with the undergraduates

The scale of the hovercraft is right for an introduction to the physics of motion. By

allowing students to experience the motions, the forces and the torques over time periods of

studied, floor surfaces can be mapped, hovercraft launchers can be designed and tested, or students can study the operation of the hovercraft itself. In a word the personal hovercraft

experiments. Use of the hovercraft in a laboratory setting stimulates creativity.. and it is

can generate quantitative data for a wide range of kinematic and dynamic physics

The personal hovercraft has been used with primary level students as well as our

and installation of batteries, blower and controller to make the hovercraft self contained and able to operate continuously for periods up to one hour; and the addition of a marking system that records the position of two hovercraft at regular, synchronized time intervals. The hovercraft were used in a cooperative lab for our introductory physics course during the spring term this year. Based on that experience and with the help and suggestions of many students we have developed a variety of useful accessories for the hovercraft. These include a vinyl floor marked with a grid of squares 10[cm] on each side, a simple push stick, and a fishing reel winch that produces constant forces. Fig. 1 shows our hovercraft

several seconds the connections between concepts are immediate and concrete. With the built in timer marking two hovercraft positions twice every second it is straight forward to analyze those marks in terms of two dimensional position vectors, velocity vectors, acceleration vectors, momentum vectors and force vectors. Speeds, potential energies and kinetic energies can also readily be obtained as functions of time or position from these marks on the floor. Students can organize themselves in groups to study motions modified by one or more constant forces, centripetal forces, or spring forces. Collisions can be

Appendix A: Instructions For Use

and several of these accessories.

termed her lab, "The greatest!"

CONCLUSION

fun.

Find the switch marked BLOWER under the seat and turn it on. When the hovercraft lifts off sit down on the seat and get your feet on board. Look down at the bubble level and adjust your position until the bubble is in the center. Now look around you. The question that usually comes to mind is, "How do I stop this thing?". When riding a hovercraft the natural state of motion is easily seen to be a constant vector velocity

both your vector velocity and your spin rate. You may grab a push stick made of foam so that you can have some control over forces that you apply to yourself. You will soon find out that in order to push yourself you need something to push against. Forces come in pairs that are oppositely directed; forces that affect your motion also affect another object. In addition to Newton's three laws of motion (as mentioned above) it will become obvious that to avoid spinning or tilting the hovercraft you must apply the forces in line with a particular point that is located a few decimeters above the center of the hovercraft. This point is the center of mass of the combination of the hovercraft and your body; it is easier to find this point experimentally, than it is to define it. When you have learned how to control the motion and attitude of the hovercraft you are ready to do quantitative experiments with one or more hovercraft. Place the hovercraft on the vinyl coordinate flooring and decide what experiment you want to do. Practice the choreography of the experiment a few times. When you are ready to do the experiment

turn on the switch under the seat labeled MARKER, this energizes the motor that moves the

marker. Then locate the hand held clock transmitter, it should be near the CLOCK TRANSMITTER sign. The marker will hit the floor twice a second for as long as the transmitter trigger is held in the squeezed position. The marker stops when either switch is

Construction of Personal Hovercraft

This is attached with screws, but no glue, so that it is easy to remove.

turned off. Enjoy!

Appendix B:

with a constant rate of rotation. A sloping floor will definitely change your velocity vector without changing your rate of rotation, but a push from a bystander will generally change

The parts shown in Fig. 2 are put together as per the instructions in Ref. 2. The parts labeled in Fig 3. represent our contribution to the hovercraft design and we treat them individually: Braces: Cut from a two by six spruce stud and attached to the base with glue and screws. Seat: Made from scraps of the 4x8 sheet of half inch plywood that was used for the bases.

Hardware \$2.97) was applied to the upper side of the 4 mil polyethylene sheet before it was attached to the base.

Noise Suppression: A spider web pattern Liquid Nails Clear Seal adhesive (Aubuchon

integrated circuit (Digikey \$0.56) is used to provide a clock signal to a 49.86[MHz] FM

Batteries: Hawker Energy, Pure Lead Tin Monoblock Type 0859-0012, 6[V] 8[A hour].

(Allied Electronics #221-0012, \$20.74) Four batteries for each hovercraft. BAT NETWORK G 682

Blower: AMETEK #11645-00 Three stage 24[V] 10[A] three phase DC blower. Output 20[ft3/min] at a pressure of 20[inch (water)]. (C&H Sales, Stock #DCB9601 \$12.50) Attached by three bolts to a brace with the output pressed into the hole through the base of

Blower Controller: W.W.Grainger Blower controller with heat sink #3HV53 \$193.59.

Screwed to the base. A hole 1[cm] by 5[cm] was made at the base of the brace when the blower was attached and the wiring harmons for the brace when t Screwed to the base. A hole 1[cm] by 5[cm] was made at the base of the brace where the blower was attached and the wiring harness from the controller to the blower was passed through this hole. through this hole. Radio Controlled Timer/Marker: A National Semiconductor LM555CN-ND timer

radio transmitter. The clock signal is picked up by an FM receiver mounted on the hovercraft that controls a motor that pushes a marker against the floor. The transmitter, receiver and motor assembly are a radio controlled car, somewhat modified. New Bright Radio Control 1:24 Scale Sports Car #346 (Ben Franklin \$12.99). Markers are Crayola Washable Markers #7809, Package of 8. (Ben Franklin \$3.99).

Vinyl Flooring: A square three yards on a side. (Countryside Carpet and Paint \$90.00).

Any plain, light colored, smooth flooring that will lie flat will do.


Push Stick: Life-Like Products, Wacky Noodle (Ben Franklin \$2.99) a 7[cm] diameter rod of high density foam 1.6[m] long.

Constant Magnitude Force Winch: A Diawa fishing reel (Outdoor Sports Model 3500 Regal-X \$63.59) is attached to the handle section of a rod (Outdoor Sports Free). Drag is adjusted by a knob with numbers on the front of the reel.

¹ For example start at http://www.hovercraftofamerica.org/ or do a search on "hovercraft".

² Ken Altshuler, "The Human Airpuck", Phys. Teach.. 27, 615-617, (Nov 1989).

Fig.1(a). Hovercraft Top View



(b). Side View





Fig. 2(a) Airpuck Bottom View



(b). Side View







Fig. 3(a). Hovercraft (Seat Removed)






Attn Ernest Provo 413-538 2357



P.01/02

220 678 8227 TO 14135382357

OCT 26 '98 16:20 FR AMETEK-RTMD



CB 960 IN UN T U (yGN EQUAUALENT 3 75 GR UNO HALL HALL YEL SENSOR + SUPPLY HAU GRW OUT HALL BLU HAU Purce HAU GRAY BROWN REO REO 1 2 0 RN BROWN

ORANGE