


Analytical Ultracentrifugation of Nanoparticles

Helmut Cölfen

Max-Planck-Institute of Colloids and Interfaces,
Colloid Chemistry, Am Mühlenberg 2, D-14424 Potsdam
Email: Coelfen@mpikg-golm.mpg.de

Why fractionating analytics in solution ?



Svedbergs Colloid Research




- T. Svedberg, J.B. Nichols, *J. Amer. Chem. Soc.* 45 (1923) 2910
- T. Svedberg, H. Rinde, *J. Amer. Chem. Soc.* 45 (1923) 943
- T. Svedberg, H. Rinde, *J. Amer. Chem. Soc.* 46 (1924) 2677
- T. Svedberg, *Kolloid-Z. Zsigmondy Festschrift, Erg.-Bd. Zu 36* (1925) 53

1925: Nobel prize for colloid work


1926: First protein work

Principle of AUC


Radial Detection



Rayleigh Interference optics



UV/VIS Absorption optics




Experiments very simple, evaluation not always

Different AUC experiments


Sedimentation velocity

High centrifugal force
Sedimentation stronger than back diffusion




Sedimentation equilibrium

Moderate centrifugal force
Sedimentation in the order of back diffusion



Density gradient

High centrifugal force for distribution of a low molecular salt or a second solvent




Lamm equation

$$\frac{\partial c}{\partial t} = \frac{1}{r} \frac{d}{dr} \left(\underbrace{rD \frac{dc}{dr}}_{\text{Diffusion term}} - \underbrace{s\omega^2 r^2 c}_{\text{Sedimentation term}} \right)$$

Experiment	Effective term in the Lamm equation	Characteristics
Sedimentation velocity	Sedimentation term much bigger than diffusion term	High rotational speed
Synthetic boundary experiment for the determination of D	Only diffusion term effective	Synthetic boundary cell, low rotational speed
Sedimentation velocity	Sedimentation and diffusion term effective / equilibrium between sedimentation and diffusion	Moderate / low rotational speed
Density gradient (special case of sedimentation equilibrium)	Sedimentation and diffusion term effective / equilibrium between sedimentation and diffusion	Moderate / high rotational speed, locally dependent solution density

Basic evaluation important for nanoparticles



1 step means
1 component

Flat baseline
indicates
purity

- Determine distance travelled in given time interval
- Calculate sedimentation coefficient s or its distribution
- Calculate particle size or distribution

$$s = \frac{u}{\omega^2 r}$$

$$s = \frac{\ln(r/r_m)}{\omega^2 t}$$

$$d_i = \sqrt{\frac{18\eta s_i}{\rho_2 - \rho}}$$

Sedimentation velocity

$$M = \frac{sRT}{D(1 - \bar{v}\rho)} \quad f = \frac{RT}{N_A D} \quad d_i = \sqrt{\frac{18\eta s_i}{\rho_2 - \rho}}$$

Sedimentation velocity depends on:

Molar mass / particle size

Density

Shape / Friction

Charge


for given exptl. parameters (temperature, solvent density & viscosity etc.)

Particle size distributions calculated on a hard sphere basis


Sedimentation velocity

Gold in
H₂O at
5000 RPM,
25 °C


a) Raw data




b) S-
evaluation



c) s-
distribution




d) d-
distribution




Common Problems

- Extremely broad s-distributions, Big aggregates or small impurities are not detected
- Colloids aggregate or grow during centrifugation (concentration dependent aggregation)
- Density of hybrid colloids is unknown to access particle size
- Electrostatic stabilization complicates analysis due to charge contributions
- Particle polydispersity in size, shape, density and hydration
- High particle density often makes density gradients or density variation methods impossible
- Often multicomponent mixtures


Fractionation of heterogeneous samples



pH 2, κ -carrageenan




+ Fe^{3+}



pH 2, Fe^{3+} crosslinking of κ -carrageenan

+ OH^-

pH > 5




Hollow space

Protected partially mobile FeOOH

Latex mixture

A. Völkel



$$d = \sqrt{\frac{18\eta s}{\rho_P - \rho_{LM}}}$$




- **Particles analyzed over colloidal range in one experiment**

given		reproduced	
d [nm]	%	d [nm]	%
60	20	60	20
103	20	100	20
160	20	164	18
347	20	389	20
2800	20	3400	22

Very small colloids



Reported to be definite cluster with 55 Au, Transition between metal and molecule





High resolution PSD

Pt colloid in MeOH / HAc


0.1 nm Baseline resolution

Problem:
Elimination of diffusion broadening by self-sharpening effects is not common

Baseline resolution > 1 Angström




TEM



Particles slightly bigger than from AUC (Density)

Pt-Colloid growth

▲ = Stabilizer ● = Pt-Particle



Critical
crystal
nucleus

Coalescence until particles
are stabilized

Further growth through
reduction of PtCl_4 at the
crystal surface

$$\frac{dr}{dt} = \text{const.}$$


Differences in the particle size after
coalescence are conserved

Growing ZnO Colloid

Zn(Ac)₂ (1 mmol/l)
+ 20 mmol/l NaOH
in water free
Isopropanol



Dilution 1 : 5 and
Heating to 65 °C

Start of heating defines
start of the reaction




Problem: AUC has a low time resolution



Bio-labeling




Differential distribution




+
BSA




Phase transfer and bio-labeling







Dopamine functionalized TiO₂



Concept of turbidimetric immunoassay




Concept of turbidimetric immunoassay




**Big particles
aggregate first,
then the small
ones**

TEM of Immunoassay

0 mg/L CRP




4.28 mg/L CRP




Statistics, drying artifacts ?

25 mg/L CRP





156 mg/L CRP



Scale bars 500 nm



SLS of Immunoassay



- Fast measurement without equilibration needs
- Kinetic measurements possible
- Ambient conditions without any pressure effects
- Low resolution in particle size
- Low statistical resolution (No fractionation)

92 wt-% small latices

AUC of Immunoassay



- AUC has equilibration time > 10 min before experiment
- Fractionation enables to determine correct particle quantities and sizes by turbidity detection (MIE correction) with speed profile even over decades of s -values
- 92 wt-% small latices in mixture correctly detected

Microgels






Figure from:
 H.G. Müller, A.
 Schmidt, D. Kranz;
 Progr. Colloid
 Polym. Sci. **86**, 70
 (1991)

$$Q = \frac{b \cdot d_h^2}{s} \cdot \frac{\rho_P - \rho_s}{18 \cdot \eta}$$



b = mass ratio of soluble parts

PNIPAM Microgels

Fringe shift of crosslinked part: 7.90 **76%**
Fringe shift of free chains : 2.37 **24%**



Swelling degree distribution





$$Q = \left(\frac{s_{\text{unswollen}}}{s_{\text{swollen}}} \right)^3$$

In principle also possible with soluble fraction b

$$Q = \frac{b \cdot d_h^2}{s} \cdot \frac{\rho_P - \rho_s}{18 \cdot \eta}$$

Solvent density variation


PS₁₂₃-P4VP₁₄₅ in toluene resp. d-toluene



Combination of 2 s-distributions in 2 chemically identical solvents yields **particle size and density distribution**

- Density in good agreement with values from density meter
- Particle size in agreement with results from DLS or viscosity measurements
- Same result for core crosslinked micelles

PS-PMAA Latex mixture 40 : 60




$$d = \sqrt{\frac{18(s_2 \eta_{s2} - s_1 \eta_{s1})}{\rho_{s1} - \rho_{s2}}}$$

$$\rho_p = \frac{s_1 \eta_{s1} \rho_{s2} - s_2 \eta_{s2} \rho_{s1}}{s_1 \eta_{s1} - s_2 \eta_{s2}}$$

Synthetic hybrid colloids

Brushite $\text{CaHPO}_4 \cdot 2\text{H}_2\text{O}$
TG 18% Mineral




A: Star like crystals/whisker, $d = 17 \text{ nm}$, $l = 184 \text{ nm}$

B: Transversal growth, destabilization of filaments, further densification of core

C: Compact aggregates with periodicity („Chrysanthemum“ structure)


Synthetic hybrid colloids



- Superposition of particle size and density distribution
- Densities can exceed 2 g/ml

Synthetic hybrid colloids

density increase



- Clear monitoring of density increase upon Ca²⁺ complexation
- Global Analysis could reveal density and molar mass

Combination of AUC with DLS or better FI-FFF is highly desirable

Particle size and density

UV : $d = 1.65 \text{ nm}$

DLS : $d = 1.6 \text{ nm}$




AUC: " $d = 1.25 \text{ nm}$ " ($s = 3.6 \text{ S}$)

calculated with $\rho_{\text{CdS}} = 4.83 \text{ g/cm}^3$

CdS



$\rho_{\text{particle}} = 3.2 \text{ g/cm}^3$
 $d_{\text{particle}} = 1.6 \text{ nm}$




$$\rho_p = \frac{18\eta s}{d_p^2} + \rho_s$$

Density core and shell are known

- $d_{\text{total}} = 1.6 \text{ nm}$
- $d_{\text{core}} = 1.1 \text{ nm}$
- $d_{\text{shell}} = 0.5 \text{ nm}$
(Lit.: $d_{\text{shell}} = 0.4 \text{ nm}$)

With independent size measurement, hybrid particle is fully characterized



Hybrid Colloids



- Iron storage protein Ferritin forms robust hollow sphere from 24 dimeric subunits
- Ferritin oligomerizes
- Different amounts of ferrihydrite inside the core (up to 4500 Fe)
- **Simultaneous particle size and density distribution**
- **Size and amount of oligomers ?**
- **Amount of ferrihydrite inside the core for the different oligomers ?**

Ferritin



FI-FFF separates only according to particle size, AUC according to particle size and density



AUC no baseline separation of oligomers (density distribution)


Global analysis

	Monomer	Dimer	Trimer
s^* at 25 °C (AUC)	59.8 S	96.9 S	127.9 S
D^* at 25 °C (FI-FFF)	$3.68 \times 10^{-7} \text{ cm}^2/\text{s}$	$2.41 \times 10^{-7} \text{ cm}^2/\text{s}$	$1.90 \times 10^{-7} \text{ cm}^2/\text{s}$
d_H (FI-FFF)	11.9 nm	18.1 nm	23.0 nm
ρ_P	1.764 g/cm³	1.531 g/cm³	1.436 g/cm³
M_w	586,700 g/mol	954,100 g/mol	1,367,400 g/mol
Oligomer amount	69.4 wt.-%	19.5 wt.-%	11.1 wt.-%
$n \text{ Fe}^{3+}$	1476	667	34



First direct observation of near critical size clusters by AFM


S.T. Yau, P.G. Vekilov „Quasi-planar nucleus structure in apoferritin crystallization“, Nature **406**, 494 - 497; 3 August 2000



D.W. Oxtoby „Catching crystals at birth“, Nature **406**, 464 - 465; 3 August 2000:


„Although the first small crystallites involved in nucleation *form in the bulk of the solution*, Yau and Vekilov can observe them only once they have fallen to the bottom of the vessel and attached themselves to the surface there. How can the authors be sure that they are observing the critical early stages of crystallization described by nucleation ?“

Synthetic boundary cell of the Vinograd Type



Reservoir


Capillaries




Reactant (■) containing solution is layered onto the second reactant via the capillary (■)

Formation of a sharp reaction boundary


Particle formation in a synthetic boundary cell



Transformation of a time distribution into a radial distribution



Synthetic boundary crystallization ultracentrifugation



All known stable CdS growth species simultaneously accessible in one experiment. Even subcritical clusters detectable in solution

Conclusion

- AUC is a very versatile instrument for colloid chemistry (Range of possible methods still not fully explored !!!)
- Fractionation allows the investigation of complex mixtures or very polydisperse samples with s ranging over decades (speed profile)
- High statistical accuracy
- Critical crystal nuclei and subcritical clusters can be resolved
- Crystal growth reactions can be investigated by transformation of time distribution into radial distribution

BUT:

- Colloids can express various unfavourable properties which makes them problematic to investigate
- Fast multidetection systems together with global analysis can potentially address even most complex mixtures with particle size and VBAR distribution but role of charge still unaddressed

Future Perspectives

- Fast speed profiles and detectors to suppress diffusion broadening
- Multi-detection systems (RI, UV-Vis, SLS and others)
- Global Analysis (promising with SEDVEL & SEDEQUIL + DLS & FI-FFF) to obtain s , M , D & VBAR distributions
- New ultracentrifugation methods could allow to access new quantities

Surface tension from phase transfer, particle charge determination via sedimentation in pH or charge gradient, Conformational changes in solvent quality gradients, new synthetic boundary techniques for membrane / crystal growth investigations or chemical reactions in the AUC, etc.

Recommended reading

. Svedberg, K.O. Pedersen, *The ultracentrifuge*, Clarendon Press, Oxford (1940)

the classical textbook about analytical ultracentrifugation still with much impact today

.K. Schachman, *Ultracentrifugation in biochemistry*, Academic Press, New York (1959)

compact and useful book covering experimental and theoretical aspects

.E. Harding, A.J. Rowe and J.C. Horton; *Analytical ultracentrifugation in biochemistry and polymer science*, The royal society of chemistry, Cambridge, ISBN 0-85186-345-0 (1992)

the most comprehensive modern book about analytical ultracentrifugation. A very good overview about methods & techniques of analytical ultracentrifugation and a valuable source of modern applications.