



# Structural Characterization of Nanoparticles

Nicola Pinna

Max Planck Institute of Colloids and Interfaces

e-mail: pinna@mpikg-golm.mpg.de - http://www.pinna.cx





- 1. Transmission Electron Microscopy
  - Selected Area Electron Diffraction
  - High Resolution Transmission Electron Microscopy
  - HRTEM Simulations
- 2. X-Ray Diffraction
  - Particle Size Broadening
  - Sherrer equation
  - Debye scattering equation



# Transmission Electron Microscopy







# Transmission Electron Microscopy

C





## Transmission Electron Microscopy







# Assembly of Nanoparticles







# **Electron** Diffraction







#### HRTEM & Power Spectrum





 $N_{pix-image}$ : N° of pixels along the quadratic frame (128, 256, 512 ...) of the image  $P_{size}$ : Pixel size in the HR image  $N_{pix}$ : N° of pixel between the origin and the reflection hkl studied



#### HCP - FCC structures\*





Crystals oriented along the [110] direction

\*The wurtzite structure is composed by two hcp networks, one occupied by the sulfur and the other by the cadmium, shifted by  $OO\frac{1}{4}$ . The zinc blende structure is composed by two fcc arrays, as for the hcp occupied by the sulfur and the zinc respectively, shifted by  $\frac{1}{4}\frac{1}{4}\frac{1}{4}$ .











Computer simulation of the HRTEM images on the basis of the structure of model particles with the multislice technique:

- 1. The construction of one or more atomic models of the nanocrystals
- 2. The calculation of the HRTEM image of these models
- 3. The calculation of the PS of the calculated HRTEM image
- 4. The comparison of the HRTEM images and the PS calculated with the data obtained from the experimental HRTEM







Electrons are assumed to scatter only in a forward direction with small diffraction angles. With this approximation the crystal can be divided in sub-slices with a thickness  $\Delta z$  perpendicular to the incident beam.

- 1. The crystal is divided in slices perpendicular to the electron beam
- 2. The electrostatic potential V(x, y) with in-plane coordinates x,y of the sliced crystal or supercell is projected for each slice of the included atoms onto its exit surface
- 3. On the basis of  $V_P(x, y)$  the amplitude of the electron wave function is calculated
- 4. Calculate the propagation of the electron wave throught all the slices



\*J.W. Cowley, A.F. Moodie, Acta Cryst. 10, 609, 1957



G





#### **HRTEM Simulations**







# **HRTEM Simulations**

G





## **HRTEM** Simulations







# Copper Nanoparticles

C





# Copper Nanoparticles























C













![](_page_23_Picture_0.jpeg)

#### **Oriented Attachement**

![](_page_23_Picture_2.jpeg)

![](_page_23_Figure_3.jpeg)

J. Polleux, N. Pinna, M. Antonietti, M. Niederberger, Adv. Mat. 2004, 16, 436

![](_page_24_Picture_0.jpeg)

![](_page_24_Picture_2.jpeg)

Inelastic scattering of an electron of the incident beam and the atomic electrons of the solid.

- Transition from an inner-shell (K, L M ...) to an unoccupied energy level (i.e. above the Fermi level)
- Transition of a valence electron across the energy gap (insulator, semiconductor) or excitation of a plasma resonance (mostly in metals metals).

![](_page_24_Figure_6.jpeg)

![](_page_25_Picture_0.jpeg)

#### Electron Energy Loss Spectroscopy

![](_page_25_Picture_2.jpeg)

![](_page_25_Picture_3.jpeg)

![](_page_25_Figure_4.jpeg)

![](_page_26_Picture_0.jpeg)

![](_page_26_Picture_2.jpeg)

### $YBa_2CU_3O_7$

#### Carbon allotropes

![](_page_26_Figure_5.jpeg)

![](_page_27_Picture_0.jpeg)

## Electron Energy Loss Spectrometry

![](_page_27_Picture_2.jpeg)

![](_page_27_Figure_3.jpeg)

![](_page_28_Picture_0.jpeg)

#### Electron Energy Loss Spectrometry

![](_page_28_Figure_2.jpeg)

![](_page_29_Picture_0.jpeg)

![](_page_29_Picture_2.jpeg)

- TEM permits the structural characterization of a collection or isolated nanoparticles
- Electron Diffraction  $\rightarrow$  structure of single or many particles
- HRTEM  $\rightarrow$  structure, orientation, crystallinity, defaults
- Image processing: fondamental tool for structural studies in electron microscopy
- Spectroscopy: EELS, EDX  $\rightarrow$  local structure, band structure, composition
- Do not permit to study structural properties of the whole sample  $\rightarrow$  Necessity to compare the results with other techniques

![](_page_30_Picture_0.jpeg)

![](_page_30_Picture_2.jpeg)

#### Particle Size Broadening

$$I = I_e F^2 * \frac{\sin^2(\pi/\lambda)(\mathbf{s} - \mathbf{s_0})N_1 \mathbf{a_1}}{\sin^2(\pi/\lambda)(\mathbf{s} - \mathbf{s_0})\mathbf{a_1}} * \frac{\cdots N_2 \mathbf{a_2}}{\cdots \mathbf{a_2}} * \frac{\cdots N_3 \mathbf{a_3}}{\cdots \mathbf{a_3}}$$
(2)

 $N_1, N_2, N_3$  Number of the unit cells along the  $a_1, a_2, a_3$  directions Normally  $N_1, N_2, N_3$  are large numbers  $\rightarrow$  the three quotients differs from zero only if the three Laue equations are closely satisfied. If  $N_1, N_2, N_3$  are small, the quotients broaden.

![](_page_30_Figure_6.jpeg)

![](_page_31_Picture_0.jpeg)

Approximations:

- 1) Cubic crystal  $N_1, N_2, N_3 = N$
- 2) Crystal free from strains and faulting  $\rightarrow$  peak broadening is only due to the small crystallite size
- 3) Each of the three quotients of equation 2 by a Gaussian function

![](_page_31_Figure_6.jpeg)

![](_page_32_Picture_0.jpeg)

![](_page_32_Picture_2.jpeg)

(4)

The intensity distribution spherical averaged over the reciprocal space is described by the Debye formula:

$$I_N(b) = \sum_{n,m\neq n}^N f_n f_m \frac{\sin(2\pi b r_{nm})}{2\pi b r_{nm}}$$

 $b = \frac{1}{d} = \frac{2sin\theta}{\lambda}$   $r_{nm}$  distance between atom n, m $f_n, f_m$  atomic scattering factors

- General equation valid for any form of matter in which there is a random orientation: gases, liquids, amorphous solids, and crystalline powders.
- No limitation on the number of different kinds of atoms in the sample.
- The number of terms increases proportional to the sixth order!

![](_page_33_Picture_0.jpeg)

![](_page_33_Picture_1.jpeg)

![](_page_33_Picture_2.jpeg)

![](_page_33_Figure_3.jpeg)

![](_page_34_Picture_0.jpeg)

# $BaTiO_3 6 nm$

![](_page_34_Picture_2.jpeg)

![](_page_34_Figure_3.jpeg)

![](_page_35_Figure_0.jpeg)

![](_page_35_Figure_1.jpeg)

![](_page_35_Picture_2.jpeg)

![](_page_35_Figure_3.jpeg)

![](_page_36_Picture_0.jpeg)

## Final Example: HfO<sub>2</sub>

![](_page_36_Picture_2.jpeg)

Ī13

Ī12

111

111

011

• Ī11

![](_page_36_Figure_3.jpeg)

![](_page_37_Picture_0.jpeg)

Final Example: HfO<sub>2</sub>

![](_page_37_Picture_2.jpeg)

![](_page_37_Figure_3.jpeg)

![](_page_38_Picture_0.jpeg)

#### Conclusion

![](_page_38_Picture_2.jpeg)

X-Ray diffraction associated to calculations is a powerfull tool to study the:

- Structure
- Crystallinity
- Particle size and size distribution
- Particle shape
- Homogeneity of the whole sample

It is the perfect tool to be associated with transmission electron microscopy study