




CNT in Microscopy




Atomic Force Microscopy is a powerful technique for imaging, nanomanipulation, as platform for sensor work, nanolithography...

Conventional silicon or tungsten tips wear out quickly.
CNT tip is robust, offers amazing resolution.



Simulated Mars dust






2 nm thick Au on Mica




NASA Ames Research Center
Ramsey Stevens, Lance Delzeit, Cattien Nguyen

MWNT Scanning Probe






280 nm line/space. Array of polymeric resist on a silicon substrate.





Conventional Si Pyramidal
Cantilever




MWNT Probe




AFM Image with a MWNT Tip 193 nm IBM Version 2 Resist







DUV Photoresist Patterns Generated by Interferometric Lithography




Nguyen et al., *App. Phys. Lett.*, 81, 5, p. 901 (2002).



AFM Imaging with Single Wall Nanotube Tips




2 nm thick Au on Mica



5 nm thick Ir on Mica





Si₃N₄ on Silicon substrate




Nguyen et al., Nanotechnology, 12, 363 (2001).

High Resolution Imaging of Biological Materials



DNA



PROTEIN