


# **LbL Technology in industrial applications**

## **Particles 2009 Berlin**


Lars Dähne

Surflay Nanotec GmbH Berlin



## The Company

- Surflay (Surface Layers) founded 2008 Berlin
- Based on Layer-by-Layer (LbL) Technology
- Spin-Off from Capsulation Nanoscience AG
- Splitting of activities in Pharma-Applications (Capsulation) and None-Pharma (Surflay)
- Licence on relevant patents of MPG and Caps.



## Interdisciplinary Team

### Administration

Dr. habil Lars Dähne (CEO)

Dr. Claudia Aldenhoven (Biology)

### Research & Development

Dr. Lars Dähne (Phys. Chemistry)


Dr. Gabriella Egri (Biotechnology)

Dr. Zhongbing Huang (Org. Chemistry)

Dipl. Ing. Jing Kang (Material Science)

Barbara Baude (Chem. Techn. Ass.)

Moritz Klickermann (Chemical Ingenieur)



## The Technology

Specialized on coating of colloidal and planar surfaces for material functionalization.

### **Key technology**


LbL-Technology with functionalized polyelectrolytes and nanoparticles.

### **Application fields**



Specific separation materials

Pigment stabilization


Sensoric and diagnostic Particles.



# LbL-Technology



Conc. 1 g/L



## Products

- **Functionalized particles (fluorescence, magnetic, biozide, enzymatic, DNA specific etc.)**
- **Encapsulated macromolecules and hydrophobic materials**
- **Ultrasound contrast agents (Microbubbles)**
- **Polyelectrolyte capsules 200 nm – 20 µm**
- **Labelled polyelectrolytes**


## Examples of Industrial developments with LbL


- Plant protection Particles against miner moths




### Possible Capsule Solution:

- Sticking Layer
- Appetizer
- Pheromones
- Insecticide (release during digestion)





## Oligonucleotid Beads




Magnetic Beads with specific Oligonucleotide sequence

Outermost layer labelled with FRET-Donor


Detection of hybridization by FRET-Signal (Acceptor-DNA)

High hybridization capacity due to nanoporosity of LbL







# Sun protection




OMC (o-Methoxycinnamate, oil) causes strong radical formation after penetrating the skin → accelerated aging, cancer


Solution: oil encapsulation in 5 μm beads





a) EP 6.1 1layer    b) EP 6.2    c) EP 6.3    d) EP 6.4




e) EP 6.5    f) EP 6.6    g) EP 6.7    h) EP 6.8




# Drug reservoirs on stent surfaces




REM-image of calcinated (600°C) Sol-Gel-sample (size 15 μm x 15 μm)




Confocal image model drug filled pores




## Reversible micro glucose sensors



4  $\mu\text{m}$  capsules filled with  
ConA-DyeA Dextrane-Dye B  
Complex + glucose release of  
dextrane



Fluorescence spectra of microsensor-  
capsules  $\lambda_{\text{exc}} = 585 \text{ nm}$



## Our partners

