



Instrumentation for Ambient Ultrafine Particle Measurement

**historic perspectives
and
recent developments**

**Susanne Hering
Aerosol Dynamics Inc.**


Ultrafine Particles Dominate Particle Number Concentrations




Rural and Urban Sites : Similar Relationship between Ultrafine Particle Number Concentrations and Total Particle Number Concentrations

Simplest Indicator for Ultrafine Particles:

Particle Number Concentration



Bricard et al, 1976



Hering et al, 2005


1889: First Measurement of Atmospheric Particle Number Concentrations

J. AITKEN ON THE NUMBER OF DUST PARTICLES IN THE ATMOSPHERE.

John Aitken,


“On the Number of Dust
Particles in the
Atmosphere”,

Transactions of the
Royal Society of
Edinburgh, 1889



How Aitken's Instrument Worked

- Humidified air sample
- Expanded adiabatically
- Droplets formed around particles, which then settled
- Counted manually
- Repeated with dilution



What Aitken Observed:


Number of Dust Particles in Air.

| Source of Air. | Number per c.c. | Number per c. in. |
|--------------------------|-----------------|-------------------|
| Outside (Raining), . . . | 32,000 | 521,000 |
| Outside (Fair), . . . | 130,000 | 119,000 |
| Room, | 1,860,000 | 30,318,000 |
| Room near ceiling, . . . | 5,420,000 | 88,346,000 |
| Bunseu flame, | 30,000,000 | 489,000,000 |

“The reason of the greater number of particles in the room than that found outside was due to the particles produced by the two gas flames burning in the room at the time.”

- **1912: Wilson Cloud Chamber**
 - **1929: Nobel Prize for particle physics**
 - **Determined precise expansion ratios for avoiding homogeneous and ion-induced nucleation of particles**
- **1930s: Scholz: Automated expansion for particle counting**
- **1950s: Vonnegut: Automated counting by recognizing particles grow to uniform size.**

1950s: Vonnegut, Automated Counting




Inventor:
 Bernard Vonnegut,
 by Paul A. Frank
 His Attorney

July 20, 1954
 Filed Nov. 23, 1949
B. VONNEGUT
 METHOD AND APPARATUS FOR MEASURING THE
 CONCENTRATION OF CONDENSATION NUCLEI
 2,684,008
 2 Sheets-Sheet 2

Automated, but not continuous
 (condensing vapor: water)

1970s: Thermally Diffusive Condensation Particle Counter (CPC)



Bricard et al, 1976

Automated & Continuous (condensing vapor: butanol or other alcohol)
Widely used, suitable as detector for size distribution instruments
Many models, sold by several companies

Can you have a Continuous, Automated Particle Counter without Butanol?

- **Challenge:** too small for direct optical detection
- **Approach:** Create region of supersaturation to activate particle growth => form droplets
- **Why Supersaturation?** Equilibrium vapor pressure over a droplet is greater than over a flat surface due to free energy associated with surface (surface tension)
- **Kelvin Relation:**


$$P_{\text{droplet}} = P_{\text{flat surface}} \exp\left(\frac{2 \sigma v}{kT\rho R} \right)$$

Surface tension

Particle radius

Thermally Diffusive CPCs: Operational Principle

- Saturate flow with vapor
- Flow into cold-walled tube
- Vapor condenses on particles
- Requires slowly diffusing vapor (e.g. butanol)



Thermal Diffusivity,
air = 0.215 cm²/s


Mass Diffusivity,
butanol = 0.081 cm²/s

Note: diffusivity of water = 0.265 cm²/s > air

Does not work well with water

2003: Water Condensation Particle Counter (WCPC)

- Cold flow enters warm wet-walled tube
- Water vapor diffuses more quickly than flow warms
- Supersaturation, particle activation and growth occurs inside of a warm, wet walled tube.



2003 WCPC
(~5 nm) :

Saturator, 20C

Condenser, 60C

Optics Head


2004:
Nano-WCPC (~2.5 nm)

Saturator, 12C

Condenser, 75C

Thermal Diffusivity,
air = 0.215 cm²/s


Mass Diffusivity,
water = 0.265 cm²/s



Comparison of Centerline Water Saturation Ratios


| |
|---|
| WCPC: Inlet flow at 20°C & 100%RH, Wetted walls at 60°C |
| <hr style="border: none; border-top: 1px dashed black;"/> |
| Traditional: Inlet flow 60°C & 100%RH Wetted walls at 20°C |

Water Condensation Particle Counter




Hering et al, 2005

First laminar-flow WCPC Response to Ambient Aerosols & Vehicle Emissions



Tunnel Measurements with Antonio Miguel, Arantza Eiguren-Fernandez, UCLA

Calculated Supersaturation Profiles within the Ultrafine Water – CPC




α_v = vapor mass diffusivity


z = axial distance

Q = volumetric flow rate

Calibration of the Ultrafine WCPC Water Residue Particles



Comparison of Water CPC and TSI Ultrafine, Queen's College, NY (Univ. at Albany, EPA Supersite)



Compare:
0.2-sec average
11-sec running average



counts/cm³ measured below nanoSMPS
(not ambient size distribution – no charging correction)

Comparison to Butanol Ultrafine CPC with Nano-DMA

Particle counts per 0.2 sec window (average over two scans)

TSI-3025: 1.1 ± 1.0 particles (0.5cc/s)

WCPC: 35 ± 6 particles (16.7 cc/s) --- much better counting statistics



In collaboration with ASRC, University at Albany




Field Comparisons Among CPCs





Total Particle Number Concentrations for Traffic Emissions

Ultrafine WCPC (TSI-3786) compared to Butanol UCPC (TSI-3025)



Ambient Sampling In Riverside, California

Ultrafine WCPC – 3786 & Butanol UCPC - 3025



2005: Micro-Environmental WCPC

Size: 7" x 7" x 5"

Weight: 5 lb

Power 12 V, 30 Watts

Features:

internal data logging


one week unattended

up to 10^6 particles/cm³


in single count mode



Response to Near-Monodisperse Aerosols



ME-WCPC Measurements in a Residential Kitchen



Summary

- **1890's**

Aitken, first explorations of airborne particles number concentration

Identified combustion as source of particles

- **1920s-1950s**

Advances on Aitken's approach: automated instruments

- **1970's**

First continuous flow condensation particle counters

widely used, especially as detectors for mobility size distributions

- **2003**

Introduction of continuous water-based condensation counters

Acknowledgements

- California Air Resources Board, who provided funding for many of the field comparisons.
- ASRC, University at Albany and the EPA Supersites Program who made possible the measurements in New York City.
- TSI Inc., who loaned equipment for our field and laboratory testing.
- Quant Technologies LLC, who provided the engineering and detailed instrument design.