

# **NANOSTRUCTURED ENVIRONMENTAL BIOCHEMICAL SENSOR FOR WATER MONITORING**

## **Nanotechnology and the Environment**

Günther Proll, Günter Gauglitz

Institute of Physical and Theoretical Chemistry  
Eberhard-Karls-University of Tübingen




## **Outline**

- **Biosensor Principles and Components**
- **Immunochemistry & Assay Design**
- **AWACSS Project: Overview & Intentions**
- **Instrument & Network**
- **Multi-analyte Measurements**
- **Environmental Assays**
- **Summary & Acknowledgements**



# Biosensor principles



# Biosensor principles

## Labeled Systems

- ELISA (heterogeneous, competitive)
- Direct test format with label:  
Liegler et al.; Krämer et al.
- Binding-inhibition-assay with label  
RIANA; **AWACSS**




# Immunochemistry

## Antigen-Antibody-Interaction:

### *Binding Inhibition Assay*

Specific antibodies which are labelled with a fluorescence dye (Cy5.5<sup>TM</sup>) recognise organic compounds in aqueous samples.



# Components & parameters

Analytical tool

Antibody

- Affinity-chromatography
- Fluorescent dye (Cy5.5, Alexa 680, QDs)

Method

Assay

- Surface chemistry
- Regeneration
- Automation

Transduction

IO-Chip

- Evanescant field
- Enhancement structure

Detection

Optics  
Electronics

- Filters / PMF
- Photodiodes
- Amplifier

Data Treatment

AD  
SW

- Drift correction
- Calibration
- Chemometrics

Fully automated biosensor

# Nano-scale

Analytical tool

Antibody  
ng

•Affinity-chromatography

•Fluorescent dye (Cy5.5, Alexa 680, QDs)  
nanoparticles)

Method

Assay

•Surface chemistry  
pL - nL

•Regeneration  
•Automation

Transduction

nanostr.  
IO-Chip

•Evanescent field  
nm

•Enhancement structure  
nm

Detection

Optics  
Electronics

•Filters / PMF

•Photodiodes  
pW - nW

•Amplifier

Data Treatment

AD  
SW

•Drift correction

•Calibration




•Chemometrics

Fully automated biosensor


LOD: ng L<sup>-1</sup>



# QDs as labels



# AWACSS-Project




**Automated  
Water  
Analyser  
Computer  
Supported  
System**

- EVK1-2000-00045
- Duration: 3 Years, RTD Project
- (March 2001 – March 2004)
- Visit: <http://barolo.ipc.uni-tuebingen.de/awacss>

*Biosensors & Bioelectronics 20 (2005) Part I & II (1499-1508 & 1509-1519)*

## Goal

Develop a cost-effective, on-line, water monitoring biosensor that will measure a variety of small organic pollutants in short-time with remote control and surveillance.





# Software & Network



Local station

## Local stations:

- Multi-analyte detection
- Data acquisition
- Data evaluation
- Local station control



Local station

## Control station:

- Global survey & early warning
- Global trend analyses
- Local station control
- Detection limit control prescribed by EU-Directives



Control station

# Target Compounds

## Ecotoxicology:

- Herbicides
- Fungicides
- Insecticides
- Antibiotics
- Blue Algae Toxins
- Endocrine Disrupting Compounds
- Suspected Carcinogens
- Industrial Wastes / Chemical Markers




# Integrated Optical Chips

Goal: Develop IO-Chips capable of simultaneously illuminating 32 analysis windows in a miniaturised format.

## Advantages:

- Single mode waveguides allow for efficient and consistent analysis window illumination.
- Compact dimensions allow for microfluidic design and reduce sample & reagent volumes.
- Spatially defined but resolved analysis windows allow efficient detection and suppress signal crosstalk.




# Integrated Optical Chips

Goal: Develop IO-Chips capable of simultaneously illuminating 32 analysis windows in a miniaturised format.


## Derivatives:

-  Estrone
-  Bisphenol A
-  Atrazine
-  Sulphonamides
-  Isoproturon
-  Propanil




18 minutes per cycle and up to 500 measurements!

# Instrument



# Multi-analyte Measurements







# Validation Parameters

Analyte	A <sub>2</sub>	p	IC <sub>50</sub>	LOD	SD <sub>ZS</sub>
Estrone	16.7	1.46	0.093	<b>0.007</b>	0.6
Bisphenol A	29.5	0.61	1.585	<b>0.008</b>	0.9
Atrazine	20.8	0.76	0.217	<b>0.010</b>	2.7
Sulphamethizole	34.3	0.61	0.419	<b>0.018</b>	2.8
Isoproturon	8.5	1.3	1.263	<b>0.020</b>	2.0
Propanil	9.6	1.08	0.232	<b>0.019</b>	1.9
Dimension	[ % ]	[ ]	[ µg L <sup>-1</sup> ]	[ µg L <sup>-1</sup> ]	[ % ]



# Drinking water

Analyt	Estron	BPA	Atrazin	SUL	IP	PRN
RR (0,05 ppb)	106,2 ±3,6	31,6 ±4,2	86,7 ±18,4	110,8 ±5,3	117,4 ±4,4	96,7 ±15,9
RR (0,10 ppb)	119,4 ±6,4	143,8 ±12,8	87,0 ±16,8	67,4 ±4,0	96,9 ±3,8	80,0 ±3,9
RR (0,15 ppb)	113,8 ±11,2	101,8 ±24,8	111,5 ±10,3	109,6 ±9,8	98,0 ±4,4	115,8 ±5,3

- RR recovery rate (70 – 120 % as recommended by the AOAC International)
- $RR = c_d / c_t \cdot 100 \%$  (d: detected, t: true)



# Collaborative Trial

## Objectives:


- **Evaluation of the AWACSS performance in comparison to:**
  - Immunoassay techniques – RIANA and ELISA
  - Conventional LC- and GC-based analytical techniques




# Collaborative Trial

- **Matrices**
  - Milli-Q water
  - River sediments (homogenized 63  $\mu\text{m}$  fraction, freeze dried)
  - Mixed industrial and municipal wastewater - Novaky Chemical Plant, Slovak Republic
- **Analytes**
  - Atrazine - pesticides, on the list of EU WFD Priority Substances
  - Bisphenol A - endocrine disrupting compound, industrial pollutant
  - Estrone - endocrine disrupting compounds, hormones, municipal WWTP
- **Spiked levels**
  - 0.1 and 1  $\mu\text{g L}^{-1}$  - water samples
  - 50 and 500  $\text{ng g}^{-1}$  - sediments


# Collaborative Trial




# Collaborative Trial



# Collaborative Trial



# Collaborative Trial






# Collaborative Trial

Compound	Assigned value	Matrix		Assigned value		
		Milli-Q water			Sediment	
	[ $\mu\text{g L}^{-1}$ ]	Mean [ $\mu\text{g L}^{-1}$ ]	RR	[ $\mu\text{g L}^{-1}$ ]	Mean [ $\mu\text{g L}^{-1}$ ]	RR
<b>Atrazine</b>						
Level 1	0.1	0.11	111	0.1	0.1	100
Level 2	1	1.15	115	1	0.77	77
<b>Bisphenol A</b>						
Level 1	0.1	0.08	80	0.1	0.133	133
Level 2	1	1.25	125	1	0.85	85
<b>Estrone</b>						
Level 1	0.1	0.08	80	0.1	0.071	71
Level 2	1	1.04	104	1	0.71	71

RR: Recovery Rate





# Field Test - Sampling System



# AWACSS Biosensor (Field Test)

Calibration: atrazine, isoproturon, estrone, bisphenol A  
Fully automated sampling and detection!



- Estrone below LOQ (AWACSS and HPLC-DAD)
- Bisphenol A contamination from tubes



# Assays: How far it can go ...

Jens Tschmelak · Guenther Proll · Guenter Gauglitz

**Immunosensor for estrone with an equal limit of detection as common analytical methods**




Table 1 Values for the logistic fit function and the calculated results for  $SD_{\text{blank}}$ , LOD and LOQ

Parameters	Values	Calculated results	Values
$A_2$	$17.16 \pm 1.54$	$SD_{\text{blank}}$	0.67%
$x_0$	$0.062 \pm 0.004$	LOD	$0.19 \text{ ng L}^{-1}$
$p$	$0.639 \pm 0.052$	LOQ	$1.39 \text{ ng L}^{-1}$


**LOD: 190 ppq**

Fig 1 Calibration for estrone from 0 to  $9 \mu\text{g L}^{-1}$  in Milli-Q water with ovalbumin from chicken eggs, phosphate buffered saline and 3 ng anti-total-estrogen in each sample, with total sample volume of 1 mL.





# Progesterone in raw milk (1:10)



## Summary

The AWACSS biosensor is a fully automated water monitoring system and capable for remote control:

- Quick!
- No pre-treatment nor pre-concentration!
- Multi-analyte measurements!
- Real samples: tap, ground and surface!
- Successful field test
- LODs below EU limits!
- LOD in ppq-range for several compounds: Progesterone, Testosterone, Propanil, Estrone, Isoproturone, Atrazine, ...



# Acknowledgements

- EU for funding AWA CSS
- All partners of the AWA CSS project for the excellent cooperation!

**Central Research Laboratories**  
Richard Nudd & Michael Jackson

**Optoelectronics Research Centre**  
James S. Wilkinson, Ping Hua and J. Patrick Hole

**Siemens AG**  
Joachim Kaiser, Peter Krämer and Luis Barzaga

**King `s College London**  
Ram Abuknesha et al.

**DVGW-Technologiezentrum Wasser**  
Frank Sacher and Jan Stien

**Institute of Chem. & Environ. Chem. IIQAB**

Damià Barceló, Sara Rodríguez-Mozaz and Maria J. López de Alda

**Environmental Institute**  
Jaroslav Slobodník, Peter Oswald, Helena Kozmenko and Eva Korenková

**Water Research Institute**  
Livia Tóthová and Zoltan Kracszenits



# Thank you for your attention!