

**WHO International Standard
Prolactin, Human.
NIBSC code: 84/500
Instructions for use
(Version 3.0, Dated 10/12/2007)**

1. INTENDED USE

The material consists of a batch of ampoules (coded 84/500) which was established as the 3rd International Standard at the 39th Meeting of the Expert Committee on Biological Standardization of WHO, Geneva^{1,2}.

2. CAUTION

This preparation is not for administration to humans.

The preparation contains material of human origin, and either the final product or the source materials, from which it is derived, have been tested and found negative for HBsAg, anti-HIV and HCV RNA. As with all materials of biological origin, this preparation should be regarded as potentially hazardous to health. It should be used and discarded according to your own laboratory's safety procedures. Such safety procedures should include the wearing of protective gloves and avoiding the generation of aerosols. Care should be exercised in opening ampoules or vials, to avoid cuts.

3. UNITAGE

Each ampoule of the International Standard contains 53 MILLI INTERNATIONAL UNITS (by definition).

4. CONTENTS

Country of origin of biological material: United Kingdom.

Each ampoule of the material contains the residue, after freeze-drying, of a solution which contained:

Human prolactin	approx 2.5 µg
Human albumin	" 1 mg
Lactose	" 5 mg
Ammonium formate	" 0.63 mg

Pure dry nitrogen at slightly less than atmospheric pressure.


5. STORAGE

Unopened ampoules should be stored at -20°C.

Please note: because of the inherent stability of lyophilized material, NIBSC may ship these materials at ambient temperature.

6. DIRECTIONS FOR OPENING

Tap the ampoule gently to collect the material at the bottom (labelled) end. Ensure ampoule is scored all round at the narrow part of the neck, with a diamond or tungsten carbide tipped glass knife file or other suitable implement before attempting to open. Place the ampoule in the ampoule opener, positioning the score at position 'A'; shown in the diagram below. Surround the ampoule with cloth or layers of tissue paper. Grip the ampoule and holder in the hand and squeeze at point 'B'. The ampoule will snap open. Take care to avoid cuts and projectile glass fragments that enter eyes. Take care that no material is lost from the ampoule and that no glass falls into the ampoule.



Side view of ampoule opening device containing an ampoule positioned ready to open. 'A' is the score mark and 'B' the point of applied pressure.

7. USE OF MATERIAL

For all practical purposes, each ampoule contains the same quantity of the substances listed above. Dissolve the total contents of the ampoule in a known volume of a suitable solvent with carrier protein where extensive dilution is required.

No attempt should be made to weigh out any portion of the freeze-dried material.

For economy of use, it is recommended that the solution be sub-divided into several small containers and stored at -40°C, or below.

The ampoules do not contain bacteriostat and solutions of them should not be assumed to be sterile.

8. PREPARATION OF AMPOULES

A quantity of highly purified extract was donated by Dr S. S Lynch and Professor W R Butt (Birmingham & Midland Hospital for Women, Birmingham, UK). The material consisted of a single batch (H-Prol SFK2) provided as a frozen solution in 0.01M ammonium acetate. The bulk solution (35.46ml; 10mg protein) was dissolved in 4.0L of a solution, pH 6.8, containing 0.01M ammonium formate, 0.1% human albumin and 0.5% lactose and sterilized by membrane filtration. The solution was filtered through a Millex HA (Millipore, Bedford, MA, USA) membrane, mean pore diameter 0.4 micrometres. Some 4000 ampoules were filled with 1.0g of the solution maintained at 4°C. The ampouled solution was then pre-frozen to -35°C, freeze-dried and further desiccated *in vacuo*. The ampoules were then filled with pure dry nitrogen and sealed by glass fusion (Annex 4, 29th ECBS Report, 1978)³. The mean weight of solution filled in 80 weighed ampoules was 1.002gm, range ± 0.2% of the mean.

9. CONTAMINANTS

Radioimmunoassay indicated contamination with human growth hormone of about 0.5% by weight (some of this may represent cross-reaction with hPRL).

10. STABILITY

Reference materials are held at NIBSC within assured, temperature-controlled storage facilities. Reference Materials should be stored on receipt as indicated on the label.

NIBSC follows the policy of WHO with respect to its reference materials. Ampoules of 84/500 which had been stored for 84 and 692 days at elevated temperatures were examined in a few assay systems. Comparison of these samples with the samples stored continuously at -20°C, suggests that in general no immunologically detectable degradation has occurred; at -20°C the detectable loss is likely to be less than 0.1% per year, although this may be dependent on the assay system used.

11. REFERENCES

- 39th Report, WHO Expert Committee on Biological Standardization (1989) WHO Tech Rep Ser No. 786.
- Schulster D, Gaines Das R E & Jeffcoate S L (1989). International standards for human prolactin: calibration by international collaborative study. J Endocrinol 121:157-166.
- Annex 4, 29th Report WHO Expert Committee on Biological Standardization (1978). WHO Tech Rep Ser No. 626.

12. FURTHER INFORMATION

Further information can be obtained as follows;

This material:

enquiries@nibsc.hpa.org.uk

WHO Biological Standards:

Http://www.who.int/biologicals/en/

JCTLM Higher order reference materials:

Http://www.bipm.org/en/committees/jc/jctlm/

Derivation of International Units:

http://www.who.int/biologicals/reference_preparations/en/

Ordering standards from NIBSC:

[Http://www.nibsc.ac.uk/products/ordering_information/frequently_asked_questions.aspx](http://www.nibsc.ac.uk/products/ordering_information/frequently_asked_questions.aspx)
NIBSC Terms & Conditions:
[Http://www.nibsc.ac.uk/terms_and_conditions.aspx](http://www.nibsc.ac.uk/terms_and_conditions.aspx)

13. CUSTOMER FEEDBACK

Customers are encouraged to provide feedback on the suitability or use of the material provided or other aspects of our service. Please send any comments to enquiries@nibsc.hpa.org.uk

14. CITATION

In all publications, including data sheets, in which this material is referenced, it is important that the preparation's title, its status, the NIBSC code number, and the name and address of NIBSC are cited and cited correctly.

15. MATERIAL SAFETY SHEET

Physical and Chemical properties	
Physical appearance: Freeze dried powder	Corrosive: No
Stable: Yes	Oxidising: No
Hygroscopic: Yes	Irritant: No
Flammable: No	Handling: See caution, Section 2
Other (specify):	Contains material of human origin
Toxicological properties	
Effects of inhalation:	Not established, avoid inhalation
Effects of ingestion:	Not established, avoid ingestion
Effects of skin absorption:	Not established, avoid contact with skin
Suggested First Aid	
Inhalation:	Seek medical advice
Ingestion:	Seek medical advice
Contact with eyes:	Wash with copious amounts of water. Seek medical advice
Contact with skin:	Wash thoroughly with water.
Action on Spillage and Method of Disposal	
Spillage of ampoule contents should be taken up with absorbent material wetted with an appropriate disinfectant. Rinse area with an appropriate disinfectant followed by water. Absorbent materials used to treat spillage should be treated as biological waste.	

16. LIABILITY AND LOSS

Information provided by the Institute is given after the exercise of all reasonable care and skill in its compilation, preparation and issue, but it is provided without liability to the Recipient in its application and use.

It is the responsibility of the Recipient to determine the appropriateness of the standards or reference materials supplied by the Institute to the Recipient ("the Goods") for the proposed application and ensure that it has the necessary technical skills to determine that they are appropriate. Results obtained from the Goods are likely to be dependant on conditions of use by the Recipient and the variability of materials beyond the control of the Institute.

All warranties are excluded to the fullest extent permitted by law, including without limitation that the Goods are free from infectious

agents or that the supply of Goods will not infringe any rights of any third party.

The Institute shall not be liable to the Recipient for any economic loss whether direct or indirect, which arise in connection with this agreement.

The total liability of the Institute in connection with this agreement, whether for negligence or breach of contract or otherwise, shall in no event exceed 120% of any price paid or payable by the Recipient for the supply of the Goods.

If any of the Goods supplied by the Institute should prove not to meet their specification when stored and used correctly (and provided that the Recipient has returned the Goods to the Institute together with written notification of such alleged defect within seven days of the time when the Recipient discovers or ought to have discovered the defect), the Institute shall either replace the Goods or, at its sole option, refund the handling charge provided that performance of either one of the above options shall constitute an entire discharge of the Institute's liability under this Condition.

17. INFORMATION FOR CUSTOMS USE ONLY

Country of origin for customs purposes*: United Kingdom
* Defined as the country where the goods have been produced and/or sufficiently processed to be classed as originating from the country of supply, for example a change of state such as freeze-drying.
Net weight: 7mg
Toxicity Statement: Non-toxic
Veterinary certificate or other statement if applicable.
Attached: No