



529-0043-00L  
Analytical Strategy  
ETH Zurich

# Analysis of pesticides in food

Thomas M. Amrein  
Head of Trace Analysis  
Coop Central Laboratory  
CH-4133 Pratteln

## Agenda

- Central laboratory of Coop
- Department of trace analysis @ Coop
- What is important in the job?
- Question 1: Analyze: what? when? targets?
- Question 2: Quantiation: Ideal vs reality
- Question 3: Confirmation: How and when?  
+ examples
- Question 4: Retro-analysis: Potential problems?
- Question 5: Glyphosat: Analytical approach?
- Question 6: Assessment of residues

# Coop Central Laboratory

- Responsible for the analytical work and related questions for the whole Coop group.
- Together with quality assurance (total ca. 60 persons)  
→ whole QM-team in the same building
- Located in Pratteln (near Basel)
- Laboratory: 35 employees, 5 departments
- ISO 17025 accreditation
- No contract lab
- Broad range of analytical methods and samples (food, feed, non-food):  
→ 815 active methods covering 2467 parameters in e.g.
  - microbiology
  - GMO
  - food composition
  - physical & visual properties
  - food additives, vitamins
  - mycotoxins
  - **trace analysis**
  - non-food testing

11.12.2012

TMA

3

# Coop Central Laboratory



11.12.2012

TMA

4

## Department of trace analysis

- Staff: 8 (6 technicians, 2 academics)
- Instrumentation:  
LC-Q-ToF, 2 LC-MS/MS, GC-MS/MS, 2 GC-MS  
GC-FID, GC-PFPD, GC-sniff  
ICP-OES, GFAAS
- Main analyses: - veterinary drugs  
- pesticides  
- heavy metals  
- plasticizers  
- polycyclic aromatic hydrocarbons  
- illegal dyes  
- taints / off-flavours
- Main samples come from our competence labels, e.g.  
- organic fruits and vegetables  
- Naturafarm & Bio meat and fish  
- Naturaline textile

11.12.2012

TMA

5

## Range of samples



11.12.2012

TMA

6

## What is important in the job?

- Legal aspects
  - Multidisciplinary approaches
  - Explaining and presenting
  - People: understand them, motivate them, lead them
  - Efficiency: a lot is about time and money!
  - Network
  - Stay curious
- 
- And of course: knowledge!


11.12.2012

TMA

7

## Question 1

### Analyse: What? When? How?



- Known pesticides
- GC-MS/MS amenable
- LC-MS/MS amenable
- Group methods
- Single methods



11.12.2012

TMA

8

# Question 1

## Analyse: What? When? How?

- Residue data, check data bases:
  - how often is the compound found?
  - on what crops?
  - country of origin? Important for supply chain?
  - MRL exceedencies?
  - Are there registered applications for the compound?
- Authorisation of use? On what crops?
- RASFF alerts? <https://webgate.ec.europa.eu/rasff-window/portal/>
- Season vs climate vs susceptibility of crop
- Amenable to multi-methods? ↔ efficiency?!
- Long-term experience / own data
- Information from your network
- Public interest / pressure
- Relevance to company / economic factors
- Supplier: performance in preceeding years; new?
- ...
- ... .. but also do the unexpected!



11.12.2012

TMA

9

## Question 1: Examples

- Lettuce in January / February
- Strawberries from cold and wet summers
- Grapes from harvests with havy rainfall
- Vegetables from Thailand

### Lettuce, France, February 2011

CS46420 PestizideQuECHERS LC		
Pestizide QuECHERS LC	nw	
Acetamiprid	0.053	mg/kg
Iprodion	0.191	mg/kg
Propamocarb	3.46	mg/kg

CS46430 PestizideQuECHGC		
Pestizide QuECHERS GC	nw	
Bifenthrin	0.107	mg/kg
Cyprodinil	0.343	mg/kg
Fludioxonil	0.461	mg/kg
lambda-Cyhalothrin	0.0459	mg/kg
Pencycuron	0.0112	mg/kg
Pyrimethanil	0.456	mg/kg

CS46460*		
Pestizide QuECHERS GC	nw	
Folpet	0.139	mg/kg

### Hot-Chili, Thailand, March 2010

Methode	Resultat	Einheit	Coop-Norm	Toleranzwert	Grenzwert oder Textnorm
<b>Messgrösse</b>					
<b>Pestizide/Begasungsmittel</b>					
CS46420 PestizideQuECHERS LC					
Pestizide QuECHERS LC	nw				
Carbendazim	0.205	mg/kg		max. 0.100	
Methomyl	0.0352	mg/kg		max. 0.2000	
Profenofos	1.21	mg/kg			max. 0.05 (EU Höchstwert)
Summe Methomyl	0.0352	mg/kg		max. 0.2000 (EU Höchstwert)	
Ergebniskommentar zu Pestizide QuECHERS LC: Das Untersuchungsspektrum umfasst 85 Wirkstoffe Bestimmungsgrenze: 0.01 mg/kg Ergebniskommentar zu Summe Methomyl: Summe von Methomyl und Thiodicarb berechnet als Methomyl					
CS46430 PestizideQuECHGC					
Pestizide QuECHERS GC	nw				
Chlorpyrifos(ethyl)	0.171	mg/kg		max. 0.050	
Cypermethrin	0.0288	mg/kg		max. 0.5000	

11.12.2012

TMA

10

## Question 1: RASFF

**RASFF Portal**

Notifications list : 23 results

Search criteria | Notified from 01/01/2012 | Notified till 30/11/2012 | Product category fruits and vegetables | Hazard category pesticides residues | Origin country THAILAND (TH)

<< First << << Previous 100 << Notifications 1 to 23 of 23 >> Next 100 >> >> Last >>

Classification	Date of case	Last change	Reference	Country	Subject
border rejection	28/11/2012	28/11/2012	2012_OON	FR	tetradifon (0.03 mg/kg - ppm) in fresh subergines from Thailand
border rejection	13/11/2012	13/11/2012	2012_GGO	FR	carbendazim (4.8 mg/kg - ppm) in aubergines from Thailand
information for attention	13/11/2012	13/11/2012	2012_1370	NL	methomyl (0.10 mg/kg - ppm) in green papaya from Thailand
border rejection	18/10/2012	07/11/2012	2012_CFB	GB	thiacloprid in fresh okra from India
border rejection	30/10/2012	30/10/2012	2012_CZM	NO	carbofuran (0.075 mg/kg - ppm) in chilled coriander leaves from Thailand
border rejection	11/10/2012	11/10/2012	2012_BXZ	DK	diclofop (1.2 mg/kg - ppm) in fresh red chili from Thailand
border rejection	08/08/2012	08/08/2012	2012_BOG	FR	indoxacarb (0.047 mg/kg - ppm) in asparagus bean from Thailand
information for attention	30/07/2012	30/07/2012	2012_1088	FR	carbendazim (4.2 mg/kg - ppm) in longan from Thailand
information for attention	19/07/2012	19/07/2012	2012_1020	NL	carbendazim (2.2 mg/kg - ppm) in durian from Thailand
border rejection	04/06/2012	04/06/2012	2012_BEE	FI	methomyl (0.074 mg/kg - ppm), omethoate and dimethoate (sum: 0.095 mg/kg - ppm) in and unauthorised genetically modified green papaya from Thailand
information for attention	23/05/2012	30/05/2012	2012_0698	CH	profenofos (266 µg/kg - ppb) and omethoate (87 µg/kg - ppb) in small eggplants from Thailand
border rejection	21/05/2012	29/05/2012	2012_BCE	DK	carbendazim (0.28 mg/kg - ppm) in chili from Thailand
information for attention	22/05/2012	23/05/2012	2012_0694	CH	unauthorised substance dinotefuran (110 µg/kg - ppb) in chilled kale shoot from Thailand
information for attention	22/05/2012	22/05/2012	2012_0692	CH	unauthorised substance chlorfluazuron (30 µg/kg - ppb) in chilled kale from Thailand
information for attention	26/04/2012	04/05/2012	2012_0993	DE	omethoate (0.69 mg/kg - ppm) in okra (Abelmoschus esculentus) from Thailand
border rejection	25/04/2012	25/04/2012	2012_AXX	FR	dimethoate (0.17 mg/kg - ppm) in eggplants from Thailand
border rejection	07/03/2012	07/03/2012	2012_AOY	DE	dimethoate (0.91 mg/kg - ppm) in water mimosa from Thailand
border rejection	17/02/2012	23/02/2012	2012_ALN	DE	omethoate and dimethoate (sum: 1.039 mg/kg - ppm) in water mimosa from Thailand


11.12.2012

TMA

11

## Question 2: Quantitation

- **Analysis by LC-MS/MS and GC-MS/MS**  
some targets can only be analyzed either by LC or GC.
- **MS/MS mode:** A parent ion (e.g.  $[M+H]^+$ ) is selected in Q1, fragmented in Q2, and two specific daughter ions are monitored in Q3 (Selected Reaction Monitoring, **SRM**).  
→ exclusion of noise/matrix; high sensitivity, high selectivity
- **Identification:** a target is identified by its retention time, two specific mass transitions (SRMs), and their relative intensity.
- **Calibration of analytes and internal standards.**
- **Several hundreds of compounds can be monitored (MRM: multiple reaction monitoring)**



11.12.2012

TMA

12

## Question 2: Quantitation


- **IDEAL:** Use isotope-labeled standards for each compound
  - ↳ correction for losses during clean-up
  - ↳ correction for suppression effects in MS
  - ↳ very high costs (some targets are only rarely detected)
  - ↳ labeled standards are **not** available for all pesticides
  - ↳ more mass-transitions have to be monitored (400 → 800 MRMs)
    - dwell time per transition has to be lowered
    - or number of data points in peak becomes (too) low
  - not efficient, not feasible
- **REAL:** Use several internal standards (e.g. 3) that elute at different times during chromatography
  - ↳ some correction for losses during clean-up
  - ↳ cheap and efficient
  - ↳ limited number of extra mass-transitions
  - ↳ no reliable correction for suppression effects
  - feasible
  - quantitation via standard addition for more accurate results

11.12.2012

TMA

13

## Question 2: Quantitation



11.12.2012

TMA

14

## Question 3: Confirmation

### Options:

- Second analysis from the very beginning
- Change chromatography:  
LC ↔ GC, other stationary phase
- Change detector, e.g. from MS/MS to ToF with accurate mass
- Monitor more mass transitions and compare to standard
- Spiking experiment: recovery? Peak shape unchanged?
- Standard addition for more accurate quantitation  
→ matrix effects are compensated for

Plausibility of residue: does the result make sense?

11.12.2012

TMA

15

## Question 3: Further investigation

- Is a legal limit exceeded?  
with/without measurement uncertainty?
- Is the substance prohibited (e.g. DDT)?
- Is the result not plausible?
- Is the sample especially important,  
e.g. organic?
- Is there health concern?
- Is the application of the pesticide illegal  
on the present sample?
- Is the sample analysed not fully  
homogenous?
- Can a mistake in samples not be  
excluded?

one YES! is enough to  
make further investigation!

False positives and false negatives must be avoided! Which one is worse?

11.12.2012


TMA

16



## Question 3: Example

Example: GC-MS/MS Analysis indicates Bupropfen in Koriander.  
Sample is analysed by LC-MS/MS using 4 MRMs.



- Only the 4<sup>th</sup> mass transition reveals that the result is false positive!
- Intensities of the transitions differ between standard and sample (② : ③!)

11.12.2012

TMA


17



## The nasty case...



Tasmanian pepper with ethoxyquin?!  
Concentration above MRL! False positive?




11.12.2012

TMA


18

# Tasmanian pepper

Fullscan measurement with LC-Q-ToF  
Ethoxyquin: exact mass of  $[M+H]^+$  218.1539



Ethoxyquin




11.12.2012

19

# Tasmanian pepper

MS-MS experiment with LC-Q-ToF  
Fragments of  $m/z$  218



Some shifts of retention time are observed but are not unusual with spices.  
Is it really positive? If yes → no purchase!



11.12.2012

TMA

20


# Tasmanian pepper

Fullscan measurement with LC-Q-ToF  
 Ethoxyquin: exact mass = 218.1539  
 Extracted ion chromatograms (EIC)




# Tasmanian pepper

Fullscan measurement with LC-Q-ToF  
 Mass spectra of the suspected Ethoxyquin-peak



## Tasmanian pepper: more SRMs and spiking




There is an interference showing all 5 (!) mass transitions (SRMs) ....  
... but is there also some real ethoxyquin?

11.12.2012

TMA

23

## Tasmanian pepper: LC-Q-ToF fragments of m/z 218




11.12.2012

TMA

24

# Tasmanian pepper: finally a closer look: no ethoxyquin!!


LC-Q-ToF: fragments of m/z 218



11.12.2012

25

# Tasmanian pepper: but what is it?



Polygodial is a natural component of tasmanian pepper

Exact mass: [M+H]<sup>+</sup> = 235.1693

During ionization polygodial undergoes a spontaneous loss of water

Its <sup>13</sup>C isotope thereby generates m/z 218.1621 (ethoxyquin: 218.1539)

This ion produces virtually the same fragment spectrum as ethoxyquin!

11.12.2012

TMA

26

## Tasmanian pepper: finally, we sell it!

- Polygodial is a natural component of tasmanian pepper, concentration is in the %-range.
- Polygodial can simulate the presence of ethoxyquin.
- Falso positive results can be avoided by careful investigations using triple quad and/or Q-ToF technologies.
- A false positive result whould have stopped the purchase
- Polygodial causes a numbing / tingling sensation on the tongue and thereby contributes an important part to the product characteristics.



11.12.2012

TMA

27

## Question 4: potential problems in retro-analysis


- Was the compound of interest covered by the sample preparation?  
If not → risk of false-negatives!
- Quantitation: is the response of the detector well comparable between actual and former measurement?  
→ risk of over-/under-estimation
- If no reference standard is available and information on retention time is lacking: how can a "positive hit" be confirmed?  
→ risk of false-positives!
- Is the sample still available for confirmation analysis, including extraction?
- Is the analyte of interest stable in the sample?

11.12.2012

TMA

28

## Question 5: Glyphosate?



- Non-selective, systemic herbicide; Number 1 pesticide in the world!  
"kills everything green"
- First marketed in 1974 by Monsanto: RoundUp
- RoundUp Ready Soja: GMO-Soy from Monsanto, resistant to glyphosate
- Inzwischen gibt es zahlreiche Anbieter von Glyphosat-haltigen Herbiziden, z.B. Syngenta, Dow, etc.
- Application:
  - viticulture
  - fruits (e.g. blackberries, apples, ...)
  - against weeds on uncultivated land/der Brache
  - roadsides
  - siccation before harvesting (lentils, wheat, soy, ...)



11.12.2012



TMA




29

## Question 5: how to analyse glyphosate?

### Option 1:

- Aqueous extract
  - Anion exchange chromatography
  - MS/MS detection
  - Isotope labeled internal standards
- ☞ simple, cheap, fast
- ☞ a lot of co-extractives  
→ some matrices ruin column with 5 injections
- ☞ LOQ may be too high (0.01 mg/kg is the goal)



11.12.2012

TMA

30



## Question 5: how to analyse glyphosate?

### Option 2:

- Aqueous extract
- Derivatisation with FMOC-Cl
- Derivate is much less polar and can be enriched on a SPE-cartridge
- Classical C18-LC-MS/MS works fine
- Isotope labeled internal standards

- rather clean extracts
- low LOQ can be achieved (0.01 mg/kg)
- more expensive (time, chemicals, work)

→ Method of choice!



11.12.2012

TMA

31

## Question 7: Assessment of residues

- FIV: Fremd- und Inhaltsstoffverordnung  
<http://www.admin.ch/ch/d/sr/8/817.021.23.de.pdf>  
Swiss Maximum Residue Levels (MRLs)
- European MRLs:  
[http://ec.europa.eu/sanco\\_pesticides/public](http://ec.europa.eu/sanco_pesticides/public)
- Swiss pesticide database:  
[www.blw.admin.ch/psm/produkte/index.html?lang=de](http://www.blw.admin.ch/psm/produkte/index.html?lang=de)
- Pesticides allowed for application on fruit for "Suisse garantie"  
[www.swissfruit.ch/m/mandanten/239/download/2012\\_Saio\\_wirkstoffe\\_liste\\_d\\_komplett.pdf](http://www.swissfruit.ch/m/mandanten/239/download/2012_Saio_wirkstoffe_liste_d_komplett.pdf)

11.12.2012

TMA

32



## Question 7: Assessment of residues

Pesticide	Concentration mg/kg	Horwitz mg/kg	MRL mg/kg	Application allowed on blackberries?	Complaint?
Buprofezin	0.09	0.021	0.1	yes	no
Bifenthrin	0.67	0.11	0.3	No!	Yes!
Spinosad	0.02	0.0057	0.5	No!	Yes!
Folpet	2.24	0.32	3	yes	no
Cyprodinil	0.01	0.0032	10	yes	no

11.12.2012

TMA

33

# Thank you for your attention!

Acknowledgement to Thomas Döring  
for LC-Q-ToF measurements and data