

For a complete list of references, please visit antibiotics.toku-e.com

Microorganism Genus, Species, and Strain (if shown)	Concentration Range (µg/ml)
<i>Bacteroides distasonis</i>	≤0.06 – 256
<i>Bacteroides ovatus</i>	2 – >128
<i>Bacteroides thetaiotaomicron</i>	0.5 – 128
<i>Bacteroides uniformis</i>	≤0.06 – 128
<i>Bacteroides vulgatus</i>	1 – 128
<i>Bifidobacterium spp.</i>	2 – >128
<i>Clostridium bifementans</i>	0.25 – 2
<i>Clostridium cadaveris</i>	0.25 – 2
<i>Clostridium difficile</i>	2 – 4
<i>Clostridium perfringens</i>	0.5 – 2
<i>Clostridium sordellii</i>	0.25 – 2
<i>Clostridium spp.</i>	0.25 – 2
<i>Clostridium tertium</i>	0.25 – 2
<i>Escherichia coli</i>	0.125 – 16
<i>Eubacterium lentum</i>	2 – >128
<i>Fusobacteria</i>	0.06 – 2
<i>Fusobacterium mortiferum</i>	0.5 – 1
<i>Fusobacterium necrophorum</i>	0.06 – 0.25
<i>Fusobacterium nucleatum</i>	0.06 – 1
<i>Fusobacterium varium</i>	0.25 – 2
<i>Lactobacillus spp.</i>	2 – >128
<i>Peptostreptococcus anaerobius</i>	0.125 – 4
<i>Peptostreptococcus asaccharolyticus</i>	0.125 – 4
<i>Peptostreptococcus magnus</i>	0.125 – 4
<i>Peptostreptococcus tetradius</i>	0.125 – 4
<i>Porphyromonas</i>	0.125 – 64
<i>Porphyromonas asaccharolytica</i>	0.125 – 64
<i>Prevotella bivia</i>	0.5 – 16
<i>Prevotella buccae</i>	0.125 – 64
<i>Prevotella corporis</i>	0.25 – 4
<i>Prevotella disiens</i>	0.5 – 64
<i>Prevotella intermedia</i>	0.5 – 16
<i>Prevotella melaninogenica</i>	1 – 16
<i>Prevotella oralis</i>	0.125 – 64
<i>Prevotella oris</i>	0.125 – 64
<i>Propionibacterium acnes</i>	0.5 – 4
<i>Propionibacterium spp.</i>	0.5 – 4
<i>Pseudomonas aeruginosa</i>	256

The data above is sourced from The Antimicrobial Index. For further assistance, please contact us at info@toku-e.com or visit www.toku-e.com.