
# Product Information



# **DL**- $\alpha$ -Lipoic Acid

Item No. 10005728

CAS Registry No.:	1077-28-7
Formal Name:	1,2-dithiolane-3-pentanoic acid
Synonym:	Thioctic acid
MF:	$C_8H_{14}O_2S_2$
FW:	206.3
Purity:	≥95%
Stability:	≥2 years at -20°C
Supplied as:	A crystalline solid
UV/Vis.:	λ <sub>max</sub> : 203, 333 nm



# Laboratory Procedures

For long term storage, we suggest that DL-α-lipoic acid be stored as supplied at -20°C. It should be stable for at least two years.

 $DL-\alpha$ -lipoic acid is supplied as a crystalline solid. A stock solution may be made by dissolving the  $DL-\alpha$ -lipoic acid in an organic solvent purged with an inert gas. DL-α-lipoic acid is soluble in organic solvents such as ethanol, DMSO, and dimethyl formamide. The solubility of DL- $\alpha$ -lipoic acid in these solvents is approximately 30 mg/ml.

DL-α-lipoic acid is sparingly soluble in aqueous buffers. For maximum solubility in aqueous buffers, DL-α-lipoic acid should first be dissolved in ethanol and then diluted with the aqueous buffer of choice. DL- $\alpha$ -lipoic acid has a solubility of approximately 0.25 mg/ml in a 1:8 solution of ethanol:PBS (pH 7.2) using this method. We do not recommend storing the aqueous solution for more than one day.

DL-α-lipoic acid is a cyclic disulfide antioxidant that interconverts with its reduced dithiol form. It is an essential cofactor for decarboxylation reactions of the citric acid cycle, and acts as a general antioxidant.<sup>1-3</sup> DL- $\alpha$ -lipoic acid can act as a direct radical scavenger, as a cofactor to regenerate reduced glutathione, and as a metal chelator.

# References

- 1. Biewenga, G.Ph., Haenen, G.R.M.M., and Bast, A. The pharmacology of the antioxidant lipoic acid. General Pharmacology: The Vascular System 29(3), 315-331 (1997).
- Randle, P.J. Regulatory interactions between lipids and carbohydrates: The glucose fatty acid cycle after 35 years. Diabetes Metab. Rev. 14(4), 263-283 (1998).

3. Biewenga, G.Ph., Haenen, G.R.M.M., and Bast, A. The role of lipoic acid in the treatment of diabetic polyneuropathy\*. Drug Metab. Rev. 29(4), 1025-1054 (1997).

# **Related Products**

For a list of related products please visit: www.caymanchem.com/catalog/10005728

### WARNING: This product is for laboratory research only: not for administration to humans. Not for human or veterinary DIAGNOSTIC OR THERAPEUTIC USE.

#### MATERIAL SAFETY DATA

This material should be considered hazardous until information to the contrary becomes available. Do not ingest, swallow, or inhale. Do not get in eyes, on skin, or on clothing. Wash thoroughly after handling. This information contains some, but not all, of the information required for the safe and proper use of this material. Before use, the user must review the complete Material Safety Data Sheet, which has been sent via email to your institution.

#### WARRANTY AND LIMITATION OF REMEDY

Cayman Chemical Company makes no warranty or guarantee of any kind, whether written or oral, expressed or implied, including without limitation, any warranty of fitness for a particular purpose, suitability and merchantability, which extends beyond the description of the chemicals hereof. Cayman warrants only to the original customer that the material will meet our specifications te time of delivery. Cayman will carry out its delivery obligations with due care and skill. Thus, in no event will Cayman have **any obligation or liability**, whether in tort (including negligence) or in contract, for

any direct, incidental or consequential damages, even if Cayman is informed about their possible existence. This limitation of liability does not apply in the case of intentional acts or negligence of Cayman, its directors or its employees. Buyer's **exclusive remedy** and Cayman's sole liability hereunder shall be limited to a <u>refund</u> of the purchase price, or at Cayman's option, the <u>replacement</u>, at no cost to Buyer, of all material that

does not meet our specifications

Said refund or replacement is conditioned on Buyer giving written notice to Cayman within thirty (30) days after arrival of the material at its destination. Failure of Buyer to give said notice within thirty (30) days shall constitute a waiver by Buyer of all claims hereunder with respect to said material.
For further details, please refer to our Warranty and Limitation of Remedy located on our website and in our catalog.
Copyright Cayman Chemical Company, 06/18/2012

# Cayman Chemical

# **Mailing address**

1180 E. Ellsworth Road Ann Arbor, MI 48108 USA

Phone (800) 364-9897 (734) 971-3335

Fax (734) 971-3640

#### E-Mail

custserv@caymanchem.com

# Web

www.caymanchem.com