

IMMPORT DATA SUBMISSION USER GUIDE



IMMPORT

BIOINFORMATICS FOR THE FUTURE OF IMMUNOLOGY

Version 2.28

Bioinformatics Integration Support Contract (BISC), Phase III

Contents

1.0 Highlights of ImmPort Versions	5
2.0 Introduction	5
2.1 Data Submission Guide Organization.....	5
2.2 What are the Data Submission options?	5
2.3 Quick Start	5
3.0 ImmPort Research Data management and Submission for New Users	7
3.1 ImmPort Research Data Model.....	8
4.0 Metadata Upload templates or Data Submission Templates	10
4.1 Spreadsheet Based Metadata upload File Features.....	12
4.1.1 Required and Optional Columns	13
4.2 ImmPort Accessions	13
4.3 Overview of ImmPort Data Submission Templates	14
5.0 Protocols.....	15
6.0 Study Design	16
7.0 Subjects.....	17
7.1 subjectsHuman.xls	18
7.2 subjectsAnimal.xls.....	18
8.0 Experiment Samples.....	18
9.0 Treatments	20
10.0 reagents.....	20
11.0 Lab Tests.....	21
12.0 Assessments and Adverse Events.....	21
13.0 Experiments.....	22
14.0 BioSamples	22
15.0 Lab Test Panels.....	23
16.0 Review of the ImmPort Metadata Files.....	23
17.0 Derived Data Templates	24
17.1 Flow Cytometry Data Report.....	24
17.2 HLA Typing Summary Report	24
17.3 MBAA Results Report.....	25
17.3.1 Control Sample	26
17.3.2 Standard Curve	27
17.3.3 Experiment Sample MBAA	27
17.3.4 MBAA Results Template.....	27
18.0 Data Submission Package	28
18.1 Save the Spreadsheet-Based File as Tab-Delimited Text.....	28

18.2	Creating an ImmPort Data Submission Package.....	29
18.2.1	ImmPort Metadata upload Tab-Delimited Format File Names (after conversion from .xls files when saving)	29
19.0	Data Submission Package Validator	30
20.0	Send the Data Submission Package	32
21.0	Data Submission History	34
APPENDIX 1	ImmPort Template Description File	1
APPENDIX 2	Saving MBAA primary Bead Level Data	1
Section 2.01	Extracting bead level data from xPONENT	1
Section 2.02	Extracting bead level data using MasterPlex CT	2
Section 2.03	Extracting bead level data using Bio-Plex Manager	3
APPENDIX 3	Affymetrix GCOS/MAS 5 Gene Expression	3
APPENDIX 4	Illumina Bead Studio Final Report	5
APPENDIX 5	Glossary.....	6
APPENDIX 6	FAQs (Frequently Asked Questions)	9

BISC Documentation Style Guide Version History

Version	Date	Description
2.28	May, 2015	<ul style="list-style-type: none"> - Combined templates for experiment samples, biosamples and experiment. - Combined templates for lab test panel and lab tests. - New template for assessments. -
2.27	December, 2014	<ul style="list-style-type: none"> - No changes to upload templates.
2.26	October, 2014	<ul style="list-style-type: none"> -New templates available for Lab Test panel and Lab Test. -Experiments must link to a study. -Enhanced support for standardized cell population name and gating definitions in derived flow cytometry template.
2.25	July, 2014	<ul style="list-style-type: none"> - New templates for Adverse Events and Substance Merging with Subject -Human Subject Population column displays http://allelefrequency.net values.
2.24	June, 2014	<ul style="list-style-type: none"> -Significant re-work of study, bioSample, experimentSample, reagent templates. - New data dictionary is available with templates to describe templates, columns, rows, lists, data types, etc.
2.19	May, 2013	<ul style="list-style-type: none"> - Changed TREATMENT to a required Experiment Sample attribute. -Only one biological sample can be linked from an experiment sample. - Added MBAA result support enhancements - Discontinued support for XML template submission.
2.18	February 2013	<ul style="list-style-type: none"> -Inclusion/exclusion criteria, publication references added to study design template -Links from experiment samples to protocols removed
2.17	December 2012	<ul style="list-style-type: none"> -Upgrade subject demographics -Support HAI and Neutralizing antibody assay results
2.13	June 2012	<ul style="list-style-type: none"> -References to arms from subjects and from bioSamples to studies required. Modified subject and sample templates. Created new assay specific experiment sample templates.
2.8	04/28/2010	<ul style="list-style-type: none"> - Study day as a required field in the Biosamples.xls - On all edit pages, the delete button removed.
2.7	02/02/2010	<ul style="list-style-type: none"> - Changed 'gender' to a required SUBJECT attribute. - Removed status from templates, changed affection_phenotype to subject_phenotype and made this a required field. - Allowed multiple subject phenotypes to be entered into the metadata template using a standard separator to be defined. - Require a SUBJECT link from a BIOSAMPLE (BS), Experiment Sample (ES) link to a Biosample (BS). - Added upld_tct_num column in treatment_info table. - Allow users to batch upload Flow text files for subsequent analysis
2.3.2	09/30/2008	<ul style="list-style-type: none"> - Added support for ELISA, ELISPOT, MBAA results per the DAIT Minimum Information Guidelines. - Modified the experiment samples and reagent templates for ELISA, ELISPOT, MBAA results. - Created new results templates for ELISA, ELISPOT and MBAA results. - Modified the user guide and tutorials to support ELISA, ELISPOT and MBAA results.

Version	Date	Description
2.3	07/23/2008	Made revisions and incorporated DAIT Minimum Information Standards.
2.2	03/19/2008	Improvements made on the Data Submission process upload to allow greater than 100MB of data. Modified the ImmPort user interface. Implemented an off-line Data QC pipeline for HLA typing results.
2.1.2	December 2007	No changes to Data Submission.
2.1.1	November 2007	Data Submission Validator Beta to check for errors before submitting packages into ImmPort.
2.1	August 2007	Added the ability to update research metadata for subjects, experiments, biological samples, experiment samples, reagents and protocols.
2.0	May 2007	Implemented generation of Linkage Pedigree (pre MAKEPED) format (.ped and .info) files from uploaded genotype results. Enabled upload of derived flow cytometry results.
1.4.1	December 2006	No changes to Data Submission.
1.4	October 2006	Created tutorials for TagSNP, Data Submission. Enabled download of result file(s).
1.3.1	June 2006	No changes to Data Submission.
1.3	June 2006	Batch submission of genotyping data from Illumina panels in two formats.
1.2.1	May 2006	XML based system for collecting descriptions of experiment results.
1.2	April 2006	No changes to Data Submission.
1.1	January 2006	No changes to Data Submission.
1.0	October 2005	Initial release of the ImmPort website.

1.0 HIGHLIGHTS OF IMMPORT VERSIONS

ImmPort releases updates on the data capture process to reflect changes in assay technology and improvements to support data sharing and re-use. These changes can include converting optional descriptive elements to required attributes, support new file formats, and adding new descriptive data concepts.

For more detailed information on each of the versions indicated, please refer to the ImmPort website templates history at <https://immport.niaid.nih.gov/immportWeb/experimental/templateHistory.do>

2.0 INTRODUCTION

2.1 DATA SUBMISSION GUIDE ORGANIZATION

Welcome to the ImmPort Data Submission User Guide. The ImmPort system supports the National Institute of Allergy and Infectious Diseases Division of Allergy, Immunology, and Transplantation mandate to facilitate storage, sharing and analysis of research data and metadata (descriptive information about data). This guide is intended as an introduction to ImmPort and specifically the data submission process. It is meant to supplement, but not replace interaction with the ImmPort staff. The user should contact the ImmPort staff for the preparation of the data for submission and using the analytical tools.

This guide is organized into the following sections

- A Quick Start Guide
- An Overview of the ImmPort Research Data Model and Data Submission process
- An Explanation of the Metadata upload File Templates
- Supplemental Information

2.2 WHAT ARE THE DATA SUBMISSION OPTIONS?

Contact the ImmPort staff for help (BISC_Helpdesk@niaid.nih.gov) with organizing, packaging, and loading of the data into ImmPort.

Read this guide then contact Helpdesk to ask questions and provide feedback and suggestions.

Submit data to ImmPort and contact the Helpdesk (BISC_Helpdesk@niaid.nih.gov) with any comments, questions, or recommendations.

2.3 QUICK START

Quick Start is intended for the experienced user who needs few hints or reminders on how to assemble and submit data to ImmPort.

If the user wants to take advantage of ImmPort's analysis, querying and sharing tools, please go to [Step 1](#) of the data submission process in ImmPort and download the metadata upload templates that you need.

Complete the metadata upload templates required for your submission. In each template enter detailed information about one or more samples used in the experiment. Please complete all data entry categories that are not shaded gray.

When using .xls files for data entry save the files as tab delimited text files before zipping into an archive file.

Use a ZIP utility (e.g. 7-Zip, Ziplt or WinZIP) to put the files and completed metadata upload files into a .zip file.

If the user wants to store results or assay data only, use a utility to put the files into a .zip archive and proceed to Step 2 of the Data Submission Process.

Go to [Step 2](#) of the Data Submission process in ImmPort (figure 2a) to install and run the Data Submission Validator.

Figure 1a. Step 2 Check the Data Submission Package with the Validator tool.






Figure 1b. ImmPort Data Submission Validator online validation



3.0 IMMPORT RESEARCH DATA MANAGEMENT AND SUBMISSION FOR NEW USERS


ImmPort is currently available to NIAID/DAIT funded researchers providing them with the ability to store, analyze, and share their research data. ImmPort also provides access to public reference data bases allowing researchers to utilize some of this information for their analyses (figure 2a).

Figure 2a. ImmPort’s capabilities



While ImmPort offers this wide array of tools and data sharing capabilities, research data initially submitted by an investigator is placed into a private work space that is only available to that investigator and any other personal they choose. This work space can contain a number of different studies and allows an investigator and their staff to take full advantage of ImmPort’s tools and shared data bases while their research data is not available to other registered users (figure 2b).

Figure 2b. ImmPort’s Management of Data Access



The first step in the data submission process is to visit ImmPort’s data submission page and download metadata template files. The template files are available in .xls format. After downloading the template files, the ImmPort data submission process can be summarized as:

- Enter your data into the metadata upload template(s)
- Generate a .zip format archive file that includes the completed metadata upload files and any other files that are to be submitted (e.g. assay results files, protocol files, assay platform annotations, etc.)
- Send the .zip file via a secure transmission method to the ImmPort system.

A more complete description of the steps involved in the data submission process follows.

3.1 IMMPORT RESEARCH DATA MODEL

The ability of ImmPort's tools to effectively annotate research data requires the capture of descriptive details at many levels in a study. To understand how this is achieved, the user should keep in mind the following definitions.

Metadata – collections of files containing detailed information and descriptions of the data

Templates – are .xls files for the structured listing of the detailed information and descriptions of the data

Transfer / Upload- submit, send or transmit data to ImmPort

Research Data / Assay Results – actual experiment data submitted by registered users

Reference data – data available from a number of public sources

ImmPort's model for handling research data is to partition information into seven metadata categories. Each category has its own .xls file or set of files for the entry of different types of information associated with an experiment. These categories are also called metadata classes in this guide.

The metadata categories and their functions are as follows:

Study: Studies provide the context and organization of a research effort. Studies organize subjects into groups (e.g. arms or cohorts) based on phenotype and/or treatment. The encounter schedule provides a guide as to the temporal relationship between samples and encounters (i.e. is the sample taken before or after a treatment).

Subjects: Subjects may be patients or animals from which samples are taken for analysis. Two .xls templates (one for human and one for animal subjects) are available for recording subject information. In these files, treatment protocols used on the subjects can also be listed as well as many other details. Subjects are assigned to a single group (arms or cohorts) within a study and maybe linked to multiple studies.

Biological Samples: Describe the types of samples taken from subjects or cell culture and processed for the experiment (i.e. organs, tissue, blood, plasma, cell culture name, etc.), when the samples were taken in the course of a study and protocols used in the sample collection, processing, and/or treatment. Samples are linked to a single study.

Experiments: Describe the type of experiment, measurement technique and links to protocols used in the experiment.

Experiment Samples: The biological samples analyzed in an experiment are linked to the assay reagent, protocol and results via the experiment sample record. Since ImmPort offers the ability to submit data for a number of immunological, proteomic and genomic analysis methods, several different template files are available for listing sample details for different experiment approaches.

Treatments: Describe the experimental conditions for specific biosamples or experiment samples. Treatments link to experiment samples as well as biosamples if needed.

Protocols: Describe the methods and procedures in studies, subject recruitment/treatment, sample collection/preparation/treatment and experiments. Protocols may be PDF files, Word documents, excel or other file types.

Reagents: Provide detailed information about the reagents used in an experiment. Since different analysis platforms employ very different reagents, the reagents.xls file contains several tabs for entering platform specific details of the reagents that were used.

One of the important properties of these metadata categories is that the details listed in these files allow archived data to be linked by common information. For example, several different biological samples can be taken from a subject and used for different assays. If the experiments for these assays have been properly formatted and submitted to ImmPort, a user can search their private and/or shared databases for

what experiments and samples came from each subject to select data from a broader number of subjects or samples to perform a specific analysis. The linkage of archived data by the different metadata categories facilitating this capability is diagrammed in figure 3.

Lab Test Panels: Provide a means to group clinical lab tests and links to study and protocols.

Lab Tests: Provide a means to describe and capture clinical lab test results.

Figure 3a. ImmPort Research Data Model Classes

How the metadata classes are organized and linked in the ImmPort Research Data Model.





Figure 3b. ImmPort Data Model Alternative View

Illustration of the relationships between research data components



4.0 METADATA UPLOAD TEMPLATES OR DATA SUBMISSION TEMPLATES

Metadata upload templates enable the user to provide detailed information about how a study was conducted in a structured format for ImmPort to parse and store for subsequent queries or analysis.

Metadata templates for entering your data can be downloaded from the ImmPort Data Submission web page at [“Step 1: Download and Fill Templates”](#) (figure 4a). It is strongly recommended to download templates prior to submitting a data package to ensure use of the most current template file versions. The order in which metadata files are completed by the user depends on the data submission goals of the research team and what information they have already submitted to ImmPort.

The file ImmPort.Upload.Templates.Description.pdf is a complete listing of all templates, their usage and content.

The templates are available in a spreadsheet-based .xls format or in a tab separated text format. **It should be noted that after entering your data into spreadsheet based .xls files, the files need to be saved as tab-delimited text files before they are zipped into a submission package and sent to ImmPort.**

Because metadata upload files are used to define and describe metadata items (e.g. samples, protocols, experiments, etc.) and define links or associations between metadata (e.g. linking a sample preparation protocol to a biological sample or experiment), the primary identifier used when defining and linking metadata is an item's User Defined ID. User Defined IDs are intended to be short names for each metadata item and are the primary means to link metadata in a specific project from different submitted experiments. As a result **the User Defined ID for a metadata item must be unique within the scope of a research project.** A Research Project generally refers to a private workspace created by the Principal Investigator (PI) to be used by an individual or by a group researchers collaborating on a given project. The Project created will be used to store study descriptions, experiments, results, and any accompanying data.

Some of the metadata upload files have several columns to list the User-Defined IDs for different items involved with linking data together (see figures 7a, 7b and 9a for examples). As a result, for subsequent data submissions with common samples, protocols or other details involved with data linkage, **the User Defined IDs (or ImmPort Accession number) and other identifiers must match exactly and are case sensitive.**

Examples and additional rules involved with User Defined IDs and data linkage will be presented in the following sections describing the metadata templates.

Figure 4a Step 1: Download and Fill in templates with your data

Data Submission / Resource / Data Submission Templates

Submit Data | Submission History | Resources -

Submit Data Main Page → **Step 1: Download and Fill Templates** → Step 2: Check Data in .zip file → Step 3: Send Data in .zip file → Step 4: Review Submission Status & Results

- Which Data Submission Templates do you need?
Please contact us by email at BISC_Helpdesk@niaid.nih.gov.
The [User Guide](#) is a reference you can use to determine which templates need to be completed.
- Complete the templates that are needed.
Note: Please save spreadsheet .xls templates as tab delimited .txt files.
- Create a .zip file that contains the files you want to submit (e.g. results, protocols, bioSamples template, experimentSamples template, etc.).
- Please check that you are using the [latest version](#) of the ImmPort data transfer templates.

ImmPort Research Data Class	Purpose	Spreadsheet Template	Required Data to Complete Metadata Form	Latest Version/ Date Available
Basic Study	Describes a study in terms of title, goals, endpoints, criteria for study participation, subject grouping (arms or cohorts), personel, planned visits or encounters and protocols using a single worksheet. A study design should be uploaded first.	basic_study_design.xls	<ul style="list-style-type: none"> Study User Defined ID Title Description Endpoints Inclusion and Exclusion Criteria Arms or Cohorts Personnel Period Planned Visit Protocol 	February 2013
Complete Study	The multiple worksheet workbook that enables a detailed description of a study (especially useful for clinical studies).	CompleteStudyDesign-V1.6.1.xls	<ul style="list-style-type: none"> Study User Defined ID Title, Description, Endpoints, Inclusion and Exclusion Criteria, Arms or Cohorts, Personnel, Period, Planned Visit, Protocol 	September 2011

4.1 SPREADSHEET BASED METADATA UPLOAD FILE FEATURES

The first three rows of the spreadsheet-based .xls metadata upload file contain file formatting structure and explanatory information for the user (e.g. comments). The column headers are used to ensure that data entered by the user is properly processed. **Do not edit or delete these column headers.**

There are drop lists with controlled vocabulary for some columns (e.g. Measurement technique in experiments.xls). These are activated by clicking on the data entry cell in the column.

Note that whenever there is a comma with text in any cell, upon conversion to a .txt file Microsoft Excel will add a double quote at the beginning of the text and a double quote at the end of the text. This could result in a failed data submission and will not allow ImmPort to properly use this data.

Upon completion of entering data into these templates, they must be saved as a tab-delimited text file with its ImmPort file name. As a result, for each data submission it is advisable to create a folder containing all the blank metadata templates you will need for a specific submission. The completed templates in this folder can then be zipped into one file for submission as described in sections 15 and 16.

4.1.1 Required and Optional Columns

Please see the ImmPort.Upload.Templates.Description.pdf document for a description of the template columns.

4.2 IMMPORT ACCESSIONS

ImmPort assigns accession numbers to key research metadata elements to ensure the global uniqueness of identifiers in ImmPort. The ImmPort accessions are an alternative means to identify User-Defined IDs to reference metadata items (e.g. reagents, protocols, or samples) already stored in ImmPort. When linking metadata information either the ImmPort accession OR the user-defined ID is used to make the link, but not both.

The naming convention for ImmPort accessions has a prefix based on the data category (Table 4a) and a numerical integer generated by ImmPort. For example, if you submit an experiment with one sample and a protocol, ImmPort accession numbers could look like; EXP3428 for the experiment, BS116067 for the biological sample and PTL2187 for the stored protocol document.

ImmPort accessions assigned to previously uploaded and stored data can be obtained from a data submission report, the experiment data management pages or from experiment data queries. More information on how to access this information can be found in the instructions on searching research data.

Table 4a. ImmPort accession prefix abbreviations

ImmPort Accession	Description
BS	Biological Sample (BioSample)
EXP	Experiment
ES	Experiment Sample (ExpSample)
PTL	Protocol
REA	Reagent
SDY	Study
SUB	Subject
TRT	Treatment
LTP	Lab Test Panel
LT	Lab Test
AE	Adverse Event
AP	Assessment Panel
AC	Assessment Component
SM	Substance Merge

4.3 OVERVIEW OF IMMPORT DATA SUBMISSION TEMPLATES

The ImmPort templates used to capture descriptions of research data types/concepts are laid out together below (Figure 4d) to provide a means of understanding how they relate to each other and how they are organized.

Figure 4d ImmPort Template Overview

subjects		Schema Version 2.28	Content Version 2.28												
Please do not delete or edit this column.															
Column Name	Subject User-Defined ID	Subject Treatment Protocol ID(s)	Study Arm Or Cohort ID	Gender	Age	Age Unit	Age Event	Age Event Specify	Subject Phenotype						
assessments		Schema Version 2.28	Content Version 2.28												
Please do not delete or edit this column.															
Column Name	Assessment Panel ID	Study ID	Assessment Name	Assessment Type	Is Standard	Status	CFR File Names	Result Separator Column	Assessment Component User-Defined ID	Planned Visit ID	Subject ID	Component Name			
labtests		Schema Version 2.28	Content Version 2.28												
Please do not delete or edit this column.															
Column Name	Biological Sample ID	Lab Panel ID	Study ID	Protocol ID(s)	Source Subject ID	Planned Visit ID	Biological Sample Type	Biological Sample Sub-Type	Biological Sample Name	Biological Sample Description	Study Time Collected	Study Time Collected Unit			
protocols		Schema Version 2.28	Content Version 2.28												
Please do not delete or edit this column.															
Column Name	Protocol User-Defined ID	Protocol File Name	Protocol Name	Protocol Summary	Protocol Type										
elisa		Schema Version 2.28	Content Version 2.28												
Please do not delete or edit this column.															
Column Name	Reagent User-Defined ID	Reagent Name	Reagent Description	Manufacturer	Catalog Number	Lot Number	Web Link	Contact	Analyte Name						
elisa		Schema Version 2.28	Content Version 2.28												
Please do not delete or edit this column.															
Column Name	Experiment Sample ID	Biological Sample ID	Experiment ID	Reagent ID(s)	Treatment ID(s)	Experiment Sample Name	Experiment Sample Description	Additional Result File Names	Study ID	Protocol ID(s)	Source Subject ID	Planned Visit ID			

Table 4b is an overview of the data formats to be provided for assay results. Please contact ImmPort staff for questions and comments about assay result formats.

Table 4b. Summary of the assay specific Experiment Samples templates and results formats specified for commonly used immunological research methods

Assay	Results Format	Controls	Compensation	Standard Curves
ELISA	ELISA_Results.txt			
ELISPOT	ELISPOT_Results.txt			
Flow Cytometry	.fcs format FCM_derived_data.txt	control samples .fcs format	compensation samples .fcs format	
Gene Expression	Affymetrix CEL Illumina GEO raw NCBI GEO accession			
Genotyping	Affymetrix CEL Illumina Bead Studio dbGAP accession			
Hemagglutination Inhibition	HAI_Results.txt			
HLA	HLA_Typing.txt			
Image Histology	Custom			
KIR	KIR_Typing_Results_Column.txt or KIR_Typing_Results_Column.txt			
Mass Spectrometry	Custom			
Multiplex Soluble Protein	MBAA_Results.txt	Control_Samples.txt		Standard_Curve.txt
qRT-PCR	PCR_Results.txt			
Virus Neutralization	Virus_Neutralization_Results.txt			

5.0 PROTOCOLS

Listing Protocols allows a user to provide detailed descriptions of methods related to subject treatment, sample collection, sample processing and other experimental methods. Protocol files may be in the form of PDFs, spreadsheets, word processor documents or other file formats.

Protocols can be linked from:

- Study
- Subjects
- Biological Samples
- Experiments

Note: Protocols MUST be linked to Studies, Biological Samples, Experiments and Subject(s). This is a requirement on DAIT minimum information standards.

The protocols.xls template allows the user to define and describe protocols that can be referenced or linked by other metadata items (e.g. biological samples)

Upon completion of entering data into this template, it should be saved as a tab-delimited text file and must be saved with its ImmPort name of protocols.

Figure 5a. protocols.xls template

	A	B	C	D	E	F
1	protocols	Version 2.24				
2	Please do not delete or edit this column					
3	Column Name	Protocol User-Defined ID	Protocol File Name	Protocol Name	Protocol Summary	Protocol Type
4						
5						
6						

The file ImmPort.Upload.Templates.Description.pdf is a complete listing of all templates, their usage and content.

6.0 STUDY DESIGN

There are two versions of the study design template: basic and complete. The basic study design template is recommended for the majority of ImmPort data upload efforts and differs from the other ImmPort templates in that it has a more vertical or "form" like organization. This is a consequence of placing several domains of study design on a single form. The row and column headers in **bold text** are required elements. There are comments for row and column headers that can be viewed by placing the mouse cursor over the column header name.

Figure 6. basic_study_design.xls template

basic_study_design		Schema Version 2.28
Please do not delete or edit this column		
Column Name		
study		
User_Def_Study_ID		
Brief_Title		
Official_Title		
Study_Type		
Brief_Description		
Detailed_Description		
Hypothesis		
Objectives		
Endpoints		
Sponsoring_Organization		
Target_Enrollment		
Condition_Studied		
Minimum_Age		
Maximum_Age		
Age_Unit		
Actual_Start_Date		
Intervention_Agent		
arm_or_cohort		
User_Def_Arm_Or_Cohort_ID	Name	
study_personnel		
User_Def_Study_Personnel_ID		
Honoric		
Last_Name		
First_Name		
Suffixes		
Organization		
Email		
Title_In_Study		
Role_In_Study		
Site_Name		
period		
User_Def_Period_ID		
Title		
Order_Number		
planned_visit		
User_Def_Planned_Visit_ID	Visit_Name	

Studies link to:
 Protocols
 and may link to:
 Files
 Publications

Studies are linked from:
 Subjects
 Biological Samples
 Experiments
 Lab Test Panel
 Assessment Panel
 Substance Merge

The file ImmPort.Upload.Templates.Description.pdf is a complete listing of all templates, their usage and content.

7.0 SUBJECTS

There are human and animal specific subject templates.

Subjects link to:
 Protocols
 Study Arms or Cohorts

Subjects are linked from:
Biological Samples
Assessment Components
Substance Merges

Note: When subjects are used they MUST link to a Protocol(s), and Study Arm or Cohort. This is a requirement of DAIT minimum information standards. If multiple protocols should be linked, a semi-colon (;) may be used to separate the user-defined of each protocol listed in a cell.

The file ImmPort.Upload.Templates.Description.pdf is a complete listing of all templates, their usage and content.

7.1 SUBJECTSHUMAN.XLS

The subject information to be completed includes basic demographic and the arm or cohort to which the subject is linked in a study (figure 7). Completed files need to be saved as tab-delimited txt files before compressing into a zip archive.

Figure 7. subjectsHuman.xls template

subjectshuman	Schema Version 2.28	Content Version 2.28								
Please do not delete or edit this column										
Column Name	Subject User-Defined ID	Subject Treatment Protocol ID(s)	Study Arm Or Cohort ID	Gender	Age	Age Unit	Age Event	Age Event Specify	Subject Phenotype	P

7.2 SUBJECTSANIMAL.XLS

This template is similar to the subjectHuman.xls but excludes human specific descriptive columns and includes laboratory animal specific columns including Strain Name and Strain Characteristics. For a detailed description of the column headers, please refer to the Glossary.

8.0 EXPERIMENT SAMPLES

An experiment sample is a biological sample that is treated in vitro in some way (including not treatment) and then assayed. The experiment sample template defines and annotates the assay results for a sample by linking samples, experiments, and results together. The template enables you to define a new experiment and/or a new biological sample. as well as a new experiment sample. The template supports multiple data annotation scenarios:

- Define a new experiment, a new biological sample and a new experiment sample
- Define a new experiment, link to a biological sample defined in ImmPort and a new experiment sample
- Link to an experiment defined in ImmPort, define a new biological sample and a new experiment sample

Whenever a new record is defined, all of the required attributes for defining a record must be completed.

Experiment Sample template files are used to name and describe the control and experimental samples used in an assay. To ensure accurate recording of information for a number of different platforms and assay methods, there are several different templates to choose from for specific applications.

For some assay methods, the assay results are also recorded in the experiment sample template by completing the set of columns for assay results (e.g. for HAI, complete Virus Strain and Titer). To enter multiple assay results for an experiment sample, simply copy and paste the set of columns needed to describe the results as many times for each assay results. For example, if an experiment sample was assayed for HAI for three viral strains, there would be three sets of HAI assay result columns. The experiment sample template for ELISA, ELISPOT, HAI, Virus Neutralization, and qRT-PCR enable capturing the assay results in the experiment samples template.

The experimentSamples assay-specific .xls templates are available for downloading on the data submission home page. The first three rows and column headers of these files are used to ensure that data entered by the user is properly processed. **Do not edit or delete these column headers.** See ImmPort.Upload.Templates.Description.pdf, for a complete description or definitions of the columns and elements.

The User-Defined IDs for the experiment sample, experiment, protocol, and reagent are meant to be simple names selected by the user for these items. As an alternative to the User Defined IDs, ImmPort accession numbers can be entered if available from previously stored data. If you are entering User defined IDs, you do not have to enter ImmPort accession numbers.

In addition to these required columns, the “Result File or Folder Name” (column K) must also be entered. This option allows ImmPort to find the necessary result file that is needed by some tools while other tools may require several files to be stored in a specific folder to be accessed by ImmPort.

Experiment Samples link to:
 Biological Samples
 Reagents
 Experiments
 Result file(s)
 Treatments

Note: Experiment Samples MUST link to an Experiment, Reagent (one or more), Biological Sample (**only one**), and Result File (one or more).

Please note that only one biological sample can be linked to an experiment sample.

The experiment sample file names help enhance the accuracy of the assay annotation. Please choose the most appropriate template for your data upload.

Figure 8. expSamples.ELISA.xls template

elisa	Schema Version 2.28	Content Version 2.28										
Please do not delete or edit this column												
Column Name	Experiment Sample ID	Biological Sample ID	Experiment ID	Reagent ID(s)	Treatment ID(s)	Experiment Sample Name	Experiment Sample Description	Additional Result File Names	Study ID	Protocol ID(s)	Source Subject ID	Planned Visit ID

Figure 12b highlights an example where the experiment sample template is tailored to capture assay specific information. Flow Cytometry assay runs include control results and compensation fcs result files. These should be linked to the fcs file from the sample(s) that were assayed in a run. All experiment

samples from the same flow cytometry run should be linked to the same set of compensation and control files.

The file `ImmPort.Upload.Templates.Description.pdf` is a complete listing of all templates, their usage and content.

9.0 TREATMENTS

Treatments are required for Experiment Samples and are optional for Biological Samples. The treatment template (figure 9) allows the user to define the amount, duration and temperature of treatments specific to individual Experiment Samples or Biological Samples. **Do not edit or delete these column headers.**

Treatments are linked from:
Experiment Samples.

Figure 9. View of an empty Treatment template

	A	B	C	D	E	F	G	H	I
1	treatments	Version 2.24							
2	Please do not delete or edit this column								
3	Column Name	Treatment User-Defined ID	Treatment Name	Treatment?	Amount Value	Amount Unit	Duration Value	Duration Unit	Temperature Value
4									
5									

The file `ImmPort.Upload.Templates.Description.pdf` is a complete listing of all templates, their usage and content.

10.0 REAGENTS

The reagents to be described to ImmPort are those used in an assay to find or characterize one or more analytes. Reagents can be microarray chips for gene expression or genotyping, fluorochrome conjugated antibodies for flow cytometry, ELISA or ELISPOT antibodies, etc.

There is a reagents template for each assay type supported by ImmPort (e.g. Figure 11a). Each template is organized for entering information specific for each assay method. This organization facilitates the capture of assay specific reagent information that allows users to search for data in ImmPort.

Figure 10a. View of an ELISA reagent template

	A	B	C	D	E	F	G	H	I	J	K	L
1	elisa	Version 2.25										
2	Please do not delete or edit this column											
3	Column Name	Reagent User-Defined ID	Reagent Name	Reagent Description	Manufacturer	Catalog Number	Lot Number	Web Link	Contact	Analyte Name	Analyte Description	Antibody Registry ID
4												

The format of the reagent templates is assay specific. As an example, the reagents-HLA worksheet (figure 11b) is used to describe reagents associated with HLA typing results and has a slightly different format than the ELISpot template. On the ImmPort web site, there are two additional files under example submission packages, the HLA_Typing_System_Allele_Ambiguity file and the HLA_Typing_System_Features file. These are intended to further annotate the HLA typing system that was employed

Figure 10b. View of reagents-HLA template

	A	B	C	D	E	F	G	H	I	J	K	L
1	hla_typing_system	Version 2.25										
2	Please do not delete or edit this column											
3	Column Name	HLA Typing System User-Defined ID	Locus Genotyped	Exons / Intron Interrogated	Typing System Manufacturer	Typing System Version	Lot Number	Web Link	Contact	Typing Method	HLA Typing Date	Last Pertinent Who Report/Update
4												

11.0 LAB TESTS

Lab Test template files are used to name and describe clinical lab tests and capture the test result. Preferred values for lab test names from CDISC (cdisc.org) and SnoMED CT (www.nlm.nih.gov/snomed/) are provided to encourage use of standard terms.

Lab Tests link to:
Biological Sample
Lab Test Panel

Note: Lab tests MUST link to a biological sample and Lab Test Panel.

labtests	Schema Version 2.28	Content Version 2.28										
Please do not delete or edit this column												
Column Name	Biological Sample ID	Lab Panel ID	Study ID	Protocol ID(s)	Source Subject ID	Planned Visit ID	Biological Sample Type	Biological Sample Sub-Type	Biological Sample Name	Biological Sample Description	Study Time Collected	Study Time Collected Unit

The file ImmPort.Upload.Templates.Description.pdf is a complete listing of all templates, their usage and content.

12.0 ASSESSMENTS AND ADVERSE EVENTS

The Assessment template captures results from Case Report Forms and similar questionnaires in a standard format. This template is intended for human subjects.

The Adverse Event template records adverse events experienced by study subjects.

The file ImmPort.Upload.Templates.Description.pdf is a complete listing of all templates, their usage and content.

13.0 EXPERIMENTS

This is a legacy template that will be supported for backward compatibility, but its function is replaced by the Experiment Sample template. Experiment template files are used to provide detailed information about the type of analysis that was done and how it was done. This allows the entry of user defined definitions of names for the experiment, protocols and other information.

The experiments.xls template file is available on the data submission page and is standard for all submissions. The first three rows and column headers are used to ensure that data entered by the user is properly processed. **Do not edit or delete these column headers.**

Experiments link to:
Protocols
Study

Experiments can be linked from:
Experiment Samples
Control Samples
Standard Curves

Note: Experiments MUST link to a single Study and one or more Protocols. This is a requirement of DAIT and ImmPort minimum information standards. If multiple protocols should be linked, a semi-colon (;) may be used to separate the user-defined of each protocol listed in a cell.

The Experiment Purpose and Measurement Technique can be selected from dropdown menus by selecting the data entry cells for this information.

Figure 8. experiments.xls template

1	experiments	Version 2.24							
2	Please do not delete or edit this column								
3	Column Name	Experiment User-Defined ID	Experiment Name	Description	Experiment Purpose	Measurement Technique	Protocol ID(s)	Hypothesis	Rationale
4									
5					CD8_T_cell_recall				
6					Cellular_Activity				
7					Cellular_Phenotype				
8					Cellular_Quantifica				
9					Cytokine_Quantifica				
10					Enzyme_Activity				
					Epitope_Mapping				
					Genotyping				

The file ImmPort.Upload.Templates.Description.pdf is a complete listing of all templates, their usage and content.

14.0 BIOSAMPLES

This is a legacy template that will be supported for backward compatibility, but its function is replaced by the Experiment Sample template. ImmPort templates for biological samples allow the user to describe the characteristics of biological material (i.e. blood, sera, tissue, cell preparation, etc.) that has undergone isolation, processing and/or treatment prior to use in an experiment. In order to effectively describe a biological sample, as much information about the sources, conditions, and treatments should be captured in the template. Links to sample sources (e.g. subjects or other biological samples) and protocols can be defined.

The bioSamples.xls template file is available on the data submission page and is standard for all submissions. The first three rows and column headers are used to ensure that data entered by the user is properly processed. **Do not edit or delete these column headers.**

Biological Samples link to:
Study
Subjects
Other Biological Samples
Protocols.

Biological Samples are linked from:
Experiment Samples
Lab Tests

Note: Biological Samples MUST link to a Protocol(s) and Study. This is a requirement of DAIT and ImmPort's minimum information standards. If multiple protocols should be linked, a semi-colon (;) may be used to separate the user-defined of each protocol listed in a cell.

Fields required on this template are in white.

The file ImmPort.Upload.Templates.Description.pdf is a complete listing of all templates, their usage and content.

15.0 LAB TEST PANELS

This is a legacy template that will be supported for backward compatibility, but its function is replaced by the Lab Test template. Lab Test Panel template files are used to name and describe groups of clinical lab tests. Preferred values for Lab Test Panel names from CDISC (cdisc.org) and SnOMED CT (www.nlm.nih.gov/snomed/) are provided to encourage use of standard terms.

Lab Test Panels link to:
Study
Protocols

Lab Test Panels are linked from:
Lab Tests

Note: Lab test panels MUST link to a Study and one or more protocols.

The file ImmPort.Upload.Templates.Description.pdf is a complete listing of all templates, their usage and content.

16.0 REVIEW OF THE IMMPORT METADATA FILES

In the preceding sections describing spreadsheet-based .xls metadata upload files, examples of completed templates were displayed.

Upon completion of entering data into these templates, they should be saved as tab-delimited text files and they must be saved with their ImmPort name (original ImmPort file names are listed in section 15.2.1). As a result, for each data submission package it is advisable to create a folder containing all the blank metadata templates you will need for a specific submission. The templates in this folder can be edited to contain the information and data needed for the submission before saving them as text files. The completed .txt versions of the templates in this folder can then be zipped into one .zip file for submission.

17.0 DERIVED DATA TEMPLATES

NIAID DAIT strongly encourages the use of ImmPort derived data templates if there are applicable templates available. These templates facilitate the annotation and organization of assay results in a consistent fashion that allow the full set of descriptive data to be associated with the assay results. The derived results templates include:

FCM_derived_data.xls
HLA_Typing.xls
KIR_Typing
MBAA_Results.xls

the following templates are legacy formats that can be used to add results to existing experiment samples.

ELISA_Results.xls
ELISPOT_Results.xls
HAI_Results
PCR_Results
Virus_Neutralization_Results

17.1 FLOW CYTOMETRY DATA REPORT

Flow Cytometry (FCM) derived data refers to the analysis results of .fcs file contents. As with the metadata upload files, the columns that must be completed are white and optional columns are shaded gray. The FCM_derived_data.xls template is available from [Step 1 of Data Submission](#).

The required columns in the file are columns B, C, D, E, F, and G. The Experiment Sample User_Defined ID is used to link the derived data with the result file and the descriptions of how the data was generated and analyzed. The additional required columns capture details about analyzed FCM results that would be useful for queries such as population name and gating combination.

Figure 17b. The ImmPort FCM analyzed results template

fcml_derived_data	Version 2.24						
Please do not delete or edit this column							
Column Name	Experiment Sample ID	Population Name	Population Definition (gating combination)	Population Cell Number	Population Cell Number Unit	Workspace File	Comments

17.2 HLA TYPING SUMMARY REPORT

The ImmPort HLA_Typing.xls template supports capturing the expert determined values for the HLA genes. The HLA_Typing.xls template format follows a commonly used structure where the HLA loci symbols are the column headers (figures 14b and 14c). The locus symbols are from the IMGT/HLA curated data repository (www.ebi.ac.uk/imgt/hla/). The commonly typed HLA loci are always displayed and the data provider can add additional loci symbols to the template by selecting from the drop down list in the "Other" column header. Note that two adjacent columns must have the same locus name. The typed values for each locus can be entered in any format as "free text". An experiment sample ID is used in the template to link the HLA data with other results and metadata. This is the only required field. The

commonly typed HLA loci symbols are included as column headers. Additional HLA locus symbols are available from the dropdown list in the “Other” columns.

Figure 17c. View of HLA_Typing.xls template

HLA Typing Results	2.24									
Please do not delete or edit										
Column Name	Experiment Sample User-Defined ID*	Population Area*	HLA-A Allele 1	HLA-A Allele 2	HLA-B Allele 1	HLA-B Allele 2	HLA-C Allele 1	HLA-C Allele 2	HLA-DPA1 Allele 1	HLA-DPA1 Allele 2

Figure 17d. Completed version of HLA_Typing.xls template


HLA Typing Results	2.24									
Please do not delete or edit										
Column Name	Experiment Sample User-Defined ID*	Population Area*	HLA-A Allele 1	HLA-A Allele 2	HLA-B Allele 1	HLA-B Allele 2	HLA-C Allele 1	HLA-C Allele 2	HLA-DPA1 Allele 1	HLA-DPA1 Allele 2
	HLA typing report 1	North America	A*0101	A*0101	B*0702	B*0702				
	HLA typing report 2	North America	A*0101	A*0101	B*0702	B*0702				

17.3 MBAA RESULTS REPORT

The Multiplex Bead Array Assays (MBAA) data capture process was developed in coordination with the Human Immunology Project Consortium (www.immuneprofiling.org) Data Standards Working Group. In addition to the experiment sample to results linkage that is common for all ImmPort results, there are two additional metadata or descriptive entities: Control Samples and Standard Curves.

Figure 17e. MBAA Data Model

This diagram highlights the relationship between experiments, samples, standard curves and MBAA results.



Figures 17f-i (below) highlight the structure of the templates used to capture MBAA results..

17.3.1 **Control Sample**

Control samples are often used for quality control and to allow comparison between assay runs or across labs. These control samples are distinct from biological samples since they are often purchased in bulk, and are not linked to particular study subjects. See Appendix 1.02 for more details.

Figure 17f. Control Samples

Control Samples	Version 2.24							
Please do not delete or edit this column								
Column Name	Control Sample User-Defined ID	Experiment User-Defined ID*	Source	Catalog ID	Dilution Factor	Assay ID	Result File Name	Assay Group ID

17.3.2 Standard Curve

This is a new template for ImmPort. A standard curve is defined for each analyte in a batch (e.g., on a single plate). See Appendix 1.07 for more details.

Figure 17g. Standard Curves

Standard C: Version 2.24								
Please do not delete or edit this column								
Column Name	Standard Curve User-Defined ID	Experiment User-Defined ID	Formula	Analyte Name	Assay ID	Upper Limit	Upper Limit Unit	Lower Li

17.3.3 Experiment Sample MBAA

The expSamples.MBAA.txt template captures details such as Assay ID and Dilution Factor that are unique for this assay.

Figure 17h. MBAA Experiment Sample

	A	B	C	D	E	F	G	H	I	J	K	L	M	N	O
1	mbaa	Schema Version 2.24	Content Version 2.24												
2	Please do not delete or edit this column														
3	Column Name	Experiment Sample User-Defined ID	Experiment Sample Name	Experiment sample Description	Experiment ID	Biological Sample ID	Reagent ID(s)	Treatment ID(s)	Result File Name	Import Template?	Assay ID	Dilution Factor	Additional Result File Names	Assay Group ID	Plate Type
4															

17.3.4 MBAA Results Template

MBAA_Results.txt template captures the individual measurements of the analyte(s). There are three sources for results: experiment sample (assay run for a biological sample), control sample, and standard curve. The results file from Luminex or MSD is the primary output from the assay platform. It may consist of various formats (image, image interpreted to text, etc). One of the "commonly" seen formats is a tabular representation of the results where plate ID and plate location are critical identifying attributes. The type of material assayed (biological sample, control sample, standard curve) and the analyte assayed are provided by annotation mapping files. Since there are three sources for results, it is essential to indicate the source type when describing the source ID.

Figure 17i. MBAA Derived or Interpreted Results

	A	B	C	D	E	F	G	H	I	J	K
1	mbaa_results	Schema Version 2.24	Content Version 2.24								
2	Please do not delete or edit this column										
3	Column Name	Source ID	Source Id Type	Assay ID	Assay Group ID	Analyte Name	Mfi	Concentration Value	Concentration Unit	Mfi Coordinate	Comment
4											

In addition to the derived results, the bead level result files should also be included and linked to experiment sample, control sample and standard curve records as appropriate. See Appendix 4 for details about saving MBAA bead level result files.


18.0 DATA SUBMISSION PACKAGE

An ImmPort data submission package may include completed metadata upload files, derived results templates and other files. The spreadsheet based templates must be saved as tab delimited text files before being compressed together into an ImmPort data submission package. For .xls file submission the files are packaged together into an archive .zip file format. When submitted to ImmPort, the contents of the package are processed to store the results and make them available for the analysis, query, or sharing of data in ImmPort. If the uploading and processing is successful or fails, the status of the upload can be viewed on the Submission History page.

18.1 SAVE THE SPREADSHEET-BASED FILE AS TAB-DELIMITED TEXT

The .xls based templates must be saved as tab-delimited text files by selecting the “Save As” option under the “File” Menu bar option (figure 15a).

Figure 18a. Selecting ‘Save As’ from the ‘File’ menu list



In the pop up window opened by Save As, open the “Save as type” drop down list, scroll down the list and choose the “Text (Tab delimited) (*.txt)” document type (figure 15b)





Figure 18b. Selecting "Save as type" in a spreadsheet as tab delimited text

After selecting the file type as .txt, press the save button. A dialog box will appear regarding the conversion of the file to a .txt format (figure 15c). Confirm that you wish to save the file in the text format by selecting "Yes" in the dialog box.

Figure 18c. Saving a spreadsheet as text dialog box



18.2 CREATING AN IMMPORT DATA SUBMISSION PACKAGE

An ImmPort data submission package is a .zip file that includes all of the files the data provider wishes to send (e.g. protocols, results, metadata, derived results, etc).


It is important that the completed metadata file names in the submission package are not modified from their original ImmPort assigned names. As a result, for each data submission package it is advisable to create a folder containing all the blank metadata templates you will need for a specific submission. The templates in this folder can be edited to contain the information and data needed for the submission before saving as .txt files. The completed .txt versions of the templates in this folder can then be compressed into one .zip file for submission.

18.2.1 ImmPort Metadata upload Tab-Delimited Format File Names (after conversion from .xls files when saving)

basic_study_design.txt
 biologicalSamples.txt
 experiments.txt
 experimentSamples.[Assay Name].txt
 protocols.txt
 reagents.[Assay Name].txt
 subjects.txt
 treatments.txt

Only one copy of each of the completed metadata files should be included in the .zip file of the data submission package (figure 15d).

Figure 18d. Example of an archive .zip file showing contents and structure



19.0 DATA SUBMISSION PACKAGE VALIDATOR

The ImmPort Data Submission Package Validator is available from Step 2 of ImmPort Data Submission and is used to check for compliance with ImmPort data submission norms and identify potential conflicts with data previously submitted to ImmPort.

Figure 19a. ImmPort Data Submission Package Validator access

IMMPORT BIOMINFORMATICS FOR THE FUTURE OF IMMUNOLOGY

[Edit Profile](#) | [Sign Out](#)

[Home](#) [About ImmPort](#) [Admin](#) [Access Data](#) [Tools](#) [My Workspace](#) [Resources](#) [News & Events](#)

Data Submission / Data Submission Package Validator

[Submit Data](#) | [Submission History](#) | [Resources](#) ▾

[Submit Data Main Page](#) → [Step 1: Download and Fill Templates](#) → [Step 2: Check Data in .zip file](#) → [Step 3: Send Data in .zip file](#) → [Step 4: Review Submission Status & Results](#)

Advantages of using the Data Submission Package Validator

- ▶ The Data Submission Package Validator checks the format and content of the files in a data upload package and reports issues.
- ▶ The Validator is a downloadable software application that is run on your computer from a web browser.
- ▶ If there are any errors in the validation report, fix the first errors that are indicated and rerun the validation process. The errors encountered first may have effects on later files that are validated.

To get started, click [here](#) to select the appropriate installer/uninstaller.

The [Data Submission Package Validator Tutorial](#) includes additional explanations.

Figure 19b. ImmPort Data Submission Package Validator


IMMPORT BIOMINFORMATICS FOR THE FUTURE OF IMMUNOLOGY

[Tools](#) ▾ [About](#)

Welcome to the ImmPort Data Submission Package Validator

The Validator is designed to check data submission packages for compliance with ImmPort data submission norms and identify potential conflicts with data already submitted to ImmPort.


Run Validation against ImmPort

 **Run Validation against ImmPort**

This option requires you to be logged into ImmPort. You also need to be connected to the Internet. When the Run Validation button is clicked, the Project Data will automatically be downloaded to your computer and the Validator Application will validate the Data Submission Package File that was submitted against the downloaded Project Data and a validation report will be generated and displayed. Also, if the checkbox for "Do you want to save the Project Data File for later use?" is checked then on the report page you will get a Download Project Data link. If you click on this link you can save the Project Data File for using it in the Run Validation against local copy mode.


1. Select the Project that you want to Validate against.
2. Select the Data Submission Package File. (Please see the [ImmPort Submission Data User Guide](#) for more details) Check the "Do you want to save the Project Data File for later use?" if you want to save the Project Data file on your computer.
3. Click the Run Validation button.
4. The Validation report will be displayed. You can download the report by clicking on the Download Validation Report File link.

View Report

 **View Report**

This option will display the last report that was generated by the Validator Application.


Run Validation against local copy

 **Run Validation against local copy**

This option does not require you to be logged into ImmPort or to be connected to the Internet. The Pre-requisite for this mode is that you need to have the Project Data File downloaded to your computer. When the Run Validation button is clicked, the Validator Application will validate the Data Submission Package File that was submitted against the Project Data File and a validation report will be generated and displayed.

1. Select the Project Data File that you want to Validate against.
2. Select the Data Submission Package File. (Please see the [ImmPort Submission Data User Guide](#) for more details)
3. Click the Run Validation button.
4. The Validation report will be displayed. You can download the report by clicking on the Download Validation Report File link.

Download Project Data File

 **Download Project Data File**

This option requires you to be logged into ImmPort. You also need to be connected to the Internet. When the Download Project Data button is clicked, the Project Data will automatically be downloaded and a Download Project Data File link will be provided so that the Project Data can be saved on your computer. This Project Data File can be used in the Run Validation against local copy option mode.

1. Select the Project that you want to download the data for.

Running the Validator can be done in two modes: via internet connection against project data in ImmPort or run against a local copy of project data. The application checks for missing files, improperly formatted

metadata files, references to other metadata within the data package, verifies references to records in the ImmPort database, confirms the integrity and structure of .zip files, database column width checks, and guidance on metadata content utilizing ImmPort analysis tools. It is highly recommended to use the Validator before sending data into ImmPort as part of the data submission process.

Details on the installation and use of the Data Submission Package Validator are available in the user guide and tutorial that accompany the Validator.

20.0 SEND THE DATA SUBMISSION PACKAGE

Once files are compressed into a .zip archive file, proceed to Step 3 of the Data Submission process “Send Data in .zip file” by clicking on Step 3 in the Data Submission schematic (figure 17a). Selecting this option in the work flow diagram at the top of the page brings you to a page where the data provider must select a research project into which the data will be submitted (figure 17b). The user can also select if the data will be submitted electronically or if it is over 1 Gigabyte, it can be shipped to ImmPort for loading (table 17a).

Figure 20a. Step 3: Choose the Research Project and select Upload option

Submit Data
Main Page

Step 1:
Download and Fill Templates

Step 2:
Check Data in .zip file

Step 3:
Send Data in .zip file

Step 4:
Review Submission Status
& Results

Research data will be stored in the private project workspace of your research project. Please select the Research Project and Contract-Grant number with which your data will be associated and enter any comments in the Notes text box. If you need assistance in understanding what a Research Project is and whether you are associated with one, please contact the Help Desk.

ImmPort accepts data submission packages in the "ZIP" format. (Please see the Data Submission User Guide for more information on how to create a .zip file.)

Please DO NOT include spaces in the .zip file name.

If your .zip file is less than or equal to 1 Gigabyte, please use the Submit on-line option.

- ▶ Check the "Upload Online" check box.
- ▶ Select the .zip file using the browse button.
- ▶ Click submit and record the registration number you receive.

If your .zip file is greater than 1 Gigabyte, please use the Submit off-line option.

- ▶ Uncheck the "Upload Online" check box.
- ▶ Click submit and record the registration number you receive.
- ▶ Contact the help desk (BISC_Helpdesk@niaid.nih.gov) with the registration ticket to arrange for file delivery.

If this is your first time entering data, you will need the following:

- ▶ A set of protocol documents delineating your SOPs
- ▶ Reagent list
- ▶ A list of types of samples (blood vs lymphocytes vs serum, etc.)
- ▶ A system for uniquely defining each sample such that sample type, date of sample, type of analysis can be discerned

Please check that you are using the [latest version](#) of the ImmPort data transfer templates.

Research Project Title: VIP004 - Vaccination Cohort

Upload Notes:

Upload Online
 Upload By Aspera
 Upload Offline

Browse... No file selected.

For electronic uploading of data less \leq 1 Gigabyte, select the Research Project and Grant/Contract number with which your data will be associated and enter any comments in the Notes text box. The data provider must select a research project to which the data will be submitted. If there is no project to select, a research project must be created by the Principal Investigator or the project manager. ImmPort staff can also help with this task as requested. If there are multiple projects listed, please contact the Principal Investigator or the project manager to decide which project to select. The user can enter any extra information to be saved in the Notes text box.

After selecting the file to upload, click the submit button to start file transfer. Make sure you record the registration number you receive for the transfer. Depending on the size of the file being transferred it can take several minutes to several hours for the transfer to be completed.

If uploading transfer packages larger than 1 Gigabyte, the data cannot be transferred by the standard uploading method. In this case uncheck the upload online box (figure 17b), select the file to transfer and hit the submit button to receive a registration number for the data transfer. The next step is to contact ImmPort staff at BISC_Helpdesk@niaid.nih.gov to arrange for your data to be uploaded.

Table 17a. Summary of the data submission process (Online and Offline)	
Files less than or equal to 1 Gigabyte Can be submitted online	Files greater than 1 Gigabyte Cannot be submitted online
Check upload online check box Select the file using the browse button Click submit Record the registration number you receive	Uncheck upload online check box Click submit Record the registration number you receive Contact the helpdesk (with registration number) to arrange for file delivery

21.0 DATA SUBMISSION HISTORY

When the data submission package .zip file has been sent to ImmPort, the web page will redirect the user to the Data Submission History page (figure 18a). The History page lists the data packages that have been sent to the research project and the status of the submission. Upon completion or termination of the transfer, the user will receive a confirmation email. The email, as well as the Data Submission History page, will indicate if the transfer was successful or not. If you have any problems or other questions about the submission process, ImmPort staff can be reached at BISC_Helpdesk@niaid.nih.gov to answer any questions or resolve any problems you encounter.

Figure 21a. Step 4: Review Submission Status & Results - Data Submission History

The screenshot shows the 'Step 4: Review Submission Status & Results' page. At the top, a workflow diagram indicates the process: Submit Data Main Page → Step 1: Download and Fill Templates → Step 2: Check Data in .zip file → Step 3: Send Data in .zip file → Step 4: Review Submission Status & Results.

Below the diagram, there is a filter section: 'Filter By Project: IP2 FCM DUL' with a 'Set Filter' button. A note says: 'Click on a Ticket Number below to view a submission's details. Ticket number in brown are marked for deletion.'

The main content is a table with 6 columns: Ticket Number, .ZIP File Name, Format, Status, Submitter, and Submit Date. The table lists 83 items, with the first 25 displayed. The first row shows a failed submission for 'kongme_20080721_2829' with file 'bioSamples_GSE9650.zip'. Other rows show various 'pat-daitFund' tickets with 'BS-2FCSfilesNoSpaceFileNamesDerived.zip' and 'reagents.zip' files, some failed and some completed.

At the bottom, there is a legend explaining the status codes: Pending, Started, Completed, and Rejected. A footer note says: 'Need help? Please feel free to contact us by email at helpdesk@immport.org'.

APPENDIX 1 IMMPORT TEMPLATE DESCRIPTION FILE

The ImmPort Template Description file describes each template, the column and rows in the template, the data type for each field in the template, database column length, preferred vocabulary terms, etc. This file supersedes the appendices describing the templates.

APPENDIX 2 SAVING MBAA PRIMARY BEAD LEVEL DATA

Section 2.01 EXTRACTING BEAD LEVEL DATA FROM xPONENT

Extracting bead level data from xPONENT

Renan Sauteraud

March 14, 2013

Contents

1 General 1

2 Saving the raw bead events 1

2.1 Documentation 1

2.2 HowTo 1

2.3 xPONENT 1.x 2

2.4 xPONENT 3.x 2

3 References 2

1 General

Luminex's own software, xPONENT, outputs raw bead level data in csv (comma separated values) files.

2 Saving the raw bead events

2.1 Documentation

From xPONENT's documentation:

One can obtain the FI for each bead within each sample from a so-called rCSV file. To configure your instrument to automatically generate the rCSV files, simply go to the Admin page/CSV Options Tab and make sure the bottom-most checkbox is checked. If checked and a run is performed, the rCSV files will be located under:

C:/docs&settings/all users/application data/Luminex/xponent/batchXXX rcsv

2.2 HowTo

To save the individual bead fluorescence:

Admin section

rsautera@fhcrc.org

1

CSV Options Tab

Check 'Automatically export results CSV file when batch is complete'

Check 'Include Advanced Statistics'

Choose the destination folder in the 'CSV File Export Folder' field.

Save

The files will be located in the destination folder, in a directory EXP NAME rcsv.

2.3 xPONENT 1.x

In the rst version of xPONENT, the filenames are Run[001-096].csv.

2.4 xPONENT 3.x

In the latest versions, the filenames use the following format: EXP NAME WELL ID.csv, with well id being A1-H12.

3 References

xPONENT 3.1 Software Manual:

<http://www.luminexcorp.com/prod/groups/public/documents/lmncorp/89-00002-00-202.pdf>

Obtaining FI of each bead within each sample:

<http://www.luminexcorp.com/blog/its-all-about-the-stats/>

Section 2.02 EXTRACTING BEAD LEVEL DATA USING MASTERPLEX CT

Extracting bead level data using MasterPlex CT

Renan Sauteraud

March 7, 2013

Contents

1 General 1

2 Saving the raw bead events 1

3 Locating the files 2

4 Summary file 2

5 References 2

1 General

Classical statistical analysis of xMap, and notably MiraiBio's own analysis software, 'MasterPlex QT', only uses the MFI. Therefore, saving the bead level information is set as an option in the acquisition software 'MasterPlex CT' and needs to be enabled.

2 Saving the raw bead events

In order to save the bead level information, some options should be checked before starting the acquisition.

In the 'Home' tab

Click the Setup 'Modify' button to open the 'Acquisition Setup' dialog box.

rsautera@fhcrc.org

1

'Setup' tab

v1.0 In the Output section, check 'Save binary file also.'

v1.2.0.7 In the Output section, check 'Save copy of output file in the output folder' and 'Save copies of the individual well run files in the output folder'

Select a 'Plate ID' and an 'Output Folder'

OK

3 Locating the files

The files are located in the selected 'Output Folder' and are based on the given 'Plate ID'.

PLATE ID.lxd is the summary file used by the vendor's analysis software. It contains the calculated MFIs as well as information on the Setup of the acquisition and the analytes used. PLATE ID WELL.lxb are the binary files containing the raw bead events required to run LumiR.

4 Summary file

Unlike with other acquisition softwares, LumiR make use of the summary file produced by MasterPlex CT. If the .lxd file is provided, the package will extract the information regarding the matching of the bead ID with the analyte name. See the LumiR's User guide for more information.

5 References

MasterPlex CT 1.0 User's Manual:

<http://www.miraibio.com/download-document/masterplex-ct-v1-user-s-manual.html>

MasterPlex CT Tutorial:

<http://www.miraibio.com/download-document/tutorial-for-masterplex-ct.html>

2

Section 2.03 EXTRACTING BEAD LEVEL DATA USING BIO-PLEX MANAGER

Extracting bead level data using Bio-Plex Manager

Renan Sauteraud

March 6, 2013

Contents

1 General 1

2 Saving the raw bead events 1

3 Exporting to XML 2

3.1 Export a single file 2

3.2 Export as a routine 2

4 References 2

1 General

This acquisition software from BioRad uses a proprietary file format that cannot be read directly by LumiR. These files with a .rbx extension should be exported in XML in order to be read by read.experiment.

2 Saving the raw bead events

The first step is to make sure that the bead level information is saved by Bio-Plex. In the 'Save As' dialog box, when saving protocols (extension .pbx) and results (.rbx), the compression mode must be disabled (default). Compression mode removes the raw bead event from the saved file in order to save space.

rsautera@fhcrc.org

1

3 Exporting to XML

For Bio-Plex data, LumiR requires a single .xml file per plate.

3.1 Export a single file

To export an existing file into XML:

Open a .rbx file

File Menu > Document Export

Choose 'Bio-Plex XML'

Choose The destination File

Export

3.2 Export as a routine

To export automatically at the end of a run:

File Menu > Document Export Properties

Choose 'Bio-Plex XML'

Check 'Use Folder' in the Destination section and Choose the desired folder.

OK

Create or load a Protocol

Advanced Settings Button

Check 'Auto save after run' and 'Auto XML export after run'

OK

4 References

Bio-Plex Manager 6.1 User Guide:

<http://www.bio-rad.com/webroot/web/pdf/lsr/literature/10022815.pdf>

2

APPENDIX 3 AFFYMETRIX GCOS/MAS 5 GENE EXPRESSION

MAS5.0 software generates six output files

- Chip File (.CHP)
- Expression Analysis File (.TXT)
- Report File (.RPT)
- Experiment Information File (.EXP)
- Cell Intensity File (.CEL)
- Data file (.DAT)

ImmPort requires a .CEL file for Affymetrix results. The other file types may be submitted as an additional result files.

A. Data File (.DAT)

This file contains raw image data of the chip.

B. Cell Intensity File (.CEL)

The .CEL file contains fluorescence intensities for each probe on the microarray. When the .CEL file is opened in either MAS 5.0 or dChip, these probe-specific intensity values are used to reconstruct the scanned image of the hybridized array. It is recommended that the investigator view the .CEL images for each sample to make sure there are no obvious chip defects. The probe-specific intensities in the .CEL file are also used in the Statistical Algorithm to calculate the probe-set-level Signals and Presence Calls recorded in the .CHP file.

C. Chip File (.CHP) and Expression Analysis File (.TXT)

The .CHP file contains Signal values and Presence Calls for each probe set on the microarray and is processed from the .CEL file. The .TXT file is a text version of the .CHP file with the same information in a tab-delimited format. The columns in the .TXT file are listed in the table below

D. Report File (.RPT)

The Report File summarizes background noise, housekeeping information, and spiked-in controls in the following format for the chip of interest. The file includes information such as:

- Filename
- Probe array type
- Algorithm
- Background
- Noise
- Total probesets and average signal
- Present probesets and average signal
- Absent probesets and average signal
- Marginal probesets and average signal
- Housekeeping controls
- Spike controls

E. Experiment Information File (.EXP)

The Experiment Information File is a text file containing sample information, fluidics settings and scanner settings. It includes information such as:

- Sample Info
 1. Chip type
 2. Sample type
 3. Description
 4. Project
- Fluidics
 1. Protocol
 2. Date of hybridization
- Scanner

1. Filter
2. Scan date
3. Scanner type
4. Number of Scans

APPENDIX 4 ILLUMINA BEAD STUDIO FINAL REPORT

The Bead Studio report format is highly customizable in terms of content (i.e. which columns of data are included) and format (e.g. tab delimited or matrix). The required Final Report columns and their headers to be included in an Illumina Bead Studio Final Report are listed below. These columns may appear in any order within the file. Additional columns may be included in the Final Report.

Sample ID SNP Name Allele1 – Top Allele2 – Top GC Score

An extract from an Illumina Final Report in tab delimited form with the required columns is presented below.

[Header]

BSGT Version 2.3.15.14130
Processing Date 5/25/2006 13:33
Content GS000XXXX-OPA
Num SNPs 1152
Total SNPs 1152
Num Samples 96
Total Samples 96

[Data]

SNP Name Sample ID Allele1 – Top Allele2 – Top GC Score
rs977008 1445888_R001_C001 G G 0.8092

Sample Identifiers

The genotype results analytical support is envisioned to allow the investigator guided creation of linkage format files and marker information files that are acceptable to algorithms such as Haploview (<http://www.broad.mit.edu/mpg/haploview>) which are part of the ImmPort genetic analysis tool set. In order to generate linkage format files, it is necessary to associate the genotype results for a sample with the affection status, gender and pedigree information for the individual or organism. To ease the task of generating these links, ImmPort supports designating whether the sample identifiers in the Sample ID column of the result file are subject identifiers or biological sample identifiers.

Annotation

The .opa manifest file describes an Illumina custom panel and serves as the annotation file. The .opa manifest file can be included in a data package and is declared in the reagents document.

An extract from an .opa manifest file is presented below for reference.

© Illumina, Inc. 2004

[Heading]

OPA, GS0006526-OPA
Test Version, 1227160
Test Version Name, 1227160
Test Name, Internal MHC OPAs
Bundle Index, 1 of 2
SNP Position, 1
Assay Format, Golden Gate
Date Manufactured, 2/1/2005
Loci Count, 1228 of 2521

Comment,
 Ilmn ID, Name, oligo1, oligo2, oligo3, IllumiCode name, IllumiCode Seq, Ilmn Strand, SNP, CHR, Ploidy,
 Species, MapInfo, TopGenomicSeq, CustomerStrand
 rs1003878-121.1_T_F_1155586, rs1003878,
 ACTTCGTCAGTAACGGACGACAGGCAGCTTACCAATAGGTCCTA,GAGTCGAGGGTCATATCGTGACA
 GGCAGCTTACCAATAGGTCCTG,
 AGTTTGCACTAAAAAGAAATTCAGTCCCAATCGAGCGGATAACATGTCTGCCTATAGTGAGTC,
 0003, CTGCCAATCGAGCGGATAACAT, top,
 [A/G],
 6, diploid, Homo sapiens, 32407800,
 TAACCCAGATACTTAAACAGTTGGAGCCAGTCTCCTTCAGACATAGTAAGAAGCCAGTAGAGATAAG
 TTGATATATACAGGCAGCTTACCAATAGGTCCT
 [A/G]
 GAGTTTGCACTAAAAAGAAATTCAAATTGGCATATTAGTACAGTATTGGAGAGTGTATTTTCACTA
 ATTTTATCCTAGAAGTGAGGCTTTGAGAGGT,
 BOT

APPENDIX 5 GLOSSARY

D-1 Glossary Table		
Term	Definition	Examples
Biological Sample	Biological material that has undergone isolation, processing and/or treatment prior to use in an experiment. As much information about the sources, conditions, and treatments of the Biological Sample as possible should be captured and parsed into the database. Biological Sample attributes include all treatments and reagents used for sample processing and EXCLUDE reagents used to measure analytes; these would be captured in Experiment Sample description as Reagents.	Blood from patient X Protein lysate from mouse brain RNA isolated from T cells purified from the spleen of patient X with disease phenotype A and placed into cell culture and treated with IL-2 for 2 hrs at 37°C RNA labeled with Cy3 for microarray experiment
Biological Sample Protocol	A protocol file submitted by the user that describes how a Biological Sample was isolated, enriched, processed, treated, or otherwise prepared for experimental use. The Biological Sample Treatment Protocol serves to capture any information which does not need to be parsed into the database and will be stored as a blob.	Purification of mRNA from mouse lymphocytes.txt Labeling RNA with Cy3 and Cy5.txt B cell enrichment from bone marrow.txt
Reagent	Something used to measure an Analyte(s). In the case of a microarray, this would be the chip type and name. In many cases, this may include information about probes: name, type, description, manufacturer, software used to select probe, etc.	Affymetrix Gene Chip HGU95 Information in the .gpl and .cdf files for custom arrays CD4-FITC, Pharmingen, clone 123 DNA oligo, sequence atgctgatccgaat
Analyte	The target compound being measured in a single Biological Sample using Reagent(s) used to detect a single target compound.	IFNG mRNA CD4 protein SNP rs12345 "A" allele

D-1 Glossary Table		
Term	Definition	Examples
Experiment Sample	Experiment Sample is distinguished from Biological Sample by the presence of Reagent(s). Reagent(s) captured for Experiment Sample are those specifically used to measure one or more analyte in the Biological Sample during the experiment. An Experiment Sample may also contain more than one Biological Sample. Thus, the combination of Biological Sample(s) + Reagent(s) used to measure an analyte(s) makes up an Experiment Sample. Replicates samples are stored as UNIQUE Experiment Samples that share the same Biological Sample, Experiment Sample attributes, and Experiment Sample Protocol, but were processed separately in the experiment.	RNA hybridized to Affy chip with 40,000 probes (one-color microarray) Control RNA labeled with Cy3 plus patient RNA labeled with Cy 5 hybridized to cDNA array with 16,000 probes (2-color microarray) T cells in one tube stained with CD4-FITC and CD8-PE (flow cytometry) Protein lysate in one anti-CD3-coated well of a 96-well plate (ELISA) 1 patient sample hybridized to one well of an Affy SNP chip (genotyping)
Experiment Sample Result	Result includes all of the measured results for a single Experiment Sample.	In the case of a two-color microarray, this would correspond to the single results file that contains the results for multiple analytes measured in multiple Biological Samples. This would also be the type of result captured in a single flow cytometry results file
Experiment Sample Results File	The file containing the minimally processed results generated from a single Experiment Sample.	Affymetrix .chp, .file ELISA experiment and contains data for all Experiment Samples
Protocol	A protocol file that describes how things are done.	Hybridization of labeled RNA to cDNA spotted array chip.txt (microarray) Staining using CD4-FITC antibody (flow cytometry)
Experiment	One or more Experiment Sample(s) containing at least one biological sample from analytes are measured, evaluated, or collected for the purposes of testing hypotheses or theories, demonstrating known facts or theories or measuring/collecting data.	A series of Affy gene expression chips A series of cDNA chips 5 tubes of T cells stained with CD4-FITC, CD8-PE, none, or both
Experiment	A label for the general type of experiment	Genotyping

D-1 Glossary Table

Term	Definition	Examples
Type	performed.	Molecular_Quantification Cellular_Quantification
Experiment Measurement Technique	The experimental system upon which an experiment is based or performed.	Array ELISPOT ELISA MBAA
Treatment	Captures descriptive information of experimental conditions for specific Experimental Samples and Biological Samples.	In cases where experiment samples may differ based on incubation temperature

APPENDIX 6 FAQs (FREQUENTLY ASKED QUESTIONS)

Q. What can you send?

A. You can send just about any file or set of files to ImmPort. There are some size limitations when using the ImmPort web site interface, but if you have large files (1 Gigabyte or more) you can request a secure ftp account or make other arrangements to have data loaded.

In addition, ImmPort can process descriptive information about your data which enables the use of analytical tools, advanced queries and sharing. ImmPort supports tab delimited files derived from spreadsheets (e.g. Excel).

Q. What should you send?

A. The question of what gets sent to ImmPort is a policy decision made by NIH staff, specifically NIAIAD DAIT program officers. The ImmPort staff is happy to help you organize and prepare your data submission packages.

Q. What would you like to send to ImmPort?

A. Although just about any file can be sent to ImmPort, what you or your team wish to do with the data once it is in ImmPort has a significant impact on how and what you send. In part that's because you can send just about any data and ImmPort will store it for you. To do more than just store it though, ImmPort has a higher threshold for how the data is described. For instance, if you want to use an ImmPort gene expression analysis tool, you'll need to tell ImmPort that you are sending gene expression data from well described samples that were assayed by a particular platform.

Q. Why are there so many files to complete and what are they for?

A. There is a metadata upload file for each metadata category. The metadata in one category can be linked to metadata in another category and to other files. By making the necessary associations between metadata and results, you can search for specific data to use with the ImmPort analysis tools. For example, a genotype assay result file from a custom panel can be submitted.

Q. What is a data submission package?

A. A data submission package is a .zip file containing a collection of metadata template and other files with many details about an experiment or study and the results.

Q. What happens to a data submission package?

A. Data submission packages are encrypted and transferred by secure ftp to the ImmPort system. The contents of the package will be processed to store the results and make them available for analysis, query or sharing in ImmPort.

Q. Which metadata upload file is required to submit?

A. Please refer to the Research Data Submission Instruction for complete details.

Q. Where do I find if my data has submitted successfully?

A. A Submission History page summarizes all the data that has been submitted within your project. You will also receive an email about the status of your submission.

Q. How do I view my submitted data?

A. By clicking on the Access Data link on the menu bar, you can see the options to search and display your data.

Q. My submission failed because an invalid value was provided for a controlled vocabulary element.

A. All controlled vocabulary values are case sensitive to ImmPort. The value specified in the completed ImmPort template file must exactly match one of the preferred vocabulary values found for that attribute.