

# Improved Methodology for Assessing Workplace Uranium Intakes Based on Accelerator Mass Spectrometric Measurements of Uranium-236 ( $^{236}\text{U}$ )



**Terry Hamilton**

**Center for Accelerator Mass Spectrometry**

**Global Security Principal Directorate**

This work was performed under the auspices of the U.S. Department of Energy by Lawrence Livermore National Laboratory under Contract DE-AC52-07NA27344.

Lawrence Livermore National Laboratory



GlobalSecurity  
Anticipate • Innovate • Deliver

# Acknowledgments

**Title:** Improved Methodology for Assessing Workplace Uranium Intakes Based on Accelerator Mass Spectrometric Measurements of Uranium-236 ( $^{236}\text{U}$ )

**Coauthors:**

T.F. Hamilton

T.A. Brown

A.R. Wood-Zika

S.J. Tumey

R.E. Martinelli

B.A. Buchholz

R.W. Williams

S.R. Kehl

Wm.G. Mansfield



## Why an interest in uranium isotopes!

- ❑ **Emerging needs in a number of different fields**
  - **Human Health**
  - **Nuclear Forensics**
  - **Environmental assessments**
  
- ❑ **Advances in measurement technologies especially in relation to the use of mass spectrometry**



# Chemistry, Materials and Life Sciences (CMELS)



GlobalSecurity  
Anticipate • Innovate • Deliver



**MC ICP-MS system for low-level actinides measurement**

Lawrence Livermore National Laboratory





## AMS Applications:


- ❑ Routine measurements of rare/stable isotope ratios, e.g., very well developed for  $^{14}\text{C}/^{12}\text{C}$
- ❑ More recently demonstrated for long-lived radionuclides such as the actinides,  $^{129}\text{I}$ , and  $^{99}\text{Tc}$
- ❑ AMS heavy-element line at LLNL was designed specifically for low-level detection of actinide elements and  $^{129}\text{I}$



# Accelerator Mass Spectrometry at LLNL

## Features of the heavy element AMS system at CAMS

- ❑ Rapid electrostatic switching between masses of interest
- ❑ High abundance sensitivity and a very wide dynamic range
- ❑ Simple chemistry and relatively high sample through-put (helps improve operational efficiency)
- ❑ AMS is a very robust measurement technique for actinide measurements with a very low susceptibility to interferences



HVEC Model FN Tandem Van de Graaff accelerator operated at voltages up to 9 MV with gas or foil stripping



# Uranium Isotopes in the Environment or Workplace

Isotope	Natural Abundance (atom %)	Half-Life (years)	Origin
$^{234}\text{U}$	0.005	$4.468 \times 10^5$	$^{238}\text{U}$ decay chain
$^{235}\text{U}$	0.72	$2.34 \times 10^7$	Primordial
$^{238}\text{U}$	99.27	$7.038 \times 10^8$	Primordial
$^{236}\text{U}$	$\sim 10^{-10}$ - $10^{-12}$	$2.446 \times 10^5$	$^{235}\text{U} + 1n \rightarrow ^{236}\text{U}$







# Initial $^{236}\text{U}$ AMS Application to Bioassay Monitoring

Global Security  
Anticipate • Innovate • Deliver

- ❑ **First sets of analyses performed in 2004**
  
- ❑ **Initial application of the technique included;**
  - **20 Reagent Blank samples**
  - **17 urine samples from the radiation worker population**
    - **9 workers from 4 programs/operations**
    - **Known or potential workplace intake(s) of uranium**
  - **8 “Baseline” urine samples**
    - **8 workers from several LLNL programs**
    - **No current or recent uranium intake potential**
    - **Range of residential locations and  $^{238}\text{U}$ -in-urine levels roughly similar to radiation worker population**







# $^{236}\text{U}$ content of reagent blank samples

## Data Table

Sample ID	$^{236}\text{U}$ , ng per sample	1- $\sigma$
Reagent Blank_1	$6.7 \times 10^{-7}$	$2.9 \times 10^{-7}$
Reagent Blank_2	$7.9 \times 10^{-7}$	$2.8 \times 10^{-7}$
Reagent Blank_3	$1.0 \times 10^{-7}$	$1.1 \times 10^{-7}$
Reagent Blank_4	$7.5 \times 10^{-7}$	$2.7 \times 10^{-7}$
Reagent Blank_5	$4.2 \times 10^{-7}$	$1.9 \times 10^{-7}$
Reagent Blank_6	$3.5 \times 10^{-7}$	$1.8 \times 10^{-7}$
Reagent Blank_7	$2.4 \times 10^{-7}$	$1.7 \times 10^{-7}$
Reagent Blank_8	$9.5 \times 10^{-7}$	$2.9 \times 10^{-7}$
Reagent Blank_9	$7.9 \times 10^{-7}$	$3.0 \times 10^{-7}$
Reagent Blank_10	$6.4 \times 10^{-7}$	$2.9 \times 10^{-7}$
Reagent Blank_11	$6.1 \times 10^{-7}$	$2.9 \times 10^{-7}$
Reagent Blank_12	$7.1 \times 10^{-7}$	$2.9 \times 10^{-7}$
Reagent Blank_13	$3.6 \times 10^{-7}$	$2.2 \times 10^{-7}$
Reagent Blank_14	$2.7 \times 10^{-7}$	$2.4 \times 10^{-7}$
Reagent Blank_15	$2.8 \times 10^{-7}$	$2.4 \times 10^{-7}$
Reagent Blank_16	$4.5 \times 10^{-7}$	$2.4 \times 10^{-7}$
Reagent Blank_17	$6.1 \times 10^{-7}$	$2.5 \times 10^{-7}$
Reagent Blank_18	$1.7 \times 10^{-7}$	$2.2 \times 10^{-7}$
Reagent Blank_19	$2.7 \times 10^{-7}$	$2.1 \times 10^{-7}$
Reagent Blank_20	$4.6 \times 10^{-7}$	$2.4 \times 10^{-7}$

## Summary Statistics

### $^{236}\text{U}$ Mean Value

**$4.9 \pm 2.4 \times 10^{-7}$  ng per sample**

AMS Limit of Detection

$\sim 7 \times 10^{-7}$  ng  $^{236}\text{U}$

c/w

MCICP-MS Limits of Detection

$\sim 0.2$  ng  $^{238}\text{U}$

$\sim 0.04$  ng  $^{235}\text{U}$



## <sup>236</sup>U Results for the Initial LLNL Baseline Cohort

Sample ID	<sup>236</sup> U, ng per sample	1- $\sigma$	<sup>236</sup> U, ng per gram	1- $\sigma$
Worker_01	$1.7 \times 10^{-7}$	$3.0 \times 10^{-7}$	$1.8 \times 10^{-9}$	$5.9 \times 10^{-9}$
Worker_02	$4.6 \times 10^{-7}$	$3.4 \times 10^{-7}$	$8.2 \times 10^{-9}$	$6.8 \times 10^{-9}$
Worker_03	$3.6 \times 10^{-7}$	$3.6 \times 10^{-7}$	$6.7 \times 10^{-9}$	$6.8 \times 10^{-9}$
Worker_04	$3.5 \times 10^{-7}$	$3.7 \times 10^{-7}$	$3.8 \times 10^{-9}$	$4.0 \times 10^{-9}$
Worker_05	$1.2 \times 10^{-7}$	$3.0 \times 10^{-7}$	$1.2 \times 10^{-9}$	$3.0 \times 10^{-9}$
Worker_06	$2.6 \times 10^{-7}$	$3.7 \times 10^{-7}$	$4.6 \times 10^{-9}$	$6.6 \times 10^{-9}$
Worker_07	$1.3 \times 10^{-7}$	$3.0 \times 10^{-7}$	$2.4 \times 10^{-9}$	$5.7 \times 10^{-9}$
Worker_08	$3.9 \times 10^{-7}$	$3.7 \times 10^{-7}$	$5.0 \times 10^{-9}$	$4.7 \times 10^{-9}$
<b>Mean Value</b>	<b><math>2.8 \pm 1.3 \times 10^{-7}</math></b>		<b><math>4.2 \pm 2.4 \times 10^{-9}</math></b>	

### Notes:

1. The associated <sup>238</sup>U concentrations ranged from ~0.3 to 20 ng per sample



# Preliminary $^{236}\text{U}$ Data on Bioassay Samples Collected from Radiation Workers at LLNL



GlobalSecurity  
Anticipate • Innovate • Deliver

Sample ID	$^{236}\text{U}$ , ng per sample	1- $\sigma$	$^{236}\text{U}$ , ng per gram	1- $\sigma$
RadWorker_01	$1.7 \times 10^{-4}$	$6.7 \times 10^{-6}$	$3.2 \times 10^{-7}$	$1.3 \times 10^{-7}$
RadWorker_02	$1.5 \times 10^{-6}$	$5.5 \times 10^{-7}$	$2.0 \times 10^{-8}$	$7.4 \times 10^{-9}$
RadWorker_03	$4.6 \times 10^{-6}$	$7.5 \times 10^{-7}$	$5.9 \times 10^{-8}$	$9.6 \times 10^{-9}$
RadWorker_04	$8.4 \times 10^{-5}$	$3.9 \times 10^{-6}$	$1.1 \times 10^{-6}$	$4.9 \times 10^{-8}$
RadWorker_05	$4.5 \times 10^{-6}$	$7.5 \times 10^{-7}$	$5.0 \times 10^{-8}$	$8.5 \times 10^{-9}$
RadWorker_06	$8.0 \times 10^{-5}$	$5.1 \times 10^{-6}$	$3.9 \times 10^{-7}$	$2.5 \times 10^{-8}$
RadWorker_07	$3.4 \times 10^{-5}$	$2.7 \times 10^{-6}$	$1.6 \times 10^{-7}$	$1.3 \times 10^{-8}$
RadWorker_08	$2.5 \times 10^{-4}$	$9.1 \times 10^{-6}$	$1.5 \times 10^{-6}$	$5.3 \times 10^{-8}$
RadWorker_09	$9.8 \times 10^{-6}$	$1.2 \times 10^{-6}$	$7.5 \times 10^{-8}$	$9.0 \times 10^{-9}$
RadWorker_10	$5.4 \times 10^{-5}$	$3.0 \times 10^{-6}$	$4.3 \times 10^{-7}$	$2.3 \times 10^{-8}$
RadWorker_11	$<7.0 \times 10^{-7}$		$<4.0 \times 10^{-8}$	
RadWorker_12	$6.3 \times 10^{-6}$	$9.5 \times 10^{-7}$	$3.1 \times 10^{-8}$	$4.6 \times 10^{-9}$
RadWorker_13	$1.6 \times 10^{-5}$	$1.5 \times 10^{-6}$	$5.8 \times 10^{-8}$	$5.6 \times 10^{-9}$
RadWorker_14	$2.7 \times 10^{-6}$	$6.7 \times 10^{-7}$	$3.2 \times 10^{-8}$	$7.9 \times 10^{-9}$
RadWorker_15	$1.3 \times 10^{-4}$	$6.2 \times 10^{-6}$	$1.6 \times 10^{-6}$	$8.0 \times 10^{-8}$
RadWorker_16	$3.8 \times 10^{-5}$	$2.6 \times 10^{-6}$	$1.9 \times 10^{-7}$	$1.3 \times 10^{-8}$
RadWorker_17	$4.9 \times 10^{-5}$	$3.0 \times 10^{-6}$	$6.1 \times 10^{-7}$	$3.7 \times 10^{-8}$
<b>Mean Value</b>	<b><math>5.8 \pm 7.1 \times 10^{-5}</math></b>		<b><math>4.1 \pm 5.3 \times 10^{-7}</math></b>	



## <sup>236</sup>U versus <sup>238</sup>U Content of Selected Bioassay Samples Collected from Radiation Workers at LLNL



GlobalSecurity  
Anticipate • Innovate • Deliver

Sample ID	<sup>238</sup> U, ng per sample (ICP-MS)	1-σ	<sup>236</sup> U, ng per sample (AMS)	1-σ
RadWorker_01	7.17	0.14	$1.7 \times 10^{-4}$	$6.7 \times 10^{-6}$
RadWorker_02	ND		$1.5 \times 10^{-6}$	$5.5 \times 10^{-7}$
RadWorker_03	0.23	0.06	$4.6 \times 10^{-6}$	$7.5 \times 10^{-7}$
RadWorker_04	3.83	0.09	$8.4 \times 10^{-5}$	$3.9 \times 10^{-6}$
RadWorker_05	0.26	0.06	$4.5 \times 10^{-6}$	$7.5 \times 10^{-7}$
RadWorker_06	5.79	0.12	$8.0 \times 10^{-5}$	$5.1 \times 10^{-6}$
RadWorker_07	3.00	0.08	$3.4 \times 10^{-5}$	$2.7 \times 10^{-6}$
RadWorker_08	9.54	0.18	$2.5 \times 10^{-4}$	$9.1 \times 10^{-6}$
RadWorker_09	2.40	0.08	$9.8 \times 10^{-6}$	$1.2 \times 10^{-6}$
RadWorker_10	3.39	0.09	$5.4 \times 10^{-5}$	$3.0 \times 10^{-6}$
RadWorker_11	0.03	0.06	ND	
RadWorker_12	16.0	0.28	$6.3 \times 10^{-6}$	$9.5 \times 10^{-7}$

ND, not detected






# What Do These Preliminary Measurements Show?

## Main Points

- **The  $^{236}\text{U}$  content of bioassay samples collected from the LLNL radiation workers appears to be clearly elevated over those of the baseline cohort**
- **Data confirms that some workers have been previously exposed to an anthropogenic source of uranium containing  $^{236}\text{U}$**



# Subsequent Applications of $^{236}\text{U}$ Bioassay at LLNL



- ❑ **2 additional AMS  $^{236}\text{U}$  runs were performed during the period between 2005 and 2006**

## 2005 Measurements

- I. **6 baseline samples (drinking water and urine)**
- II. **54 radiation worker samples (30 workers, 3 programs)**


## 2006 Measurements

- I. **6 baseline samples (drinking water and urine)**
- II. **74 radiation worker samples (32 workers, 3 programs)**



# $^{236}\text{U}$ Measurement Data from 2005 and 2006

Example 2006 AMS U236 Run Results: Low to High



Legend: ■ No Intake Potential ■ Low Intake Potential ■ Medium Intake Potential ■ High Intake Potential ■ ORNL U238 Spike - - - Det. Limit







GlobalSecurity  
Anticipate • Innovate • Deliver

## Applications in $^{236}\text{U}$ Bioassay - Findings & Uses

- ❑ In general, the level of  $^{236}\text{U}$  in bioassay samples increased with the relative exposure potential predicted by a workplace review in advance of the AMS analysis (one exception from 2006)





GlobalSecurity  
Anticipate • Innovate • Deliver

## Applications in $^{236}\text{U}$ Bioassay - Findings & Uses

- ❑ The  $^{236}\text{U}$  content of bioassay samples from the LLNL radiation worker cohort were useful for discerning relatively low-level occupational intakes in the presence of fluctuating environmental uranium excretion rates

### Examples;

- $^{238}\text{U}$  in urine levels from the radiation worker cohort ranged from 0.02 to 1.2 ug/L within same year, while the associated  $^{234}\text{U}/^{238}\text{U}$  ratios were within the expected environmental range
- The  $^{236}\text{U}$  content of such bioassay samples were often comparable across co-workers performing similar operation(s) with the same relative exposure potential (even in the presence of large variations in the total uranium content of samples)





GlobalSecurity  
Anticipate • Innovate • Deliver

## Applications in $^{236}\text{U}$ Bioassay - Findings & Uses

- ❑ The  $^{236}\text{U}$  excretion rates/patterns and  $^{236}\text{U}:^{238}\text{U}$  ratios (for operations involving the same and/or very similar uranium source terms), combined with worker  $^{238}\text{U}$  excretion history and workplace information, can be used to identify:
  - **The presence of multi-intake scenarios;**
  - **Previously unrecognized intake potential within specific operations and worker handling practices; and**
  - **Material bio-solubility information**



## Summary Points

- ❑ AMS offers a robust measurement technique for low-level detection of  $^{236}\text{U}$  in bioassay samples in the presence of fluctuating amounts of environmental uranium
- ❑ The application of  $^{236}\text{U}$  in bioassay studies appears to offer new opportunities for improving the standard of occupational safety and risk management at LLNL or elsewhere around the DOE complex
- ❑  $^{236}\text{U}$  is a potentially a useful ‘fingerprint’ for assessing the presence of anthropogenic sources of uranium (to include reprocessed, DU or enriched uranium) either inside the human body or in the environment



**Thank you for your attention**

**Terry Hamilton  
Center for Accelerator Mass Spectrometry  
Lawrence Livermore National Laboratory  
PO Box 808, L-397  
Livermore, CA 94551**

**telephone  
email**

**(925) 422 6621  
hamilton18@llnl.gov**