


Portable Oxygen Concentrators (POC) Performance Variables that Affect Therapy

Robert McCoy BS RRT FAARC

### **Old LTOT Patients and Products**



### New LTOT Patients and Products

You've come a long way baby













# LTOT Market Trends

- Early diagnosis
- Early treatment
- Maintain active (normal) life
- More options and more choices
- Patients involved in therapy

# New LTOT Products

- Can we see the forest for the trees?
  - Products are tools for clinicians to use. (trees)
  - Clinically effective therapy with positive patient outcomes is the objective (forest)






### Variables Related to Oxygenation

### Activity / Environment

- □ Rest
- □ Exercise
- □Sleep
- Altitude
- Oxygen %
- Patient / physiology / disease

### Dose vs. Saturation

The numbers on the dial are reference pointsObtain the right saturation and then look at the dial



# **Oxygen Therapy Prescription**

- 2 lpm continuous flow ?
- What is missing?
  - Titration with device the patient will be using at activity level
  - Understanding device operation
  - Respiratory rate and tidal volume
  - Follow up




Figure 7. Arterial blood gas response in 17 patients with chronic airway obstruction in whom progressive Pco<sub>2</sub> rise did not develop while they were receiving low flow oxygen. (Cross-hatched area répresents normal Pao<sub>2</sub> at Denver.)

### Recommendations from LTOT Consensus Conferences

- Titration LTOT with activity has been specified or implied since the 2nd Consensus Conference.
- Recommendation 8 from 2nd LTOT Consensus Conference (1987) "Clinical evaluation should include regular assessments of patients compliance with prescribed therapy, potential complications, potential hazards and the need for continued education. Patients receiving LTOT share responsibility with the prescribing physicians for remaining in communication with their physician in order to assure continued appropriate care for their condition"

#### Calculated Oxygen Concentrations Delivered by a Demand Valve with Pure Oxygen at Rest, Lessthan-pure Oxygen at Rest and Less-than-pure Oxygen with Reduced Volume during Exercise

	A 20 b/min, rest			<b>B</b> 20 b/min rest			C 30 b/min, exercise		
L/min flow	Vox ml	Cox %	Cf %	Vox ml	Cox %	Cf %	Vox ml	Cox %	Vox %
"1"	10	100	22.6	10	95	22.5	6	95	21.9
"2"	20	100	24.2	20	90	23.8	12	90	22.7
"3"	30	100	25.7	30	85	24.8	18	85	23.3
"4"	40	100	27.3	40	85	26.1	24	85	24.1
"5"	50	100	28.9	50	80	26.9	32	80	24.8

Editorial: Dr. Shigeoka & Linda Gallegos RRT Respir. Care Jan. 2006 Vol 51 No1

# Does 93% vs. 99% Make a Difference?

 Bolton CE, Annandale JA, Ebden P. Comparison of an oxygen concentrator and wall oxygen in the assessment of patients undergoing long term oxygen therapy assessment. Chron Respir Dis. 2006;3(1):49-51.

?

# Normobaric Hypoxia Inhalation Test vs. Response to Airline Flight in Healthy Passengers

Aviation, Space, and Environmental Medicine. Vol. 77, No 11 Nov. 2006






Fig. 1. Group  $S_{PO_2}$  response (mean  $\pm$  SD) to the NHIT.





# Altitude Adventure

- 13 LTOT patients and 5 clinicians
- Bus trip to 10,500 ft
- Continuous oxygen monitoring
- 12 patients maintained > 90% sat during the trip



My sat @ 84% at >10,000 ft with exercise

### **Check Saturation Levels**



# **Portable Concentrator**

- Improved operational range with access to ac or battery
- FAA approved
- Issues
  - Max. oxygen generation
  - ☐ Max. oxygen dose
  - Sleeping
  - 24 hour use












### Oxygen Dose Per Setting 15, 20, 25, 30 Breath per Minute



Breathing	Patterns	in	COPE	)
-----------	----------	----	------	---

<b>3</b>				Mean and (Standard Deviation)					
Investigator	Severity	n	Re F, min <sup>-1</sup>	est Vt,cc	Sub Ma F	x Work <u>Vt</u>	Max F	Work <u>Vt</u>	Ti/ <u>∏</u> t
Bliss/NOTT '99 ('80)	FEV1 32%	99	24 (6.8)		<b>26.7 (</b> 8.1)	825 (411)	<b>32.1</b> (8.3)		
Schanning '76 *	FEV1 45%	14	<b>17</b> (4)	650 (170)	<b>28 (</b> 5)	<b>1250</b> (270)	32 (6)	<b>1480</b> (290)	
Schanning '78 *	FEV1 66%	18	<b>16.7</b> (3.4)	<b>920</b> (180)	<b>25.2</b> (5.4)	<b>1460</b> (320)	<b>33.7</b> (6.0)	<b>1910</b> (530)	0.44
Tobin '84 PCO2<44 PCO2>44	FEV1 35% FEV1 32%	16 12	<b>20.4</b> (4.1) <b>23.3</b> (3.3)	<b>447</b> (139) <b>476</b> (158)					0.35
Sorlie '78 Supine PCO2 38 mean PCO2 50 mean Seated PCO2 38 mean PCO2 50 mean	FEV1 38% FEV1 22% FEV1 38% FEV1 22%	8 7 8 7	<b>16.9</b> (3.9) <b>24.1</b> (9.1)	930 (120) 640 (180) 710 (90) 560 (70)					0.39 0.41
Bradley '79		20	<b>21.6</b> (1.6)	<b>580</b> (40)					0.34
Hagarty '97 * Note that in the two	FEV1 1.0L	15			27 (6)	<b>1000 (</b> 400)	<b>30 (</b> 5)	<b>1300 (</b> 700)	

\* Note that in the two studies by Schanning, the subjects were not as severely obstructed than in other studies.

In these trials, the diseased patients maintained tidal volumes almost as large as the normal, control group at rest and at sub-maximal exercise.

#### References

- Bradley GW, Fleetam, Anthonisen N. Ventilatory Control in Hypoxic COPD. Am Review Resp Dis 1979; 120:21-30
- Sorli J, Grassino A, Lorange G, Milic-Emili J. Control of Breathing in Patients with Chronic Obstructive Lung Disease. *Clinical Science and Molecular Medicine 1978*; 54: 295-304
- Hagarty EM, Skrodin MS, Langbein WE, Hultman CI, Jessen JA, Maki KC. Comparison of Three Oxygen Delivery Systems during Exercise in Hypoxemic Patients with Chronic Obstructive Pulmonary Disease. Am J Respir Crit Care Med 1997;155: 893-898
- Tobin MG, Tejvir SC, Jenouri G, Birch ST, Gazeroglu HB, Sackner MA. Breathing Patterns 2. Diseased Subjects. *Chest 1983*; 84:286-294
- Schanning J. Respiratory Cycle Time Duration during Exercise in Patients with Chronic Obstructive Lung Disease. Scand J Resp Dis. 1978;59:313-318
- Schanning J. Ventilatory and Heart Rate Adjustments during Submaximal and Maximal Exercise in Patients with Chronic Obstructive Lung Disease. Scand J Resp Dis. 1976;57:63-72
- Bliss P. Analysis of Nocturnal Oxygen Therapy Trial Data Set. Unpublished. Valley Inspired Products 1999

## **Breath Rate Tracking**


Visual counting of breath rate: count for 15 seconds and multiply by 4, most recorded breath rates = 16 (4x4)



Dr. Tiep prototype

#### **VIP** Prototype






Duel Lumen Cannula


# Average Respiratory Rate from Two Locations, using a Breath Rate Counter

#### Milwaukee WI




Average Invacare 23.18 Average Spirit 20.54 Combined 22 Breath per minute

#### Albuquerque NM




Average Invacare 19.79 Average Spirit 22.25 Combined 21 Breath per minute

# Oxygen Delivered in 1 breath




# Oxygen Delivered in 1 breath




# Issues Differentiating OCDs

- Dose per Setting
- Maximum O2 Delivery
- Oxygen Purity
- Response to Rate Increases
- Sensitivity

# Pulse Volume per 1 lpm Setting




# Maximum FIO2 (20 bpm)




## FIO2 vs. Pulse Volume

Respiratory rate had an impact on 99% oxygen systems




### **Testing Pulse Products**

Device   Setting	1	2	3	4	5	6
Chad OM 301/2400	*	*	*	35	1	-
Chad OM 401/302	*	*	*	40	50	60
Chad OM 411/311	10	20	30	40	50	60
AirSep Impulse Select A Mode	*	*	*	35	44	52
AirSep Impulse Select B mode	17	33	50	66	83	99
Caire Spirit	15	30	45	60	75	-
Penox Escort	16	32	48	64	80	96



PulseMeter

### % of Patients Served



### POC Sleep Study Demographics

Patient	Age	Gender	Diagnosis	AHI	Years on Oxygen
1	73	М	Emphysema	1.0	3
2	65	F	Emphysema	2.6	6
3	58	М	Pulmonary Fibrosis	0.8	3
4	73	М	Emphysema	9.0	3
5	77	М	Emphysema	4.4	.4
6	71	F	Emphysema	5.0	6
7	74	F	Emphysema	5.0	2.5
8	76	F	Emphysema	3.8	7
9	72	F	Emphysema	Not Available	1.5
10	64	F	Emphysema	1.6	5

#### Results

			<b>DD00D</b>		of data	М	ean SpO₂(%)		Mean	HR (beats	s/min)
Patient	CF	PDOCD Sensitivity	PDOCD	CF	PDOC D	CF	PDOCD	Δ	CF	PDOC D	Δ
1	0.75	S	1	7.4	9.2	94.0	93.2	-0.8	75	73.4	-1.6
2	2	S	3	7.1	6.6	96.0	95.8	-0.2	89.6	89.8	0.2
3	2	D	5	8.1	9.0	90.1	90.3	0.2	74.1	82.7	8.6
4	3	S	3	9.5	9.0	97.5	96.1	-1.4	79.7	69.6	-10.1
5	2	S	2	7.2	8.6	96.9	94.6	-2.3	64.3	63.9	-0.4
6	2	S	3	7.2	9.6	96.5	97.2	0.7	64.3	69.6	5.3
7	2	S	3	6.6	6.6	97.2	93.3	-3.9	69.3	69.6	0.3
8	2	D	2.5	7.2	9.0	96.9	86.3	-10.6	79.9	84.3	4.4
9	2	S	3.5	5.1	7.2	96.5	94.0	-2.5	80	80.5	0.5
10	2.5	D	3	5.6	6.5	95.0	91.6	-3.4	97	96.0	-1.0
			Mean	7.1	8.1	95.7	93.2	-2.4	77.3	77.9	0.6
			SD	1.2	1.2	2.2	3.2	3.3	10.4	10.3	5.0

# **Conserving Devices**

- Multiple choices
   How do they work

   Capabilities
   Limitations
- Titrate for dose
- Understand the objective (patient oxygenation)



### Oxygen Saturation and Heart Rate

Patient	Spirit	Helios	Spirit	Helios
-	<u>SpO2</u>	<u>SpO2</u>	<u>HR</u>	HR
1	91.8	86.0	76.0	78.8
2	96.1	92.0	85.8	94.1
3	91.9	89.0	109.0	114.8
4	86.9	84.0	106.5	103.4
5	98.9	94.3	99.5	98.8
6	91.0	89.3	103.4	106.1
7	89.5	86.9	102.6	104.3
8	85.6	85.6	104.9	114.3
9	92.9	91.0	101.4	108.4
10	97.5	93.3	98.0	115.8
Average	92.2	89.1	98.7	103.9

## **Exercise and POCs**

- Dose volume and maximum dose is an important factor
- Respiratory rate has an impact on oxygen dose and purity
- Patient should be titrated on the POC at the activity level they will be using the POC

### NOTT Patients Matched for Ambulation





Figure 8. Re-examination of the Nocturnal Oxygen Therapy Trial by pretreatment ambulation status. Origins of the 80 matched patients.



# DEAL <sup>(f)</sup>NO DEAL



- Evidence based research - when findings conflict
  - □ Accept what you see
    - Deal!
  - Repeat the study yourself
    - No Deal!

# Thank You