

## BRUCELLOSIS (INFECTION WITH *B. ABORTUS*, *B. MELITENSIS* AND *B. SUIS*)

---

### SUMMARY

**Description of the disease:** *Brucellosis* is the generic name used for the animal and human infections caused by several species of the genus *Brucella*, mainly *Brucella abortus*, *B. melitensis* and *B. suis*. Infection with *Brucella* in cattle is usually caused by *B. abortus*, less frequently by *B. melitensis*, and occasionally by *B. suis*. *Brucella melitensis* is the main causative agent of infection with *Brucella* in sheep and goats. Infection with *Brucella* in pigs is due to *B. suis* biovars 1–3, but the disease caused by biovar 2 differs in its host range, its limited geographical distribution and its pathogenicity. In some areas, *B. suis* infection has become established in wild pigs. Clinically, infection with *Brucella* in animals is characterised by one or more of the following signs: abortion, retained placenta, orchitis, epididymitis and, rarely, arthritis, with excretion of the organisms in uterine discharges and in milk. Unequivocal diagnosis depends on the isolation of *Brucella* from abortion material, udder secretions or from tissues removed at post-mortem. *Brucella abortus*, *B. melitensis* and *B. suis* are highly pathogenic for humans, and potentially contaminated tissues, cultures and materials must be handled under appropriate containment conditions.

**Identification of the agent:** Evidence of *Brucella* is provided by the demonstration of *Brucella*-like organisms in abortion material or vaginal discharge using modified acid-fast staining, and is considered presumptive, especially if supported by serological tests. Polymerase chain reaction (PCR) methods are additional means for detection of the presence of *Brucella* DNA in a sample. Whenever possible, *Brucella* spp. should be isolated by culturing samples from uterine discharges, aborted fetuses, udder secretions or selected tissues, such as lymph nodes and male and female reproductive organs. Species and biovars should be identified by phage lysis, and by cultural, biochemical and serological tests. PCR can provide a complementary identification and typing method based on specific genomic sequences.

**Serological and cellular immunity tests:** The buffered *Brucella* antigen tests (rose bengal test and buffered plate agglutination test), the complement fixation test, the enzyme-linked immunosorbent assays (ELISA) or the fluorescence polarisation assay, are suitable tests for screening of herds/flocks and individual small ruminants, camelids and bovines (cattle and buffaloes). However, no single serological test is appropriate in each animal species and all epidemiological situations, and some of these tests are not adequate for diagnosing brucellosis in pigs. Therefore, the reactivity of samples that are positive in screening tests should be assessed using an established confirmatory or complementary strategy. The indirect ELISA or milk ring test performed on bulk milk samples is effective for screening and monitoring dairy cattle. The brucellin skin test can be used in both unvaccinated ruminants and swine as either a screening or a confirmatory herd test when positive serological reactors occur in the absence of obvious risk factors.

**Requirements for vaccines and diagnostic biologicals:** *Brucella abortus* strain 19 and *B. melitensis* strain Rev.1 remain the reference vaccines for the control of *Brucella* infections in cattle and in sheep and goats, respectively, with which any other vaccines should be compared. Both should be prepared from adequately derived seed cultures. The rough *B. abortus* strain RB51 vaccine has also become the official vaccine for prevention of *B. abortus* infection in cattle in some countries. No suitable vaccines exist for the control of *Brucella* infection in swine. Brucellin preparations must be free of smooth lipopolysaccharide, and antigens for serological tests must be prepared from smooth *B. abortus* strain 1119-3 or 99 and, in the case of indirect ELISA, from

smooth *B. melitensis* strain 16M as well. Vaccines and brucellin preparations must comply with relevant standards.

## A. INTRODUCTION

Brucellosis is the generic name used for the animal and human infections caused by several species of the genus *Brucella*, mainly *Brucella abortus*, *B. melitensis* and *B. suis*. Infection of sheep with *B. ovis* is described separately in Chapter 3.7.7 *Ovine epididymitis* (*Brucella ovis*).

**Causal pathogens:** Genetic and immunological evidence indicates that all members of the *Brucella* genus are closely related. Nevertheless, based on relevant differences in host preference and epidemiology displayed by the major variants, as well as molecular evidence of genomic variation, the International Committee on Systematics of Prokaryotes, Subcommittee on the Taxonomy of *Brucella* took a clear position in 2005 on a return to pre-1986 *Brucella* taxonomic opinion; the consequences of this statement imply the re-approval of the six classical *Brucella* nomenclatures with their corresponding recognised biovars. The classical names related to the six *Brucella* nomenclatures are validly published in the Approved Lists of Bacterial Names, 1980, and the designated type strains are attached to these validly published names: *B. abortus*, *B. melitensis*, *B. suis*, *B. neotomae*, *B. ovis* and *B. canis*<sup>1</sup>. The first three of these are subdivided into biovars based on cultural and serological properties (see Tables 2 and 3). Strains of *Brucella* have been isolated in the last decade from marine mammals but these strains cannot be ascribed to any of the above-recognised species. Investigations are continuing to establish their proper position in the taxonomy of the genus and it has been proposed that they be classified into two new species: *B. ceti* and *B. pinnipedialis* (Foster *et al.*, 2007). A new species, named *B. microti*, was also isolated from the common vole (*Microtus arvalis*) as well as from foxes and soil in Central Europe (Scholz *et al.*, 2008). Novel isolates from human breast implant and lung infections (strains BO1 and BO2) and from baboons that had delivered stillborn offspring have also been described, although the natural reservoir of these isolates remains unclear. While only two isolates of each new type have been described, they have been formally published as the tenth and eleventh *Brucella* species, *B. inopinata* and *B. papionis* respectively (Scholz *et al.*, 2010; Whatmore *et al.*, 2014). Finally, strains isolated from rodents, foxes and frogs were characterised as atypical *Brucella* strains distinct from the currently described species. They have not yet been approved as new *Brucella* species.

*Brucella* is a member of the *Brucellaceae* family, in the order Rhizobiales, class Alphaproteobacteria. It shows close genetic relatedness to some plant pathogens and symbionts of the genera *Agrobacterium* and *Rhizobium*, as well as animal pathogens (*Bartonella*) and opportunistic or soil bacteria (*Ochrobactrum*).

### 1. Description of the disease

#### 1.1. Infection with *Brucella* in cattle

Infection with *Brucella* in cattle is usually caused by biovars (bv.) of *Brucella abortus*. In some countries, particularly in southern Europe, Africa and western Asia, where cattle are kept in close association with sheep or goats, infection can also be caused by *B. melitensis* (Verger, 1985). Occasionally, *B. suis* may cause a chronic infection of the mammary gland of cattle, but it has not been reported to cause abortion or spread to other animals. Infection with *Brucella* in cattle is widespread globally but several countries in northern and central Europe, Canada, Japan, Australia and New Zealand are believed to be free from both *B. abortus* and *B. melitensis*.

The disease is usually asymptomatic in young animals and non-pregnant females. Following infection with *B. abortus* or *B. melitensis*, pregnant adult females develop a placentitis usually resulting in abortion between the fifth and ninth month of pregnancy. Even in the absence of abortion, profuse excretion of the organism occurs in the placenta, fetal fluids and vaginal discharges. The mammary gland and associated lymph nodes may also be infected, and organisms may be excreted in the milk. Subsequent pregnancies are usually carried to term, but uterine and mammary infection recurs, with reduced numbers of organisms in afterbirth products and milk. In acute infections, the organism is present in most major body lymph nodes. Adult male cattle may develop orchitis/epididymitis and brucellosis may be a cause of infertility in both sexes. Hygromas, usually involving leg joints, are a common manifestation of brucellosis in some tropical countries and may be the only obvious indicator of infection; the hygroma fluid is often infected with *Brucella*.

---

1 <http://www.the-icsp.org/subcoms/Brucella.htm>

## 1.2. Infection with *Brucella* in sheep and goats

Infection with *Brucella* in sheep and goats (excluding *B. ovis* infection) is primarily caused by one of the three biovars of *B. melitensis*. Sporadic infections caused by *B. abortus* or *B. suis* have been observed in sheep and goats, but such cases are extremely rare. Infection with *Brucella* in sheep and goats is endemic in the Mediterranean region, but infection is widespread world-wide. North America (except Mexico) is believed to be free from the agent, as are northern and central Europe, south-east Asia, Australia and New Zealand. Pathologically and epidemiologically, *B. melitensis* infection in sheep and goats is very similar to *B. abortus* infection in cattle. In most circumstances, the primary routes of transmission of *Brucella* are the placenta, fetal fluids and vaginal discharges expelled by infected ewes and goats when they abort or have full-term parturition. Shedding of *Brucella* is also common in udder secretions and semen, and *Brucella* may be isolated from various tissues, such as lymph nodes from the head, spleen and organs associated with reproduction (uterus, epididymides and testes), and from arthritic lesions (Alton *et al.*, 1988).

## 1.3. Infection with *Brucella* in pigs

Infection with *Brucella* in pigs is primarily caused by biovars 1, 2 or 3 of *B. suis*. Sporadic infections caused by *B. abortus* or *B. melitensis* have been also observed in pigs, but such cases are rare. The disease occurs in many countries where pigs are raised. Generally, the prevalence is low, but in some regions, such as South America and south-east Asia, the prevalence may be much higher. Porcine brucellosis may be a serious, but presently unrecognised, problem in some countries. *Brucella suis* bv. 1 infection has been reported from feral pigs in some of the southern states of the United States of America (USA), and in Queensland (Australia) and several other countries in Oceania. In these countries, a number of human infections have been reported from people who hunt and handle material taken from feral pigs. The disease is generally transmitted by consumption of feed contaminated by birth or abortion products and uterine discharges. Pigs will readily eat aborted fetuses and placental membranes. Transmission during copulation also occurs frequently, and *B. suis* excretion in semen has implications for those practising artificial insemination. In pigs, as in ruminants, after the initial bacteraemia, *B. suis* colonises cells of the reproductive tract of either sex. In females, placentas and fetuses are invaded, while in males, invasion occurs in one or more of the following: testes, prostate, epididymides, seminal vesicles or bulbo-urethral glands. In males the lesions, which are most often unilateral, start with a hyperplasia that may progress to abscess formation; the final stage is characterised by sclerosis and atrophy. Arthritis may occur in various joints, and sometimes spondylitis occurs. The most common manifestation of brucellosis in female pigs is abortion, occurring at any time during pregnancy, but most frequently between day 50 and 110 of gestation. Vaginal discharge is not often evident, and, in chronically infected herds, infertility rather than abortion is the most relevant clinical sign of the disease. In males, brucellosis is more likely to be persistent, with lesions in the genital tract often leading to interference with sexual activity, which can be temporary or permanent. The boar may excrete *Brucella* in the semen without any apparent abnormality in the sex organs or interference with sexual activity. In both sexes, there may be swollen joints and tendon sheaths, lameness and, occasionally, posterior paralysis. A significant proportion of both male and female pigs will recover from the infection, often within 6 months, but many will remain permanently infected (Olsen *et al.*, 2012).

Infection caused by *B. suis* bv. 2 differs from infection caused by bv. 1 and bv. 3 in its host range, distribution, and in pathogenicity. Historically, the geographical distribution of *B. suis* bv. 2 has been in a broad range between Scandinavia and the Balkans. The prevalence in wild boars appears to be high throughout continental Europe (EFSA, 2009). In outbreaks in Europe, wild boars were implicated as the source of transmission of bv. 2 to outdoor reared pigs, and are considered as the main wild reservoir of this infection (EFSA, 2009). *Brucella suis* bv. 2 causes miliary lesions, particularly in reproductive tissues, that often become purulent. To date, bv. 2 has rarely been reported as the cause of human brucellosis. However, bv. 2 infections have been reported in immuno-compromised hunters, who had been extensively exposed through gutting or skinning boars or hares. Moreover, rare cases of *B. suis* bv. 2 infection without clinical signs have been reported in Europe in cattle or sheep exposed to infected wild boars.

## 1.4. Infection with *Brucella* in other domestic, captive-wild or wild species

Infection with *B. abortus* or *B. melitensis* has been reported in the one-humped camel (*Camelus dromedarius*) and the two-humped camel (*C. bactrianus*), as well as in the South American camelids: llama (*Lama glama*), alpaca (*Lama pacos*), guanaco (*Lama guanicoe*), and vicuña (*Vicugna vicugna*), and is related to contact with large and small ruminants infected with *B. abortus* or *B. melitensis*.

In addition, brucellosis has been observed in the domestic buffalo (*Bubalus bubalis*), American and European bison (*Bison bison* and *B. bonasus*, respectively), yak (*Bos grunniens*), elk/wapiti (*Cervus elaphus*), African buffalo (*Syncerus caffer*) and various African antelope species. The clinical manifestations of brucellosis in these animals are similar to those seen in cattle or sheep and goats.

*Brucella melitensis* infection in wild ruminants may occur when these species are in close contact with sheep and goats in enzootic areas. The manifestations of brucellosis in these animals are similar to those in cattle or sheep and goats. However, in several wild ruminant species (e.g. chamois [*Rupicapra rupicapra*], Alpine ibex [*Capra ibex*] and the Iberian wild goat [*Capra pyrenaica*]), purulent or calcified arthritis and orchitis as well as uveitis and neurological troubles have been reported. These species are considered as dead-end carriers, and the disease usually disappears naturally as soon as *Brucella* infection has been eradicated from domestic livestock, unless anthropogenic effects take place.

There are two different types of epidemiological situation with regard to *B. suis* infection in other non-porcine species. In the first case, *B. suis* infection occurs in animals that are not the natural host of the particular infection through the ingestion of contaminated materials or by co-habitation with infected natural hosts. For example, Arctic foxes and wolves may contract *B. suis* bv. 4 from reindeer; dogs and rodents, such as rats and mice, may acquire other *B. suis* biovars by cohabitation with infected hosts; cattle and horses may become infected by cohabitation or interaction with infected swine. The infecting bacteria are invariably the well-defined biovars of the natural host species. In the second case, wildlife species are natural hosts for *B. suis* or *B. suis*-like infections. One example is the so-called murine brucellosis of the Commonwealth of Independent States (CIS) and the Baltic countries, where small rodents are infected with *B. suis* bv. 5. Another situation has been reported from Australia where strains isolated from rodents resembling *B. suis* but with different characteristics were involved; they were finally considered to be different from *B. suis* on genetic.

In addition to wild boar, the European hare (*Lepus europaeus*) is also considered to be a reservoir for *B. suis* bv. 2 and has been implicated as a possible source of transmission to domestic livestock. The disease in the European hare is characterised by the formation of nodules, varying in size from that of a millet seed to a cherry or even larger; these often become purulent. Such nodules may occur in almost any location, sometimes subcutaneously or intramuscularly, in the spleen, liver or lung and in the reproductive organs of either sex. The bodily condition of the hare may be surprisingly unaffected. Other species may also become infected by cohabitation with *B. suis* bv. 2 infected swine, wild boars or hares. Gutting or skinning wild boars in cattle sheds could be a route of transmission to cattle.

*Brucella suis* bv. 4 causes a serious zoonotic disease in wild or domesticated reindeer or caribou (*Rangifer tarandus* and its various subspecies) throughout the Arctic region, including Siberia, Canada and Alaska. *Rangifer tarandus* is very susceptible to *B. suis* infection, which causes fever, depression and various local signs, such as abortion, retained placentas, metritis, sometimes with blood-stained discharge, mastitis, bursitis and orchitis. Transmission to humans may be by direct contact or through consumption of raw milk and other inadequately heated products from reindeer, bone marrow in particular.

### 1.5. Zoonotic risk and biosafety requirements

*Brucella* infection is readily transmissible to humans, causing acute febrile illness – undulant fever – which may progress to a more chronic form and can also produce serious complications affecting the musculo-skeletal, cardiovascular, and central nervous systems. Precautions should be taken to prevent human infection. Infection is essentially acquired by the oral, respiratory, or conjunctival routes, but ingestion of raw milk products constitutes the main risk to the general public where the disease is endemic. There is an occupational risk to veterinarians, abattoir workers and farmers who handle infected animals/carcasses and aborted fetuses or placentas. Brucellosis is also one of the most easily acquired laboratory infections, and all laboratory manipulations with live cultures or potentially infected/contaminated material must be performed at an appropriate biosafety and containment level determined by biorisk analysis (see Chapter 1.1.4 *Biosafety and biosecurity: Standard for managing biological risk in the veterinary laboratory and animal facilities*). Specific recommendations have been made for the biosafety precautions to be observed with *Brucella*-infected materials (for further details see Alton *et al.*, 1988; Joint FAO/WHO Expert Committee on Brucellosis, 1986; WHO, 1953; WHO, 2004; Chapter 1.1.3 *Transport of biological materials*).

## B. DIAGNOSTIC TECHNIQUES

**Table 1.** Test methods available for the diagnosis of infection with *Brucella abortus*, *melitensis* or *suis*

Method	Purpose					
	Population freedom from infection	Individual animal freedom from infection <sup>a</sup>	Contribute to eradication policies <sup>b</sup>	Confirmation of suspect or clinical cases <sup>c</sup>	Herd/flock prevalence of infection – surveillance	Immune status in individual animals or populations post-vaccination
<b>Identification of the agent</b>						
Staining methods	–	–	–	+	–	–
Culture	–	–	–	+++	–	–
PCR <sup>d</sup>	–	–	–	+ / ++	–	–
<b>Detection of immune response</b>						
BBAT (RBT or BPAT)	+++	++	+++	+	+++	–
FPA	++	++	+	++	++	–
CFT	++	++	+++	++	+++	–
I-ELISA	+++	++	+++	++	+++	–
C-ELISA	++	+	+	+	++	–
BST	++	–	+	+++	++	–
SAT	++	+	+	–	+	–
NH and cytosol protein-based tests <sup>e</sup>	–	–	+	++	–	–
Bulk milk tests <sup>f</sup> Milk I-ELISA or Milk ring-test	+++	–	+++	+	+++	–

Key: +++ = recommended for this purpose; ++ recommended but has limitations;

+ = suitable in very limited circumstances; – = not appropriate for this purpose.

PCR = polymerase chain reaction; BBAT = buffered *Brucella* antigen tests (i.e. RBT [rose bengal test] and BPAT [buffered plate agglutination test]); FPA = fluorescence polarisation assay; CFT = complement fixation test; I- or C-ELISA = indirect/competitive enzyme-linked immunosorbent assay; BST = brucellin skin test; SAT = serum agglutination test; NH = native hapten

<sup>a</sup>This applies only to herds/flocks, countries or zones free from infection with *Brucella*.

<sup>b</sup>To increase the efficiency of eradication policies in infected herds/flocks, it is recommended to associate tests in parallel so as to increase the sensitivity of the diagnosis, i.e. two serological tests at least, e.g. BBAT or FPA and CFT or I-ELISA. The sensitivity is further increased by parallel testing by both serology and BST.

<sup>c</sup>In low-prevalence or almost-free zones, the predictive value of positive results to serological tests may be very low. In such situations, agent identification is usually needed to confirm clinical cases.

In infected herds/flocks, a positive result to any serological test may be considered as confirmation of a clinical case. Any reactor in any serological test should be considered to be infected even in the absence of clinical signs.

In low-prevalence or almost-free zones, singleton serological reactors may be confirmed by culture (or PCR) or BST.

In free countries or zones, suspect animals are those positive to both a screening and a confirmatory serological test (tests in series) and may be confirmed by culture (or PCR) and/or BST.

<sup>d</sup>False-positive results may occur.

<sup>e</sup>In zones where subcutaneous S19 or Rev.1 vaccination is practised, this test may help in differentiating antibodies due to vaccination from those due to infection.

<sup>f</sup>Dairy cattle only.

All cases of abortion as well as orchitis in cattle, sheep and goats, and pigs, should be considered as suspected brucellosis and should be investigated through the herd/flock history and submission of specimens for laboratory testing. The clinical signs are not pathognomonic and unequivocal diagnosis of *Brucella* infections can be made only by the isolation and identification of *Brucella*, but in situations where bacteriological examination is not practicable, diagnosis must be based on molecular or immunological methods.

## 1. Identification of the agent

There is no single test by which a bacterium can be identified unequivocally as *Brucella*. Accordingly, for a definitive identification, a combination of growth characteristics, serological, bacteriological or molecular methods is required (Alton *et al.*, 1988; Joint FAO/WHO Expert Committee on Brucellosis, 1986).

All samples from suspect cases should be cooled immediately after they are taken, and transported to the laboratory by the most rapid means. If they are to spend more than 12 hours in transit, all samples apart from vaginal swabs, should be frozen. On arrival at the laboratory, samples that are not to be cultured immediately should be frozen (Alton *et al.*, 1988). In all cases, the shorter the shipment and storage time, the higher is the probability of *Brucella* isolation, all the more so as the initial amount of *Brucella* is low in the sample. No specific transport medium has been demonstrated to improve *Brucella* survival in animal samples.

### 1.1. Staining methods

*Brucella* are coccobacilli or short rods measuring from 0.6 to 1.5 µm long and from 0.5 to 0.7 µm wide. They are usually arranged singly, and less frequently in pairs or small groups. The morphology of *Brucella* is fairly constant, except in old cultures where pleomorphic forms may be evident. *Brucella* are nonmotile. They do not form spores, and flagella, pili, or true capsules are not produced. *Brucella* are Gram negative and usually do not show bipolar staining. They are resistant to decolourisation by weak acids and thus stain red by the Stamp's modification of the Ziehl–Neelsen's method (Alton *et al.*, 1988). With this method, in smears of organs or biological fluids previously fixed with heat or ethanol, *Brucella* organisms stain red against a blue background. A fluorochrome or peroxidase-labelled antibody conjugate-based technique could also be used. The presence of intracellular, weakly acid-fast organisms of *Brucella* morphology or immuno-specifically stained organisms is a presumptive evidence of brucellosis. However, these methods are not feasible or have a low sensitivity in milk and dairy products where *Brucella* are often present in small numbers, and interpretation is frequently impeded by the presence of fat globules. Care must be taken as well in the interpretation of positive results in the Stamp's method because other organisms that cause abortions, e.g. *Chlamydia abortus* or *Coxiella burnetii*, may be difficult to differentiate from *Brucella* organisms in these preparations. The results, whether positive or negative, should be confirmed by culture.

Direct DNA probes or polymerase chain reaction (PCR) methods can also be used to demonstrate the agent in various biological samples (Bricker, 2002; Whatmore & Gopaul, 2011), but for the moment, the sensitivity and specificity of these approaches remain low with respect to the classical bacteriology. However, some of these molecular assays facilitate definitive identification and typing of the *Brucella* isolates (see Section B.1.4).

### 1.2. Collection of samples and culture

Bacteriological isolation is slow, expensive and cumbersome, but it should be performed whenever possible to confirm the disease and to determine the *Brucella* species/biovars involved. Although often considered not sensitive, it can be very effective when the type and number of samples, their adequate storage, amount seeded and culture media used are optimised.

#### 1.2.1. Basal media

Direct isolation and culture of *Brucella* are usually performed on solid media. This is generally the most satisfactory method as it enables the developing colonies to be isolated and recognised clearly. Such media also limit the establishment of non-smooth mutants and excessive development of contaminants. However, the use of liquid media may be recommended for voluminous samples or for the purpose of enrichment. A wide range of commercial dehydrated basal media is available, e.g. *Brucella* medium base, tryptose (or trypticase)–soy agar (TSA). The addition of 2–5% bovine or equine serum is necessary for the growth of strains such as *B. abortus* bv. 2, and many laboratories systematically add serum to basal media, such as blood agar base or Columbia agar, with excellent results. Other satisfactory media, such as serum–dextrose agar (SDA) or glycerol–dextrose agar, can be used (Alton *et al.*, 1988). SDA is usually preferred for observation of colonial morphology. A non-selective, biphasic medium, known as Castañeda's medium, is recommended for the isolation of *Brucella* from blood and other body fluids or milk, where enrichment culture is advised. Castañeda's medium is used because brucellae tend to dissociate in broth medium, and this interferes with biotyping by conventional bacteriological techniques.

### 1.2.2 Selective media

All the basal media mentioned above can be used for the preparation of selective media. Appropriate antibiotics are added to suppress the growth of organisms other than *Brucella*.

The most widely used selective medium is the modified Farrell's medium (FM) (Stack *et al.*, 2002), added to 1 litre of agar: polymyxin B sulphate (5000 units = 5 mg); bacitracin (25,000 units = 25 mg); natamycin (50 mg); nalidixic acid (5 mg); nystatin (100,000 units); vancomycin (20 mg). A corresponding freeze-dried antibiotic supplement is available commercially. However, nalidixic acid and bacitracin, at the concentration used in FM, have inhibitory effects on some *B. abortus*, *B. melitensis* and *B. suis* strains. Accordingly, the simultaneous use of FM and the less selective Thayer–Martin's modified (mTM) culture media has been considered the strategy of choice for *Brucella* primary isolation from field veterinary samples. However, the mTM is not translucent because of the haemoglobin contained as a basal component, being thus unsuitable for the direct observation of colonial morphology, probably the most practical procedure for the presumptive identification of *Brucella* (Alton *et al.*, 1988). A new selective and translucent culture medium (named CITA) has been recently formulated (De Miguel *et al.*, 2011). For its preparation, blood agar base is used as a basal component, supplemented with 5% sterile calf serum and containing vancomycin (20 mg/litre), colistin methanesulfonate (7.5 mg/litre), nitrofurantoin (10 mg/litre), nystatin (100,000 International Units [IU]/litre), and amphotericin B (4 mg/litre). This antibiotic mixture can be prepared as follows: weigh vancomycin, colistin and nystatin in the same 50 ml sterile container, then rehydrate the mixture with 10 ml of a 1:1 solution of absolute methanol in sterile purified water. Weigh then nitrofurantoin in a sterile tube and dissolve it with 1 ml of 0.1 M NaOH solution (sterilised previously by filtration through a 0.22 µm filter). Finally, weigh 10 mg of amphotericin B in a 20 ml sterile container and dissolve with 1 ml dimethyl sulphoxide. Once fully dissolved (5–10 minutes are required), add 9 ml of 10 mM sterile phosphate-buffered saline (PBS) (pH=7.2 ± 0.2). The final concentration of amphotericin B would be 1 mg/ml; a total of 4 ml of this solution are required for 1 litre of medium. The remaining Amphotericin B suspension can be kept at 5°C ± 3°C for several days for further uses. This new CITA medium inhibits most contaminant microorganisms but allows simultaneously the growth of all *Brucella* species and is more sensitive than both mTM and Farrell's media for isolating all smooth *Brucella* species from field samples, being thus the selective medium of choice for overall *Brucella* isolation, although the maximal diagnostic sensitivity is obtained using both FM and CITA simultaneously (De Miguel *et al.*, 2011).

Contrary to the situation with several *B. abortus* biovars as well as *B. ovis*, the growth of *B. melitensis* or *B. suis* is not dependent on an incubating atmosphere containing 5–10% CO<sub>2</sub> (Table 2), but such a CO<sub>2</sub> enriched-atmosphere is optimal for the culture of all *Brucella*.

As the number of *Brucella* organisms is likely to be lower in milk, colostrum and some tissue samples than in abortion material, enrichment can be advisable. In the case of milk, results can be improved by centrifugation and culture from both the cream and the pellet, but strict safety measures should be implemented in this case to avoid aerosols. A more practical way to increase the sensitivity of milk culture while avoiding the risks of centrifugation is increasing the number of both FM and CITA culture plates per each milk sample tested (two plates per udder quarter should be a minimum), each plate being inoculated with ca. 0.5 ml of milk. Enrichment can be carried out in liquid medium consisting of serum–dextrose broth, tryptose broth (or trypticase)–soy broth (TSA) or *Brucella* broth supplemented with an antibiotic mixture of at least amphotericin B (1 µg/ml), and vancomycin (20 µg/ml) (all final concentrations). The enrichment medium should be incubated at 37°C ± 2°C in air supplemented with 5–10% (v/v) CO<sub>2</sub> for up to 6 weeks, with weekly subcultures on to solid FM and CITA selective media. If preferred, a biphasic system of solid and liquid selective medium in the same bottle (Castañeda's method) may be used to minimise subculture. A selective biphasic medium composed of the basal Castañeda's medium with the addition of the following antibiotics to the liquid phase, is sometimes recommended for isolation of *Brucella* in milk (quantities are per litre of medium): polymyxin B (sulphate) (6000 units = 6 mg); bacitracin (25,000 units = 25 mg); natamycin (50 mg); nalidixic acid (5 mg); amphotericin B (1 mg); vancomycin (20 mg); D-cycloserine (100 mg).

All culture media used should be subjected to quality control with the reference strains to show that it performs properly. The use of a small inoculum of fastidious strains, such as *B. abortus* bv. 2, *B. ovis* or *B. suis* bv. 2, is preferred.

On suitable solid media, *Brucella* colonies can be clearly visible after a 3- to 4-day incubation period. After 4-days' incubation, *Brucella* colonies are round, 1–2 mm in diameter, with smooth margins. They are translucent and a pale honey colour when plates are viewed in the daylight

through a transparent medium. When viewed from above, colonies appear convex and pearly white. Later, colonies become larger and slightly darker. Smooth (S) *Brucella* cultures have a tendency to undergo variation during growth, especially with subcultures, and to dissociate to rough (R) forms. Colonies are then much less transparent, have a more granular, dull surface, and range in colour from matt white to brown in reflected or transmitted light. Checking for dissociation is easily tested by crystal violet staining: rough colonies stain red/violet and smooth colonies do not uptake dye or stain pale yellow. If the colonies are smooth, they should be checked against antiserum to smooth *Brucella*, or, if available, against anti-A and -M monospecific sera. In the case of non-smooth colonies, isolates should be checked with antiserum to *Brucella* R antigen. Changes in the colonial morphology are generally associated with changes in virulence, serological properties or phage sensitivity. Typical colonial morphology and positive agglutination with specific *Brucella* antiserum, followed by the oxidase and urease tests (see Tables 2 and 3), allow preliminary identification of the isolate as *Brucella*. However, it is recommended that subsequent confirmation and typing is performed by a reference laboratory.

### 1.2.3 Collection and culture of samples

For the diagnosis of animal brucellosis by cultural examination, the choice of samples usually depends on the clinical signs observed. The most valuable samples include vaginal secretions (swabs), aborted fetuses (stomach contents, spleen and lung), fetal membranes, and milk, semen and arthritis or hygroma fluids. From animal carcasses, the preferred tissues for culture are those of the reticulo-endothelial system (i.e. head, mammary and genital lymph nodes and spleen), the pregnant or early post-parturient uterus, and the udder. Growth normally appears after 3–4 days, but cultures should not be discarded as negative until 7–10 days have elapsed.

#### 1.2.3.1. Tissues

Samples are removed aseptically with sterile instruments. The tissue samples are prepared by removal of extraneous material (e.g. fat), cut into small pieces, and macerated using a 'Stomacher' or tissue grinder with a small amount of sterile PBS, before being inoculated on to solid media.

#### 1.2.3.2. Vaginal discharge

A vaginal swab taken after abortion or parturition is an excellent source for the recovery of *Brucella* and far less risky for the personnel than abortion material. The swab is then streaked directly onto solid media.

#### 1.2.3.3. Milk

Samples of milk must be collected cleanly after washing and drying the whole udder and disinfecting the teats. It is essential that samples should contain milk from all quarters, and 10–20 ml of milk should be taken from each teat, changing or disinfecting the gloves from one animal to the next to avoid cross-contamination of the samples. The first streams are discarded and the sample is milked directly into a sterile vessel or container. Care must be taken to avoid contact between the milk and the milker's hands. The milk can be centrifuged and the cream and deposit are spread on solid selective medium, either separately or mixed or streaked directly as indicated above. If *Brucella* are present in bulk milk samples, their numbers are usually low, and isolation from such samples is very unlikely.

#### 1.2.3.4. Dairy products

Dairy products, such as cheese, should be cultured on the media described above. As these materials are likely to contain small numbers of organisms, enrichment culture is advised. Samples need to be carefully homogenised before culture, after they have been ground in a tissue grinder or macerated and pounded in a 'Stomacher' or an electric blender with an appropriate volume (avoiding over-dilution) of sterile PBS. Superficial strata (rind and underlying parts) and the core of the product should be cultured. As brucellae grow, survive or disappear quite rapidly, their distribution throughout the different parts of the product varies according to the local physico-chemical conditions linked to specific process technologies.

#### 1.2.3.5. Arthritis/hygroma fluids – abscesses content

Such samples must be collected aseptically and spread directly on solid selective media.

All the above samples should be cooled (4–10°C) immediately after sampling and transported to the laboratory in the fastest way. Otherwise, the samples should be frozen to avoid viability losses. On arrival at the laboratory, milk and tissue samples and other biological liquids should be frozen if they are not to be cultured immediately.

Use of laboratory animals should be avoided unless absolutely necessary, but may sometimes provide the only means of detecting the presence of *Brucella*, especially when samples have been shown to be heavily contaminated or are likely to contain a low number of *Brucella* organisms. Animal inoculation may be intravenously or intraperitoneally in mice or intramuscularly, subcutaneously or intraperitoneally in guinea-pigs. This work must be carried out under appropriate biosafety conditions as outlined in chapter 1.1.4. The spleens of inoculated animals are cultured at 7 days (mice) or 3–6 weeks (guinea-pigs) after inoculation. Serum samples can be collected by intra-cardiac puncture before necropsy from guinea-pigs and subjected to buffered *Brucella* antigen tests (BBAT); a positive serological result is highly suggestive of brucellosis (Alton *et al.*, 1988).

### 1.3. Identification and typing

Any colonies showing the characteristic *Brucella* morphology should be examined using a Gram-stained-smear. As the serological properties, dyes and phage sensitivity are usually altered in the non-smooth phases, attention to the colonial morphology is essential in the typing tests described below. The recommended methods for observing colonial morphology are Henry's method by obliquely reflected light, the acriflavine test described by Braun & Bonestell, or White & Wilson's crystal violet method of staining colonies (Alton *et al.*, 1988).

Identification of *Brucella* organisms can be carried out by a combination of the following tests: organism morphology after Gram or Stamp's staining, direct observation of colonial morphology, growth characteristics, urease and oxidase tests, and the slide agglutination test with a polyclonal anti-*Brucella* serum. Species and biovar identification requires elaborate tests (such as phage lysis and agglutination with anti-A, -M or -R monospecific sera), the performance of which should be left to reference laboratories with accredited expertise in these methods. The simultaneous use of several phages e.g. Tbilisi (Tb), Weybridge (Wb), Izatnagar<sub>1</sub> (Iz<sub>1</sub>) and R/C provides a phage-typing system that, in experienced hands, allows a practical identification of the *Brucella* species. However, several characteristics, for example added CO<sub>2</sub> requirement for growth, production of H<sub>2</sub>S (detected by lead acetate papers), and growth in the presence of basic fuchsin and thionin, are revealed by routine tests that can be performed in moderately equipped non-specialised laboratories (see Tables 2 and 3).

For the maintenance of *B. abortus*, *B. melitensis* or *B. suis* strains as well as for sending them to a reference laboratory for typing, it is essential that only smooth colonies be selected. Cultures may be maintained for short periods at 5°C ± 3°C, but for longer periods they should be lyophilised or stored in a screw-capped tube at a temperature ≤ –16°C in tryptose broth with 15% (v/v) glycerol. For shipment, cultures should be lyophilised and sealed in ampoules packed in screw-capped canisters or subcultured onto appropriate nutrient agar slopes contained in screw-capped bottles. The strains could also be sent suspended in transport media (e.g. Amies), but this could cause dissociation.

For transporting *Brucella* cultures, the caps of the bottles or canisters should be screwed tightly down and sealed with PVC (polyvinyl chloride) tapes. Bottles should be wrapped in absorbent paper or cotton wool, sealed in polyethylene bags and packed into a rigid container (triple packaging) in accordance with the requirements of the International Air Transport Association (IATA) for shipping dangerous goods (IATA, 2013). These regulations are summarised in Chapter 1.1.2 *Collection, submission and storage of diagnostic specimens* and Chapter 1.1.3 *Transport of biological materials*, and they must be followed.

### 1.4. Nucleic acid recognition methods

The PCR, including the real-time format, provides an additional means of detection and identification of *Brucella* sp. (Bricker, 2002; Lopez-Goñi *et al.*, 2011; Ocampo-Sosa *et al.*, 2005; Whatmore & Gopaul, 2011). Despite the high degree of DNA homology within the genus *Brucella*, several molecular methods including PCR, PCR restriction fragment length polymorphism (RFLP) and Southern blot, have been developed allowing, to a certain extent, the differentiation of *Brucella* species and some of their biovars (for a review see Bricker, 2002; Moreno *et al.*, 2002; Whatmore & Gopaul, 2011). Pulse-field gel electrophoresis has been developed that allows the differentiation of several *Brucella* species. PCR can identify satisfactorily *Brucella* species and distinguish vaccine strains but there has been limited validation of the PCR for direct diagnosis.

The first species-specific multiplex PCR assay for the differentiation of *Brucella* was described by Bricker & Halling. The assay, named AMOS-PCR, was based on the polymorphism arising from species-specific localisation of the insertion sequence IS711 in the *Brucella* chromosome, and comprised five oligonucleotide primers that could identify without differentiating *B. abortus* bv. 1, 2 and 4 but could not identify *B. abortus* bv. 3, 5, 6, and 9. Modifications to the assay have been introduced over time to improve performance, and additional strain-specific primers were incorporated for identification of the *B. abortus* vaccine strains, and other biovars and species (Ocampo-Sosa *et al.*, 2005). A new multiplex PCR assay (Bruce-ladder) has been proposed for rapid and simple one-step identification of *Brucella*. The major advantage of this assay over previously described PCRs is that it can identify and differentiate in a single step most *Brucella* species as well as the vaccine strains *B. abortus* strain 19 (S19), *B. abortus* RB51 and *B. melitensis* Rev.1. In contrast to other PCRs, Bruce-ladder is also able to detect DNA from *B. neotomae*, *B. pinnipedialis* and *B. ceti*. The reference *B. abortus* bv. 3, 5, 6, 9, and *B. suis* bv. 2, 3, 4, 5 can be distinguished by this new multiplex PCR. An update to the original Bruce-ladder PCR protocol has been described. This new version (Bruce-ladder v2.0), that has been validated in several laboratories, is also able to discriminate between *B. suis* and *B. canis*, and allows the differentiation of *B. microti*. Similarly, a new multiplex PCR assay (Suis-ladder), has been developed for fast and accurate identification of *B. suis* strains at the biovar level (Lopez-Goñi *et al.*, 2011).

Another advanced multiplex PCR is also able to discriminate between *B. suis* and *B. canis* and between *B. suis* and *B. microti* in just one step, and between the vaccine strains *B. abortus* S19, *B. abortus* RB51 and *B. melitensis* Rev.1 (Kang *et al.*, 2011). This test could also allow the differentiation of the two marine mammal species, but this deserves further validation on field strains.

Other tests such as *omp25*, 2a and 2b PCR/RFLP are available and may be useful to identify some *Brucella* species.

Alternative approaches allowing identification of all *Brucella* species, *B. suis* biovars and vaccine strains based on single nucleotide polymorphism (SNP) discrimination by either primer extension or real-time PCR or the ligase-chain-reaction have been described. These tests are rapid, simple, unambiguous, and based on a robust population genetic analysis that helps ensure the species/biovar specificity of markers used (Whatmore & Gopaul, 2011).

A number of other methods adding useful epidemiological information have also been described. These include a multilocus sequencing scheme (Whatmore & Gopaul, 2011) and several typing schemes based on the use of MLVA (multiple locus variable number of tandem repeats analysis) (Bricker *et al.*, 2003; Le Flèche *et al.*, 2006; Whatmore & Gopaul, 2011). Depending on the particular markers chosen, these methods allow isolates to be identified at species level and provide epidemiological information at the subspecies level.

### 1.5. Identification of vaccine strains

Vaccine strains *B. abortus* S19, *B. melitensis* Rev.1 and *B. abortus* RB51 may be identified using specific PCRs (Kang *et al.*, 2011; Lopez-Goñi *et al.*, 2011), or by their growth characteristics in culture.

*Brucella abortus* S19 has the typical properties of a bv. 1 strain of *B. abortus*, but does not require CO<sub>2</sub>, does not grow in the presence of benzyl-penicillin (3 µg/ml = 5 IU/ml), thionin blue (2 µg/ml), or i-erythritol (1 mg/ml) (all final concentrations), and presents a high L-glutamate use (Alton *et al.*, 1988).

*Brucella melitensis* strain Rev.1 has the typical properties of a bv. 1 strain of *B. melitensis*, but develops smaller colonies on solid media, does not grow in the presence of basic fuchsin, thionin (both at 20 µg/ml) or benzyl-penicillin (3 µg/ml), but does grow in the presence of streptomycin at 2.5 or 5 µg/ml (5 IU/ml) (Alton *et al.*, 1988).

*Brucella abortus* strain RB51 can be distinguished from its *B. abortus* biovar 1 smooth counterparts by its rough morphology and growth in presence of rifampicin (250 µg per ml of media).

**Table 2.** Differential characteristics of species of the genus *Brucella*

Species	Colony morphology <sup>b</sup>	Serum requirement	Lysis by phages <sup>a</sup>					Oxidase	Urease activity	Preferred host
			Tb		Wb	Iz <sub>1</sub>	R/C			
			RTD <sup>c</sup>	10 <sup>4</sup> RTD	RTD	RTD	RTD			
<i>B. abortus</i>	S	– <sup>d</sup>	+	+	+	+	–	(+) <sup>e</sup>	(+) <sup>f</sup>	Cattle and other <i>Bovidae</i>
<i>B. melitensis</i>	S		–	–	(–) <sup>g</sup>	+	–	+	+ <sup>h</sup>	Sheep and goats
<i>B. suis</i>	S	–	–	+	(+) <sup>i</sup>	(+) <sup>i</sup>	–	+	+ <sup>j</sup>	Bv. 1: swine Bv. 2: swine, hare Bv. 3: swine Bv. 4: reindeer Bv. 5: rodents
<i>B. neotomae</i>	S	–	– <sup>k</sup>	+	+	+	–	–	+ <sup>j</sup>	Desert wood rat <sup>l</sup>
<i>B. ovis</i>	R	+	–	–	–	–	+	–	–	Sheep
<i>B. canis</i>	R	–	–	–	–	–	+	+	+ <sup>j</sup>	Dogs
<i>B. ceti</i>	S	ND	(–)		(+)	(+)	–	(+)	+	Cetaceans
<i>B. pinnipedialis</i>	S	ND	(–)		(+)	(+)	–	(+)	+ <sup>h</sup>	Pinnipeds
<i>B. microti</i>	S	–	–	+	+	+	ND	+	+ <sup>h</sup>	Common voles
<i>B. inopinata</i>	S	ND	–	PL <sup>m</sup>	ND	ND	ND	ND	+ <sup>j</sup>	Unknown
<i>B. papionis</i>	S		PL <sup>m</sup>	PL <sup>m</sup>	+	ND	–	–	+ <sup>j</sup>	Unknown

From Alton *et al.* (1988), Joint FAO/WHO Expert Committee on Brucellosis (1986), Whatmore (2009), and Whatmore *et al.* (2014).

(+)/(–) Most isolates positive/negative

a Phages: Tbilisi (Tb), Weybridge (Wb), Izatnagar1 (Iz<sub>1</sub>) and R/C

b Normally occurring phase: S: smooth, R: rough

c RTD: routine test dilution

d *B. abortus* bv. 2 generally requires serum for growth on primary isolation

e Some African isolates of *B. abortus* bv. 3 are negative

f Intermediate rate, except strain 544 and some field strains that are negative

g Some isolates are lysed by Wb

h Slow rate, except some strains that are rapid

i Some isolates of *B. suis* bv. 2 are not or only partially lysed by phage Wb or Iz<sub>1</sub>

j Rapid rate

k Minute plaques

l *Neotoma lepida*

m Partial lysis

ND Not determined

**Table 3.** Differential characteristics of the biovars of *Brucella* species

Species	Biovar	CO <sub>2</sub> requirement	H <sub>2</sub> S production	Growth on dyes <sup>a</sup>		Agglutination with monospecific sera			Reference strain <sup>2</sup>		
				Thionin	Basic Fuchsin	A	M	R	Strain	ATCC	NCTC
<i>B. abortus</i>	1	(+) <sup>b</sup>	+	-	+	+	-	-	544	23448	10093
	2	(+) <sup>b</sup>	+	-	-	+	-	-	86/8/59	23449	10501
	3 <sup>c</sup>	(+) <sup>b</sup>	+	+	+	+	-	-	Tulya	23450	10502
	4	(+) <sup>b</sup>	+	-	(+)	-	+	-	292	23451	10503
	5	-	-	+	+	-	+	-	B3196	23452	10504
	6 <sup>c</sup>	-	(-)	+	+	+	-	-	870	23453	10505
	9	+/-	+	+	+	-	+	-	C68	23455	10507
<i>B. melitensis</i>	1	-	-	+	+	-	+	-	16M	23456	10094
	2	-	-	+	+	+	-	-	63/9	23457	10508
	3	-	-	+	+	+	+	-	Ether	23458	10509
<i>B. suis</i>	1	-	+	+	(-)	+	-	-	1330	23444	10316
	2	-	-	+	-	+	-	-	Thomsen	23445	10510
	3	-	-	+	+	+	-	-	686	23446	10511
	4	-	-	+	(-)	+	+	-	40	23447	11364
	5	-	-	+	-	-	+	-	513	ND	11996
<i>B. neotomae</i>	-	+	- <sup>d</sup>	-	+	-	-	5K33	23459	10084	
<i>B. ovis</i>	+	-	+	(-)	-	-	+	63/290	25840	10512	
<i>B. canis</i>	-	-	+	(-)	-	-	+	RM6/66	23365	10854	
<i>B. ceti</i>	(-)	-	(+)	(+)	+	(-)	-	B1/94 BCCN 94-74	ND	12891	
<i>B. pinnipedialis</i>	(+)	-	+	(+)	(+)	(-)	-	B2/94 BCCN 94-73	ND	12890	
<i>B. microti</i>	-	-	+	+	-	+	-	CCM4915 BCCN 07-01 CAPM 6434	ND	ND	
<i>B. inopinata</i>	-	+	+	+	-	+ <sup>e</sup>	-	BO1 BCCN 09-01 CAPM 6436	ND	ND	
<i>B. papionis</i>	-	-	-	-	+	-	-	F8/08-60 CIRMBP 0958	ND	13660	

From Alton *et al.* (1988), Joint FAO/WHO Expert Committee on Brucellosis (1986), Whatmore (2009), and Whatmore *et al.* (2014).

(+)/(-) Most isolates positive/negative

a Dye concentration in serum dextrose agar: 20 µg/ml

b Usually positive on primary isolation

c For more certain differentiation of bv. 3 and 6, thionin at 40 µg/ml is used in addition: bv. 3 = +, bv. 6 = -

d Growth at a concentration of 10 µg/ml thionin

e Weak agglutination

ND Not determined

<sup>2</sup> <http://www.the-icsp.org/taxa/Brucellalist.htm>

## 2. Serological tests

No single serological test is appropriate in all epidemiological situations and all animal species; all tests have limitations especially when screening individual animals. Consideration should be given to all factors that impact on the relevance of the test method and test results to a specific diagnostic interpretation or application. In epidemiological units where vaccination with smooth *Brucella* is practised, and depending on the vaccination method (dose/route) used, positive serological reactions may be expected among the vaccinated animals because of antibodies cross-reacting with wild strain infection. Moreover, a number of bacteria, in particular *Yersinia enterocolitica* O:9, may induce antibody responses that cause false positive serological reactions (FPSR) in brucellosis tests, impeding accurate serological diagnosis. These FPSR may occur in all animal species at variable rates according to the time and the region.

The serum agglutination test (SAT) is generally regarded as being unsatisfactory for the purposes of international trade. The complement fixation test (CFT) is more specific than the SAT, and has also a standardised system of unitage. The diagnostic performance characteristics of some enzyme-linked immunosorbent assays (ELISAs) and the fluorescence polarisation assay (FPA) are comparable with or better than that of the CFT, and as they are technically simpler to perform and more robust, their use may be preferred. The diagnostic performance of these tests has been compared in cattle, small ruminants and swine.

For the control of brucellosis at the national or local level, BBATs (the rose bengal test [RBT] and the buffered plate agglutination test [BPAT]), ELISA and FPA, are considered as suitable screening tests. Depending on the purpose of testing, positive reactors could be retested using a suitable confirmatory or complementary method.

In other species, for example, buffaloes (*Bubalus bubalis*), American and European bison (*Bison bison*, *Bison bonasus*), yak (*Bos grunniens*), elk/wapiti (*Cervus elaphus*), camels (*Camelus bactrianus* and *C. dromedarius*), and South American camelids, *Brucella* sp. infection follows a course similar to that in cattle. The same serological procedures may be used for these animals, but each test should be validated for its fitness in the corresponding animal species.

### 2.1. Reference sera

The OIE reference standards are those against which all other standards are compared and standardised. These reference standards are all available to national reference laboratories and should be used to establish secondary or national standards against which working standards can be prepared and used in the diagnostic laboratory for daily routine use.

These sera have been developed and designated by the OIE as International Standard Sera<sup>3</sup>. The use of these reagents promotes international harmonisation of diagnostic testing and antigen standardisation:

- i) For RBT, CFT, SAT and milk ring test (MRT), the OIE International Standard Serum (OIEISS, previously named the WHO Second International standard anti-*Brucella abortus* Serum; WHO, 1953) is used. This serum is of bovine origin and contains 1000 IU (SAT) and 1000 ICFTU (international complement fixation test units).
- ii) For indirect ELISA (I-ELISA), competitive or blocking ELISA (C-ELISA) and FPA in cattle, three OIE ELISA Standard Sera are available for use. These are also of bovine origin and consist of a strong positive (OIEELISA<sub>SPSS</sub>), a weak positive (OIEELISA<sub>WPSS</sub>) and a negative (OIEELISA<sub>NSS</sub>) standard.
- iii) For I-ELISA, C-ELISA and FPA in sheep and goats, the International standard anti-*Brucella melitensis* Serum (ISaBmS) is used (McGiven *et al.*, 2011).
- iv) For I-ELISA, C-ELISA and FPA in pigs, there is no international standard serum available at present.

### 2.2. Production of antigens

*Brucella abortus* strain 99 (Weybridge) (S99)<sup>3</sup> or *B. abortus* strain 1119-3 (USDA) (S1119-3)<sup>4</sup> should always be used for the production of antigens for the BBATs, SAT, CFT and FPA. These *B. abortus* strains can be also used as a source of soluble antigen extracts (smooth lipopolysaccharide [S-LPS] or O-polysaccharide [OPS]) for the ELISAs or the Native Hapten tests, but *B. melitensis* strain 16M is also suitable for such a purpose. It should be emphasised that antigen made with any of the two *B. abortus* or *B. melitensis* 16M strains is used to test for any infections due to smooth *Brucella* species.

3 Obtainable from the OIE Reference Laboratory for Brucellosis in the United Kingdom.

4 Obtainable from the United States Department of Agriculture (USDA), National Veterinary Services Laboratories (NVSL) 1800 Dayton Road, Ames, Iowa, United States of America.

The strains must be completely smooth and should not auto-agglutinate in saline and 0.1% (w/v) acriflavine. They must be pure cultures and conform to the characteristics of CO<sub>2</sub>-independent strains of *B. abortus* bv. 1 or *B. melitensis* bv. 1. The original seed cultures should be propagated to produce a seed lot that must conform to the properties of these strains, and should be preserved by lyophilisation or by freezing in liquid nitrogen.

For antigen production, the seed culture is used to inoculate a number of potato-infusion agar slopes that are then incubated at 37°C ± 2°C for 48 hours. SDA and TSA, to which 5% equine or new-born calf serum or 0.1% yeast extract may be added, are satisfactory solid media provided a suitable seed is used as recommended above. The growth is checked for purity, resuspended in sterile PBS, pH 6.4, and used to seed layers of potato-infusion agar or glycerol-dextrose agar in Roux flasks. These are then incubated at 37°C ± 2°C for 72 hours with the inoculated surface facing down. Each flask is checked for purity by Gram staining samples of the growth, and the organisms are harvested by adding 50–60 ml phenol saline (0.5% phenol in 0.85% sodium chloride solution) to each flask. The flasks are gently agitated, the suspension is decanted, and the organisms are killed by heating at 80°C for 90 minutes. Following a viability check, the antigen is stored at 5°C ± 3°C.

Alternatively, the cells may be produced by batch or continuous culture in a fermenter, using a liquid medium containing (per litre of purified water) D-glucose (30 g), a high-grade peptone (30 g), yeast extract (Difco) (10 g), sodium dihydrogen phosphate (9 g) and disodium hydrogen phosphate (3.3 g). The initial pH is 6.6 (± 0.2), but this tends to rise to pH 7.2 (± 0.2) during the growth cycle. Care should be taken to check batches of peptone and yeast extract for capacity to produce good growth without formation of abnormal or dissociated cells. Vigorous aeration and stirring is required during growth, and adjustment to pH 7.2 (± 0.2) by the addition of sterile 0.1 M HCl may be necessary. The seed inoculum is prepared as described above. The culture is incubated at 37°C ± 2°C for 48 hours. Continuous culture runs can be operated for much longer periods, but more skill is required to maintain them. In-process checks should be made on the growth from either solid or liquid medium to ensure purity, an adequate viable count and freedom from dissociation to rough forms. Cells for use in the preparation of all antigens should be checked for purity and smoothness at the harvesting stage.

The culture is harvested by centrifugation to deposit the organisms, which are resuspended in phenol saline. The organisms are killed by heating at 80°C for 90 minutes and are stored at 5°C ± 3°C. They must form stable suspensions in physiological saline solutions and show no evidence of auto-agglutination. A viability check must be performed on the suspensions and no growth must be evident after 10 days of incubation at 37°C ± 2°C. The packed cell volume (PCV) of the killed suspensions can be determined by centrifuging 1 ml volumes in Wintrobe tubes at 3000 **g** for 75 minutes.

### 2.3. Buffered *Brucella* antigen tests (BBAT)

#### 2.3.1. Rose bengal test

This test is a simple spot agglutination test using antigen stained with rose bengal and buffered to a low pH, 3.65 ± 0.05 (Morgan *et al.*, 1969).

##### 2.3.1.1. Antigen production

Antigen for the RBT is prepared by depositing killed *B. abortus* S99 or S1119-3 cells by centrifugation at 23,000 **g** for 10 minutes at 5°C ± 3°C, and uniformly resuspending in sterile phenol saline (0.5%) at the rate of 1 g to 22.5 ml. (Note: if sodium carboxymethyl cellulose is used as the sedimenting agent during preparation of the cell concentrate, insoluble residues must be removed by filtering the suspension through an AMF-CUNO Zeta-plus prefilter [Type CPR 01A] before staining.) To every 35 ml of this suspension, 1 ml of 1% (w/v) rose bengal (CI No. 45440) in sterile purified water is added, and the mixture is stirred for 2 hours at room temperature. The mixture is filtered through sterile cotton wool, and centrifuged at 10,000 **g** to deposit the stained cells, which are then uniformly resuspended at the rate of 1 g cells to 7 ml of diluent (21.1 g of NaOH dissolved in 353 ml of sterile phenol saline, followed by 95 ml of lactic acid, and adjusted to 1056 ml with sterile phenol saline). The colour of this suspension should be an intense pink and the supernatant of a centrifuged sample should be free of stain; the pH should be 3.65 ± 0.05. After filtration through cotton wool, the suspension is filtered twice through a Sartorius No. 13430 glass fibre prefilter, adjusted to a PCV of approximately 8%, pending final standardisation against serum standardised against the OIEISS, and stored at 5°C ± 3°C in the dark. The antigen should be stored as recommended by the manufacturer. It should not be frozen.

### 2.3.1.2. Antigen standardisation

When used in the standard test procedure, the RBT antigen should give a clearly positive reaction with 1/45 dilution, but not 1/55 dilution, of the OIEISS diluted in 0.5% phenol saline or normal saline.

Additional checks may be performed with the ISaBmS. The highest dilution (in negative goat serum) of this standard that must give a positive result and the lowest dilution (in negative goat serum) that must simultaneously give a negative result have been established at 1/16 and 1/200, respectively (McGiven *et al.*, 2011).

It is also be advisable to compare the reactivity of new and previously standardised batches of antigen using a panel of well-defined reference sera.

However the above standardisation against the OIEISS is probably a cause of the reduced sensitivity of some RB antigen batches for diagnosing *B. melitensis* infection in small ruminants and of the discrepancies with the CFT (Blasco *et al.*, 1994a). When testing small ruminants, the discrepancies with the CFT can be minimised by using three volumes of serum and one volume of antigen (e.g. 75 µl and 25 µl, respectively) in place of an equal volume of each as mentioned in the standard test procedure. However, this modification of the RBT should not be recommended for testing cattle and pig sera.

### 2.3.1.3. Test procedure

- i) Bring the serum samples and antigen to room temperature (22°C ± 4°C); only sufficient antigen for the day's tests should be removed from the refrigerator.
- ii) Place 25–30 µl of each serum sample on a white tile, enamel or plastic plate, or in a WHO haemagglutination plate.
- iii) Shake the antigen bottle well, but gently, and place an equal volume of antigen near each serum spot.
- iv) Immediately after the last drop of antigen has been added to the plate, mix the serum and antigen thoroughly (using a clean glass or plastic rod for each test) to produce a circular or oval zone approximately 2 cm in diameter.
- v) The mixture is agitated gently for 4 minutes at room temperature (22°C ± 4°C) on a rocker or three-directional agitator (if the reaction zone is oval or round, respectively).
- vi) Read for agglutination immediately after the 4-minute period is completed. Any visible reaction is considered to be positive. A control serum that gives a minimum positive reaction should be tested before each day's tests are begun to verify the sensitivity of test conditions.

The RBT is very sensitive. However, like all other serological tests, it could sometimes give a positive result in cattle because of *B. abortus* S19 vaccination or of FPSR. The same phenomenon occurs in small ruminants or pigs affected by FPSR and in small ruminants vaccinated with *B. melitensis* Rev.1. Therefore positive reactions should be investigated using suitable confirmatory or complementary strategies (including epidemiological investigation). Conversely, false-negative reactions occur rarely. Nevertheless RBT appears to be adequate as a screening test for detecting infected herds or to guarantee the absence of infection in brucellosis-free herds or flocks.

## 2.3.2. Buffered plate agglutination test

### 2.3.2.1. Antigen production

Antigen for the BPAT is prepared from *B. abortus* S1119-3 according to the procedure described by Angus & Barton (1984).

Two staining solutions are required: brilliant green (2 g/100 ml) and crystal violet (1 g/100 ml) both certified stains dissolved in purified water. Once prepared, the two solutions should be stored separately for a period of 24 hours, and then mixed together in equal volumes in a dark bottle and stored in a refrigerator for a period of not less than 6 months before use. The mixed stain may only be used between 6 and 12 months after initial preparation.

Buffered diluent is prepared by slowly dissolving sodium hydroxide (150 g) in 3–4 litres of sterile phenol saline. Lactic acid (675 ml) is added to this solution, and the final volume is adjusted to 6 litres by adding sterile phenol saline. The pH of the solution should be 3.65 ± 0.05.

*Brucella abortus* S1119-3 packed cells are diluted to a concentration of 250 g/litre in phenol saline; 6 ml of stain is added per litre of cell suspension, and the mixture is shaken thoroughly before being filtered through sterile absorbent cotton. The cells are centrifuged at 10,000 **g** at 5°C ± 3°C, and the packed cells are then resuspended at a concentration of 50 g/100 ml in buffered diluent (as described above). This mixture is shaken thoroughly for 2 hours, and is then further diluted by the addition of 300 ml of buffered diluent per 100 ml of suspended cells (i.e. final concentration of 50 g packed cells/400 ml buffered diluent). The mixture is stirred at room temperature for 20–24 hours before the cell concentration is adjusted to 11% (w/v) in buffered diluent. This suspension is stirred overnight before testing. Pending final quality control tests, the antigen is stored at 5°C ± 3°C until required for use. The antigen should not be frozen.

The pH of the buffered plate antigen should be 3.70 ± 0.03 and the pH of a serum–antigen mixture at a ratio of 8:3 should be 4.02 ± 0.04. The 11% stained-cell suspension should appear blue–green. Each batch of buffered plate antigen should be checked by testing at least 10 weakly reactive sera and comparing the results with one or more previous batches of antigen. If possible, the antigen batches should be compared with the standard antigen prepared by the NVSL, USDA (see footnote 4 for address). There is, however, no international standardisation procedure established for use with either the OIEISS or with the ISaBmS.

#### 2.3.2.2. Test procedure

- i) Bring the serum samples and antigen to room temperature (22°C ± 4°C); only sufficient antigen for the day's tests should be removed from the refrigerator.
- ii) Shake the sample well. Place 80 µl of each serum sample on a glass plate marked in 4 × 4 cm squares
- iii) Shake the antigen bottle well, but gently, and place 30 µl of antigen near each serum spot.
- iv) Immediately after the last drop of antigen has been added to the plate, mix the serum and antigen.
- v) Thoroughly (using a clean glass or plastic rod for each test) to produce a circular zone approximately 3 cm in diameter.
- vi) After the initial mixing, the plate should be rotated three times in a tilting motion to ensure even dispersion of the reagents, and then incubated for 4 minutes in a humid chamber at ambient temperature.
- vii) The plate should be removed and rotated as above, and then returned for a second 4-minute incubation.
- viii) Read for agglutination immediately after the 8-minute period is completed. Any visible reaction is considered to be positive. A control serum that gives a minimum positive reaction should be tested before each day's tests are begun to verify the sensitivity of test conditions.

Like the RBT, the test is very sensitive in cattle, especially for detection of vaccine-induced antibody, and positive samples should be retested using a confirmatory or complementary test(s). False-negative reactions may occur, usually due to prozoning, which may be overcome by diluting the serum or retesting after a given time. While the BPAT has been extensively used with apparent good results in small ruminants and pigs in some countries, its diagnostic value in these species has not been reported at international level.

#### 2.4. Complement fixation test

The CFT is widely used but it is complex to perform, and requires good laboratory facilities and adequately trained staff to accurately titrate and maintain the reagents. There are numerous variations of the CFT in use, but this test is most conveniently carried out in a microtitre format. Warm or cold fixation may be used for the incubation of serum, antigen and complement: either 37°C ± 2°C for 30 minutes or 5°C ± 3°C for 14–18 hours. A number of factors affect the choice of the method: anti-complementary activity in serum samples of poor quality is more evident with cold fixation, while fixation at 37°C ± 2°C increases the frequency and intensity of prozones, and a number of dilutions must be tested for each sample.

Several methods have been proposed for the CFT using different concentrations of fresh or preserved sheep red blood cells (SRBCs) (a 2%, 2.5% or 3% suspension is usually recommended) sensitised with an equal volume of rabbit anti-SRBC serum diluted to contain several times (usually from two to

five times) the minimum concentration required to produce 100% lysis of SRBCs in the presence of a titrated solution of guinea-pig complement. The latter is independently titrated (in the presence or absence of antigen according to the method) to determine the amount of complement required to produce either 50% or 100% lysis of sensitised SRBCs in a unit volume of a standardised suspension; these are defined as the 50% or 100% haemolytic unit of complement/minimum haemolytic dose (C'H or MHD<sub>50</sub> or C'H or MHD<sub>100</sub>), respectively. It is generally recommended to titrate the complement before each set of tests, a macromethod being preferred for an optimal determination of C'H<sub>50</sub>. Usually, 1.25–2 C'H<sub>100</sub> or 5–6 C'H<sub>50</sub> are used in the test.

Barbital (veronal) buffered saline is the standard diluent for the CFT. This is prepared from tablets available commercially; otherwise it may be prepared from a stock solution of sodium chloride (42.5 g), barbituric acid (2.875 g), sodium diethyl barbiturate (1.875 g), magnesium sulphate (1.018 g), and calcium chloride (0.147 g) in 1 litre of purified water, diluted by the addition of four volumes of 0.04% gelatine solution before use. However, this buffer contains barbituric derivatives that are no longer available in several countries. Satisfactory results may be also obtained with a barbituric-free solution of sodium chloride 0.85% containing calcium and magnesium, prepared by adding 1 ml of a stock solution of 1 M magnesium chloride and 0.3 M calcium chloride (anhydrous MgCl<sub>2</sub>: 9.5 g CaCl<sub>2</sub>: 3.7 g; purified water: up to 100 ml) (stored in small amounts at 5°C ± 3°C) to 1 litre of saline solution (Alton *et al.*, 1988). The pH is critical and must be strictly adjusted to 7.35 (± 0.05). The replacement of the veronal buffer by this barbituric-free buffer has been validated in the OIE Brucellosis Reference Laboratory in France<sup>5</sup>.

#### 2.4.1. Antigen production

Numerous variations of the test exist but, whichever procedure is selected, the test must use an antigen that has been prepared from an approved smooth strain of *B. abortus*, such as S99 or S1119-3, and standardised against the OIEISS. Antigen for the CFT can be prepared following specialised procedures (Alton *et al.*, 1988) or a whole cell antigen can be used after diluting the stock suspension such that the PCV of the concentrated antigen suspension for CFT is approximately 2% before standardisation against the OIEISS.

#### 2.4.2. Antigen standardisation

The antigen should be standardised to give 50% fixation at a dilution of 1/200 of the OIEISS and must also show complete fixation at the lower serum dilutions, because too weak (or too strong) a concentration of antigen may not produce 100% fixation at the lower dilutions of serum. When two dilutions of antigen are suitable, the more concentrated antigen suspension must be chosen in order to avoid prozone occurrence. The appearance of the antigen, when diluted 1/10 must be that of a uniform, dense, white suspension with no visible aggregation or deposit after incubation at 37°C ± 2°C for 18 hours. It must not produce anti-complementary effects at the working strength for the test. The antigen is stored at 5°C ± 3°C and should not be frozen.

#### 2.4.3. Test procedure (example)

The undiluted test sera and appropriate working standards should be inactivated for 30 minutes in a water bath at 60°C ± 2°C. If previously diluted with an equal volume of veronal buffered saline, these sera could be inactivated at 58°C ± 2°C for 50 minutes. Usually, only one serum dilution is tested routinely (generally 1/4 or 1/5 depending on the CF procedure chosen), but serial dilutions are recommended for trade purposes and when clinical signs have been reported in order to detect prozone.

Using standard 96-well microtitre plates with round (U) bottoms, the technique is usually performed as follows:

- i) Volumes of 25 µl of diluted inactivated test serum are placed in the well of the first, second and third rows. The first row is an anti-complementary control for each serum. Volumes of 25 µl of CFT buffer are added to the wells of the first row (anti-complementary controls) to compensate for lack of antigen. Volumes of 25 µl of CFT buffer are added to all other wells except those of the second row. Serial doubling dilutions are then made by transferring 25 µl volumes of serum from the third row onwards; 25 µl of the resulting mixture in the last row are discarded.
- ii) Volumes of 25 µl of antigen, diluted to working strength, are added to each well except in the first row.

---

5 Obtainable from the OIE Reference Laboratory for Brucellosis in France.

- iii) Volumes of 25 µl of complement, diluted to the number of units required, are added to each well.
- iv) Control wells containing:
  - a) diluent only,
  - b) complement + diluent,
  - c) antigen + complement + diluent,are set up to contain 75 µl total volume in each case.

A control serum that gives a minimum positive reaction should be tested in each set of tests to verify the sensitivity of test conditions.

- v) The plates are incubated at 37°C ± 2°C for 30 minutes or at 5°C ± 3°C overnight, and a volume (25 or 50 µl according to the technique) of sensitised SRBCs is added to each well. The plates are re-incubated at 37°C ± 2°C for 30 minutes.
- vi) The results are read after the plates have been centrifuged at 1000 *g* for 10 minutes at 5°C ± 3°C or left to stand at 5°C ± 3°C for 2–3 hours to allow unlysed cells to settle. The degree of haemolysis is compared with standards corresponding to 0, 25, 50, 75 and 100% lysis. The absence of anti-complementary activity is checked for each serum in the first row.
- vii) Standardisation of results of the CFT

A unit system that is based on the OIEISS exists for the standardisation of results. This serum contains 1000 ICFTU per ml. If this serum is tested in a given method and gives a titre of, for example 200 (50% haemolysis), then the factor for an unknown serum tested by that method can be found from the formula:  $1000 \times 1/200 \times \text{titre of test serum} = \text{number of ICFTU of antibody in the test serum per ml}$ . The OIEISS contains specific IgG; national standard sera should also depend on this isotype for their specific complement-fixing activity. Difficulties in standardisation arise because different techniques selectively favour CF by different immunoglobulin isotypes. It is recommended that any country using the CFT on a national scale should obtain agreement among the different laboratories performing the test to use the same method in order to obtain the same level of sensitivity. To facilitate comparison between countries, results should always be expressed in ICFTUs, calculated in relation to those obtained in a parallel titration with a standard serum, which in turn may be standardised against the OIEISS.

- viii) Interpretation of the results: Sera giving a titre equivalent to 20 ICFTU/ml or more are considered to be positive.

Animals that have been vaccinated with *B. abortus* S19 or *B. melitensis* Rev.1 between 3 and 6 months are usually considered to be infected if the sera give positive fixation at a titre of 30 or greater ICFTU/ml when the animals are tested at an age of 18 months or older.

This procedure is an example, other volumes and quantities of reagents could be chosen provided that the test is standardised against the OIEISS as described above and the results expressed in ICFTU/ml.

The CFT is usually very specific but less sensitive than RBT and ELISA, particularly in the case of swine, as swine complement interacts with guinea-pig complement to produce a pro-complementary activity that reduces the sensitivity. Thus, the CFT has a reduced sensitivity for diagnosing *B. suis* infection, is not capable of eliminating the FPSR problem, and can be recommended only as a complementary test in swine. Moreover, like most serological tests, the CFT can result positive in ruminants after *B. abortus* S19 or *B. melitensis* Rev.1 vaccination and it is not specific enough in presence of FPSR. Therefore, CFT results should be investigated using suitable confirmatory or complementary strategies.

## 2.5. Enzyme-linked immunosorbent assays

### 2.5.1. Indirect ELISA

Numerous variations of the I-ELISA have been described for cattle, small ruminants and pigs employing different antigen preparations, antiglobulin-enzyme conjugates, and substrate/

chromogens. *Brucella abortus* strain 99 (Weybridge) (S99)<sup>3</sup> or *B. abortus* strain 1119-3 (USDA) (S1119-3)<sup>4</sup> should be used for production of these antigens, but *B. melitensis* strain 16M can also be suitable for such a purpose. Several commercial I-ELISAs using whole cell, S-LPS or the OPS as antigens that have been validated in extensive field trials are available and are in wide use. Nevertheless, the technique used and the interpretation of results must have been validated in accordance with the principles laid down in Chapter 1.1.6 *Principles and methods of validation of diagnostic assays for infectious diseases*.

#### 2.5.1.1. Test standardisation of the I-ELISA (EU, 2008; McGiven *et al.*, 2011)

##### 2.5.1.1.1. Infection with *Brucella* in cattle

- i) A 1/2 pre-dilution of the OIEELISA<sub>WPSS</sub> or a 1/16 pre-dilution of the OIEELISA<sub>SPSS</sub> made up in a negative bovine serum (or in a negative pool of bovine sera) must give a positive reaction;  
and
- ii) A 1/8 pre-dilution of the OIEELISA<sub>WPSS</sub> or a 1/64 pre-dilution of the OIEELISA<sub>SPSS</sub> made up in a negative bovine serum (or in a negative pool of bovine sera) must give a negative reaction;  
and
- iii) The OIEELISA<sub>NSS</sub> must always give a negative reaction.

##### 2.5.1.1.2. Infection with *Brucella* in sheep and goats:

- i) A 1/64 pre-dilution of the ISaBmS made up in a negative goat serum (or in a negative pool of goat sera) must give a positive reaction;  
and
- i) A 1/750 pre-dilution of the ISaBmS made up in a negative goat serum (or in a negative pool of goat sera) must give a negative reaction;  
and
- i) The above-mentioned negative goat serum (or negative pool of goat sera) must always give a negative reaction.

##### 2.5.1.1.3. Infection with *Brucella* in pigs:

- i) In the absence of international standard serum for porcine brucellosis the test should be duly validated and the cut-off established in the test population with appropriate validation techniques (see chapter 1.1.6).

The I-ELISAs that use S-LPS or OPS as antigens are highly sensitive for the detection of anti-*Brucella* antibodies in cattle, small ruminants and pigs, but are not capable of fully resolving the problem of differentiating between antibodies resulting from *B. abortus* S19 and *B. melitensis* Rev.1 vaccination. The *B. abortus* RB51 vaccine may also interfere in S-LPS-based I-ELISAs.

Positive reactions should be investigated using suitable confirmatory or complementary strategies as for BBAT and CFT.

Using I-ELISA standardised against the OIE standard sera described above, the diagnostic sensitivity should be equal to or greater than that of the BBATs (RBT/BPAT) or the CFT in the testing of infected cattle, small ruminants and pigs. However, the specificity would usually be lower (Praud *et al.*, 2012).

The problem of FPSR may be reduced but not fully resolved, in pigs in particular, by performing I-ELISAs using extracts from rough strains of *Brucella*. Most FPSR are a result of cross reaction with the OPS portion of the S-LPS molecule; cross-reaction among core epitopes is less frequent but does exist. The use of chaotropic I-ELISAs or procedures using heterologous extracts or *Brucella* cytosolic proteins as antigens does not resolve this problem, at least in cattle (Muñoz *et al.*, 2005). Moreover, in the context of FPSR, the most specific diagnostic procedure remains the brucellin skin test (see Section B.3.1 below).

Monoclonal, polyclonal antiglobulin or protein G or AG enzyme conjugates may be used depending on availability and performance requirements. A monoclonal antibody (MAb) specific for the heavy chain of bovine IgG<sub>1</sub> may provide some improvement in specificity at the possible

cost of some loss of sensitivity while a protein G or AG enzyme conjugate may provide a reagent useful for testing a variety of mammalian species.

Several commercial I-ELISAs are available. Some protocols are less sensitive or less specific than others; therefore results obtained from different assays are not always comparable. I-ELISA for diagnosing anti-*Brucella* antibodies in small ruminants and pigs is essentially the same as that described for cattle, but the cut-off should have been properly established for these species using the appropriate validation techniques (see chapter 1.1.6), and, moreover, I-ELISA for sheep and goats should be standardised against the ISaBmS (McGiven *et al.*, 2011).

Whatever the I-ELISA format used:

- i) A positive and a negative control are included in each plate. OD (optical density) ranges to be obtained with these two controls must be established to define the criteria for validating each plate results. The OD of the positive control is the one with which the OD of each test serum is compared to establish the final result (negative or positive).
- ii) An additional positive serum (internal control) must be included in each plate to validate the repeatability of the test from plate to plate and from day to day.

## 2.5.2. Competitive ELISA

Several variations of the C-ELISA, using S-LPS or OPS as antigens, have been described for cattle, small ruminants and pigs employing different antiglobulin-enzyme conjugates, substrate or chromogens and antigens prepared from different smooth *Brucella* strains. Nevertheless, the technique used and the interpretation of results must have been validated in accordance with the principles laid down in chapter 1.1.6.

### 2.5.2.1. Test standardisation of the C-ELISA (EU, 2008; McGiven *et al.*, 2011)

#### 2.5.2.1.1. Infection with *Brucella* in cattle:

- i) A 1/2 pre-dilution of the OIEELISA<sub>WPSS</sub> or a 1/16 pre-dilution of the OIEELISA<sub>SPSS</sub> made up in a negative bovine serum (or in a negative pool of bovine sera) must give a positive reaction;  
and
- ii) A 1/8 pre-dilution of the OIEELISA<sub>WPSS</sub> or a 1/64 pre-dilution of the OIEELISA<sub>SPSS</sub> made up in a negative bovine serum (or in a negative pool of bovine sera) must give a negative reaction;  
and
- iii) The OIEELISA<sub>NSS</sub> must always give a negative reaction.

#### 2.5.2.1.2. Infection with *Brucella* in sheep and goats:

- i) A 1/8 pre-dilution of the ISaBmS made up in a negative goat serum (or in a negative pool of goat sera) must give a positive reaction;  
and
- ii) a 1/300 pre-dilution of the ISaBmS made up in a negative goat serum (or in a negative pool of goat sera) must give a negative reaction;  
and
- iii) the above-mentioned negative goat serum (or negative pool of goat sera) must always give a negative reaction.

#### 2.5.2.1.3. Infection with *Brucella* in pigs:

- i) In the absence of an international standard serum for porcine brucellosis, the test should be duly validated and the cut-off established in the test population with appropriate validation techniques (see chapter 1.1.6).

The C-ELISA using an MAb specific for one of the epitopes of the *Brucella* sp. OPS has been shown in cattle, sheep and swine to have usually, but not always, higher specificity but lower sensitivity than the BBAT or I-ELISA (Muñoz *et al.*, 2005; 2012; Nielsen *et al.*, 1995; Praud *et al.*, 2012; Stack *et al.*, 1999).

It has been shown however that the C-ELISA eliminates some but not all FPSR caused by cross-reacting bacteria in cattle (Muñoz *et al.*, 2005) and swine (Praud *et al.*, 2012). In some cases, in ruminants or pigs, FPSR may be observed in C-ELISA while not in other S-LPS-based tests, including I-ELISA. Moreover, the C-ELISA reduces but not fully eliminates the reactions caused by antibodies produced in response to vaccination.

It is highly probable that much of the specificity improvement is due to a reduction in sensitivity of the C-ELISA compared with BBAT and I-ELISA. Therefore, positive C-ELISA reactions should be investigated using suitable confirmatory or complementary strategies as for BBAT, CFT and I-ELISA.

The choice of MAb and its more or less high and unique specificity and affinity will have a distinct influence on the diagnostic performance characteristics of the assay. As with any MAb-based assays, the universal availability of the MAb or the hybridoma must also be considered with respect to international acceptance and widespread use.

Several commercial C-ELISA are available. Some protocols are less sensitive or less specific than others; therefore results obtained from different assays are not always comparable. C-ELISA for diagnosing anti-*Brucella* antibodies in small ruminants and swine is essentially the same as that described for cattle, but the cut-off should have been properly established for these species using the appropriate validation techniques (see chapter 1.1.6; moreover, C-ELISA for sheep and goats should be standardised against the ISaBmS (McGiven *et al.*, 2011).

Whatever the C-ELISA format used:

- i) A positive and a negative control are included in each plate. OD ranges to be obtained with these two controls must be established to define the criteria for validating each plate results. The OD of the positive control is the one with which the OD of each test serum is compared to establish the final result (negative or positive).
- ii) An additional positive serum (internal control) must be included in each plate to validate the repeatability of the test from plate to plate and from day to day.

## 2.6. Fluorescence polarisation assay

The FPA is a simple technique for measuring antigen/antibody interaction. It is a homogeneous assay in which analytes are not separated and it is therefore very rapid. However, unlike another homogeneous method (e.g. RBT) a blank/background read is required for each sample before adding the antigen. Thus it is a two-step assay.

For the diagnosis of brucellosis, a small molecular weight fragment (average 22 kD) of the OPS of *B. abortus* strain *B. abortus* 1119-3 S-LPS is labelled with fluorescein isothiocyanate (FITC) and used as the antigen. Once the blank/background read is performed (2–3 minutes), this antigen is added to diluted serum and a measure of the antibody content is obtained in about 2 minutes after the addition of antigen using a fluorescence polarisation analyser (FPM) (Nielsen *et al.*, 1996).

### 2.6.1. Antigen production (example)

OPS from 5 g dry weight (or 50 g wet weight) of *B. abortus* S1119-3 is prepared by adding 400 ml of 2% (v/v) acetic acid, autoclaving the suspension for 15 minutes at 121°C and removing the cellular debris by centrifugation at 10,000 *g* for 10 minutes at 5°C ± 3°C. The supernatant solution is then treated with 20 g of trichloroacetic acid to precipitate any proteins and nucleic acids. The precipitate is again removed by centrifugation at 10,000 *g* for 10 minutes at 5°C ± 3°C. The supernatant fluid is dialysed against at least 100 volumes of purified water and freeze dried; 3 mg of OPS are dissolved in 0.6 ml of 0.1 M sodium hydroxide (4 g NaOH/litre) and incubated at 37°C ± 2°C for 1 hour, followed by the addition of 0.3 ml of FITC isomer 1 at a concentration of 100 mg/ml in dimethyl sulphoxide and a further incubation at 37°C ± 2°C for 1 hours. The conjugated OPS is applied to a 1 × 10 cm column packed with DEAE (diethylaminoethyl) Sephadex A 25 equilibrated in 0.01 M phosphate buffer, pH 7.4 ± 0.2. The first fraction (after 10–15 ml of buffer) is bright green, after which the buffer is switched to 0.1 M phosphate, pH 7.4 ± 0.2. This results in the elution of 10–15 ml of buffer followed by 25–40 ml of green fluorescent material. The latter material is the antigen used in the FPA. Antigen preparation may be scaled up proportionally.

The amount of antigen used per test is determined by diluting the material derived above until a total fluorescence intensity of 250,000–300,000 is achieved using the FPM.

The antigen can be stored as a liquid for several years at 5°C ± 3°C in a dark bottle or it may be freeze-dried in dark bottles. Labelled antigen may be obtained from a limited number of commercial sources.

## 2.6.2. Test standardisation of the FPA (EU, 2008; McGiven *et al.*, 2011)

### 2.6.2.1. Infection with *Brucella* in cattle

- i) the OIEELISA<sub>WPSS</sub> and OIEELISA<sub>SPSS</sub> consistently give a positive reaction;  
and
- ii) a 1/8 pre-dilution of the OIEELISA<sub>WPSS</sub> or a 1/64 pre-dilution of the OIEELISA<sub>SPSS</sub> made up in a negative bovine serum (or in a negative pool of bovine sera) must give a negative reaction;  
and
- iii) the OIEELISA<sub>NSS</sub> must always give a negative reaction.

### 2.6.2.2. Infection with smooth *Brucella* in sheep and goats

- i) a 1/16 pre-dilution of the ISaBmS made up in a negative goat serum (or in a negative pool of goat sera) must give a positive reaction;  
and
- ii) a 1/200 pre-dilution of the ISaBmS made up in a negative goat serum (or in a negative pool of goat sera) must give a negative reaction;  
and
- iii) the above-mentioned negative goat serum (or negative pool of goat sera) must always give a negative reaction.

### 2.6.2.3. Infection with *Brucella* in pigs

- i) In the absence of international standard serum for porcine brucellosis the test should be duly validated and the cut-off established in the test population with appropriate validation techniques (see chapter 1.1.6).

The diagnostic sensitivity and specificity of the FPA for bovine brucellosis are almost identical to those of the C-ELISA. The FPA for diagnosing anti-*Brucella* antibodies in small ruminants and swine is essentially the same as that described for cattle, but the cut-off should have been properly established for these species using the appropriate validation techniques (see chapter 1.1.6), and the test should be standardised against the corresponding international Standards, as mentioned above.

The FPA is capable of reducing but not fully eliminating the reactions due to residual antibody produced in response to vaccination (Nielsen *et al.*, 1996). Moreover, the specificity of FPA in FPSR conditions is currently unknown in cattle and small ruminants, but it has been clearly shown that it does not resolve the FPSR problem in swine (Praud *et al.*, 2012).

Accordingly, like all other serological tests, positive reactions should be investigated using suitable confirmatory and/or complementary strategies.

## 2.6.3. Test procedure (Nielsen *et al.*, 1996)

The FPA can be performed in glass tubes or a 96-well plate format.

Bovine sera are diluted 1/10 for the plate test or 1/100 for the tube test.

Sheep and goat as well as porcine sera are diluted 1/10 for the plate test or 1/25 (goat and porcine) and 1/40 (sheep) for the tube test.

The diluent used is 0.01 M Tris (1.21 g), containing 0.15 M sodium chloride (8.5 g), 0.05% Igepal CA630 (500 µl) (formerly NP40), 10 mM EDTA (3.73 g) per litre of purified water, pH 7.2 ± 0.2 (Tris buffer).

An initial reading to assess light scatter is obtained with the FPM after mixing. Suitably labelled titrated antigen is added, mixed and a second reading is obtained in the FPM about 2 minutes later.

A reading (in millipolarisation units, mP) over the established threshold level is indicative of a positive reaction.

A typical threshold level is 90–100 mP units, however, the test should be standardised locally against the corresponding OIE reference standard sera (as mentioned above). A strong positive, a weak positive and a negative working standard serum (standardised against the above-mentioned OIE reference standard sera) should be included.

#### 2.6.3.1. Example for bovine sera

- i) 1 ml of diluent is added to a 10 × 75 mm borosilicate glass tube followed by 10 µl of serum or 20 µl of EDTA-treated blood. For the 96-well format, 20 µl of serum is added to 180 µl of buffer. It is important to mix well. A reading is obtained on the FPM to determine light scatter.
- ii) A volume of antigen, which results in a total fluorescence intensity of 250–300 × 10<sup>3</sup>, is added to the tube and mixed well. This volume will vary from batch to batch, but is generally in the range of about 10 µl. A second reading is obtained on the FPM after incubation at room temperature (22°C ± 4°C) for approximately 2 minutes for serum and 15 seconds for EDTA-treated blood.
- iii) A reading above the predetermined threshold is indicative of a positive reaction.

### 2.7. Serum agglutination test (cattle only)

The SAT has been used with success for many years in surveillance and control programmes for bovine brucellosis, particularly in northern Europe. Its specificity is significantly improved with the addition of EDTA to the antigen (MacMillan & Cockrem, 1985).

The antigen represents a bacterial suspension of *B. abortus* strain 99 or *B. abortus* strain 1119-3 in phenol saline (NaCl 0.85% [w/v] and phenol at 0.5% [v/v]). Formaldehyde must not be used. Antigens may be delivered in the concentrated state provided the dilution factor to be used is indicated on the bottle label. EDTA may be added to the antigen suspension to 5 mM final test dilution to reduce the level of false-positive results. Subsequently the pH of 7.2 ± 0.2 must be readjusted in the antigen suspension.

The OIEISS contains 1000 IUs of agglutination. The antigen should be prepared without reference to the cell concentration, but its sensitivity must be standardised in relation to the OIEISS in such a way that the antigen produces either 50% agglutination with a final serum dilution of 1/600 to 1/1000 or 75% agglutination with a final serum dilution of 1/500 to 1/750. It may also be advisable to compare the reactivity of new and previously standardised batches of antigen using a panel of defined sera.

The test is performed either in tubes or in microplates. The mixture of antigen and serum dilutions should be incubated for 16–24 hours at 37°C ± 2°C. If the test is carried out in microplates, the incubation time can be shortened to 6 hours. At least three dilutions must be prepared for each serum in order to refute prozone negative responders. Dilutions of suspect serum must be made in such a way that the reading of the reaction at the positivity limit is made in the median tube (or well for the microplate method).

*Interpretation of results:* The degree of *Brucella* agglutination in a serum must be expressed in IU per ml. A serum containing 30 or more IU per ml is considered to be positive.

### 2.8. Native hapten and cytosol protein-based tests (ruminants only)

In cattle, native hapten tests<sup>6</sup> are highly specific in *B. abortus* S19 vaccination contexts, and have been used successfully in combination with the RBT as a screening test. The optimal sensitivity (close to that of CFT but significantly lower than that of RBT and S-LPS based I-ELISAs) is obtained in a reverse radial immunodiffusion (RID) system in which the serum diffuses into a hypertonic gel containing the polysaccharide, but the double gel diffusion assay is also useful (Muñoz *et al.*, 2005). Calves

---

<sup>6</sup> The detailed procedure can be obtained from the Brucellosis Laboratory, Centro de Investigación y Tecnología Agroalimentaria / Gobierno de Aragón, Avenida Montañana 930, 50059 Zaragoza, Spain.

vaccinated subcutaneously with the standard dose of *B. abortus* S19 at 3-5 months of age are negative usually by 2 months after vaccination, and adult cattle vaccinated subcutaneously 4–5 months previously with the reduced dose of *B. abortus* S19 do not give positive reactions unless the animals become infected and shed the vaccine strain in their milk. The conjunctival vaccination (both in young and adults) reduces significantly the time to obtain a negative response in native hapten tests. A characteristic of the native hapten tests is that a positive result correlates with *Brucella* shedding as shown in experimentally and in naturally infected cattle (Jones *et al.*, 1980). In case of FPSR caused by *Yersinia enterocolitica* O:9 and FPSR of unknown origin in cattle, gel precipitation tests using native hapten or *Brucella* cytosol proteins usually give negative results (Muñoz *et al.*, 2005).

These native hapten tests are also of interest in sheep and goats as they are very specific for discriminating the serological responses of infected animals (positive) from those induced in *B. melitensis* Rev.1 vaccinated animals (usually negative after a given time post-vaccination). The optimal diagnostic sensitivity (around 90%) is obtained in the double gel diffusion or RID tests for sheep and goats, respectively.

## 2.9. Milk tests

An efficient means of screening dairy herds is by testing milk from the bulk tank. When a positive test result is obtained, all cows contributing milk should be blood tested. The milk I-ELISA is a sensitive and specific test, and is particularly valuable for testing large herds. The MRT is a suitable alternative if the milk I-ELISA is not available. However, the MRT is not suitable in milk from small ruminants.

### 2.9.1. Milk I-ELISA (cattle and sheep and goats only)

As with the serum I-ELISA, many variations of the milk I-ELISA are possible. Several commercial I-ELISAs are available that have been validated in extensive field trials and are in wide use. In the interests of international harmonisation, the three OIE ELISA Standard Sera should be used by national reference laboratories to check or standardise a particular test for use in cattle.

#### 2.9.1.1. Test standardisation of the Milk I-ELISA in dairy cattle

The milk I-ELISA for cattle should be standardised such that the OIE ELISA strong positive standard when diluted 1/125 in negative serum and further diluted 1/10 in negative milk consistently tests positive (EU, 2008).

Bulk milk samples are generally tested at much lower dilutions than sera, i.e. undiluted to 1/2 to 1/10 in diluent buffer, with the remainder of the assay being similar to that described for serum.

FPSR may be observed in milk I-ELISA, but usually less frequently than with blood tests.

The I-ELISA for diagnosing anti-*Brucella* antibodies in sheep or goat milk is essentially the same as that described for cattle, but the cut-off should be properly established for these species using the appropriate validation techniques (see chapter 1.1.6). However there are no international standardisation recommendations of the milk I-ELISA against the corresponding ISaBmS.

### 2.9.2. Milk ring test (cattle only)

In lactating cattle, the MRT can be used for screening herds for brucellosis.

In large herds (> 100 lactating cows), the sensitivity of the test becomes less reliable. The MRT may be adjusted to compensate for the dilution factor from bulk milk samples from large herds. The samples are adjusted according to the following formula: herd size < 150 animals use, 1 ml bulk milk; 150–450 animals, use 2 ml milk sample, 451–700 animals, use 3 ml milk sample.

False-positive reactions may occur in cattle vaccinated less than 4 months prior to testing, in samples containing abnormal milk (such as colostrum) or in cases of mastitis. Therefore, it is not recommended to use this test in very small farms where these problems have a greater impact on the test results.

### 2.9.2.1. Antigen production

MRT antigen is prepared from concentrated, killed *B. abortus* S99 or S1119-3 cell suspension, grown as described previously. It is centrifuged at, for example, 23,000 *g* for 10 minutes at 5°C ± 3°C, followed by resuspension in haematoxylin-staining solution. Various satisfactory methods are in use; one example is as follows: 100 ml of 4% (w/v) haematoxylin (CI No. 75290) dissolved in 95% ethanol is added to a solution of ammonium aluminium sulphate (5 g) in 100 ml of purified water and 48 ml of glycerol; 2 ml of freshly prepared 10% (w/v) sodium iodate is added to the solution. After standing for 30 minutes at room temperature (22°C ± 4°C), the deep purple solution is added to 940 ml of 10% (w/v) ammonium aluminium sulphate in purified water. The pH of this mixture is adjusted to 3.1 ± 0.2, and the solution must be aged by storage at room temperature in the dark for 45–90 days.

Before use, the staining solution is shaken and filtered through cotton wool. The packed cells are suspended in the staining solution at the rate of 1 g per 30 ml stain, and held at room temperature (22°C ± 4°C) for 48 hours (some laboratories prefer to heat at 80°C for 10 minutes instead). The stained cells are then deposited by centrifugation, and washed three times in a solution of sodium chloride (6.4 g), 85% lactic acid (1.5 ml) and 10% sodium hydroxide (4.4 ml) in 1.6 litres of purified water, final pH 3.0 ± 0.2. The washed cells are resuspended at the rate of 1 g in 27 ml of a diluent consisting of 0.5% phenol saline, adjusted to pH 4.0 ± 0.2 by the addition of 0.1 M citric acid (approximately 2.5 ml) and 0.5 M disodium hydrogen phosphate (approximately 1 ml) and maintained at 5°C ± 3°C for 24 hours. The mixture is filtered through cotton wool, the pH is checked, and the PCV is determined and adjusted to approximately 4%.

### 2.9.2.2. Antigen standardisation

The antigen should be standardised against the OIEISS so that a 1/500 dilution is positive and 1/1000 dilution is negative. The sensitivity of the new batch should be compared as well with a previously standardised batch using a panel of samples of varying degrees of reaction prepared by diluting a positive serum in milk.

The antigen should be stored as recommended by the manufacturer but usually at 5°C ± 3°C.

The pH of the antigen should be 3.5 (± 0.2) and its colour should be dark blue. A little free stain in the supernatant of a centrifuged sample is permissible. When diluted in milk from a brucellosis-free animal, the antigen must produce a uniform coloration of the milk layer with no deposit and no coloration of the cream layer.

### 2.9.2.3. Test procedure

The test is performed on bulk tank milk samples. If necessary, samples could be pre-treated with preservative (0.1% formalin or 0.02% bronopol) for 2–3 days at 5°C ± 3°C prior to use.

- i) Bring the milk samples and antigen to room temperature (22°C ± 4°C); only sufficient antigen for the day's tests should be removed from the refrigerator.
- ii) Gently shake the antigen bottle well.
- iii) The test is performed by adding 30–50 µl of antigen to a 1–2 ml volume of whole milk (the volume of milk may be increased for bulk samples from larger herds – see above).
- iv) The height of the milk column in the tube must be at least 25 mm. The milk samples must not have been frozen, heated, subjected to violent shaking or stored for more than 72 hours.
- v) The milk/antigen mixtures are normally incubated at 37°C ± 2°C for 1 hour, together with positive and negative working standards. However, overnight incubation at 5°C ± 3°C increases the sensitivity of the test and allows for easier reading.
- vi) A strongly positive reaction is indicated by formation of a dark blue ring above a white milk column. Any blue layer at the interface of milk and cream should be considered to be positive as it might be significant, especially in large herds.
- vii) The test is considered to be negative if the colour of the underlying milk exceeds that of the cream layer.

- viii) When the MRT is adjusted for large herd sizes (2 or 3 ml of milk used), 0.1 ml of pooled negative cream is added to the test tube and is followed by 30–50 µl of the ring test antigen. After mixing, the test is incubated and read in the same manner as the unadjusted MRT. The negative pooled cream is collected from the separation of composite, unpasteurised milk from a brucellosis negative herd of 25 or more cows.

## 2.10. Serological tests in wildlife

Serological investigations in wild species are usually carried out for screening purposes. In these particular circumstances, adequate specificity is of paramount importance. The RBT can be recommended as a general purpose diagnostic test in all wildlife species. The CFT could also be recommended for such purpose, but the selection of the complement inactivation temperature and the cut-off titres have not been properly documented in all wildlife species. Both tests require the use of high quality serum samples that are not easy to obtain in wildlife studies. When poor quality serum samples are tested in both RBT and CFT, the results are frequently uninterpretable. The I- and C-ELISAs appear to be useful for epidemiological sero-surveys in wild animals as both are generally more reliable than both RBT and CFT, and, moreover, can be used with poor quality and haemolysed sera (Stack *et al.*, 1999). Another advantage of the ELISAs is that if serum is not available, it is possible to test meat juice samples. Attention must be paid to the conjugate used in the I-ELISA that must have a satisfactory affinity for the corresponding antibody isotypes of the wild species under study. However, in wild species, the interpretation of ELISA results may be problematic, due to the lack of validation studies. Whenever possible, the cut-off of ELISAs should be properly established for the particular species using the appropriate validation techniques (see chapter 1.1.6). Nevertheless, where positive or doubtful serological results are found, a bacteriological investigation should be conducted, when possible, to clarify the diagnosis.

## 3. Tests for cellular immunity

### 3.1. Brucellin skin test

An alternative immunological assay is the brucellin skin test, which can be used for screening unvaccinated herds, provided that a purified (free of S-LPS) and standardised antigen preparation is used. The brucellin skin test has a very high specificity, such that serologically negative unvaccinated cattle that are positive reactors to the brucellin test should be regarded as infected animals (Pouillot *et al.*, 1997). The brucellin skin test also has a high sensitivity for the diagnosis of *B. melitensis* infection in small ruminants and, in the absence of vaccination, is considered one of the most specific diagnostic tests.

Brucellin was developed for use in ruminants, but is also effective for confirming the disease at the herd level in pigs. Field trials have also shown its good sensitivity in *Brucella*-infected pigs (Dieste-Pérez *et al.*, 2014). Also, results of this test may aid the interpretation of serological reactions thought to be FPSR due to infection with cross-reacting bacteria, as FPSR affected animals always give negative results in the skin test (Dieste-Pérez *et al.*, 2014; Pouillot *et al.*, 1997).

Animals vaccinated with *B. melitensis* Rev.1, *B. abortus* S19 or RB51 can give positive results in this test for years (Pouillot *et al.*, 1997; De Massis *et al.*, 2005). Therefore this test cannot be recommended either as the sole diagnostic test or for the purposes of international trade in areas where *Brucella* vaccines are used. Moreover, not all infected animals react, therefore this test alone cannot be recommended as an individual diagnostic test or for the purposes of international trade. However, due to its high specificity and its adequate sensitivity at the herd or flock level, it can be recommended for herd/flock surveillance in brucellosis-free areas.

It is essential to use a standardised, defined brucellin preparation that does not contain S-LPS antigen, as this may provoke nonspecific inflammatory reactions or interfere with subsequent serological tests. One such preparation is brucellin INRA prepared from a rough strain of *B. melitensis* (strain B115) that is commercially available, but suitable alternative preparations from rough *Brucella* mutants could be also available in the near future (Dieste-Pérez *et al.*, 2014).

Although the brucellin test is probably the most specific indirect assay for diagnosing brucellosis (in unvaccinated animals), the final diagnosis should not be made solely on the basis of positive intradermal reactions given by a few animals in the herd, and should be supported by a complementary diagnostic test. The intradermal inoculation of brucellin induces a temporary anergy in the cellular immune response, at least in some animal species (Blasco *et al.*, 1994b). Therefore an interval of 6 weeks is generally recommended between two tests on the same animal.

### 3.1.1. Test procedure

#### 3.1.1.1. Infection with *Brucella* in cattle

- i) A volume of 0.1 ml of brucellin (2000 Units/ml) is injected intra-dermally into the caudal fold, the skin of the flank, or the side of the neck.
- ii) The test is read after 48–72 hours.
- iii) The skin thickness at the injection site is measured with vernier callipers before injection and at re-examination.
- iv) A strong positive reaction is easily recognised by local swelling and induration. However, borderline reactions require careful interpretation. Skin thickening of 1.5–2 mm would be considered as a positive reaction.

#### 3.1.1.2. Infection with *Brucella* in sheep and goats

- i) A volume of 0.1 ml of brucellin (2000 Units/ml) is injected intra-dermally into the lower eyelid.
- ii) The test is read after 48 hours.
- iii) Any visible or palpable reaction of hypersensitivity, such as an oedematous reaction leading to an elevation of the skin or thickening of the eyelid ( $\geq 2$  mm), should be interpreted as a positive reaction.

#### 3.1.1.3. Infection with *Brucella* in pigs

As a diagnostic agent in pigs, 0.1 ml of the allergen suspension (2000 Units/ml) is injected intra-dermally into the skin at the base of the ear or preferably, next to the base of the tail. The latter appears more practical and less hazardous. The reaction is assessed by visual inspection and palpation of the inoculated area after 48 hours and a positive reaction is characterised by erythema of non-pigmented skin and an oedematous swelling. In some cases, there may also be some haemorrhage/necrosis.

### 3.2. Interferon gamma release assay

The interferon gamma release assay involves stimulation of lymphocytes in whole blood with a suitable antigen such as brucellin. The resulting gamma interferon (IFN- $\gamma$ ) production is detected through a capture ELISA. This test could be useful in the discrimination of brucellosis from FPSR but more specific antigens are still needed and the protocol needs to be standardised and properly validated in the different animal species and epidemiological conditions. For the moment, a fully validated protocol is not available.

## C. REQUIREMENTS FOR VACCINES AND DIAGNOSTIC BIOLOGICALS

### 1. Vaccines

As mentioned previously, brucellosis is one of the most easily acquired laboratory infections, and strict safety precautions should be observed. Laboratory manipulation of live cultures of *Brucella*, including vaccine strains, is hazardous and must be performed at an appropriate biosafety and containment level determined by biorisk analysis (see Chapter 1.1.4). The S19, RB51 and Rev.1 vaccines have some virulence for humans, and a hazard warning should be included on the label of the final containers. Medical advice should be sought in the event of accidental inoculation or exposure (see section C.1.2.3.2.3 Precautions) (Ashford *et al.*, 2004; Joint FAO/WHO Expert Committee on Brucellosis, 1986; USDA, 2003).

#### 1.1. Background

##### 1.1.1. *Brucella abortus* strain 19 vaccine

A widely used vaccine for the prevention of brucellosis in cattle is *B. abortus* S19, which remains the reference vaccine with which any other vaccines must be compared. It is used as a live vaccine and is normally given to female calves aged between 3 and 6 months as a single subcutaneous dose of  $5\text{--}8 \times 10^{10}$  viable organisms. A reduced dose of from  $3 \times 10^8$  to  $5 \times 10^9$  organisms can be administered subcutaneously to adult cattle, but some animals can develop persistent antibody titres and may abort and excrete the vaccine strain in the milk. Alternatively, the vaccine can be administered to cattle of any age as either one or two doses of  $5 \times 10^9$

viable organisms, given by the conjunctival route. This vaccination procedure induces protection against both *B. abortus* (Nicoletti *et al.*, 1978) and *B. melitensis* (Jiménez de Bagües *et al.*, 1991) without a persistent antibody response and reduces the risks of abortion and excretion in milk when vaccinating adult cattle.

*Brucella abortus* S19 vaccine induces good immunity to moderate challenge by virulent *B. abortus* or *B. melitensis* organisms. The vaccine must be prepared from USDA-derived seed (see footnote 4 for address) and each batch must be checked for purity (absence of extraneous microorganisms), viability (live bacteria per dose) and smoothness (determination of dissociation phase). Seed lots for *B. abortus* S19 vaccine production should be regularly tested for residual virulence and immunogenicity in mice.

Control procedures for this vaccine are detailed in Section C.1.2.2.3 *In-process controls*.

### 1.1.2. *Brucella abortus* strain RB51 vaccine

Since 1996, *B. abortus* strain RB51 has become the official vaccine for prevention of brucellosis in cattle in several countries. However there is a disagreement over the protective performance of strain RB51 compared with strain S19 in cattle (Moriyón *et al.*, 2004). Each country uses slightly different methods to apply this vaccine. In the USA (a country that was almost free of bovine brucellosis before RB51 was introduced), calves are vaccinated subcutaneously between the ages of 4 and 12 months with  $1\text{--}3.4 \times 10^{10}$  viable organisms. Vaccination of cattle over 12 months of age is carried out only under authorisation from the State or Federal Animal Health Officials, and the recommended dose is  $1\text{--}3 \times 10^9$  viable organisms (USDA, 2003). In other countries, it is recommended to vaccinate cattle as calves (4–12 months of age) with a  $1\text{--}3.4 \times 10^{10}$  dose, with revaccination from 12 months of age onwards with a similar dose to elicit a booster effect and increase immunity.

However, it has been reported that full doses of *B. abortus* strain RB51 when administered intravenously in cattle induce severe placentitis and placental infection in most vaccinated cattle, and that there is excretion in milk in a relevant number of vaccinated animals. Field experience also indicates that it can induce abortion and increased perinatal mortality if applied to pregnant cattle. These observations have led to the recommendation to avoid vaccination of pregnant cattle with *B. abortus* RB51. One way to reduce the side effects of *B. abortus* RB51 is to reduce the dose. When using the reduced dose of this vaccine ( $1 \times 10^9$  colony-forming units [CFU]), on late pregnant cattle, no abortions or placentitis lesions have been reported, but the vaccine strain can be shed by a significant proportion of vaccinated animals. However, this reduced dose does not protect against *B. abortus* when used as a vaccine in calves, and the protection against *B. abortus* is only moderate when used as an adult vaccine. The protection conferred by *B. abortus* RB51 against *B. melitensis* infection in cattle is unknown.

Control procedures for this vaccine are detailed in Section C.1.2.2.3.

### 1.1.3. *Brucella melitensis* strain Rev.1 vaccine

It is not infrequent to isolate *B. melitensis* in cattle in countries with a high prevalence of this infection in small ruminants (Verger, 1985). There has been some debate on the protective efficacy of *B. abortus* S19 against *B. melitensis* infection in cattle, but there is published evidence proving that this vaccine is able to control *B. melitensis* in cattle. It has been hypothesised that *B. melitensis* Rev.1 should be a more effective vaccine than *B. abortus* S19 in these conditions. However there is very little information related to this issue (Joint FAO/WHO Expert Committee on Brucellosis, 1986), and no experiments have been reported showing the efficacy of Rev.1 against *B. melitensis* infection in cows. Moreover, the safety of Rev.1 vaccine is practically unknown in cattle. Accordingly, until the safety of Rev.1 in cattle of different physiological status and efficacy studies against *B. melitensis* under strictly controlled conditions are performed, this vaccine should not be recommended for use in cattle.

*Brucella melitensis* Rev.1 is the most widely used vaccine for the prevention of brucellosis in sheep and goats, and, despite its drawbacks, remains the reference vaccine with which any other vaccines should be compared. By contrast, the rough *B. abortus* RB51 vaccine is not effective against *B. melitensis* infection in sheep. The Rev.1 vaccine is used as a freeze-dried suspension of live *B. melitensis* Rev.1 strain for the immunisation of sheep and goats. It should be given to lambs and kids aged between 3 and 5 months as a single subcutaneous or conjunctival inoculation, 5 months being the upper time limit to minimise the antibody response to make this vaccination compatible with further serological testing. No matter the inoculation route, the standard dose must be between  $0.5 \times 10^9$  and  $2.0 \times 10^9$  viable organisms. The reduced doses confer a significantly lower protection than the standard doses, and should not be recommended for vaccinating sheep and goats. The subcutaneous vaccination induces long-

lasting serological responses, causing strong interferences in serological tests and should not be recommended for use in combined eradication programmes. However, when this vaccine is administered conjunctivally at the standard dose, it produces a similar protection without inducing a persistent antibody response, thus facilitating the application of eradication programmes combined with vaccination. Care must be taken when using *B. melitensis* Rev.1 vaccine to avoid the risk of contaminating the environment or causing human infection. In many developing countries and endemic areas, vaccination of the whole population has to be considered as the best option for the control of the disease (Blasco, 1997). However, Rev.1 vaccine is known to often cause abortion and excretion in milk when animals are vaccinated during pregnancy, either with a full or reduced dose (Blasco, 1997). These side-effects are considerably reduced when adult animals are vaccinated conjunctivally (full dose) during lambing/kidding, lactation or before mating. Therefore, when mass vaccination is the only means of controlling the disease, a vaccination campaign should be recommended using the standard dose of Rev.1 administered by the conjunctival route when the animals are not pregnant or during the late lambing/kidding and pre-breeding season (Blasco, 1997).

The subcutaneous vaccination of young animals and the vaccination of adult animals, even at reduced doses, may lead to long-term persistence of vaccinal antibodies in a significant proportion of Rev.1 vaccinated animals that creates serious interferences in the serological diagnosis of brucellosis. As indicated above, conjunctival vaccination minimises these problems (particularly when the upper limit of age for vaccination is 5 months) and thus it is the method of choice for combined eradication programmes. Therefore, the serological diagnosis of brucellosis should take into account the vaccinal state of the flock and the overall frequency distribution of antibody titres detected in the group of animals tested.

Control procedures for this vaccine are detailed in Section C.1.2.2.3.

#### 1.1.4. Vaccination in pigs

Attempts have been made to develop a suitable vaccine to immunise pigs against *B. suis*, but none has been found fully effective. Only one vaccine, *Brucella suis* strain 2 (S2), has been reported effective after extensive field use in China (People's Rep. of), but efficacy data against *B. suis* infection under strictly controlled conditions are not available. It is also reported that *B. abortus* strain RB51 vaccine is ineffective for the protection of swine against exposure to *B. suis* (Stoffregen *et al.*, 2006).

### 1.2. Outline of production and requirements

#### 1.2.1. Characteristics of the seed

##### 1.2.1.1. Biological characteristics of the master seed

*Brucella abortus* S19 original seed for vaccine production must be obtained from the USDA (see footnote 4 for address), and used to produce a seed lot that is preserved by lyophilisation or by freezing at liquid nitrogen temperature. The properties of this seed lot must conform to those of a pure culture of a CO<sub>2</sub>-independent *B. abortus* bv. 1 that is also sensitive to benzyl-penicillin (3 µg [5 IU]/ml), thionin blue (2 µg/ml) and i-erythritol (1 µg/ml), and that displays minimal pathogenicity for guinea-pigs.

*Brucella abortus* RB51 original seed for vaccine production is available commercially. It is also obtainable from the USDA (see footnote 4 for address). *Brucella abortus* RB51 has the normal properties of a bv. 1 strain of *B. abortus*, but is 100% in the rough phase and does grow in the presence of rifampicin (250 µg/ml).

*Brucella melitensis* strain Rev.1 original seed for vaccine production can be obtained commercially. A European reference Rev.1 strain that possesses the characteristics of the Rev.1 original seed is also obtainable from the OIE Reference Laboratory for Brucellosis in France (see footnote 5 for address). Strain Rev.1 must conform to the characteristics of *B. melitensis* bv. 1, except that it should grow more slowly. Additionally, when incubated in air (atmospheres containing CO<sub>2</sub> alter the results) at 37°C ± 2°C, it should grow on agar containing streptomycin (2.5 µg/ml), and it should be inhibited by the addition to a suitable culture medium of sodium benzyl-penicillin (3 µg [5 IU]/ml), thionin (20 µg/ml) or basic fuchsin (20 µg/ml).

PCR and molecular techniques have been used to further characterise the S19, RB51 or Rev.1 vaccines (see Section B.1.4).

The specific requirements for S19 and Rev.1 vaccine production recommend that each seed lot (i.e. the culture used to inoculate medium for vaccine production) should be no more than three passages removed from an original seed culture and that the harvest of a vaccine lot should be no more than three passages from a seed lot or an original seed. The original seed culture should always be checked for the absence of dissociation before use. The recommended method for preparing seed material is given in Alton *et al.* (1988).

#### 1.2.1.2. Quality criteria

*Brucella abortus* S19 and RB51 as well as *B. melitensis* Rev.1 master seeds should be checked for purity, identity and, where appropriate, smoothness or roughness. S19 and Rev.1 seed lots must also conform to the characteristics of residual virulence and immunogenicity in mice of the original seed.

##### 1.2.1.2.1. Purity

Tests for purity and freedom from contamination of biological materials may be found in Chapter 1.1.9 *Tests for sterility and freedom from contamination of biological materials intended for veterinary use*.

##### 1.2.1.2.2. Safety

The S19 and Rev.1 vaccines show reduced virulence, but should keep a minimal virulence to be efficient (see Section C.1.2.1.2.3 *Potency*). However a safety test is not routinely done. If desired, when a new manufacturing process is started and when a modification in the innocuousness of the S19 and Rev.1 vaccine preparations is expected, it may be performed on cattle (S19) and sheep and goats (Rev.1). This control should be done as follows: the test uses 12 female calves or sheep/goats respectively, aged 4–6 months. Six young females are injected with one or three recommended doses. Each lot of six young females are kept separately. All animals are observed for 21 days. No significant local or systemic reaction should occur. If, for a given dose and route of administration, this test gives good results on a representative batch of the vaccine, it does not have to be repeated routinely on seed lots or vaccine lots prepared with the same original seed and with the same manufacturing process. A safety test on S19/Rev.1 vaccine may also be performed in guinea-pigs. Groups of at least ten animals are given intramuscular injections of doses of vaccine diluted in PBS, pH 7.2 ± 0.2, to contain  $5 \times 10^9$  viable organisms. The animals should show no obvious adverse effects and there must be no mortality.

If this safety test has been performed with good results on a representative seed lot or batch of the test vaccine, it does not have to be repeated routinely on other vaccine lots prepared from the same seed lot and using the same manufacturing process.

If a safety test for RB51 is desired, 8- to 10-week-old female Balb/c mice can be injected intraperitoneally with  $1 \times 10^8$  CFUs and the spleens cultured at 6 weeks post-inoculation. Spleens should be free from RB51 and the mice should not develop anti-OPS antibodies.

##### 1.2.1.2.3. Potency

###### i) S19 vaccine

An S19 vaccine is effective if it possesses the characteristics of the S19 original strain, i.e. if it is satisfactory with respect to identity and smoothness. Moreover, it should have been produced with a given seed lot with adequate immunogenicity and residual virulence (Grilló *et al.*, 2000).

###### a) Identity

*Brucella abortus* present in the vaccine is identified by suitable morphological, serological and biochemical tests and by culture: *B. abortus* S19 has the normal properties of a *B. abortus* bv. 1 strain of, but does not require CO<sub>2</sub> for growth, does not grow in the presence of benzyl-penicillin (3 µg/ml = 5 IU/ml), thionin blue (2 µg/ml), and i-erythritol (1 mg/ml) (all final concentrations). PCR and molecular techniques have been described to identify the S19 vaccine strain (see Section B.1.4).

###### b) Smoothness (determination of dissociation phase)

The S19 vaccine reconstituted in purified water is streaked across six agar plates (serum–dextrose agar or trypticase–soy agar (TSA) with added serum

5% [v/v] or yeast extract 0.1 % [w/v] in such a manner that the colonies will be close together in certain areas, while semi-separated and separated in others. Slight differences in appearance are more obvious in adjacent than widely separated colonies. Plates are incubated at 37°C ± 2°C for 5 days and examined by obliquely reflected light (Henry's method) before and after staining (three plates) with crystal violet (White & Wilson's staining method).

*Appearance of colonies before staining:* S colonies appear round, glistening and blue to blue-green in colour. R colonies have a dry, granular appearance and are dull yellowish-white in colour. Muroid colonies (M) are transparent and greyish in colour and can be distinguished by their slimy consistency when touched with a loop. Intermediate colonies (I), which are the most difficult to classify, have an appearance intermediate between S and R forms: they are slightly opaque and more granular than S colonies.

*Appearance of colonies after staining with crystal violet:* S colonies do not take up the dye. Dissociated colonies (I, M, or R) are stained various shades of red and purple and the surface may show radial cracks. Sometimes a stained surface film slips off a dissociated colony and is seen adjacent to it.

The colony phase can be confirmed by the acriflavine agglutination test (Alton *et al.*, 1988). S colonies remain in suspension, whereas R colonies are agglutinated immediately and, if muroid, will form threads. Intermediate colonies may remain in suspension or a very fine agglutination may occur.

At least 99% of cells in seed lots should be in the smooth phase.

- c) Residual virulence (50% persistence time or 50% recovery time) (Grilló *et al.*, 2000; Pouillot *et al.*, 2004)
- 1) Prepare adequate suspensions of either the *B. abortus* S19 seed lot to be tested (test vaccine) and the S19 original seed culture (as a reference strain). For this, harvest 24–48 hours' growth of each strain in sterile buffered saline solution (BSS: NaCl 8.5 g; KH<sub>2</sub>PO<sub>4</sub> 1.0 g; K<sub>2</sub>HPO<sub>4</sub> 2.0 g; purified water 1000 ml; pH 6.8 ± 0.2) and adjust the suspension in BSS to 10<sup>9</sup> CFU/ml using a spectrophotometer (0.170 OD when read at 600 nm). The exact number of CFU/ml should be checked afterwards by plating serial tenfold dilutions on to adequate culture medium (blood agar base or TSA are recommended).
  - 2) Inject subcutaneously 0.1 ml (10<sup>8</sup> CFU/mouse) of the suspension containing the test vaccine into each of 32 female CD1 mice, aged 5–6 weeks. Carry out, in parallel, a similar inoculation in another 32 mice using the suspension containing the S19 reference strain. The original seed S19 strain, which has been shown satisfactory with respect to immunogenicity and/or residual virulence, can be obtained from USDA (see footnote 4 for address).
  - 3) Kill the mice by cervical dislocation, in groups of eight selected at random 3, 6, 9 and 12 weeks later.
  - 4) Remove the spleens and homogenise individually and aseptically with a glass grinder (or in adequate sterile bags with a Stomacher) in 1 ml of sterile BSS.
  - 5) Spread each whole spleen suspension *in toto* onto several plates containing a suitable culture medium and incubate in standard *Brucella* conditions for 5–7 days (lower limit of detection: 1 bacterium per spleen). An animal is considered infected when at least 1 CFU is isolated from the spleen.
  - 6) Calculate the 50% persistence time or 50% recovery time (RT<sub>50</sub>) by the SAS® statistical method specifically developed for RT<sub>50</sub> calculations<sup>7</sup>. For this, determine the number of cured mice (no colonies isolated in the spleen) at each slaughtering point time (eight mice per point) and calculate the percentage of cured accumulated mice over time, by the Reed and Muench method (see ref. cited in Grilló *et al.*, 2000). The function of

---

7 To obtain the specific SAS® file contact the OIE Reference Laboratory for Brucellosis in France.

distribution of this percentage describes a sigmoid curve, which must be linearised for calculating the RT<sub>50</sub> values, using the computerised PROBIT procedure of the SAS® statistical package.

- 7) Compare statistically the parallelism (intercept and slope) between the distribution lines obtained for both tested and reference S19 strains using the SAS® file specifically designed for this purpose. Two RT<sub>50</sub> values can be statistically compared exclusively when they come from parallel distribution lines. If parallelism does not exist, the residual virulence of the tested strain should be considered inadequate, and discarded for vaccine production.
- 8) If the parallelism is confirmed, compare statistically the RT<sub>50</sub> values obtained for both tested and reference S19 strains using a SAS® file specifically designed for this purpose. To be accepted for vaccine production, the RT<sub>50</sub> obtained with the tested strain should not differ significantly from that obtained with the reference S19 strain (RT<sub>50</sub> and confidence limits are usually around 7.0 ± 1.3 weeks).

The underlying basis of the statistical procedure for performing the above residual virulence calculations have been described in detail (see ref. cited in Grilló *et al.*, 2000). Alternatively, the statistical calculations described in steps 6 to 8 above can be avoided by an easy-to-use specific HTML-JAVA script program (Rev2) available free (see footnote 5 for address) (Pouillot *et al.*, 2004).

If this test has been performed with good results on a representative seed lot, it does not have to be repeated routinely on other seed lots and vaccine batches prepared from the same seed strain and using the same manufacturing process.

d) Immunogenicity in mice (Grilló *et al.*, 2000)

This test uses three groups of six female CD1 mice, aged 5–7 weeks that have been selected at random.

- 1) Prepare and adjust spectrophotometrically the vaccine suspensions as indicated above.
- 2) Inject subcutaneously a suspension containing 10<sup>5</sup> CFU (in a volume of 0.1 ml/mouse) of the vaccine to be examined (test vaccine) into each of six mice of the first group.
- 3) Inject subcutaneously a suspension containing 10<sup>5</sup> CFU of live bacteria of a reference S19 vaccine into each of six mice of the second group. The third group will serve as the unvaccinated control group and should be inoculated subcutaneously with 0.1 ml of BSS.
- 4) The exact number of CFU inoculated should be checked afterwards by plating serial tenfold dilutions onto adequate culture medium (blood agar base or TSA are recommended).
- 5) All the mice are challenged 30 days after vaccination (and immediately following 16 hours starvation), intraperitoneally with a suspension (0.1 ml/mouse) containing 2 × 10<sup>5</sup> CFU of *B. abortus* strain 544 (CO<sub>2</sub>-dependent), prepared, adjusted and retrospectively checked as above.
- 6) Kill the mice by cervical dislocation 15 days later.
- 7) Each spleen is excised aseptically, the fat is removed, and the spleen is weighed and homogenised.
- 8) Alternatively, the spleens can be frozen and kept at ≤ -16°C for from 24 hours to 7 weeks.
- 9) Each spleen is homogenised aseptically with a glass grinder (or in adequate sterile bags in Stomacher) in nine times its weight of BSS, pH 6.8 ± 0.2 and three serial tenfold dilutions (1/10, 1/100 and 1/1000) of each homogenate made in the same diluent. Spread 0.2 ml of each dilution by quadruplicate in agar plates and incubate two of the plates in a 10% CO<sub>2</sub> atmosphere (allows the growth of both vaccine and challenge

strains) and the other two plates in air (inhibits the growth of the *B. abortus* 544 CO<sub>2</sub>-dependent challenge strain), both at 37°C± 2°C for 5 days.

- 10) Colonies of *Brucella* (*B. abortus* 544) should be enumerated on the dilutions corresponding to plates showing fewer than 300 CFU. When no colony is seen in the plates corresponding to the 1/10 dilution, the spleen is considered to be infected with five bacteria. These numbers of *Brucella* per spleen are first recorded as X and expressed as Y, after the following transformation:  $Y = \log(X/\log X)$ . Mean and standard deviation, which are the response of each group of six mice, are then calculated.
- 11) The conditions of the control experiment are satisfactory when: i) the response of unvaccinated mice (mean of Y) is at least of 4.5; ii) the response of mice vaccinated with the reference S19 vaccine is lower than 2.5; and iii) the standard deviation calculated on each lot of six mice is lower than 0.8.
- 12) Carry out the statistical comparisons (the least significant differences test is recommended) of the immunogenicity values obtained in mice vaccinated with the S19 strain to be tested with respect to those obtained in mice vaccinated with the reference vaccine and in the unvaccinated control group. The test vaccine would be satisfactory if the immunogenicity value obtained in mice vaccinated with this vaccine is significantly lower than that obtained in the unvaccinated controls and, moreover, does not differ significantly from that obtained in mice vaccinated with the reference vaccine (for detailed information on this procedure, footnote 5 gives the contact address).

If this test has been performed with good results on a representative seed lot, it does not have to be repeated routinely on other seed lots and vaccine batches prepared from the same seed strain and with the same manufacturing process.

ii) RB51 vaccine

An RB51 vaccine seed lot must possess the characteristics of the RB51 original strain, i.e. if it is satisfactory with respect to identity, roughness and potency.

a) Identity

The reconstituted RB51 vaccine should not contain extraneous microorganisms. *Brucella abortus* present in the vaccine is identified by suitable morphological, serological and biochemical tests and by culture: *B. abortus* RB51 has the normal properties of a *B. abortus* bv. 1 strain, but is 100% in the rough phase and does grow in the presence of rifampicin (250 µg/ml). PCR and molecular techniques have been described to further characterise the RB51 vaccine strain (see Section B.1.4).

b) Roughness (determination of dissociation phase)

The same technical procedures indicated for S19 vaccine (see Section C. 1.2.1.2.3.i *S19 vaccine* above) have to be applied for RB51. 100% of the RB51 cells must be in the rough phase. Additionally, for RB51, all colonies should be negative on dot-blot assays with MAbs specific for the OPS antigen.

c) Potency

As dosage (CFU) of the master seed was correlated to protection as part of registration of RB51 for use cattle in the USA, *in-vivo* potency tests are not routinely conducted for serials of the RB51 vaccine. In the USA, plate counts of viable organisms have been approved and used as a measure of potency (this approach is identical to the potency test for S19 vaccine in the USA (USDA, 2003). Rough vaccines for brucellosis have been discussed in some detail (Moriyón *et al.*, 2004).

iii) Rev.1 vaccine

A Rev.1 vaccine is efficient if it possesses the characteristics of the Rev.1 original strain, i.e. if it is satisfactory with respect to identity and smoothness. Moreover, it should have been produced with a given seed lot with adequate immunogenicity, and residual virulence (Grilló *et al.*, 2000).

a) Identity

*Brucella melitensis* present in the vaccine is identified by suitable morphological, serological and biochemical tests and by culture: when incubated in air at  $37^{\circ}\text{C} \pm 2^{\circ}\text{C}$ , Rev.1 strain is inhibited by addition to the suitable culture medium of  $3 \mu\text{g}$  (5 IU) per ml of sodium benzyl-penicillin, thionin ( $20 \mu\text{g/ml}$ ) or basic fuchsin ( $20 \mu\text{g/ml}$ ); the strain grows on agar containing  $2.5 \mu\text{g}$  per ml of streptomycin. PCR and molecular techniques have been described to further characterise the Rev.1 vaccine strain (see Section B.1.4).

b) Smoothness (determination of dissociation phase)

The same technical procedures indicated for S19 vaccine (Section C. 1.2.1.2.3.i above) have to be applied for Rev.1.

Sometimes, slight and difficult to observe differences can be seen in the size of Rev.1 colonies. The small colonies (1–1.2 mm in diameter) are typical for Rev.1, but larger Rev.1 colonies can appear depending on the medium used, the amount of residual moisture in the incubator atmosphere, and the presence or absence of  $\text{CO}_2$ . The frequency of variation in colony size occurs normally at a ratio of 1 large to  $10^3$  small colonies. Both Rev.1 variants are of the S (smooth) type. To avoid an increase in this colony size variation along successive passages, it is important to always select small colonies for preparation of seed lots.

At least 99% of cells in seed lots should be in the smooth phase.

c) Residual virulence (50% persistence time or 50% recovery time) (Grilló *et al.*, 2000)

The same technical procedures indicated for  $\text{RT}_{50}$  calculation of S19 vaccine (see above) have to be applied for Rev.1, except that *B. abortus* S19 seed lot to be tested (test vaccine) and the S19 original seed culture (used as a reference strain), respectively, are replaced by the corresponding *B. melitensis* Rev.1 seed lot to be tested (test vaccine) and the *B. melitensis* Rev.1 original seed culture as the reference strain. For the reference original Rev.1 strain,  $\text{RT}_{50}$  and confidence limits are around  $7.9 \pm 1.2$  weeks. A given Rev.1 vaccine seed lot or batch should keep similar residual virulence to be acceptable.

If this test has been done with good results on a representative seed lot, it does not have to be repeated routinely on seed lots and vaccine batches prepared from the same seed strain and with the same manufacturing process.

d) Immunogenicity in mice

The same technical procedures indicated for immunogenicity calculation of S19 vaccine (see above) have to be applied for Rev.1, except that *B. abortus* S19 seed lot to be tested (test vaccine) and the *B. abortus* S19 original seed culture (used as a reference strain), respectively, are replaced by the corresponding *B. melitensis* Rev.1 seed lot to be tested (test vaccine) and the *B. melitensis* Rev.1 original seed culture as the reference strain.

Conditions of the control experiment are satisfactory when: i) the response in unvaccinated mice (mean of Y) is at least of 4.5; ii) the response in mice vaccinated with the reference Rev.1 vaccine is lower than 2.5; and iii) the standard deviation calculated on each lot of six mice is lower than 0.8.

If this test has been done with good results on a representative seed lot, it does not have to be repeated routinely on seed lots and vaccine batches prepared from the same seed strain and with the same manufacturing process.

#### 1.2.1.3. Validation as a vaccine

Numerous independent studies have confirmed the value of S19 as a vaccine for protecting cattle from brucellosis, and it has been the vaccine used (in combination with serological testing and culling) to eradicate bovine brucellosis in most currently free countries. The organism behaves as an attenuated strain when given to sexually immature cattle. In rare cases, when standard doses are applied subcutaneously, it may produce

localised infection in the genital tract particularly in males. For this reason, vaccination of males is counter-indicated. Antibody responses persisting for 6 months or longer are likely to occur in a substantial proportion of cattle that have been vaccinated subcutaneously with the standard dose as adults. Some of the cattle vaccinated subcutaneously as calves may later develop arthropathy, particularly of the femoro-tibial joints. The vaccine is safe for most animals if administered to calves between 3 and 6 months of age. It may also be used in adult animals at a reduced dose, with the advantage of reducing the diagnostic interferences. It produces lasting immunity to moderate challenge with virulent *B. abortus* strains, but the precise duration of this immunity is not well known. The length of protection of this vaccine against *B. melitensis* infection in cattle is also unknown. The vaccine strain is stable and reversion to virulence is extremely rare. It has been associated with the emergence of i-erythritol-using strains when inadvertently administered to pregnant animals. The organism behaves as an attenuated strain in mice, and even large inocula are rapidly cleared from the mice tissues.

Reports from both experimental challenge studies and field studies remain controversial as far as the value of *B. abortus* strain RB51 in protecting cattle from brucellosis is concerned (see above). The organism is attenuated in calves but can cause safety problems in adults. *Brucella abortus* strain RB51 contains minimally expressed OPS but there is no serological conversion in both standard RBT and CFT in vaccinated animals. In addition, it has also been claimed that RB51 does not induce detectable antibodies, using current OPS-based testing procedures (USDA, 2003). However, the presence of common core epitopes in both smooth and rough *Brucella* does not allow always the antibody response to RB51 to be distinguished from that induced by field smooth *Brucella* strains, no matter which S-LPS or OPS-based ELISA is used. The efficacy of RB51 against *Brucella* infection in cattle is controversial (Moriyón *et al.*, 2004), but it is claimed that it protects against moderate challenge with virulent *B. abortus*, although the precise duration of this protection is unknown. The efficacy of this vaccine against *B. melitensis* infection in cattle is also unknown. The vaccine is very stable and no reversion to smoothness has been described *in vivo* or *in vitro*. The organism behaves as an attenuated strain in a variety of animals including mice where it is rapidly cleared from the tissues.

Numerous independent studies have confirmed the value of *B. melitensis* strain Rev.1 as a vaccine for protecting sheep and goats from brucellosis. Its virulence is unchanged after passage through pregnant sheep and goats. However, abortions and excretion in milk may result when the Rev.1 vaccine is inoculated into pregnant ewes and goats. The vaccine-induced abortions are not avoided using reduced doses and doses as low as  $10^6$ , used either subcutaneously or conjunctivally in pregnant animals, have been proven to induce abortions and milk excretion of the vaccine strain (Blasco, 1997).

## 1.2.2. Method of manufacture

### 1.2.2.1. Procedure

Production of *Brucella* live vaccines is based on the seed-lot system described above (Section B.2.2) for BBAT and CFT antigens.

For the production of S19 vaccine, the procedures described above can be used, except that the cells are collected in PBS, pH  $6.3 \pm 0.2$ , and deposited by centrifugation or by the addition of sodium carboxymethyl cellulose at a final concentration of 1.5 g/litre. The yield from one fermenter run or the pooled cells from a batch of Roux flask cultures that have been inoculated at the same time from the same seed lot constitutes a single harvest. More than one single harvest may be pooled to form a final bulk, which is used to fill the final containers of a batch of vaccine. Before pooling, each single harvest must be checked for purity, cell concentration, dissociation and identity. A similar range of tests must be done on the final bulk, which should have a viable count of between 8 and  $24 \times 10^9$  CFU/ml. Adjustments in concentration are made by the addition of PBS for vaccine to be dispensed in liquid form, or by the addition of stabiliser for lyophilised vaccine. If stabiliser is to be used, loss of viability on lyophilisation should be taken into account, and should not be in excess of 50%. The final dried product should not be exposed to a temperature exceeding 35°C during drying, and the residual moisture content should be 1–2%. The contents must be sealed under vacuum or dry nitrogen immediately after drying, and stored at  $5^\circ\text{C} \pm 3^\circ\text{C}$ .

The production process for *B. abortus* strain RB51 is very similar to the one used for S19.

For the production of *B. melitensis* strain Rev.1 vaccine, the procedures described above for antigens (Alton *et al.*, 1988) can be used except that the cells are collected in a freeze-drying stabiliser and deposited by centrifugation. The yield from one fermenter run or the pooled cells from a batch of Roux flask cultures inoculated on the same occasion from the same seed lot constitutes a single harvest. More than one single harvest may be pooled to form the final bulk that is used to fill the final containers of a batch of vaccine. Before pooling, each single harvest must be checked for purity, cell concentration, dissociation and identity. The volume of the final bulk is adjusted by adding sufficient stabiliser so that a dose contains an appropriate number of viable organisms. After adjusting the cell concentration of the final bulk, tests for identity, dissociation and absence of contaminating organisms are conducted (see below).

#### 1.2.2.2. Requirements for ingredients

Strains should be cultured in a suitable medium.

*Brucella abortus* S19 and *B. melitensis* strain Rev.1 for vaccine production is grown on medium free from serum or other animal products, under conditions similar to those described above for *B. abortus* S99 or S1119-3 (Alton *et al.*, 1988). The phenol saline is replaced by a freeze-drying stabiliser.

*Brucella abortus* strain RB51 follows similar culture methods.

Serum–dextrose agar, and trypticase–soy agar, to which 5% serum or 0.1% yeast extract may be added, are among the solid media that have been found to be satisfactory for propagating the Rev.1 strain (Alton *et al.*, 1988). However, Rev.1 strain does not grow well on potato agar and generally needs 3–5 days to grow.

For all vaccines, the organisms are not killed but are stored at 5°C ± 3°C while quality control examinations are carried out as described below.

For preparation of the lyophilised vaccine, a stabiliser containing 2.5% casein digest, e.g. Tryptone (Oxoid), 5% sucrose and 1% sodium glutamate, dissolved in purified water and sterilised by filtration is recommended.

Antimicrobial preservatives must not be used in live *B. abortus* strain S19 or RB51 and *B. melitensis* strain Rev.1 vaccines.

#### 1.2.2.3. In-process controls

*Brucella abortus* S19 and RB51 as well as *B. melitensis* Rev.1 vaccines should be checked for purity, identity and, where appropriate, smoothness or roughness during preparation of the single harvests. The cell concentration of the bulks should also be checked. This can be done by opacity measurement, but a viable count must be performed on the final filling lots. The identity of these should also be checked by agglutination tests with antiserum to either *Brucella*-A, -R or -M antigen respectively.

The viable count of the final containers should not be less than the recommended doses (see above).

For S19 and Rev.1, at least 99% of cells in seed lots and 95% of cells in final lots should be in the smooth phase, while 100% of the RB51 cells must be in the rough phase. Additionally, for RB51, all colonies should be negative on dot-blot assays with monoclonal antibodies specific for the OPS antigen.

For S19 and Rev.1, the immunogenicity and the residual virulence (50% persistence time or 50% recovery time) in the mice model should also be determined on representative seed lots. If these tests have been done with good results on a representative seed lot, it does not have to be repeated routinely on seed lots and vaccine batches prepared from the same seed strain and with the same manufacturing process.

#### 1.2.2.4. Final product batch tests

With freeze-dried vaccine, the control tests should be conducted on the vaccine reconstituted in the form in which it will be used (same diluent).

1.2.2.4.1. Purity

Tests for purity and freedom from contamination of biological materials may be found in chapter 1.1.9.

1.2.2.4.2. Identity

See Section C.1.2.1.2.3 *Potency*.

1.2.2.4.3. Safety

See Section C.1.2.1.2.2 *Safety*.

1.2.2.4.4. Batch potency

i) Potency

For S19 and Rev.1 vaccines, potency can also be determined on the final lyophilised product. The procedure is as described above (identity; smoothness; residual virulence and immunogenicity checks; see Section C.1.2.1.2.3.). If residual virulence and immunogenicity checks have been performed with good results on a representative batch of the test vaccine, they do not have to be repeated routinely on other vaccine lots prepared from the same seed lot and using the same manufacturing process.

ii) Enumeration of live bacteria

Batches should also be checked for the number of viable organisms. The same procedure may be applied for the S19, Rev.1 and RB51 vaccine batches. Inoculate each of at least five plates of tryptose, serum-dextrose or other suitable agar medium with 0.1 ml of adequate dilutions of the vaccine spread with a sterile glass, wire or plastic spreader. CFU per vaccine volume unit are enumerated.

Suitable CFU counts are the following:

S19:

- a)  $0.5-1 \times 10^{11}$  CFU (standard dose; subcutaneous route);
- b)  $0.5-5 \times 10^9$  CFU (reduced dose; subcutaneous route);
- c)  $5 \times 10^9$  CFU (reduced dose; conjunctival route).

Rev.1:

- a)  $0.5-2 \times 10^9$  CFU (standard dose, subcutaneous or conjunctival route).

RB51:

- a)  $1-3.4 \times 10^{10}$  CFU (standard dose; subcutaneous route).

### 1.2.3. Requirements for authorisation/registration/licensing

#### 1.2.3.1. Manufacturing process

For registration of the vaccine, all relevant details concerning manufacture of the vaccine and quality control testing (see above) should be submitted to the authorities. This information should be provided from three consecutive vaccine batches with a volume of not less than 1/3 of the typical industrial batch volume.

#### 1.2.3.2. Safety requirements

##### 1.2.3.2.1. Target and non-target animal safety

For the potential side-effects of the *Brucella* vaccines according to the status of the animals, see Section C.1.1. *Background*.

#### 1.2.3.2.2. Reversion-to-virulence

*Brucella abortus* S19 and *B. melitensis* Rev.1 vaccines prepared from seed stock from appropriate sources are stable in characteristics, provided that the in-process and batch control requirements described above are fulfilled, and show no tendency to reversion to virulence.

*Brucella abortus* strain RB51 has shown no tendency to revert to virulent smooth organisms after many passages *in vitro* or *in vivo*. This is probably due to the nature and place of the mutations found in this strain. *Brucella abortus* strain RB51, despite carrying, among other unknown mutations, an IS711-disrupted *wboA* (a putative glycosyl-transferase gene), accumulates low amounts of cytoplasmic M-like OPS.

#### 1.2.3.2.3. Precautions

*Brucella abortus* S19 and RB51, as well as *B. melitensis* Rev.1, although attenuated strains, are still capable of causing disease in humans. Accordingly cell cultures and suspensions must be handled under appropriate conditions of biohazard containment (see chapter 1.1.4). Reconstitution and subsequent handling of the reconstituted vaccine should be done with care to avoid accidental injection or eye or skin contamination. Vaccine residues and injection equipment should be decontaminated with a suitable disinfectant (phenolic, iodophor or aldehyde formulation) at recommended concentration. Medical advice should be sought in the event of accidental exposure. The efficacy of the antibiotic treatment of infections caused by S19, RB51 or Rev.1 in humans has not been fully established. However, it must be reiterated that, while the S19 strain carries no particular antibiotic resistance compared with other *Brucella* field strains, Rev.1 and RB51 strains are respectively streptomycin- and rifampicin-resistant.

Vaccine should be identified as pathogenic for vaccinators. Manufacturers should provide adequate warnings that medical advice should be sought in the case of self-injection or exposure to vaccine (including aerosols) with warnings included on the product label/leaflet so that the vaccinator is aware of any danger.

#### 1.2.3.3. Efficacy requirements

Potency can also be determined on the final batch, but if safety/efficacy tests have been performed with good results on a representative seed lot or a batch of the test vaccine, it does not have to be repeated routinely on other vaccine lots prepared from the same seed lot and using the same manufacturing process.

#### 1.2.3.4. Duration of immunity

Vaccinating calves with a full dose of S19 vaccine is considered to give long-lasting immunity, and subsequent doses are not recommended. However, there is scanty evidence for this and revaccination within 6–12 months could be advisable in endemic areas.

The duration of immunity induced by RB51 vaccine in cattle is unknown, whatever the dose applied and the age at vaccination. Revaccination within 6–12 months has been proposed for boosting immunity in endemic areas.

It is accepted that subcutaneous or conjunctival vaccination with standard doses of Rev.1 confers a solid and durable immunity in sheep and goats. However, growing field evidence shows that the immunity conferred declines with time, and revaccination within 6–12 months could be advisable in endemic areas.

#### 1.2.3.5. Stability

*Brucella abortus* S19 and *B. melitensis* Rev.1 vaccines prepared from seed stock from appropriate sources are stable in characteristics, provided that the in-process and batch control requirements described above are fulfilled, and show no tendency to reversion to virulence. The lyophilised vaccine shows a gradual loss of viable count, but should retain its potency for the recommended shelf life. Allowance for this phenomenon is normally made by ensuring that the viable count immediately following lyophilisation is well in excess of the minimum requirement. Maintenance of a cold chain during distribution of the vaccine will ensure its viability.

*Brucella abortus* strain RB51 has shown no tendency to revert to virulent smooth organisms after many passages *in vitro* or *in vivo*. This is probably due to the nature and place of the mutations found in this strain. *Brucella abortus* strain RB51, despite carrying, among other unknown mutations, an IS711-disrupted *wboA* (a putative glycosyl-transferase gene), accumulates low amounts of cytoplasmic M-like OPS.

## 2. Diagnostic biologicals: brucellin

### 2.1. Background

Brucellin-INRA is an LPS-free extract from rough *B. melitensis* B115, and a single inoculation of this preparation does not provoke formation of antibodies reactive in BBAT, CFT or ELISAs. However, as this rough strain contains *Brucella* OPS sugars in the cytosol extract, the repeated inoculation of brucellin could elicit antibodies, interfering with other diagnostic tests. For this reason, cytosolic protein extracts have been obtained from rough *B. abortus* mutants, defective in genes strictly necessary to synthesise perosamine, then unable to generate OPS antibody response in sheep. These cytosolic protein extracts have been successfully tested in comparative studies with brucellin for the diagnosis of brucellosis in pigs (Dieste-Pérez *et al.*, 2014). These extracts could advantageously replace classical brucellin in the near future.

### 2.2. Outline of production and requirements

#### 2.2.1. Characteristics of the seed

##### 2.2.1.1. Biological characteristics of the master seed and quality criteria

Production of brucellin-INRA is based on a seed-lot system as described for antigens and vaccines. The original seed *B. melitensis* strain B115 for brucellin production<sup>5</sup> should be propagated to produce a seed lot, which should be preserved by lyophilisation or freezing at liquid nitrogen temperature. It should conform to the properties of a pure culture of a rough strain of *B. melitensis* and must not produce smooth *Brucella* LPS. It should produce reasonable yields of a mixture of protein antigens reactive with antisera to smooth and rough *Brucella* strains.

##### 2.2.1.2. Quality criteria

*Brucella melitensis* B115 seed should be checked for purity.

Tests for purity and freedom from contamination of biological materials may be found in chapter 1.1.9.

##### 2.2.1.3. Validation as an *in-vivo* diagnostic reagent

Laboratory and field studies in France have confirmed that brucellin-INRA is safe, non-toxic and specific in action. The preparation contains 50–75% proteins, mainly of low molecular weight, and 15–30% carbohydrate. It does not contain S-LPS antigens, does not provoke inflammatory responses in unsensitised animals, and it is not in itself a sensitising agent. After a single inoculation, it does not induce detectable antibodies in the standard serological tests for brucellosis. More than 90% of small ruminants infected with *B. melitensis* manifest delayed hypersensitivity to brucellin-INRA at some stage. The preparation is not recommended as a diagnostic agent for individual animals, but can be useful when used for screening herds. It is given to small ruminants in 100-µg doses by the intradermal route, and provokes a local delayed hypersensitivity reaction visible at 48–72 hours in sensitised animals. Positive reactions can be given by vaccinated as well as by infected animals (Pouillot *et al.*, 1997).

#### 2.2.2. Method of manufacture (Alton *et al.*, 1988)

##### 2.2.2.1. Procedure and requirements for ingredients

Brucellin is produced from *B. melitensis* strain B115 according to Alton *et al.*, 1988.

##### 2.2.2.2. In-process control

The crude brucellin extract should be checked for sterility after acetone extraction, to ensure killing of *Brucella* cells, and again at the end of the process to check possible contamination. The pH and protein concentration should be determined, and control tests

(as described in Section C.2.2.2.3 below) should be performed on the bulk material before filling the final containers.

### 2.2.2.3. Final product batch tests

i) Sterility

Brucellin preparations should be checked for sterility as described in chapter 1.1.9.

ii) Safety

Brucellin preparations should also be checked for abnormal toxicity. Doses equivalent to 20 cattle doses (2 ml) should be injected intraperitoneally into a pair of normal guinea-pigs that have not been exposed previously to *Brucella* organisms or their antigens. Five normal mice are also inoculated subcutaneously with 0.5 ml of the brucellin (2000 U/ml) to be examined. Animals are observed for 7 days, and there should be no local or generalised reaction to the injection.

Dermo-necrotic capacity is examined by intradermal inoculation of 0.1 ml of the product to be tested into the previously shaved and disinfected flank of three normal albino guinea-pigs that have not been exposed previously to *Brucella* organisms or their antigens. No cutaneous reaction should be observed.

Absence of allergic and serological sensitisation is checked by intradermal inoculation of three normal albino guinea-pigs, three times every 5 days, with 0.1 ml of a 1/10 dilution of the preparation to be examined. A fourth similar injection is given, 15 days later, to the same three animals and to a control lot of three guinea-pigs of the same weight that have not been injected previously. The animals should not become seropositive to the standard tests for brucellosis (RBT, CFT) when sampled 24 hours after the last injection, and should not develop delayed hypersensitivity responses.

iii) Batch potency

The potency of brucellin preparations is determined by intradermal injection of graded doses of brucellin into guinea-pigs that have been sensitised by subcutaneous inoculation of 0.5 ml of reference brucellin<sup>8</sup> in Freund's complete adjuvant from 1 to 6 months previously (the use of a live *Brucella* strain, for example Rev.1 strain, is possible provided that it produces the same level of sensitisation). The erythematous reactions are read and measured at 24 hours and the titre is calculated by comparison with a reference brucellin<sup>9</sup>. This method is only valid for comparing brucellin preparations made according to the same protocol as the sensitising allergen. Initial standardisation of a batch of allergen and the sensitisation and titration in ruminants is described (Alton *et al.*, 1988).

## 2.2.3. Requirements for authorisation/registration/licensing

### 2.2.3.1. Manufacturing process

For registration of the brucellin, all relevant details concerning manufacture of the product and quality control testing (see above) should be submitted to the regulatory authorities in accordance with their requirements. This information should be provided from three consecutive product batches with a volume of not less than 1/3 of the typical industrial batch volume.

### 2.2.3.2. Safety requirements

i) Target and non-target animal safety

No side-effects of the brucellin have ever been reported in animals.

ii) Reversion-to-virulence

Not applicable.

iii) Precautions

Brucellin is not toxic. Nevertheless it may provoke severe hypersensitivity reactions in sensitised individuals who are accidentally exposed to it. Care should be taken to avoid

---

8 An EU reference brucellin (2000 Units/ml) is obtainable from the OIE Reference Laboratory for Brucellosis in France.

9 The statistical procedure can be obtained from the OIE Reference Laboratory for Brucellosis in France.

accidental injection or mucosal contamination. Used containers and injection equipment should be carefully decontaminated or disposed of by incineration in a suitable disposable container.

Brucellin should be identified as potentially harmless for practitioners. Manufacturers should provide adequate warnings that medical advice should be sought in the case of self-injection or mucosal exposure with warnings included on the product label/leaflet so that the practitioner is aware of any danger.

#### 2.2.3.3. Efficacy requirements

Potency must be determined on the final product. The procedure is as described above.

#### 2.2.3.4. Duration of sensitivity

Duration of sensitivity is uncertain. Individual animals vary considerably in the degree of hypersensitivity manifested to brucellin. Animals in the very early stages of infection, or with long-standing infection, may not manifest hypersensitivity to intradermal injection.

#### 2.2.3.5. Stability

The freeze-dried preparation retains full potency for several years. The liquid commercial preparation should retain potency for the recommended shelf-life.

## REFERENCES

- ALTON G.G., JONES L.M., ANGUS R.D. & VERGER J.M. (1988). Techniques for the Brucellosis Laboratory. Institut National de la Recherche Agronomique, Paris, France.
- ANGUS R.D. & BARTON C.E. (1984). The production and evaluation of a buffered plate antigen for use in a presumptive test for brucellosis. *Dev. Biol. Stand.*, **56**, 349–356.
- ASHFORD D.A., DI PIETRA J., LINGAPPA J., WOODS C., NOLL H., NEVILLE B., WEYANT R., BRAGG S.L., SPIEGEL R.A., TAPPERRRO J. & PERKINS B.A. (2004). Adverse events in humans associated with accidental exposure to the livestock brucellosis vaccine RB51. *Vaccine*, **22**, 3435–3439.
- BLASCO J.M. (1997). A review on the use of *B. melitensis* Rev 1 vaccine in adult sheep and goats. *Prev. Vet. Med.*, **31**, 275–283.
- BLASCO J.M., GARIN-BASTUJI B., MARÍN C.M., GERBIER G., FANLO J., JIMÉNEZ DE BAGÜES M.P. & CAU C. (1994a). Efficacy of different rose bengal and complement fixation antigens for the diagnosis of *Brucella melitensis* infection in sheep and goats. *Vet. Rec.*, **134**, 415–420.
- BLASCO J.M., MARÍN C.M., JIMÉNEZ DE BAGÜES M.P., BARBERAN M., HERNANDEZ A., MOLINA L., VELASCO J., DIAZ R. & MORIYÓN I. (1994b). Evaluation of allergic and serological tests for diagnosing *Brucella melitensis* infection of sheep. *J. Clin. Microbiol.*, **32**, 1835–1840.
- BRICKER B.J. (2002). PCR as a diagnostic tool for brucellosis. *Vet. Microbiol.*, **90**, 435–446.
- BRICKER B.J., EWALT D.R. & HALLING S.M. (2003). *Brucella* 'HOOF-Prints' strain typing by multi-locus analysis of variable number tandem repeats (VNTRs). *BMC Microbiol.*, **3**, 15.
- DE MASSIS F., GIOVANNINI A., DI EMIDIO B., RONCHI G.F., TITTARELLI M., DI VENTURA M., NANNINI D. & CAPORALE V. (2005). Use of the complement fixation and brucellin skin tests to identify cattle vaccinated with *Brucella abortus* strain RB51. *Vet. Ital.*, **41** (4), 291–299.
- DE MIGUEL M.J., MARÍN C.M., MUÑOZ P.M., DIESTE L., GRILLÓ M.J. & BLASCO J.M. (2011). Development of a selective culture medium for primary isolation of the main *Brucella* species. *J. Clin. Microbiol.*, **49**, 1458–1463.
- DIESTE-PÉREZ L., BLASCO J.M., DE MIGUEL M.J., MARÍN C.M., BARBERÁN M., CONDE-ÁLVAREZ R., MORIYÓN I., MUÑOZ P.M. (2014). Performance of skin tests with allergens from *B. melitensis* B115 and rough *B. abortus* mutants for diagnosing swine brucellosis. *Vet. Microbiol.*, **168**, 161–168.
- EUROPEAN FOOD SAFETY AGENCY (EFSA) (2009). Scientific Opinion of the Panel on Animal Health and Welfare (AHAW) on a request from the Commission on porcine brucellosis (*Brucella suis*). *EFSA J.*, **1144**, 1–112.

- EUROPEAN UNION (2008). Commission Decision of 10 December 2008 amending Annex C to Council Directive 64/432/EEC and Decision 2004/226/EC as regards diagnostic tests for bovine brucellosis (notified under document number C [2008] 7642) (Text with EEA relevance) (2008/984/EC). *Official Journal of the European Union*, L 352/38-45.
- FOSTER G., OSTERMAN B.S., GODFROID J., JACQUES I. & CLOECKAERT A. (2007). *Brucella ceti* sp. nov. and *Brucella pinnipedialis* sp. nov. for *Brucella* strains with cetaceans and seals as their preferred hosts. *Int. J. Syst. Evol. Microbiol.*, **57**, 2688–2693.
- GRILLÓ M.J., BOSSERAY N. & BLASCO J.M. (2000). *In vitro* markers and biological activity in mice of seed lot strains and commercial *Brucella melitensis* Rev.1 and *Brucella abortus* B19 vaccines. *Biologicals*, **28**, 119–127.
- INTERNATIONAL AIR TRANSPORT ASSOCIATION (IATA) (2013). Dangerous Goods Regulations, 54<sup>th</sup> Edition. International Air Transport Association, 800 Place Victoria, P.O. Box 113; Montreal, Quebec H4Z 1M1, Canada, 815 pp <http://www.iata.org/publications/dgr/Pages/manuals.aspx>.
- JIMÉNEZ DE BAGÜES M.P., MARÍN C. & BLASCO J.M. (1991). Effect of antibiotic therapy and strain 19 vaccination on the spread of *Brucella melitensis* within an infected dairy herd. *Prev. Vet. Med.*, **11**, 17–24.
- JOINT FOOD AND AGRICULTURE ORGANIZATION OF THE UNITED NATIONS (FAO)/WORLD HEALTH ORGANIZATION (WHO) EXPERT COMMITTEE ON BRUCELLOSIS (1986). Technical Report Series 740, Sixth Report. WHO, Geneva, Switzerland.
- JONES L.M., BERMAN D.T., MORENO E., DEYOE B.L., GILSDORF M.J., HUBER J.D. & NICOLETTI P.L. (1980). Evaluation of a radial immunodiffusion test with polysaccharide B antigen for diagnosis of bovine brucellosis. *J. Clin. Microbiol.*, **12**, 753–760.
- KANG S.I., HER M., KIM J.W., KIM J.Y., KO K.Y., HA Y.M. & JUNG S.C. (2011). Advanced multiplex PCR assay for differentiation of *Brucella* species. *Appl. Environ. Microbiol.*, **77**, 6726–6728.
- LE FLÈCHE P., JACQUES I., GRAYON M., AL DAHOUK S., BOUCHON P., DENOEUDE F., NÖCKLER K., NEUBAUER H., GUILLOTEAU L.A. & VERGNAUD G. (2006). Evaluation and selection of tandem repeat loci for a *Brucella* MLVA typing assay. *BMC Microbiol.*, **6**, 9.
- LÓPEZ-GOÑI I., GARCÍA-YOLDI D., MARÍN C.M., DE MIGUEL M.J., BARQUERO-CALVO E., GUZMÁN-VERRI C., ALBERT D., GARIN-BASTUJI B. (2011). New Bruce-ladder multiplex PCR assay for the biovar typing of *Brucella suis* and the discrimination of *Brucella suis* and *Brucella canis*. *Vet Microbiol.*, **154**, 152–155.
- MCGIVEN J., TAYLOR A., DUNCOMBE L., SAYERS R., ALBERT D., BANAI M., BLASCO J.M., ELENA S., FRETIN D., GARIN-BASTUJI B., MELZER F., MUÑOZ P.M., NIELSEN K., NICOLA A., SCACCHIA M., TITTARELLI M., TRAVASSOS DIAS I., WALRAVENS K. & STACK J. (2011). The first International Standard anti-*Brucella melitensis* Serum. *Rev. sci. tech. Off. int. Epiz.*, **30**, 809–819.
- MACMILLAN A.P. & COCKREM D.S. (1985). Reduction of non-specific reactions to the *Brucella abortus* serum agglutination test by the addition of EDTA. *Res. Vet. Sci.*, **38**, 288–291.
- MORENO E., CLOECKAERT A. & MORIYON I. (2002). *Brucella* evolution and taxonomy. *Vet. Microbiol.*, **90**, 209–227.
- MORGAN W.J.B., MACKINNON D.J., LAWSON J.R. & CULLEN G.A. (1969). The rose bengal plate agglutination test in the diagnosis of brucellosis. *Vet. Rec.*, **85**, 636–641.
- MORIYÓN I., GRILLÓ M.J., MONREAL D., GONZALEZ D., MARÍN C.M., LÓPEZ-GOÑI I., MAINAR-JAIME R.C., MORENO E. & BLASCO J.M. (2004). Rough vaccines in animal brucellosis: structural and genetic basis and present status. *Vet. Res.*, **35**, 1–38.
- MUÑOZ P.M., BLASCO J.M., ENGEL B., DE MIGUEL M.J., MARÍN C.M., DIESTE L., MAINAR-JAIME R.C. (2012). Assessment of performance of selected serological tests for diagnosing brucellosis in pigs. *Vet. Immunol. Immunopathol.*, **146**, 150–158
- MUÑOZ P.M., MARIN C.M., MONREAL D., GONZALES D., GARIN-BASTUJI B., DIAZ R., MAINAR-JAIME R.C., MORIYON I. & BLASCO J. (2005). Efficacy of several serological tests and antigens for the diagnosis of bovine brucellosis in the presence of false positive serological results due to *Yersinia enterocolitica* O:9. *Clin. Diagn. Lab. Immunol.*, **12**, 141–151.

- NICOLETTI P, JONES L M & BERMAN D.T. (1978). Comparison of the subcutaneous and conjunctival route of vaccination with *Brucella abortus* strain 19 vaccine in adult cattle. *J. Am. Vet. Med. Assoc.*, **173**, 1450–1456.
- NIELSEN K., GALL D., JOLLEY M., LEISHMAN G., BALSEVICIUS S., SMITH P., NICOLETTI P. & THOMAS F. (1996). A homogenous fluorescence polarisation assay for detection of antibody to *Brucella abortus*. *J. Immunol. Methods*, **195**, 161–168.
- NIELSEN K., KELLY L., GALL D., NICOLETTI P. & KELLY W. (1995). Improved competitive enzyme immunoassay for the diagnosis of bovine brucellosis. *Vet. Immunol. Immunopathol.*, **46**, 285–291.
- OCAMPO–SOSA A.A., AGÜERO–BALBÍN J. & GARCÍA–LOBO J.M. (2005). Development of a new PCR assay to identify *Brucella abortus* biovars 5, 6 and 9 and the new subgroup 3b of biovar 3. *Vet. Microbiol.*, **110**, 41–51.
- OLSEN S.C., GARIN-BASTUJI B., BLASCO J.M., NICOLA A.M., SAMARTINO L. (2012). Swine Brucellosis. *In: Diseases of Swine*, Zimmerman J.J., Karriker A.L., Ramirez A., Schwartz K.J. & Stevenson G.W. eds, 10<sup>th</sup> Edition, John Wiley & Sons, USA, 697–708.
- POUILLOT R., GARIN–BASTUJI B., GERBIER G., COCHE Y., CAU C., DUFOUR B. & MOUTOU F. (1997). The brucellin skin test as a tool to differentiate false positive serological reactions in bovine brucellosis. *Vet. Res.*, **28**, 365–374.
- POUILLOT R., GRILLÓ M.J., ALABART J.L., GARIN–BASTUJI B. & BLASCO J.M. (2004). Statistical procedures for calculating the residual virulence of *Brucella abortus* strain 19 (S19) and *Brucella melitensis* strain Rev.1 vaccines in mice: theoretical basis and practical applications. *Rev. sci. tech. Off. int. Epiz.*, **22**, 1051–1063.
- PRAUD A., GIMENEZ O., ZANELLA G., DUFOUR B., POZZI N., ANTRAS V., MEYER L. & GARIN-BASTUJI B. (2012). Estimation of sensitivity and specificity of five serological tests for the diagnosis of porcine brucellosis. *Prev. Vet. Med.*, **104**, 94–100.
- SCHOLZ H.C., HUBALEK Z., SEDLACEK I., VERGNAUD G., TOMASO H., AL DAHOUK S., MELZER F., KÄMPFER P., NEUBAUER H., CLOECKAERT A., MAQUART M., ZYGMUNT M.S., WHATMORE A.M., FALSEN E., BAHN P., GÖLLNER C., PFEFFER M., HUBER B., BUSSE H.J. & NÖCKLER K. (2008). *Brucella microti* sp. nov., isolated from the common vole *Microtus arvalis*. *Int. J. Syst. Evol. Microbiol.*, **58**, 375–382.
- SCHOLZ H.C., NÖCKLER K., GÖLLNER C., BAHN P., VERGNAUD G., TOMASO H., AL DAHOUK S., KÄMPFER P., CLOECKAERT A., MAQUART M., ZYGMUNT M.S., WHATMORE A.M., PFEFFER M., HUBER B., BUSSE H.J. & DE B.K. (2010). *Brucella inopinata* sp. nov., isolated from a breast implant infection. *Int. J. Syst. Evol. Microbiol.*, **60**, 801–808.
- STACK J.A., HARRISON M. & PERRETT L.L. (2002). Evaluation of a selective medium for *Brucella* isolation using natamycin. *J. Appl. Microbiol.*, **92**, 724–728.
- STACK J.A., PERRETT L.L., BREW S.D. & MACMILLAN A.P. (1999). C–ELISA for bovine brucellosis suitable for testing poor quality samples. *Vet. Rec.*, **145**, 735–736.
- STOFFREGEN W.C., OLSEN S.C., BRICKER B.J. (2006). Parenteral vaccination of domestic pigs with *Brucella abortus* strain RB51. *Am. J. Vet. Res.*, **67**, 1802–1808.
- UNITED STATES DEPARTMENT OF AGRICULTURE (USDA), ANIMAL AND PLANT HEALTH INSPECTION SERVICES (APHIS) (2003). Availability of an Environmental Assessment for Licensing of *Brucella abortus* Vaccine, Strain RB–51, Live Culture. *Federal Register*, 18 Feb 2003, **68**, 7761.
- VERGER J.M. (1985). *B. melitensis* infection in cattle. *In: Brucella melitensis*, Plommet M. & Verger J.M., eds. Martinus Nijhoff Publ., Dordrecht, Netherlands, 197–203.
- WHATMORE A.M. (2009). Current understanding of the genetic diversity of *Brucella*, an expanding genus of zoonotic pathogens. *Infect. Genet. Evol.*, **9**, 1168–1184.
- WHATMORE A.M. & GOPAUL K.K. (2011). Recent advances in molecular approaches to *Brucella* diagnostics and epidemiology. *In: Brucella: Molecular Microbiology and Genomics*, López-Goñi I. & O’Callaghan D., eds, Caister Academic Press, Norfolk, UK, 57–88.
- WHATMORE A.M., DAVISON N., CLOECKAERT A., AL DAHOUK S., ZYGMUNT M.S., BREW S.D., PERRETT L.L., KOYLASS M.S., VERGNAUD G., QUANCE C., SCHOLZ H.C., DICK E.J. Jr, HUBBARD G. & SCHLABRITZ–LOUTSEVITCH N.E. (2014). *Brucella papionis* sp. nov. isolated from baboons (*Papio* spp.). *Int. J. Syst. Evol. Microbiol.*, **64**, 4120–4128.

WORLD HEALTH ORGANIZATION (1953). WHO Technical Report Series No. 68. Sixth report of the WHO Expert Committee on Biological Standardization. WHO, Geneva, Switzerland.

WORLD HEALTH ORGANIZATION (2004). WHO Laboratory Biosafety Manual, Third Edition. WHO, Geneva, Switzerland.

WORLD HEALTH ORGANIZATION (2012). Guidance on regulations for the Transport of Infectious Substances 2013-2014, WHO/HSE/GCR/2012.12, WHO, Geneva, Switzerland. <http://apps.who.int/iris/handle/10665/78075>

\*

\* \*

**NB:** There are OIE Reference Laboratories for infection with *Brucella abortus*, *B. melitensis* and *B. suis* (please consult the OIE Web site:

<https://www.oie.int/en/what-we-offer/expertise-network/reference-laboratories/#ui-id-3> ).

Please contact the OIE Reference Laboratories for any further information on diagnostic tests, reagents and vaccines for these agents

**NB:** BOVINE BRUCELLOSIS FIRST ADOPTED IN 1990; BRUCELLOSIS IN SHEEP, GOATS AND SWINE FIRST ADOPTED IN 1991. CHAPTER FIRST ADOPTED WITH CURRENT TITLE IN 2016.